



Surgical Technique

Wise-Lock "T" Wrist Volar Plate

www.auxein.com

about us

Auxein Medical is an integrated, research based, orthopaedic Implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in india and serves customers in over 75 countries worldwide.

Our Achievements



INTRODUCTION

The AUXEIN MEDICAL'S Wise-Lock "T" Wrist Volar Plate is used for fixation of fractures and osteotomies involving the distal radius. This plate is single use implantable device for long term duration (intended for continuous used for more than 30 days) contacting radius bone and its surrounding tissues.





INDICATIONS:

The Distal Radius Implants are intended for:

- Displaced extra-articular and intra-articular distal radius fractures.
- Dorsally displaced fractures
- Extra-articular fractures with metaphyseal defect (AO classification 23-A3)
- Open joint reconstruction (AO classification 23-C1, C2, C3)
- Corrective osteotomies

CONTRAINDICATIONS:

Contraindications include but are not limited to:

- Infection, local to operative site
- Signs of local inflammation
- Morbid Obesity
- Severely comminuted fractures in which bone fractures are too small or numerous to fix or maintain a reduced position
- Metal sensitivity or intolerance
- Alcohol or drug addict
- Symptomatic Arthritis
- Malignant Tumors, pimary or metastatic
- Poor bone quality
- Foreign body sensibility
- Vascular compromise





Patient Positioning: Patient is placed in supine position on an operating table. Forearm is rested on hand table in supinated position so that palm faces upwards

Skin Incision: An incision of 8cm is made over the course of flexor carpi radialis (FCR) tendon. For better access and visualization, a zigzag incision is made across the wrist flexion.

Plate Positioning: The plate is positioned in such a way that the plate conforms to the fracture line and sits on the volar surface of the radius.



Initial Plate fixation:

Firstly, the 3.5mm cortical screw is inserted into the proximal oblong hole which will allow adjusting the plate for initial fixation. Follow the below steps for 3.5mm Cortical screw insertion:

- **Predrilling Screw Hole:** The 2.5/3.5mm drill guide (1472-046) is placed along the plate hole and the 2.5mm drill bit (2100-2.5-112) is inserted through it for predrilling the hole.
- **Precaution:** Selection of correct drill bit should be carefully chosen.for the gliding hole.

The hole can be drilled at neutral or compression position according to requirement of the reduction as described below:

- 1. **Predrilling in neutral position:** The upper part of the inner sleeve of drill guide is depressed and is positioned in plate hole remote from the fracture as indicated by 1 is adjoining figure.
- 2. **Predrilling in compression position:** The upper part of the inner sleeve of drill guide is not depressed and is positioned in plate hole remote from the fracture as indicated by 2 in adjoining figure.







Determining screw length: The desired screw length is measured with the help of depth gauge (3443.37)

Precaution: Care should be taken to choose appropriate length using depth gauge. Selection of excessive or too short screw length can create complications in bone healing as excessive length can lead to soft tissue irritation or damage.

Screw Insertion: The 2.5mm tip Hexagonal screwdriver (3406.03) is used to pick and insert the 3.5mm cortical screw of selected length.

Precaution: Care should be taken while screw insertion. Avoid excessive tightening as it can lead to screw breakage and recess deformation. And also improper tightening can lead to loosening of implant and irritation or soft tissue damage.

40 50 60

30



Distal Plate Fixation: Final fracture reduction is performed before distal plate fixation. The head of the plate accepts both the 2.5mm volar screws as well as 2.0mm peg screws.

The 2.0mm peg screws provide better hold on the dorsal comminuted fragments and a strong peg to plate interface.

The 2.5mm volar screw is available in both full thread and partial thread. The volar screw with partial thread helps in securing the fragments in coronal plane.

The following steps are used for screw insertion:

• The k-wire is passed through the k- wire hole in plate head to secure the bone fragment to the plate.

The Drill sleeve Ø2.0 (1312-14-2.0) is attached to the proximal plate threaded head hole and drill bit Ø2.0 (2100-2.0-100) is used to drill through the bone. The marking on the drill bit coinciding with the one on the drill sleeve indicates the required screw length for insertion.



.



- The Star head Screwdriver shaft is assembled to 0.8Nm Torque limited screwdriver (1312-21) and then attached into the quick coupling handle (1312-22). This assembly is then used to insert the 2.5mm volar screws or the 2.0mm peg screws as per requirement into the predrilled holes.
- After filling the proximal head holes, all the distal head holes are filled with the volar/peg screws.

Precaution: Care should be taken while screw insertion. Avoid excessive tightening or loosening. Choose appropriate size torque limiter.



Proximal Plate Fixation: The locking shaft holes of the plate accept 3.5mm Wise lock screws.

• The 3.5mm threaded drill guide (3441.18) is aligned with the trajectory of the threaded plate hole and the guide is approximately given a quarter counter clockwise turn until the thread engages. Finally the guide is advanced clockwise until firmly seated. The 2.8mm drill bit (2103-2.8-165) is then inserted into the guide to drill the bone for screw insertion.





• The depth gauge (3443.37) is then used to measure the required screw length for insertion.





• The 1.5Nm Torque Screwdriver (TQ-2.5) is attached to the Torque Screw Driver Handle (1472-064). The Star head Screwdriver shaft is further coupled into the 1.5Nm Torque screwdriver. Now this screwdriver assembly is used to insert the 3.5mm Wise lock into the predrilled hole.





2100-1.8-110	Drill Bit Ø1.8		
2100-2.4-100	Drill Bit Ø2.0		
2100-2.0-100	Drill Bit Ø2.4		.82.4
2100-2.7-100	Drill Bit Ø2.7		ø27
2105-01	Guide Pin Ø1.5		
2105-02	Threaded Guide Pir	nØ1.5	
3767-2.4	Double Drill Guide	Ø2.4/1.8mm	















1312-20-2.7	Screw holding sleeve for 2.7mm cortical sc	2.7 «шиц.)
1312-20-2.4	Screw holding sleeve for 2.4mm cortical sc	rew 2.1 mmm
1312-30-2.4	Screw holding sleeve for 2.4mm wise lock s	screw 2.1 amment
2104-02	Bone Tap Quick Coupling for Cortical Screv	vs, Ø3.5mm
2100-2.5-112	Drill Bit Ouick Couplina.Ø2.5mm x 112mm	OCE p 2 0mm x 112mm 200-20-112 2000
2103-2.8-165	Drill Bit Quick Coupling with Stopper, Ø 2.8	Bmm x Length 165mm
1472-046	Drill Guide, Ø2.5/3.5mm	A CE 147 DAS 27000 ALLER SA CE INTO ALLE
3441.18	Threaded Drill Guide, Ø3.5mm for drill bit 2.8mm	
TQ-2.5	Torque Screw Driver - 2.5mm Tip, 1.5Nm	@CE 10-25-2700



Wise-Lock "T" Wrist Volar Plate

	Left Direction	Right Direction
Holes	Titanium	Titanium
3	6137/06	6136/06
4	6137/01	6136/01

2.5mm Volar Screws, Self Tapping, Partial Thread



2.5mm Volar Screws, Self Tapping, Full Thread

Length (mm)	Titanium	
10	6152/11	1
12	6152/01	
14	6152/02	
16	6152/03	
18	6152/04	
20	6152/05	
22	6152/06	
24	6152/07	
26	6152/08	a a a a a a a a a a a a a a a a a a a
28	6152/09	
30	6152/10	

Titanium
6151/10
6151/01
6151/02
6151/03
6151/04
6151/05
6151/06
6151/07
6151/08
6151/09
6151/11

2.0mm Peg Screws

Length (mm)	Titanium
10	6153/09
12	6153/10
14	6153/11
16	6153/01
18	6153/02
20	6153/03
22	6153/04
24	6153/05
26	6153/06
28	6153/07
30	6153/08

3.5mm Cortical Screws

Length (mm)	Titanium
10	TI-104.010
12	TI-104.012
14	TI-104.014
16	TI-104.016
18	TI-104.018
20	TI-104.020
22	TI-104.022
24	TI-104.024
26	TI-104.026
28	TI-104.028
30	TI-104.030
32	TI-104.032
34	TI-104.034
36	TI-104.036
38	TI-104.038
40	TI-104.040
42	TI-104.042
44	TI-104.044
46	TI-104.046
48	TI-104.048
50	TI-104.050



3.5mm Wise-Lock Screws, Self Tapping

Length (mm)	Titanium
10	TI-117.010
12	TI-117.012
14	TI-117.014
16	TI-117.016
18	TI-117.018
20	TI-117.020
22	TI-117.022
24	TI-117.024
26	TI-117.026
28	TI-117.028
30	TI-117.030
32	TI-117.032
34	TI-117.034
36	TI-117.036
38	TI-117.038
40	TI-117.040
42	TI-117.042
44	TI-117.044
46	TI-117.046
48	TI-117.048
50	TI-117.050
52	TI-117.052
54	TI-117.054
56	TI-117.056
58	TI-117.058
60	TI-117.060

and the second division of the second divisio	
and the second se	
and the second s	
and the second se	
the second se	
and the second se	
the second se	
and the second se	
The second s	
No. of Concession, Name	
and the second se	
of the owner of the	
the second se	
the second se	
and the second se	
and the second second	
and the second se	
and the second se	
and the second se	
and the second se	
and the second se	
and the second se	
and the second se	
and the second se	
the second se	
and the second se	
the second se	
and the second se	
and the second se	
the second se	
and the second se	
and the second se	
and the second se	
the second se	
and the second s	
the second se	
and the second s	
Concession of the local division of the loca	
and the second s	
Concession of the local	
and the second s	
and the second s	
The second se	
The second se	
and the second second	
1000	



USA

Auxein Inc. 1500 Nw 89th Court, Suite 107-108 Doral, Florida 33172 Tel: +1 305 395 6062 E Fax: +1 305 395 6262 Email: USoffice@auxein.com

www.auxein.com

MEXICO

Auxein México, S.A. de C.V. Tepic 139 int 801, Colonia Roma Sur, Alcaldía Cuauhtémoc, CDMX, México, C.P. 06760 Tel: +521 55 7261 0318 Email: info@auxein.mx

INDIA

Auxein Medical Pvt. Ltd. Plot No. 168-169-170, Phase-4, Kundli Industrial Area, HSIIDC, Sector-57, Sonepat - 131028, Haryana Tel: +91 99106 43638 | Fax: +91 86077 70197 Email: info@auxein.com