

auxilock® 


SMALL JOINT Arthroscopy System

Hand & Wrist




Foot & Ankle

SMALL JOINT
Arthroscopy System

AUXILOCK[®]
AUXSUTURE
ANCHOR



Ordering Information

Code	Product Description	Suture/Tape Colour	Needle
6-040-01	AUXILOCK [®] 1.2mm AuxSuture Anchor With One #0 BioBraid: White/Blue		
6-040-02	AUXILOCK [®] 1.2mm AuxSuture Anchor With One #0 BioBraid: White/Blue, With Needles		

Recommended Instrumentation







Code	Product Description	Compatible Implants
7-139-02	8 Point Slotted Drill Guide, Length 50mm	6-040-01, 6-040-02
7-139-52	Ø1.2mm X 110mm Drill Bit	6-040-01, 6-040-02

STERILE 

AUXILOCK[®] TITANIUM SCREW-IN SUTURE ANCHOR



Ordering Information


Code	Product Description	Suture/Tape Colour	Needle
6-005-16	AUXILOCK [®] 2.0mm Titanium Screw-In Suture Anchor With One #0 BioBraid: White/Blue		
6-005-11	AUXILOCK [®] 2.0mm Titanium Screw-In Suture Anchor With One #0 BioBraid: White/Blue, With Needles		
6-005-15	AUXILOCK [®] 2.5mm Titanium Screw-In Suture Anchor With One #1 BioBraid: White/Blue		
6-005-12	AUXILOCK [®] 2.5mm Titanium Screw-In Suture Anchor With One #1 BioBraid: White/Blue, With Needles: MO-6		

Recommended Instrumentation

Code	Product Description	Compatible Implants
7-139-01	Ø1.2mm X 100mm Drill Bit	6-005-11, 6-005-16
7-139-02	8 Point Slotted Drill Guide, Length 50mm	6-005-11, 6-005-12, 6-005-16, 6-005-15
7-139-66	Ø1.8mm X 100mm Drill Bit	6-005-15, 6-005-12

STERILE EO

AUXILOCK[®] PEEK CF INTERFERENCE SCREW



Ordering Information


Code	Product Description
6-009-27	AUXILOCK [®] 3mm x 8mm PEEK CF Interference Screw (With Driver)
6-009-21	AUXILOCK [®] 4mm x 10mm PEEK CF Interference Screw
6-009-22	AUXILOCK [®] 4mm x 15mm PEEK CF Interference Screw
6-009-28	AUXILOCK [®] 5mm x 10mm PEEK CF Interference Screw
6-009-23	AUXILOCK [®] 5mm x 15mm PEEK CF Interference Screw
6-009-24	AUXILOCK [®] 5mm x 20mm PEEK CF Interference Screw
6-009-29	AUXILOCK [®] 6mm x 15mm PEEK CF Interference Screw
6-009-25	AUXILOCK [®] 6mm x 20mm PEEK CF Interference Screw
6-009-26	AUXILOCK [®] 6mm x 25mm PEEK CF Interference Screw

Related Instrumentation

Code	Product Description
7-139-54	Ø1.2mm X 150mm Threaded Guide wire
7-139-55	Ø1.2mm X 180mm Blunt tip Guide wire
7-139-68	Ø3.0mm X 100mm Cannulated Reamer
7-139-67	Ø3.2mm X 100mm Cannulated Reamer
7-139-56	Ø3.5mm X 100mm Cannulated Reamer
7-139-57	Ø4.0mm X 100mm Cannulated Reamer
7-139-58	Ø4.5mm X 100mm Cannulated Reamer
7-139-59	Ø5.0mm X 100mm Cannulated Reamer
7-139-60	Ø5.5mm X 100mm Cannulated Reamer
7-139-61	Ø6.0mm X 100mm Cannulated Reamer
7-139-62	Ø6.5mm X 100mm Cannulated Reamer
7-139-63	Graft Sizer, 4.0mm to 7.0mm (0.5mm Steps)
7-139-65	Screwdriver For Ø4mm, Ø5.0mm & Ø6.0 mm Polymer Interference Screw

STERILE EO

AUXILOCK[®] PEEK OPTIMA INTERFERENCE SCREW



Ordering Information

Code	Product Description
6-014-27	AUXILOCK [®] 3mm x 8mm PEEK OPTIMA Interference Screw (With Driver)
6-014-21	AUXILOCK [®] 4mm x 10mm PEEK OPTIMA Interference Screw
6-014-22	AUXILOCK [®] 4mm x 15mm PEEK OPTIMA Interference Screw
6-014-28	AUXILOCK [®] 5mm x 10mm PEEK OPTIMA Interference Screw
6-014-23	AUXILOCK [®] 5mm x 15mm PEEK OPTIMA Interference Screw
6-014-24	AUXILOCK [®] 5mm x 20mm PEEK OPTIMA Interference Screw
6-014-29	AUXILOCK [®] 6mm x 15mm PEEK OPTIMA Interference Screw
6-014-25	AUXILOCK [®] 6mm x 20mm PEEK OPTIMA Interference Screw
6-014-26	AUXILOCK [®] 6mm x 25mm PEEK OPTIMA Interference Screw

Related Instrumentation


Code	Product Description
7-139-54	Ø1.2mm X 150mm Threaded Guide wire
7-139-55	Ø1.2mm X 180mm Blunt tip Guide wire
7-139-68	Ø3.0mm X 100mm Cannulated Reamer
7-139-67	Ø3.2mm X 100mm Cannulated Reamer
7-139-56	Ø3.5mm X 100mm Cannulated Reamer
7-139-57	Ø4.0mm X 100mm Cannulated Reamer
7-139-58	Ø4.5mm X 100mm Cannulated Reamer
7-139-59	Ø5.0mm X 100mm Cannulated Reamer
7-139-60	Ø5.5mm X 100mm Cannulated Reamer
7-139-61	Ø6.0mm X 100mm Cannulated Reamer
7-139-62	Ø6.5mm X 100mm Cannulated Reamer
7-139-63	Graft Sizer, 4.0mm to 7.0mm (0.5mm Steps)
7-139-65	Screwdriver For Ø4mm, Ø5.0mm & Ø6.0 mm Polymer Interference Screw

STERILE EO

AUXILOCK[®]

ROTADOR PEEK OPTIMA

SCREW-IN ANCHOR



Ordering Information

Code	Product Description
6-020-01	AUXILOCK [®] Rotador 4.75mm x 15mm PEEK OPTIMA Screw-In Anchor


Recommended Instrumentation

Code	Product Description	Compatible Implants
7-139-69	Tap for 4.75 mm Rotador Screw-in Anchor	6-020-01

STERILE

EO

AUXILOCK[®] SYNDESMOSIS ULTIMATE KIT



AUXILOCK[®] Syndesmosis Ultimate Kit can be used for the fixation of Ankle Syndesmosis joint. The versatile design of the kit is ideal for use with or without a fracture plate. The kit comprises of adjustable loop with round shaped button (Dia. 5.5mm) and oblong shaped button (L 12mm X W 3.2mm X H 1.5mm). The round button has a tapered slot design which can easily fix on any fracture plate. The kit is passed through 3.5mm drilling tunnel from lateral to medial provides the fixation of oblong shaped button on the medial side and round button on the lateral side of the syndesmosis joint. The syndesmosis kit provides a double locking mechanism which eliminates the need for knot tying.

Ordering Information


Code	Product Description
6-041-01	AUXILOCK [®] Syndesmosis Ultimate Kit

Related Instrumentation (Include in Kit)

Code	Product Description
7-147-01	Cannulated Drill Bit with Quick Coupling End, Ø3.5mm x Length 175mm
7-147-02	Ø1.5mm K-Wire
7-147-03	Ø2.4mm Passing Pin
7-147-04	Suture Pull Handle

STERILE EO

AUXILOCK[®]
ADVANCE SYNDESMOSIS
ULTIMATE KIT



The AUXILOCK[®] Advance Syndesmosis implant system features a unique delivery mechanism that allows the user to insert a Syndesmosis implant without pulling a needle through the medial skin. Advance Syndesmosis Handle system gives the user control to flip the medial oblong button below the skin for less soft-tissue interposition.

Ergonomic Advance Syndesmosis implant tensioning handles have also been added to the system. The handle is designed tension the syndesmosis implant without hurting the surgeon fingers. The Avance Syndesmosis implant system also features a new trocar style drill bit that reduces skiving off the fibula when drilling in the transmalleolar plane 30° anterior to the coronal plane.

Ordering Information

Code	Product Description
6-051-01	AUXILOCK [®] Advance Syndesmosis Ultimate Kit

Related Instrumentation (Include in Kit)

Code	Product Description
7-147-01	Cannulated Drill Bit with Quick Coupling End, Ø3.5mm x Length 175mm
7-147-02	Ø1.5mm K-Wire

STERILE EO

Instruments

7-139-01
Ø1.2mm X 100mm Drill Bit
(For 2.0mm Ti Anchor)




7-139-02
8 point Slotted Drill Guide, Length 50mm



7-139-52
Ø1.2mm X 110mm Drill Bit
(For 1.2mm AuxSuture Anchor)



7-139-66
Ø1.8mm X 100mm Drill Bit
(For 2.5mm Ti Anchor)



7-139-54
Ø1.2mm X 150mm Threaded Guide wire



7-139-55
Ø1.2mm X 180mm Blunt tip Guide wire




Cannulated Reamer

Code	Diameter
7-139-68	Ø3.0mm X 100mm Cannulated Reamer
7-139-67	Ø3.2mm X 100mm Cannulated Reamer
7-139-56	Ø3.5mm X 100mm Cannulated Reamer
7-139-57	Ø4.0mm X 100mm Cannulated Reamer
7-139-58	Ø4.5mm X 100mm Cannulated Reamer
7-139-59	Ø5.0mm X 100mm Cannulated Reamer
7-139-60	Ø5.5mm X 100mm Cannulated Reamer
7-139-61	Ø6.0mm X 100mm Cannulated Reamer
7-139-62	Ø6.5mm X 100mm Cannulated Reamer



7-139-63 Graft Sizer, 4.0mm to 7.0mm (0.5mm Steps)



7-139-65 Screwdriver For Ø4.0mm, Ø5.0mm & Ø6.0mm Polymer Interference Screw



7-139-69 Tap for 4.75 mm Rotador Screw-in Anchor



7-139-70
3.5mm Drill Guide, Length 20mm



7-139-71
Drill Guide Sleeve for Ø1.5mm K-Wire



OPTIONAL INSTRUMENTS

7-139-05
Exploring Probe, Graduated, Curved



Mini Rasp

Code	Shape
7-139-06	Mini Rasp, Concave Jaw
7-139-07	Mini Rasp, Convex Jaw



7-139-09
Mini Knife With V - Knotch



7-139-11
Mini Bone Curette



7-139-21
Scalpel Knife



7-139-22
Banana Knife



7-139-23
Hook Knife, Reverse Cutting



7-139-27
Mini Micro Fracture Awl, 45°



7-139-28
Mini Micro Fracture Awl, 60°



7-139-31
Mini Punch



7-139-32
Mini Scissors, Upwards Cutting



7-139-33
Mini Grasping Forceps



7-139-43
Blunt Nose Mini Basket Punch, Straight Jaws, Straight Shaft



7-139-48
Shark Mini Grasping Forceps, Straight Jaws, Straight Shaft



ORDER INFORMATION

7-139 Small Joint Arthroscopy Instrument Set

Code	Set Consisting of	Unit
7-139-01	Ø1.2mm X 100mm Drill Bit	1
7-139-02	8 Point Slotted Drill Guide, Length 50mm	1
7-139-52	Ø1.2mm X 110mm Drill Bit	1
7-139-66	Ø1.8mm X 100mm Drill Bit	1
7-139-54	Ø1.2mm X 150mm Threaded Guide wire	1
7-139-55	Ø1.2mm X 180mm Blunt tip Guide wire	1
7-139-68	Ø3.0mm X 100mm Cannulated Reamer	1
7-139-67	Ø3.2mm X 100mm Cannulated Reamer	1
7-139-56	Ø3.5mm X 100mm Cannulated Reamer	1
7-139-57	Ø4.0mm X 100mm Cannulated Reamer	1
7-139-58	Ø4.5mm X 100mm Cannulated Reamer	1
7-139-59	Ø5.0mm X 100mm Cannulated Reamer	1
7-139-60	Ø5.5mm X 100mm Cannulated Reamer	1
7-139-61	Ø6.0mm X 100mm Cannulated Reamer	1
7-139-62	Ø6.5mm X 100mm Cannulated Reamer	1
7-139-63	Graft Sizer, 4.0mm to 7.0mm (0.5mm Steps)	1
7-139-65	Screwdriver For Ø4.0mm, Ø5.0mm & Ø6.0mm Polymer Interference Screw	1
7-139-69	Tap for 4.75mm Rotador Screw-In Anchor	1
7-139-70	3.5mm Drill Guide, Length 20mm	1
7-139-71	Drill Guide Sleeve for Ø1.5mm K-Wire	1
7-147-01	Cannulated Drill Bit with Quick Coupling End, Ø3.5mm X Length 175mm	1
7-147-02	Ø1.5mm K-Wire	1
7-147-03	Ø2.4mm Passing Pin	1

OPTIONAL INSTRUMENTS

7-139-05	Exploring Probe, Graduated, Curved	1
7-139-06	Mini Rasp, Concave Jaw	1
7-139-07	Mini Rasp, Convex Jaw	1
7-139-09	Mini Knife With V - Knotch	1
7-139-11	Mini Bone Curette	1
7-139-21	Scalpel Knife	1
7-139-22	Banana Knife	1
7-139-23	Hook Knife, Reverse Cutting	1
7-139-27	Mini Micro Fracture Awl, 45°	1
7-139-28	Mini Micro Fracture Awl, 60°	1
7-139-31	Mini Punch	1
7-139-32	Mini Scissors, Upwards Cutting	1
7-139-33	Mini Grasping Forceps	1
7-139-43	Blunt Nose Mini Basket Punch, Straight Jaws, Straight Shaft	1
7-139-48	Shark Mini Grasping Forceps, Straight Jaws, Straight Shaft	1



USA

Auxein Inc.

1500 Nw 89th Court, Suite 107-108
Doral, Florida 33172

Tel: +1 305 395 6062

E Fax: +1 305 395 6262

Email: USoffice@auxein.com

MEXICO

Auxein México, S.A. de C.V.

Tepic 139 int 801, Colonia Roma Sur,
Alcaldía Cuauhtémoc, CDMX,
México, C.P. 06760

Tel: +521 55 7261 0318

Email: info@auxein.mx

INDIA

Auxein Medical Pvt. Ltd.

Plot No. 168-169-170, Phase-4,
Kundli Industrial Area,
HSIIDC, Sector-57, Sonapat - 131028, Haryana

Tel: +91 99106 43638 | Fax: +91 86077 70197

Email: info@auxein.com

www.auxein.com