

Bone Screw




Lag screw

Kirschner Screw

IMF Screw

Anchor Screw



Neuro
Facial
Mandible
Orthognathic

Screw System	Matched Drill Bit	Pin Emergency Screw	Name	Indication
1.5mm 2.0mm 2.4mm	1.2mm 1.7mm 2.0mm	None	Kirschner Wire, Lag Screw	Applicable for fracture of ankle bump, fracture fixation in hard fixation cases by bone plates.
1.5mm 2.0mm 2.4mm	Self Tapping & Drilling Screw	None	Intermaxillary Fixation Screw, Anchorage Screw, Traction Screw	Keep normal occlusion relationship, applicable for fracture of maxillary, mandible and orthognathic surgery. Features: Screw self-tapping & self-drilling directly, ligating by titanium wire.
1.3mm 1.9mm 2.4mm	1.3mm 1.9mm 2.4mm	None	Anchor Screw	Applicable for connection between bone and tendon in head of joint because of injury, diseases, sport injury.



Anchor Screws

Length	Code
Ø1.3mmx3.5mm	200.01
Ø1.3mmx4.5mm	200.02
Ø1.9mmx4.5mm	200.03
Ø1.9mmx5.5mm	200.04
Ø2.4mmx8.5mm	200.05
Ø2.4mmx9.5mm	200.06


1.5mm Kirschner Screws


Length	Code
18mm	189.18
20mm	189.20
22mm	189.22
24mm	189.24


1.9mm Emergency Screws


Length	Code
18mm	190.18
20mm	190.20
22mm	190.22
24mm	190.24


2.0mm Kirschner Screws


Length	Code
18mm	191.18
20mm	191.20
22mm	191.22
24mm	191.24


2.3mm Emergency Screws


Length	Code
18mm	192.18
20mm	192.20
22mm	192.22
24mm	192.24


2.0mm Lag Screws


Length	Code
18mm	193.18
20mm	193.20
22mm	193.22
24mm	193.24

**1.5mm Self Tapping & Self Drilling IMF Screws,
Type-I**


Length	Code
5mm	194.05
7mm	194.07
9mm	194.09
11mm	194.11
13mm	194.13

**2.0mm Self Tapping & Self Drilling IMF Screws,
Type-I**


Length	Code
5mm	197.05
7mm	197.07
9mm	197.09
11mm	197.11
13mm	197.13

**1.5mm Self Tapping & Self Drilling IMF Screws,
Type-II**


Length	Code
5mm	195.05
7mm	195.07
9mm	195.09
11mm	195.11
13mm	195.13

**2.0mm Self Tapping & Self Drilling IMF Screws,
Type-II**


Length	Code
5mm	198.05
7mm	198.07
9mm	198.09
11mm	198.11
13mm	198.13

**1.5mm Self Tapping & Self Drilling IMF Screws,
Type-III**


Length	Code
5mm	196.05
7mm	196.07
9mm	196.09
11mm	196.11
13mm	196.13

**2.0mm Self Tapping & Self Drilling IMF Screws,
Type-III**


Length	Code
5mm	199.05
7mm	199.07
9mm	199.09
11mm	199.11
13mm	199.13