

**HERBERT SCREW (HEX HEAD),
3.0/3.9mm Cannulated Compression Screw**

HERBERT SCREW (HEX HEAD)
3.0/3.9mm Cannulated Compression Screw
Hand & Foot Surgery



about us

Auxein Medical is an integrated, research based, orthopaedic Implants & Instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and innovator products that exceed the highest expectations of our customers in term of quality and safety. The company has world-class manufacturing unit established in India and serves customers in over 60 countries worldwide.

Our Achievements


- US FDA 510K APPROVED •**
(Arthroscopy, Spine, Nailing, Locking Plates & Screws)
- EC Certificate - Class IIa & IIb •
- ISO 13485:2016 •
- ISO 9001:2015 •
- India FDA Approved •
- India MDR 2017 •

Orthopaedic Implants

DOUBLE THREADING FOR OPTIMAL COMPRESSION

Design Rationale and Main Features

The Compression Screw is easy to insert (over a guide wire) and provides efficient compression (through two separate threadings with different pitches, and an intermediate unthreaded section), thus ensuring quick, dependable internal fixation.



INDICATIONS

- Distal and proximal metatarsal osteotomies • SCARF osteotomy
- Uni and biocortical internal fixation (ex.: scaphoid) • Small bone fusion



**HERBERT SCREW (HEX HEAD),
3.0/3.9mm Cannulated Compression Screw**

Code	Length
TI-HC3.0/3.9.10	10mm, Titanium
TI-HC3.0/3.9.12	12mm, Titanium
TI-HC3.0/3.9.14	14mm, Titanium
TI-HC3.0/3.9.16	16mm, Titanium
TI-HC3.0/3.9.18	18mm, Titanium
TI-HC3.0/3.9.20	20mm, Titanium
TI-HC3.0/3.9.22	22mm, Titanium
TI-HC3.0/3.9.24	24mm, Titanium
TI-HC3.0/3.9.26	26mm, Titanium
TI-HC3.0/3.9.28	28mm, Titanium
TI-HC3.0/3.9.30	30mm, Titanium
TI-HC3.0/3.9.32	32mm, Titanium
TI-HC3.0/3.9.34	34mm, Titanium
TI-HC3.0/3.9.36	36mm, Titanium
TI-HC3.0/3.9.38	38mm, Titanium
TI-HC3.0/3.9.40	40mm, Titanium

7-042-01

Stylet Ø1.0mm



7-042-02

Protection Sleeve Ø2.2mm



7-042-03

Wire Sleeve Ø2.1/1.1mm




7-042-04

Guide Wire with Threaded Tip, Ø1.0mm x Thread Length 10mm x Length 120mm



7-042-05

Guide Wire, Ø1.1mm x Length 120mm




7-042-06

Cannulated Drill Bit Stop, Ø2.2mm x Length 150mm




7-042-07 Cannulated Drill Bit with Stopper, Ø2.2mm x Length 150mm



7-042-08 Allen Key, Hex 2.5mm for Herbert Hex Screw



7-042-09 Screw Holding Forcep for Herbert Hex Screw



7-042-10 Measuring sleeve for Guide Wire upto 60mm



7-042-11 Cannulated Screwdriver, Hex 2.0mm



7-042-12 Screwdriver, Hex 2.0mm for Hebert Hex Screw



7-042 Herbert Hex Screw Implant & Instrument Set


Code	Set Consisting of	Units
7-042-01	Stylet Ø1.0mm	1
7-042-02	Protection Sleeve Ø2.2mm	1
7-042-03	Wire Sleeve Ø2.1/1.1mm	1
7-042-04	Guide Wire with Threaded Tip, Ø1.0mm x Thread Length 10mm x Length 120mm	3
7-042-05	Guide Wire, Ø1.1mm x Length 120mm	3
7-042-06	Cannulated Drill Bit Stop, Ø2.2mm x Length 150mm	1
7-042-07	Cannulated Drill Bit with Stopper, Ø2.2mm x Length 150mm	1
7-042-08	Allen Key, Hex 2.5mm for Herbert Hex Screw	1
7-042-09	Screw Holding Forcep for Herbert Hex Screw	1
7-042-10	Measuring sleeve for Guide Wire upto 60mm	1
7-042-11	Cannulated Screwdriver, Hex 2.0mm	1
7-042-12	Screwdriver, Hex 2.0mm for Hebert Hex Screw	1
7-042-13	Container for Herbert Hex Screw Implant & Instrument Set	1



USA

Auxein Inc.
1500 Nw 89th Court, Suite 107-108
Doral, Florida 33172
Tel: +1 305 395 6062
E Fax: +1 305 395 6262
Email: USoffice@auxein.com

MEXICO

Auxein México, S.A. de C.V.
Tepic 139 int 801, Colonia Roma Sur,
Alcaldía Cuauhtémoc, CDMX,
México, C.P. 06760
Tel: +521 55 7261 0318
Email: info@auxein.mx

INDIA

Auxein Medical Pvt. Ltd.
Plot No. 168-169-170, Phase-4,
Kundli Industrial Area,
HSIIDC, Sector-57, Sonapat - 131028, Haryana
Tel: +91 99106 43638 | Fax: +91 86077 70197
Email: info@auxein.com