

# Bone Plate & Screw




**1.0mm**

1.5mm

2.0mm

2.4mm




**Neuro  
Facial  
Mandible  
Orthognathic**

Screw System	Matched Drill Bit	Pin Emergency Screw	Bone Plate (Thickness)	Indication
1.0mm	0.8mm	1.3mm	0.4mm Mini Plate	Repair and fixation of fracture of cranial base, orbital floor, nasofrontal region and maxillofacial zone of children
1.5mm	1.2mm	1.9mm	0.6mm Mini Plate	Fixation of fracture of zygoma and zygomatic arch, eye pit, maxillary
1.5mm	1.2mm	1.9mm	1.0mm Micro Plate	Correction and fixation of osteotomy of maxillary and mandible
2.0mm	1.7mm	2.3mm	0.6mm Low-incisura Plate	Fracture of Midface, correction and fixation of maxillary
2.0mm	1.7mm	2.3mm	0.8mm Micro Plate	Fracture of maxillary, fixation of orthogatha
2.0mm	1.7mm	2.3mm	1.0mm Micro Plate	Fracture of maxillary and mandible, fixation of orthogatha
2.0mm Locking	1.7mm	2.3mm	1.0mm Locking Plate	Fracture of maxillary and mandible, fixation of orthogatha and osteotomy
2.0mm Locking	1.7mm	2.3mm	1.2mm Locking Plate	Fracture of mandible, fixation of orthogatha and osteotomy
2.0mm	1.7mm	2.3mm	1.5mm One-Way Compression Plate	Fracture of mandible, fixation of osteotomy
2.4mm	2.0mm	2.7mm	2.0mm General Reconstruction Plate	Fracture of mandible, fixation of osteotomy
2.4mm Locking	2.0mm	2.7mm	2.0mm Locking Reconstruction Plate	Reconstruction and fixation of defect of mandible
2.4mm Locking	2.0mm	2.7mm	2.0mm Two-Way Compression & Locking Reconstruction Plate	Reconstruction and fixation of defect of mandible
2.4mm	2.0mm	2.7mm	2.5mm General Reconstruction Plate	Fracture of mandible, fixation of osteotomy
2.4mm Locking	2.0mm	2.7mm	2.5mm Locking Reconstruction Plate	Reconstruction and fixation of defect of mandible
2.4mm Locking	2.0mm	2.7mm	2.5mm Two-Way Compression & Locking Reconstruction Plate	Reconstruction and fixation of defect of mandible



1.0mm Straight Plate

Holes	Codes
4	01.04
6	01.06



1.0mm L-Plate, 90°

Direction	Codes
Left	02.03L
Right	02.03R



1.0mm Double Y Plate

Holes	Codes
6	03.06




1.0mm Y-Plate with Bridge

Holes	Codes
4	04.04



1.0mm Y-Plate

Holes	Codes
5	05.05




1.0mm T-Plate

Holes	Codes
6	06.06




1.0mm Square Plate

Holes	Codes
2x2	07.01
2x3	07.02



1.0mm Rectangular Grill Plate, Type-I

Holes	Codes
2x4	08.01




1.0mm Rectangular Grill Plate, Type-II

Holes	Codes
4x6	09.01




1.0mm Rectangular Grill Plate, Type-III

Holes	Codes
2x2	10.01
2x3	10.02
2x4	10.03
4x5	10.04



1.0mm Orbital Floor Plate, Type-I  
Codes 11.01




1.0mm Orbital Floor Plate, Type-II  
Codes 12.01



1.0mm Orbital Floor Plate, Type-III  
Codes 13.01



1.0mm Orbital Floor Plate, Type-IV  
Codes 13.11



1.0mm Orbital Floor Plate, Type-V  
Codes 14.01




1.0mm Orbital Floor Plate, Type-VI

Direction	Codes
Left	15.01L
Right	15.01R



1.0mm Orbital Floor Plate, Type-VII

Direction	Codes
Left	16.01L
Right	16.01R



1.0mm Orbital Floor Plate, Type-VIII

Direction	Codes
Left	17.01L
Right	17.01R


**1.0mm Self Tapping Screws**

Length	Codes
2mm	18.02
3mm	18.03
4mm	18.04
5mm	18.05
6mm	18.06
7mm	18.07


**1.0mm Self Tapping & Drilling Screws**

Length	Codes
2mm	19.02
3mm	19.03
4mm	19.04
5mm	19.05
6mm	19.06
7mm	19.07


**1.3mm Emergency Screws, Self Tapping**

Length	Codes
2mm	20.02
3mm	20.03
4mm	20.04
5mm	20.05
6mm	20.06
7mm	20.07

**7-134-11**      **Implant and Instrument Set for 1.0mm CMF System**



<b>Code</b>	<b>Item</b>	<b>Qty.</b>
<b>01.04</b>	1.0mm Straight Plate, 4 Holes	1
<b>01.06</b>	1.0mm Straight Plate, 6 Holes	1
<b>02.03L</b>	1.0mm L-Plate, 90° - Left	1
<b>02.03R</b>	1.0mm L-Plate, 90° - Right	1
<b>03.06</b>	1.0mm Double Y Plate, 6 Holes	1
<b>04.04</b>	1.0mm Y-Plate with Bridge, 4 Holes	1
<b>05.05</b>	1.0mm Y-Plate , 5 Holes	1
<b>06.06</b>	1.0mm T-Plate, 6 Holes	1
<b>07.01</b>	1.0mm Square Plate, 2 x 2 Holes	1
<b>07.02</b>	1.0mm Square Plate, 2 x 3 Holes	1
<b>08.01</b>	1.0mm Rectangular Grill Plate, Type-I, 2 x 4 Holes	1
<b>09.01</b>	1.0mm Rectangular Grill Plate, Type-II, 4 x 6 Holes	1
<b>10.01</b>	1.0mm Rectangular Grill Plate, Type-III, 2 x 2 Holes	1
<b>10.02</b>	1.0mm Rectangular Grill Plate, Type-III, 2 x 3 Holes	1
<b>10.03</b>	1.0mm Rectangular Grill Plate, Type-III, 2 x 4 Holes	1
<b>10.04</b>	1.0mm Rectangular Grill Plate, Type-III, 4 x 5 Holes	1
<b>11.01</b>	1.0mm Orbital Floor Plate, Type-I	1
<b>12.01</b>	1.0mm Orbital Floor Plate, Type-II	1
<b>13.01</b>	1.0mm Orbital Floor Plate, Type-III	1
<b>13.11</b>	1.0mm Orbital Floor Plate, Type-IV	1
<b>14.01</b>	1.0mm Orbital Floor Plate, Type-V	1
<b>15.01L</b>	1.0mm Orbital Floor Plate, Type-VI - Left	1
<b>15.01R</b>	1.0mm Orbital Floor Plate, Type-VI - Right	1
<b>16.01L</b>	1.0mm Orbital Floor Plate, Type-VII - Left	1
<b>16.01R</b>	1.0mm Orbital Floor Plate, Type-VII - Right	1
<b>17.01L</b>	1.0mm Orbital Floor Plate, Type-VIII - Left	1
<b>17.01R</b>	1.0mm Orbital Floor Plate, Type-VIII - Right	1
<b>Code</b>	<b>1.0mm Bone Screw System</b>	<b>Qty.</b>
<b>18.02</b>	1.0mm Self Tapping Screws, Length 2mm	10
<b>18.03</b>	1.0mm Self Tapping Screws, Length 3mm	10
<b>18.04</b>	1.0mm Self Tapping Screws, Length 4mm	10
<b>18.05</b>	1.0mm Self Tapping Screws, Length 5mm	10
<b>18.06</b>	1.0mm Self Tapping Screws, Length 6mm	10
<b>18.07</b>	1.0mm Self Tapping Screws, Length 7mm	10
<b>19.02</b>	1.0mm Self Tapping & Drilling Screws, Length 2mm	20
<b>19.03</b>	1.0mm Self Tapping & Drilling Screws, Length 3mm	20
<b>19.04</b>	1.0mm Self Tapping & Drilling Screws, Length 4mm	20
<b>19.05</b>	1.0mm Self Tapping & Drilling Screws, Length 5mm	20
<b>19.06</b>	1.0mm Self Tapping & Drilling Screws, Length 6mm	20
<b>19.07</b>	1.0mm Self Tapping & Drilling Screws, Length 7mm	20
<b>20.02</b>	1.3mm Emergency Screws, Self Tapping, Length 2mm	10
<b>20.03</b>	1.3mm Emergency Screws, Self Tapping, Length 3mm	10
<b>20.04</b>	1.3mm Emergency Screws, Self Tapping, Length 4mm	10
<b>20.05</b>	1.3mm Emergency Screws, Self Tapping, Length 5mm	10
<b>20.06</b>	1.3mm Emergency Screws, Self Tapping, Length 6mm	10
<b>20.07</b>	1.3mm Emergency Screws, Self Tapping, Length 7mm	10

**Implant and Instrument Set for 1.0mm CMF System**

**7-134-12** Drill Guide Ø0.8mm for CMF System



**7-134-13** Drill Bit, Ø0.8mm x Length 50mm for CMF System



**7-134-14** Drill Bit, Ø0.8mm x Length 4mm for CMF System



**7-134-15** Drill Bit, Ø0.8mm x Length 5mm for CMF System



**7-134-16** Drill Bit, Ø0.8mm x Length 6mm for CMF System





**Implant and Instrument Set for 1.0mm CMF System**

**7-134-17** Drill Bit, Ø0.8mm x Length 7mm for CMF System




**7-134-18** Self Retaining Crosshead Screwdriver Shaft for 1.0mm Screw



**7-134-19** 1.0mm Screw Caddy for CMF System



**7-134-20** Container with Tray & Lid for 1.0mm CMF Set



**Implant and Instrument Set for 1.0mm CMF System**

**7-134-11** Implant and Instrument Set for 1.0mm CMF System

Code	Item	Qty.
7-134-11	Implant and Instrument Set for 1.0mm CMF System	1
7-134-12	Drill Guide Ø0.8mm for CMF System	1
7-134-13	Drill Bit, Ø0.8mm x Length 50mm for CMF System	2
7-134-14	Drill Bit, Ø0.8mm x Length 4mm for CMF System	1
7-134-15	Drill Bit, Ø0.8mm x Length 5mm for CMF System	1
7-134-16	Drill Bit, Ø0.8mm x Length 6mm for CMF System	1
7-134-17	Drill Bit, Ø0.8mm x Length 7mm for CMF System	1
7-134-18	Self Retaining Crosshead Screwdriver Shaft for 1.0mm Screw	2
7-134-19	1.0mm Screw Caddy for CMF System	1
7-134-20	Container with Tray & Lid for 1.0mm CMF Set	1