

LOW PROFILE HAND SYSTEM

1.3mm



about us

Auxein Medical is an integrated, research based, orthopaedic Implants & Instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and innovator products that exceed the highest expectations of our customers in term of quality and safety. The company has world-class manufacturing unit established in India and serves customers in over 60 countries worldwide.

Our Achievements

US FDA 510K APPROVED •

(Arthroscopy, Spine, Nailing, Locking Plates & Screws)

EC Certificate - Class IIa & IIb •

ISO 13485:2016 •

ISO 9001:2015 •

India FDA Approved •

India MDR 2017 •

Orthopaedic Implants

1.3mm Wise-Lock Low profile Metacarpal Neck Plate

Code	Holes	
12-009-06LTI	6	Titanium
12-009-06RTI	6	Titanium



1.3mm Wise-Lock Low profile Compression Plate

Code	Holes	
12-010-04TI	4	Titanium
12-010-05TI	5	Titanium
12-010-06TI	6	Titanium



1.3mm Wise-Lock Low profile Straight Plate

Code	Holes	
12-011-03TI	3	Titanium
12-011-04TI	4	Titanium
12-011-05TI	5	Titanium
12-011-06TI	6	Titanium
12-011-07TI	7	Titanium
12-011-10TI	10	Titanium



1.3mm Wise-Lock Low profile T-Plate

Code	Holes	
12-012-04TI	4	Titanium
12-012-05TI	5	Titanium
12-012-06TI	6	Titanium
12-012-07TI	7	Titanium
12-012-08TI	8	Titanium
12-012-11TI	11	Titanium



1.3mm Wise-Lock Low Profile Rolando Fracture Hook Plate

Code	Holes	
12-013-07TI	7	Titanium



1.3mm Wise-Lock Low Profile Rotational Correction Plate

Code	Holes	
12-014-06TI	6	Titanium



1.5mm Wise-Lock Low profile Screw, Star Head

Code	Length (mm)	
12-015-05TI	5	Titanium
12-015-06TI	6	Titanium
12-015-07TI	7	Titanium
12-015-08TI	8	Titanium
12-015-09TI	9	Titanium
12-015-10TI	10	Titanium
12-015-11TI	11	Titanium
12-015-12TI	12	Titanium
12-015-13TI	13	Titanium
12-015-14TI	14	Titanium
12-015-16TI	16	Titanium
12-015-18TI	18	Titanium
12-015-20TI	20	Titanium



Star Head

1.5mm Low profile Screw, Star Head

Code	Length (mm)	
12-016-05TI	5	Titanium
12-016-06TI	6	Titanium
12-016-07TI	7	Titanium
12-016-08TI	8	Titanium
12-016-09TI	9	Titanium
12-016-10TI	10	Titanium
12-016-11TI	11	Titanium
12-016-12TI	12	Titanium
12-016-13TI	13	Titanium
12-016-14TI	14	Titanium
12-016-16TI	16	Titanium
12-016-18TI	18	Titanium
12-016-20TI	20	Titanium



Star Head

2.3mm Wise-Lock Low profile Screw, Star Head

Code	Length (mm)	
12-017-05TI	5	Titanium
12-017-06TI	6	Titanium
12-017-07TI	7	Titanium
12-017-08TI	8	Titanium
12-017-09TI	9	Titanium
12-017-10TI	10	Titanium
12-017-11TI	11	Titanium
12-017-12TI	12	Titanium
12-017-13TI	13	Titanium
12-017-14TI	14	Titanium
12-017-16TI	16	Titanium
12-017-18TI	18	Titanium
12-017-20TI	20	Titanium



Star Head

2.3mm Low profile Screw, Star Head

Code	Length (mm)	
12-018-05TI	5	Titanium
12-018-06TI	6	Titanium
12-018-07TI	7	Titanium
12-018-08TI	8	Titanium
12-018-09TI	9	Titanium
12-018-10TI	10	Titanium
12-018-11TI	11	Titanium
12-018-12TI	12	Titanium
12-018-13TI	13	Titanium
12-018-14TI	14	Titanium
12-018-16TI	16	Titanium
12-018-18TI	18	Titanium
12-018-20TI	20	Titanium



Star Head

- 7-114-01** Drill Bit with Quick Coupling Ø2.0mm x Length 90mm



- 7-114-02** Drill Bit with AO Coupling Ø2.0mm x Length 90mm



- 7-114-03** Threaded Drill Guide Ø2.0mm



- 7-095-22** Stardrive Shaft, T6



- 7-095-23** Plate Positioner 0.8/1.3mm



- 7-114-04** Depth Gauge for Ø2.3mm Screw



7-114 1.3 mm Metacarpal Implant and Instrument Case



7-114 1.3 mm Metacarpal Implant and Instrument Case
Plates

Code	Set Consisting of	Holes	Units
12-009-06LTI	1.3mm Wise-Lock Low profile Metacarpal Neck Plate, Left	6	1
12-009-06RTI	1.3mm Wise-Lock Low profile Metacarpal Neck Plate, Right	6	1
12-010-04TI	1.3mm Wise-Lock Low profile Compression Plate	4	1
12-010-05TI	1.3mm Wise-Lock Low profile Compression Plate	5	1
12-010-06TI	1.3mm Wise-Lock Low profile Compression Plate	6	1
12-011-03TI	1.3mm Wise-Lock Low profile Straight Plate	3	1
12-011-04TI	1.3mm Wise-Lock Low profile Straight Plate	4	1
12-011-05TI	1.3mm Wise-Lock Low profile Straight Plate	5	1
12-011-06TI	1.3mm Wise-Lock Low profile Straight Plate	6	1
12-011-07TI	1.3mm Wise-Lock Low profile Straight Plate	7	1
12-011-10TI	1.3mm Wise-Lock Low profile Straight Plate	10	1
12-012-04TI	1.3mm Wise-Lock Low profile T-Plate	4	1
12-012-05TI	1.3mm Wise-Lock Low profile T-Plate	5	1
12-012-06TI	1.3mm Wise-Lock Low profile T-Plate	6	1
12-012-07TI	1.3mm Wise-Lock Low profile T-Plate	7	1
12-012-08TI	1.3mm Wise-Lock Low profile T-Plate	8	1
12-012-11TI	1.3mm Wise-Lock Low profile T-Plate	11	1
12-013-07TI	1.3mm Wise-Lock Low profile Rolando Fracture Hook Plate	7	1
12-014-06TI	1.3mm Wise-Lock Low profile Rotational Correction Plate	6	1

Screws

Set Consisting of	Length	Units
1.5mm Wise-Lock Low profile screw	5mm to 20mm	6
1.5mm Low profile screw	5mm to 20mm	4
2.3mm Wise-Lock Low profile screw	5mm to 20mm	6
2.3mm Low profile screw	5mm to 20mm	4

Instruments

Code	Set Consisting of	Units
7-114-01	Ø2.0mm Drill Bit x Length 90mm with Quick Coupling	1
7-114-02	Ø2.0mm Drill Bit x Length 90mm with AO Coupling	2
7-114-03	Ø2.0mm Drill Guide	1
7-095-22	T6 Stardrive Shaft	2
7-095-23	0.8/1.3mm Plate Positioner	1
7-114-04	Depth Gauge for Ø2.3mm Screw	1

7-095-01

Bending Pliers for Low Profile Hand System

**7-095-02**


Double Periosteal Elevator, Flat, for Low Profile Hand System

**7-095-03**

Countersink for Low Profile Hand System

**7-095-04**

Position Pin for Low Profile Hand System

**7-095-05**

Retractor for Low Profile Hand System

**7-095-06**

Cleaning Hook for Low Profile Hand System



7-095-07 Bone Forcep 1.1/2.0 for Low Profile Hand System



7-095-08 Reduction Forcep with Points, for Low Profile Hand System



7-095-09 Reduction Forcep for Low Profile Hand System



7-095-10 Osteotomy Guide for Low Profile Hand System



7-095-11 Distractor for Low Profile Hand System



7-095-12 Kirschner Wire with Threaded Trocar Tip, Ø1.0mm x Thread Length 19mm x Length 150mm



7-095-13

Kirschner Wire with Threaded Trocar Tip, Ø1.2mm x Thread Length 19mm x Length 150mm

**7-095-14**

Kirschner Wire with Threaded Trocar Tip, Ø1.6mm x Thread Length 19mm x Length 150mm

**7-095-15**

Kirschner Wire, Ø1.0mm x Length 150mm

**7-095-16**

Kirschner Wire, Ø1.2mm x Length 150mm

**7-095-17**

Kirschner Wire, Ø1.6mm x Length 150mm

**7-095-18**

Handle with Mini Quick Coupling for Low Profile Hand System




7-095-19

Mini Quick Coupling (AO) for Low Profile Hand System




7-095-20


Instrument Trays for 0.8-1.3mm Low Profile Hand Instrument Set



7-095-21

Container for 0.8-1.3mm Low Profile Hand Instrument Set



7-095 0.8-1.3mm Low Profile Hand System Instrument Set


Item Code	Item Name	Unit
7-095-01	Bending Pliers for Low Profile Hand System	2
7-095-02	Double Periosteal Elevator, Flat, for Low Profile Hand System	1
7-095-03	Countersink for Low Profile Hand System	2
7-095-04	Position Pin for Low Profile Hand System	4
7-095-05	Retractor for Low Profile Hand System	2
7-095-06	Cleaning Hook for Low Profile Hand System	1
7-095-07	Bone Forcep 1.1/2.0 for Low Profile Hand System	1
7-095-08	Reduction Forcep with Points, for Low Profile Hand System	1
7-095-09	Reduction Forcep for Low Profile Hand System	1
7-095-10	Osteotomy Guide for Low Profile Hand System	1
7-095-11	Distractor for Low Profile Hand System	1
7-095-12	Kirschner Wire with Threaded Trocar Tip, Ø1.0mm x Thread Length 19mm x Length 150mm	2
7-095-13	Kirschner Wire with Threaded Trocar Tip, Ø1.2mm x Thread Length 19mm x Length 150mm	2
7-095-14	Kirschner Wire with Threaded Trocar Tip, Ø1.6mm x Thread Length 19mm x Length 150mm	2
7-095-15	Kirschner Wire, Ø1.0mm x Length 150mm	2
7-095-16	Kirschner Wire, Ø1.2mm x Length 150mm	2
7-095-17	Kirschner Wire, Ø1.6mm x Length 150mm	2
7-095-18	Handle with Mini Quick Coupling for Low Profile Hand System	1
7-095-19	Mini Quick Coupling (AO) for Low Profile Hand System	1
7-095-20	Instrument Trays for 0.8-1.3mm Low Profile Hand Instrument Set	1
7-095-21	Container for 0.8-1.3mm Low Profile Hand Instrument Set	1

**USA**

Auxein Inc.
1500 Nw 89th Court, Suite 107-108
Doral, Florida 33172
Tel: +1 305 395 6062
E Fax: +1 305 395 6262
Email: USoffice@auxein.com

MEXICO

Auxein México, S.A. de C.V.
Tepic 139 int 801, Colonia Roma Sur,
Alcaldía Cuauhtémoc, CDMX,
México, C.P. 06760
Tel: +521 55 7261 0318
Email: info@auxein.mx

INDIA

Auxein Medical Pvt. Ltd.
Plot No. 168-169-170, Phase-4,
Kundli Industrial Area,
HSIIDC, Sector-57, Sonapat - 131028, Haryana
Tel: +91 99106 43638 | Fax: +91 86077 70197
Email: info@auxein.com