

2.0mm WISE-LOCK MODULAR HAND & FOOT SYSTEM



about us

Auxein Medical is an integrated, research based, orthopaedic Implants & Instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and innovator products that exceed the highest expectations of our customers in term of quality and safety. The company has world-class manufacturing unit established in India and serves customers in over 60 countries worldwide.

Our Achievements

- US FDA 510K APPROVED •**
(Arthroscopy, Spine, Nailing, Locking Plates & Screws)
- EC Certificate - Class IIa & IIb •
- ISO 13485:2016 •
- ISO 9001:2015 •
- India FDA Approved •
- India MDR 2017 •

Orthopaedic Implants



2.0mm Wise-Lock Distal Ulna Hook Plate

Titanium	Stainless Steel	Holes
10-063-07TI	10-063-07SS	7

2.0mm Cortical Screw, Self-Tapping, (Star Head)

Code	Length
TI-102.206	6mm
TI-102.207	7mm
TI-102.208	8mm
TI-102.209	9mm
TI-102.210	10mm
TI-102.211	11mm
TI-102.212	12mm
TI-102.213	13mm
TI-102.214	14mm
TI-102.215	15mm
TI-102.216	16mm
TI-102.218	18mm
TI-102.220	20mm
TI-102.222	22mm
TI-102.224	24mm
TI-102.226	26mm
TI-102.228	28mm
TI-102.230	30mm
TI-102.232	32mm
TI-102.234	34mm
TI-102.236	36mm
TI-102.238	38mm
TI-102.240	40mm


2.0mm Wise-Lock Screw, Self-Tapping, (Star Head)

Code	Length
TI-1131.06	6mm
TI-1131.07	7mm
TI-1131.08	8mm
TI-1131.09	9mm
TI-1131.10	10mm
TI-1131.11	11mm
TI-1131.12	12mm
TI-1131.13	13mm
TI-1131.14	14mm
TI-1131.15	15mm
TI-1131.16	16mm
TI-1131.18	18mm
TI-1131.20	20mm
TI-1131.22	22mm
TI-1131.24	24mm
TI-1131.26	26mm
TI-1131.28	28mm
TI-1131.30	30mm
TI-1131.32	32mm
TI-1131.34	34mm
TI-1131.36	36mm
TI-1131.38	38mm
TI-1131.40	40mm



1136-01 Self-Centering Double Drill Guide $\varnothing 2.0/1.5\text{mm}$



1136-02 Wise-Lock Drill Sleeve 2.0mm, For Drill Bit, $\varnothing 1.5\text{mm}$



1136-03 Bending Pin with thread For 2.0mm Wise-Lock Plates



1136-04 Drill Bit Quick Coupling $\varnothing 1.5\text{mm}$ x Length 96/82mm, 2-Flute



1136-05 Drill Bit Quick Coupling $\varnothing 2.0\text{mm}$ x Length 67/55mm, 2-Flute



1136-06 Countersink For 2.0/2.4/2.7mm Screws



1136-10 Screwdriver Shaft, T6 for 2.0mm Screws



1136-08 Screw Holding Sleeve For 2.0mm Cortical Screws




1136-09 Screw Holding Sleeve For 2.0mm Wise-Lock Screws



1138-09 Container with Tray and Lid For 2.0mm Implant & Instrument Set



1138-09 Container with Tray and Lid For 2.0mm Implant & Instrument Set

Instruments for 2.0mm Wise-Lock Modular Hand System

Code	Items	Qty
1136-01	Self-Centering Double Drill Guide Ø2.0/1.5mm	1
1136-02	Wise-Lock Drill Sleeve 2.0mm, For Drill Bit, Ø1.5mm	2
1136-03	Bending Pin with thread For 2.0mm Wise-Lock Plates	2
1136-04	Drill Bit Quick Coupling Ø1.5mm x Length 96/82mm, 2-Flute	2
1136-05	Drill Bit Quick Coupling Ø2.0mm x Length 67/55mm, 2-Flute	2
1136-06	Countersink For 2.0/2.4/2.7mm Screws	1
1136-10	Screwdriver Shaft, T6 for 2.0mm Screws	1
1136-08	Screw Holding Sleeve For 2.0mm Cortical Screws	1
1136-09	Screw Holding Sleeve For 2.0mm Wise-Lock Screws	1
1138-09	Container with Tray and Lid For 2.0mm Implant & Instrument Set	1

1138-09 Container with Tray and Lid For 2.0mm Implant & Instrument Set

Implants For 2.0mm Wise-lock Modular Hand System				
	Code	Items	Qty	
Plates	TI-505.204	2.0mm Wise-Lock Plate, Straight, 4 Holes	1	
	TI-505.205	2.0mm Wise-Lock Plate, Straight, 5 Holes	1	
	TI-505.206	2.0mm Wise-Lock Plate, Straight, 6 Holes	1	
	TI-505.208	2.0mm Wise-Lock Plate, Straight, 8 Holes	1	
	TI-821.212	2.0mm Adaption Plate, Straight, 12 Holes	1	
	TI-559.207	2.0mm Wise-Lock Condylar Plate, (2 Head Holes), 7 Holes	1	
	TI-555.207T	2.0mm Wise-Lock T-Adaption Plate, (2 Head Holes), 7 Holes	1	
Screws	TI-102.206	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 6mm, Titanium	4	
	TI-102.207	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 7mm, Titanium	4	
	TI-102.208	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 8mm, Titanium	4	
	TI-102.209	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 9mm, Titanium	4	
	TI-102.210	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 10mm, Titanium	4	
	TI-102.211	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 11mm, Titanium	4	
	TI-102.212	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 12mm, Titanium	4	
	TI-102.213	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 13mm, Titanium	4	
	TI-102.214	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 14mm, Titanium	4	
	TI-102.216	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 16mm, Titanium	4	
	TI-102.218	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 18mm, Titanium	4	
	TI-102.220	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 20mm, Titanium	4	
	TI-102.222	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 22mm, Titanium	4	
	TI-102.224	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 24mm, Titanium	4	
	TI-102.226	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 26mm, Titanium	4	
	TI-102.228	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 28mm, Titanium	4	
	TI-102.230	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 30mm, Titanium	4	
	TI-102.232	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 32mm, Titanium	2	
	TI-102.234	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 34mm, Titanium	2	
	TI-102.236	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 36mm, Titanium	2	
	TI-102.238	2.0mm Cortical Screws, Self-Tapping, (Star Head), Length 38mm, Titanium	2	
	Screws	TI-1131.06	2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 6mm, Titanium	4
		TI-1131.07	2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 7mm, Titanium	4
TI-1131.08		2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 8mm, Titanium	4	
TI-1131.09		2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 9mm, Titanium	4	
TI-1131.10		2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 10mm, Titanium	4	
TI-1131.11		2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 11mm, Titanium	4	
TI-1131.12		2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 12mm, Titanium	4	
TI-1131.13		2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 13mm, Titanium	4	
TI-1131.14		2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 14mm, Titanium	4	
TI-1131.16		2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 16mm, Titanium	4	
TI-1131.18		2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 18mm, Titanium	4	
TI-1131.20		2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 20mm, Titanium	4	
TI-1131.22		2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 22mm, Titanium	4	
TI-1131.24		2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 24mm, Titanium	4	
TI-1131.26		2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 26mm, Titanium	4	
TI-1131.28		2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 28mm, Titanium	4	
TI-1131.30		2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 30mm, Titanium	4	
TI-1131.32		2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 32mm, Titanium	2	
TI-1131.34		2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 34mm, Titanium	2	
TI-1131.36		2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 36mm, Titanium	2	
TI-1131.38	2.0mm Wise-Lock Screws, Self-Tapping, (Star Head), Length 38mm, Titanium	2		



USA

Auxein Inc.
1500 Nw 89th Court, Suite 107-108
Doral, Florida 33172
Tel: +1 305 395 6062
E Fax: +1 305 395 6262
Email: USoffice@auxein.com

MEXICO

Auxein México, S.A. de C.V.
Tepic 139 int 801, Colonia Roma Sur,
Alcaldía Cuauhtémoc, CDMX,
México, C.P. 06760
Tel: +521 55 7261 0318
Email: info@auxein.mx

INDIA

Auxein Medical Pvt. Ltd.
Plot No. 168-169-170, Phase-4,
Kundli Industrial Area,
HSIIDC, Sector-57, Sonapat - 131028, Haryana
Tel: +91 99106 43638 | Fax: +91 86077 70197
Email: info@auxein.com