

AV-WISELOCK

AV-WISELOCK™

Low Profile Plates



edition
2023

 **Auxein®**

about us

Auxein Medical is an integrated, research based, orthopaedic Implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in india and serves customers in over 75 countries worldwide.

Our Achievements



AV-WISELOCK™

2.4mm AV-Wiselock Distal Radius System

New Generation Low Profile Radius Plates complete System for fracture Specific Treatment.

Indication

Distal Radius plates:

- Intra- and extraarticular fractures of the distal radius
- Correction osteotomies of the distal radius

Distal Ulna plates:

- Intra- and extraarticular fractures of the distal ulna


Plate Features

- Low-Profile plates with better anatomical Fit
- AV-Wiselock Variable angle of $\pm 15^\circ$ in all directions in each screw hole
- Pre-angled Variable Angle holes for oriented screw placement
- Rounded edges and a smooth surface for soft tissue protection
- Oblong hole for variable positioning of the plate
- Universal Variable angle drill guides available for easy angulation of screws
- K-wire holes to assist with temporary plate fixation



Anatomical Contour

First distal row can be bent individually to match the anatomy



Variable angle locking


Holes allow up to 15° off-axis screw angulation in all directions.

Oblong Hole
For Plate positioning

Plates Available in

Titanium as per ISO 5832-3 and Stainless Steel as per ISO 5832-1


2.4mm AV-Wiselock Distal Radius Small Fragment M Plate, Curved



| Titanium | | | |
|----------|-----------|----------------|--------------|
| Holes | Direction | Sterile | Non-Sterile |
| 5 | Left | 11-058-01TIL-S | 11-058-01TIL |
| 5 | Right | 11-058-01TIR-S | 11-058-01TIR |



| Stainless Steel | | | |
|-----------------|-----------|----------------|--------------|
| Holes | Direction | Sterile | Non-Sterile |
| 5 | Left | 11-058-01SSL-S | 11-058-01SSL |
| 5 | Right | 11-058-01SSR-S | 11-058-01SSR |

2.4mm AV-Wiselock T-Distal Radius Small Fragment M Plate



| Titanium | | | |
|----------|-----------|----------------|--------------|
| Holes | Direction | Sterile | Non-Sterile |
| 7 | Left | 11-058-02TIL-S | 11-058-02TIL |
| 7 | Right | 11-058-02TIR-S | 11-058-02TIR |


| Stainless Steel | | | |
|-----------------|-----------|----------------|--------------|
| Holes | Direction | Sterile | Non-Sterile |
| 7 | Left | 11-058-02SSL-S | 11-058-02SSL |
| 7 | Right | 11-058-02SSR-S | 11-058-02SSR |

2.4mm AV-Wiselock L-Distal Radius Small Fragment M Plate


| Titanium | | | | Stainless Steel | | | |
|----------|-----------|----------------|--------------|-----------------|-----------|----------------|--------------|
| Holes | Direction | Sterile | Non-Sterile | Holes | Direction | Sterile | Non-Sterile |
| 8 | Left | 11-058-03TIL-S | 11-058-03TIL | 8 | Left | 11-058-03SSL-S | 11-058-03SSL |
| 8 | Right | 11-058-03TIR-S | 11-058-03TIR | 8 | Right | 11-058-03SSR-S | 11-058-03SSR |

2.4mm AV-Wiselock Lateral Distal Radius Small Fragment M Plate


| Titanium | | | Stainless Steel | | |
|----------|---------------|-------------|-----------------|---------------|-------------|
| Holes | Sterile | Non-Sterile | Holes | Sterile | Non-Sterile |
| 6 | 11-058-04TI-S | 11-058-04TI | 6 | 11-058-04SS-S | 11-058-04SS |


2.4mm AV-Wiselock Distal Radius Volar Rim M Plate


Titanium

| Holes | Direction | Sterile | Non-Sterile |
|-------|-----------|----------------|--------------|
| 13 | Left | 11-058-05TIL-S | 11-058-05TIL |
| 13 | Right | 11-058-05TIR-S | 11-058-05TIR |

Stainless Steel

| Holes | Direction | Sterile | Non-Sterile |
|-------|-----------|----------------|--------------|
| 13 | Left | 11-058-05SSL-S | 11-058-05SSL |
| 13 | Right | 11-058-05SSR-S | 11-058-05SSR |

2.4mm AV-Wiselock Volar Lunate Facet M Plate



Titanium

| Holes | Direction | Sterile | Non-Sterile |
|-------|-----------|----------------|--------------|
| 7 | Left | 11-058-06TIL-S | 11-058-06TIL |
| 7 | Right | 11-058-06TIR-S | 11-058-06TIR |

Stainless Steel

| Holes | Direction | Sterile | Non-Sterile |
|-------|-----------|----------------|--------------|
| 7 | Left | 11-058-06SSL-S | 11-058-06SSL |
| 7 | Right | 11-058-06SSR-S | 11-058-06SSR |

2.4mm AV-Wiselock FPL Volar Distal Radius M Plate



FPL Volar Distal Radius M Plate

| Titanium | | | |
|----------|-----------|----------------|--------------|
| Holes | Direction | Sterile | Non-Sterile |
| 10 | Left | 11-058-21TIL-S | 11-058-21TIL |
| 10 | Right | 11-058-22TIR-S | 11-058-22TIR |


| Stainless Steel | | | |
|-----------------|-----------|----------------|--------------|
| Holes | Direction | Sterile | Non-Sterile |
| 10 | Left | 11-058-21SSL-S | 11-058-21SSL |
| 10 | Right | 11-058-22SSR-S | 11-058-22SSR |

FPL Volar Distal Radius M Plate, Long

| Holes | Direction | Sterile | Non-Sterile |
|-------|-----------|----------------|--------------|
| 12 | Left | 11-058-23TIL-S | 11-058-23TIL |
| 12 | Right | 11-058-24TIR-S | 11-058-24TIR |


| Holes | Direction | Sterile | Non-Sterile |
|-------|-----------|----------------|--------------|
| 12 | Left | 11-058-23SSL-S | 11-058-23SSL |
| 12 | Right | 11-058-24SSR-S | 11-058-24SSR |

2.4mm AV-Wiselock Dorsal Frame Distal Radius M Plate, Narrow





| Titanium | | | |
|----------|-----------|----------------|--------------|
| Holes | Direction | Sterile | Non-Sterile |
| 18 | Left | 11-058-08TIL-S | 11-058-08TIL |
| 18 | Right | 11-058-08TIR-S | 11-058-08TIR |

| Stainless Steel | | | |
|-----------------|-----------|----------------|--------------|
| Holes | Direction | Sterile | Non-Sterile |
| 18 | Left | 11-058-08SSL-S | 11-058-08SSL |
| 18 | Right | 11-058-08SSR-S | 11-058-08SSR |

2.4mm AV-Wiselock Dorsal Frame Distal Radius M Plate


| Titanium | | | |
|----------|-----------|----------------|--------------|
| Holes | Direction | Sterile | Non-Sterile |
| 20 | Left | 11-058-09TIL-S | 11-058-09TIL |
| 20 | Right | 11-058-09TIR-S | 11-058-09TIR |


| Stainless Steel | | | |
|-----------------|-----------|----------------|--------------|
| Holes | Direction | Sterile | Non-Sterile |
| 20 | Left | 11-058-09SSL-S | 11-058-09SSL |
| 20 | Right | 11-058-09SSR-S | 11-058-09SSR |

2.4mm AV-Wiselock Volar Distal Radius M Plate, XL


| Titanium | | | |
|----------|-----------|----------------|--------------|
| Holes | Direction | Sterile | Non-Sterile |
| 17 | Left | 11-058-10TIL-S | 11-058-10TIL |
| 17 | Right | 11-058-10TIR-S | 11-058-10TIR |
| 20 | Left | 11-058-11TIL-S | 11-058-11TIL |
| 20 | Right | 11-058-11TIR-S | 11-058-11TIR |
| 24 | Left | 11-058-12TIL-S | 11-058-12TIL |
| 24 | Right | 11-058-12TIR-S | 11-058-12TIR |

| Stainless Steel | | | |
|-----------------|-----------|----------------|--------------|
| Holes | Direction | Sterile | Non-Sterile |
| 17 | Left | 11-058-10SSL-S | 11-058-10SSL |
| 17 | Right | 11-058-10SSR-S | 11-058-10SSR |
| 20 | Left | 11-058-11SSL-S | 11-058-11SSL |
| 20 | Right | 11-058-11SSR-S | 11-058-11SSR |
| 24 | Left | 11-058-12SSL-S | 11-058-12SSL |
| 24 | Right | 11-058-12SSR-S | 11-058-12SSR |

2.4mm AV-Wiselock Volar Distal Radius Fracture M Plate




Titanium

| Holes | Direction | Sterile | Non-Sterile |
|-------|-----------|----------------|--------------|
| 11 | Left | 11-058-17TIL-S | 11-058-17TIL |
| 11 | Right | 11-058-17TIR-S | 11-058-17TIR |

Stainless Steel

| Holes | Direction | Sterile | Non-Sterile |
|-------|-----------|----------------|--------------|
| 11 | Left | 11-058-17SSL-S | 11-058-17SSL |
| 11 | Right | 11-058-17SSR-S | 11-058-17SSR |

2.4mm AV-Wiselock Volar Distal Radius Fracture M Plate, Narrow, Short



Titanium

| Holes | Direction | Sterile | Non-Sterile |
|-------|-----------|----------------|--------------|
| 10 | Left | 11-058-18TIL-S | 11-058-18TIL |
| 10 | Right | 11-058-18TIR-S | 11-058-18TIR |

Stainless Steel

| Holes | Direction | Sterile | Non-Sterile |
|-------|-----------|----------------|--------------|
| 10 | Left | 11-058-18SSL-S | 11-058-18SSL |
| 10 | Right | 11-058-18SSR-S | 11-058-18SSR |


2.4mm AV-Wiselock Distal Ulna M Plate


Titanium

| Holes | Sterile | Non-Sterile |
|-------|---------------|-------------|
| 7 | 11-058-19TI-S | 11-058-19TI |
| 10 | 11-058-20TI-S | 11-058-20TI |

Stainless Steel

| Holes | Sterile | Non-Sterile |
|-------|---------------|-------------|
| 7 | 11-058-19SS-S | 11-058-19SS |
| 10 | 11-058-20SS-S | 11-058-20SS |

2.4mm AV-Wiselock Dorsal H-Distal Radius M Plate


Titanium


| Holes | Direction | Sterile | Non-Sterile |
|-------|-----------|----------------|--------------|
| 12 | Left | 11-058-07TIL-S | 11-058-07TIL |
| 12 | Right | 11-058-07TIR-S | 11-058-07TIR |

Stainless Steel

| Holes | Direction | Sterile | Non-Sterile |
|-------|-----------|----------------|--------------|
| 12 | Left | 11-058-07SSL-S | 11-058-07SSL |
| 12 | Right | 11-058-07SSR-S | 11-058-07SSR |


2.4mm Cortical Screw, Self Tapping, (Star Head)

| Dia | Length | Titanium | | Stainless Steel | |
|------|--------|--------------|-------------|-----------------|-------------|
| | | Sterile | Non-Sterile | Sterile | Non-Sterile |
| Ø2.4 | 8mm | TI-1144.08-S | TI-1144.08 | 1144.08-S | 1144.08 |
| Ø2.4 | 10mm | TI-1144.10-S | TI-1144.10 | 1144.10-S | 1144.10 |
| Ø2.4 | 12mm | TI-1144.12-S | TI-1144.12 | 1144.12-S | 1144.12 |
| Ø2.4 | 14mm | TI-1144.14-S | TI-1144.14 | 1144.14-S | 1144.14 |
| Ø2.4 | 16mm | TI-1144.16-S | TI-1144.16 | 1144.16-S | 1144.16 |
| Ø2.4 | 18mm | TI-1144.18-S | TI-1144.18 | 1144.18-S | 1144.18 |
| Ø2.4 | 20mm | TI-1144.20-S | TI-1144.20 | 1144.20-S | 1144.20 |
| Ø2.4 | 22mm | TI-1144.22-S | TI-1144.22 | 1144.22-S | 1144.22 |
| Ø2.4 | 24mm | TI-1144.24-S | TI-1144.24 | 1144.24-S | 1144.24 |
| Ø2.4 | 26mm | TI-1144.26-S | TI-1144.26 | 1144.26-S | 1144.26 |
| Ø2.4 | 28mm | TI-1144.28-S | TI-1144.28 | 1144.28-S | 1144.28 |
| Ø2.4 | 30mm | TI-1144.30-S | TI-1144.30 | 1144.30-S | 1144.30 |
| Ø2.4 | 32mm | TI-1144.32-S | TI-1144.32 | 1144.32-S | 1144.32 |
| Ø2.4 | 34mm | TI-1144.34-S | TI-1144.34 | 1144.34-S | 1144.34 |




2.4mm Variable Angle Screws, Self Tapping, (Star Head)

| Dia | Length | Titanium | | Stainless Steel | |
|------|--------|--------------|-------------|-----------------|-------------|
| | | Sterile | Non-Sterile | Sterile | Non-Sterile |
| Ø2.4 | 8mm | TI-549.008-S | TI-549.008 | 549.008-S | 549.008 |
| Ø2.4 | 10mm | TI-549.010-S | TI-549.010 | 549.010-S | 549.010 |
| Ø2.4 | 12mm | TI-549.012-S | TI-549.012 | 549.012-S | 549.012 |
| Ø2.4 | 14mm | TI-549.014-S | TI-549.014 | 549.014-S | 549.014 |
| Ø2.4 | 16mm | TI-549.016-S | TI-549.016 | 549.016-S | 549.016 |
| Ø2.4 | 18mm | TI-549.018-S | TI-549.018 | 549.018-S | 549.018 |
| Ø2.4 | 20mm | TI-549.020-S | TI-549.020 | 549.020-S | 549.020 |
| Ø2.4 | 22mm | TI-549.022-S | TI-549.022 | 549.022-S | 549.022 |
| Ø2.4 | 24mm | TI-549.024-S | TI-549.024 | 549.024-S | 549.024 |
| Ø2.4 | 26mm | TI-549.026-S | TI-549.026 | 549.026-S | 549.026 |
| Ø2.4 | 28mm | TI-549.028-S | TI-549.028 | 549.028-S | 549.028 |
| Ø2.4 | 30mm | TI-549.030-S | TI-549.030 | 549.030-S | 549.030 |
| Ø2.4 | 32mm | TI-549.032-S | TI-549.032 | 549.032-S | 549.032 |
| Ø2.4 | 34mm | TI-549.034-S | TI-549.034 | 549.034-S | 549.034 |



2.4mm Wise-Lock Screw, Self-Tapping, (Star Head)

| Dia | Length | Titanium | | Stainless Steel | |
|------|--------|--------------|-------------|-----------------|-------------|
| | | Sterile | Non-Sterile | Sterile | Non-Sterile |
| Ø2.4 | 8mm | TI-1145.08-S | TI-1145.08 | 1145.08-S | 1145.08 |
| Ø2.4 | 10mm | TI-1145.10-S | TI-1145.10 | 1145.10-S | 1145.10 |
| Ø2.4 | 12mm | TI-1145.12-S | TI-1145.12 | 1145.12-S | 1145.12 |
| Ø2.4 | 14mm | TI-1145.14-S | TI-1145.14 | 1145.14-S | 1145.14 |
| Ø2.4 | 16mm | TI-1145.16-S | TI-1145.16 | 1145.16-S | 1145.16 |
| Ø2.4 | 18mm | TI-1145.18-S | TI-1145.18 | 1145.18-S | 1145.18 |
| Ø2.4 | 20mm | TI-1145.20-S | TI-1145.20 | 1145.20-S | 1145.20 |
| Ø2.4 | 22mm | TI-1145.22-S | TI-1145.22 | 1145.22-S | 1145.22 |
| Ø2.4 | 24mm | TI-1145.24-S | TI-1145.24 | 1145.24-S | 1145.24 |
| Ø2.4 | 26mm | TI-1145.26-S | TI-1145.26 | 1145.26-S | 1145.26 |
| Ø2.4 | 28mm | TI-1145.28-S | TI-1145.28 | 1145.28-S | 1145.28 |
| Ø2.4 | 30mm | TI-1145.30-S | TI-1145.30 | 1145.30-S | 1145.30 |



1.8mm Buttress Screws, (Star Head)

| Dia | Length | Titanium | | Stainless Steel | |
|------|--------|---------------|-------------|-----------------|-------------|
| | | Sterile | Non-Sterile | Sterile | Non-Sterile |
| Ø2.4 | 12mm | TI-1103.112-S | TI-1103.112 | 1103.112-S | 1103.112 |
| Ø2.4 | 14mm | TI-1103.114-S | TI-1103.114 | 1103.114-S | 1103.114 |
| Ø2.4 | 16mm | TI-1103.116-S | TI-1103.116 | 1103.116-S | 1103.116 |
| Ø2.4 | 18mm | TI-1103.118-S | TI-1103.118 | 1103.118-S | 1103.118 |
| Ø2.4 | 20mm | TI-1103.120-S | TI-1103.120 | 1103.120-S | 1103.120 |
| Ø2.4 | 22mm | TI-1103.122-S | TI-1103.122 | 1103.122-S | 1103.122 |
| Ø2.4 | 24mm | TI-1103.124-S | TI-1103.124 | 1103.124-S | 1103.124 |



1312-01 Depth Gauge measuring upto 60mm



1312-02 T-Handle With Quick Coupling



1312-03-15 Retractor, 15mm



1312-04-1.8 Drill Bit with Quick Coupling End, Ø1.8mm x Length 140mm



1312-05-2.0 Drill Bit with Quick Coupling End, Ø2.0mm x Length 140mm




1312-06-2.4 Drill Bit with Quick Coupling End, Ø2.4mm x Length 140mm



1312-07-2.7 Drill Bit with Quick Coupling End, Ø2.7mm x Length 140mm



1312-08-1.2 Guide Wire, Ø1.2mm x Length 150mm



1312-09-1.2 Guide Wire With Threaded Tip, Ø1.2mm x Thread Length 5mm x Length 150mm



1312-10-2.4 Tap, HA 2.4mm



1312-11-2.7 Tap, HA 2.7mm



1312-13-1.8 Threaded Drill Sleeve, Ø1.8/2.4mm (Double Start Thread)



1312-15 Self-Centering Double Drill Guide Φ 1.8/2.4mm



1312-16 Self-Centering Double Drill Guide Φ 2.0/2.7mm



1312-17 Depth Gauge with Protector measuring upto 50mm



1312-18 Star Screwdriver Shaft, T8



1312-21 Torque Limiting Attachment, 0.8Nm



7-019-03 Handle for Torque Limiting Attachment



1312-22-6.5

Retractor, 6.5mm



1312-23

Double Periosteal Elevator, Round



1312-24

Double Periosteal Elevator, Flat



1312-25

Reduction Forcep with Points, Small



1312-26

Reduction Forcep with Points, Large



1312-27

Reduction Forcep with Points, Serrated



1312-28

Bending Forcep



1312-30

Universal Variable Angle Drill Guide, Ø1.8mm



1312-20-2.7

Screw Holding Sleeve for 2.7mm Cortical Screw



1312-20-2.4

Screw Holding Sleeve for 2.4mm Cortical Screw



1312-30-2.4

Screw Holding Sleeve for 2.4mm Wise-lock Screw



1312-22

Quick Coupling Handle




1312-29 Sharp Hook




1312-31 Screwdriver Star Head






1312-35 Instrument trays for 2.4mm AV-Wiselock Distal Radius Instrument Set



1312-34 Container for 2.4mm AV-Wiselock Distal Radius Instrument



1312-000 2.4mm AV-Wiselock Distal Radius Instrument Set



| Code | Set Consisting of | Units |
|--------------------|---|-------|
| 1312-01 | Depth Gauge measuring upto 60mm | 1 |
| 1312-02 | T-Handle With Quick Coupling | 1 |
| 1312-03-15 | Retractor, 15mm | 1 |
| 1312-04-1.8 | Drill Bit with Quick Coupling End, Ø1.8mm x Length 140mm | 3 |
| 1312-05-2.0 | Drill Bit with Quick Coupling End, Ø2.0mm x Length 140mm | 2 |
| 1312-06-2.4 | Drill Bit with Quick Coupling End, Ø2.4mm x Length 140mm | 1 |
| 1312-07-2.7 | Drill Bit with Quick Coupling End, Ø2.7mm x Length 140mm | 1 |
| 1312-08-1.2 | Guide Wire, Ø1.2mm x Length 150mm | 5 |
| 1312-09-1.2 | Guide Wire With Threaded Tip, Ø1.2mm x Thread Length 5mm x Length 150mm | 2 |
| 1312-10-2.4 | Tap, HA 2.4mm | 1 |
| 1312-11-2.7 | Tap, HA 2.7mm | 1 |
| 1312-13-1.8 | Threaded Drill Sleeve, Ø1.8/2.4mm (Double Start Thread) | 3 |
| 1312-15 | Self-Centering Double Drill Guide Ø1.8/2.4mm | 1 |
| 1312-16 | Self-Centering Double Drill Guide Ø2.0/2.7mm | 1 |
| 1312-17 | Depth Gauge with Protector measuring upto 50mm | 1 |
| 1312-18 | Star Screwdriver Shaft, T8 | 2 |
| 1312-21 | Torque Limiting Attachment, 0.8Nm | 1 |
| 7-019-03 | Handle for Torque Limiting Attachment | 1 |
| 1312-22-6.5 | Retractor, 6.5mm | 1 |
| 1312-23 | Double Periosteal Elevator, Round | 1 |
| 1312-24 | Double Periosteal Elevator, Flat | 1 |
| 1312-25 | Reduction Forcep with Points, Small | 1 |
| 1312-26 | Reduction Forcep with Points, Large | 1 |
| 1312-27 | Reduction Forcep with Points, Serrated | 1 |
| 1312-28 | Bending Forcep | 2 |
| 1312-30 | Universal Variable Angle Drill Guide, Ø1.8mm | 1 |
| 1312-20-2.7 | Screw Holding Sleeve for 2.7mm Cortical Screw | 1 |
| 1312-20-2.4 | Screw Holding Sleeve for 2.4mm Cortical Screw | 1 |
| 1312-30-2.4 | Screw Holding Sleeve for 2.4mm Wise-lock Screw | 1 |
| 1312-22 | Quick Coupling Handle | 1 |
| 1312-29 | Sharp Hook | 1 |
| 1312-31 | Screwdriver Star Head | 1 |
| 1312-35 | Instrument trays for 2.4mm AV-Wiselock Distal Radius Instrument Set | 2 |
| 1312-34 | Container for 2.4mm AV-Wiselock Distal Radius Instrument | 1 |

Implant Box Optional

| | | |
|----------------|--|---|
| 1312-33 | Implant Box with cover for 2.4mm Variable Angle Distal Radius system | 1 |
| 10-050 | Implant Box with cover for 2.4mm Variable Angle Volar Rim Distal Radius System | 1 |
| 11-051 | Implant Box with cover for 2.4mm AV-Wiselock Distal Radius Dorsal & Extra-Articular System | 1 |
| 11-077 | Implant Box with cover for 2.4mm AV-Wiselock Hand and Wrist - Arthrodesis System | 1 |



USA

Auxein Inc.
1500 Nw 89th Court, Suite 107-108
Doral, Florida 33172
Tel: +1 305 395 6062
E Fax: +1 305 395 6262
Email: USoffice@auxein.com

MEXICO

Auxein México, S.A. de C.V.
Tepic 139 int 801, Colonia Roma Sur,
Alcaldía Cuauhtémoc, CDMX,
México, C.P. 06760
Tel: +521 55 7261 0318
Email: info@auxein.mx

INDIA

Auxein Medical Pvt. Ltd.
Plot No. 168-169-170, Phase-4,
Kundli Industrial Area,
HSIIDC, Sector-57, Sonapat - 131028, Haryana
Tel: +91 99106 43638 | Fax: +91 86077 70197
Email: info@auxein.com