


3.0mm Cannulated Screws




3.0mm Cannulated Cancellous Screw, Short Thread

Length (mm)	Stainless Steel	Titanium
10	123.010	TI-123.010
12	123.012	TI-123.012
14	123.014	TI-123.014
16	123.016	TI-123.016
18	123.018	TI-123.018
20	123.020	TI-123.020
22	123.022	TI-123.022
24	123.024	TI-123.024
26	123.026	TI-123.026
28	123.028	TI-123.028
30	123.030	TI-123.030
32	123.032	TI-123.032
34	123.034	TI-123.034
36	123.036	TI-123.036
38	123.038	TI-123.038
40	123.040	TI-123.040
42	123.042	TI-123.042
44	123.044	TI-123.044
46	123.046	TI-123.046
48	123.048	TI-123.048
50	123.050	TI-123.050




3.0mm Cannulated Cancellous Screw, Long Thread

Length (mm)	Stainless Steel	Titanium
14	123.214	TI-123.214
16	123.216	TI-123.216
18	123.218	TI-123.218
20	123.220	TI-123.220
22	123.222	TI-123.222
24	123.224	TI-123.224
26	123.226	TI-123.226
28	123.228	TI-123.228
30	123.230	TI-123.230
32	123.232	TI-123.232
34	123.234	TI-123.234
36	123.236	TI-123.236
38	123.238	TI-123.238
40	123.240	TI-123.240
42	123.242	TI-123.242
44	123.244	TI-123.244
46	123.246	TI-123.246
48	123.248	TI-123.248
50	123.250	TI-123.250



Washer For Small Fragment Screw

	Stainless Steel	Titanium
Code	100.110	TI-100.110




STERILE
3.0mm Cannulated Cancellous Screw, Short Thread


Length (mm)	Stainless Steel	Titanium
10	123.010-S	TI-123.010-S
12	123.012-S	TI-123.012-S
14	123.014-S	TI-123.014-S
16	123.016-S	TI-123.016-S
18	123.018-S	TI-123.018-S
20	123.020-S	TI-123.020-S
22	123.022-S	TI-123.022-S
24	123.024-S	TI-123.024-S
26	123.026-S	TI-123.026-S
28	123.028-S	TI-123.028-S
30	123.030-S	TI-123.030-S
32	123.032-S	TI-123.032-S
34	123.034-S	TI-123.034-S
36	123.036-S	TI-123.036-S
38	123.038-S	TI-123.038-S
40	123.040-S	TI-123.040-S
42	123.042-S	TI-123.042-S
44	123.044-S	TI-123.044-S
46	123.046-S	TI-123.046-S
48	123.048-S	TI-123.048-S
50	123.050-S	TI-123.050-S


3.0mm Cannulated Cancellous Screw, Long Thread

Length (mm)	Stainless Steel	Titanium
14	123.214-S	TI-123.214-S
16	123.216-S	TI-123.216-S
18	123.218-S	TI-123.218-S
20	123.220-S	TI-123.220-S
22	123.222-S	TI-123.222-S
24	123.224-S	TI-123.224-S
26	123.226-S	TI-123.226-S
28	123.228-S	TI-123.228-S
30	123.230-S	TI-123.230-S
32	123.232-S	TI-123.232-S
34	123.234-S	TI-123.234-S
36	123.236-S	TI-123.236-S
38	123.238-S	TI-123.238-S
40	123.240-S	TI-123.240-S
42	123.242-S	TI-123.242-S
44	123.244-S	TI-123.244-S
46	123.246-S	TI-123.246-S
48	123.248-S	TI-123.248-S
50	123.250-S	TI-123.250-S


Washer For Small Fragment Screw

	Stainless Steel	Titanium
Code	100.110-S	TI-100.110-S



2190-2.0-150 Cannulated Drill Bit with Quick Coupling End, Ø2.0mm x Length 150mm



3767-1.1 Double Drill Guide, Ø2.0/1.1mm



3433-01 Cannulated Screwdriver, Hex 2.5mm



3446-000 Fiber Handle with Quick Coupling, Cannulated



3430-03 Cannulated Countersink with Quick Coupling for Ø3.0mm Screw




3436-09 Stylet Ø1.1mm




3426-03 Guide Wire, Ø1.1mm x Length 150mm




3438-17 Direct Measuring Device for Cannulated Screw Ø3.0mm



3425-000 Holding Clip for Washer



3415-000 Screw Holding Forcep



7-120-01 Instrument Tray for 3.0mm Cannulated Cancellous Screw Instrument Set



7-120-02 Container for 3.0mm Cannulated Cancellous Screw Instrument Set



2335-000 3.0mm Cannulated Cancellous Screw Instrument Set


Item Code	Item Name	Unit
2190-2.0-150	Cannulated Drill Bit with Quick Coupling End, Ø2.0mm x Length 150mm	2
3767-1.1	Double Drill Guide, Ø2.0/1.1mm	1
3433-01	Cannulated Screwdriver, Hex 2.5mm	1
3446-000	Fiber Handle with Quick Coupling, Cannulated	1
3430-03	Cannulated Countersink with Quick Coupling for Ø3.0mm Screw	1
3436-09	Stylet Ø1.1mm	1
3426-03	Guide Wire, Ø1.1mm x Length 150mm	10
3438-17	Direct Measuring Device for Cannulated Screw Ø3.0mm	1
3425-000	Holding Clip for Washer	1
3415-000	Screw Holding Forcep	1
7-120-01	Instrument Tray for 3.0mm Cannulated Cancellous Screw Instrument Set	1
7-120-02	Container for 3.0mm Cannulated Cancellous Screw Instrument Set	1