


5.0mm

Wise-Lock Large Fragment
LISS System



5.0mm Wise-Lock Distal Femur Plate(LISS) System

5.0mm Wise-Lock Distal Femur Plate(LISS)



Compression and
Locking Hole



Right



Left



Holes	Left Direction	Right Direction
5	10-016-05LTI	10-016-05RTI
7	10-016-07LTI	10-016-07RTI
9	10-016-09LTI	10-016-09RTI
11	10-016-11LTI	10-016-11RTI
13	10-016-13LTI	10-016-13RTI

5.0mm Wise-Lock Screw, Self-Tapping, (Hex Head)

Length (mm)	Stainless Steel	Titanium
12	119.012	TI-119.012
14	119.014	TI-119.014
16	119.016	TI-119.016
18	119.018	TI-119.018
20	119.020	TI-119.020
22	119.022	TI-119.022
24	119.024	TI-119.024
26	119.026	TI-119.026
28	119.028	TI-119.028
30	119.030	TI-119.030
32	119.032	TI-119.032
34	119.034	TI-119.034
36	119.036	TI-119.036
38	119.038	TI-119.038
40	119.040	TI-119.040
42	119.042	TI-119.042
44	119.044	TI-119.044
46	119.046	TI-119.046
48	119.048	TI-119.048
50	119.050	TI-119.050
52	119.052	TI-119.052
54	119.054	TI-119.054
56	119.056	TI-119.056
58	119.058	TI-119.058
60	119.060	TI-119.060
65	119.065	TI-119.065
70	119.070	TI-119.070
75	119.075	TI-119.075
80	119.080	TI-119.080
85	119.085	TI-119.085
90	119.090	TI-119.090



7-028-16 Screw Caddy for 5.0mm Wise-Lock LISS System



**5.0mm Wise-Lock Screw, Self-Tapping,
(Hex Head)**

Length (mm)	Stainless Steel	Titanium	Unit
26	119.026	TI-119.026	4
28	119.028	TI-119.028	5
30	119.030	TI-119.030	5
32	119.032	TI-119.032	5
34	119.034	TI-119.034	5
36	119.036	TI-119.036	5
38	119.038	TI-119.038	5
40	119.040	TI-119.040	5
42	119.042	TI-119.042	5
44	119.044	TI-119.044	5
46	119.046	TI-119.046	5
48	119.048	TI-119.048	5
50	119.050	TI-119.050	5
52	119.052	TI-119.052	5
54	119.054	TI-119.054	5
56	119.056	TI-119.056	5
58	119.058	TI-119.058	5
60	119.060	TI-119.060	5
65	119.065	TI-119.065	5
70	119.070	TI-119.070	5
75	119.075	TI-119.075	5
80	119.080	TI-119.080	6
85	119.085	TI-119.085	6
90	119.090	TI-119.090	6

**5.0mm Wise-Lock Screw,
Self-Tapping & Self-Drilling,(Hex Head)**

Length (mm)	Stainless Steel	Titanium	Unit
26	119.226	TI-119.226	5
28	119.228	TI-119.228	5
30	119.230	TI-119.230	5
32	119.232	TI-119.232	5
34	119.234	TI-119.234	5
36	119.236	TI-119.236	5
38	119.238	TI-119.238	5
40	119.240	TI-119.240	5
42	119.242	TI-119.242	6
44	119.244	TI-119.244	5
46	119.246	TI-119.246	5
48	119.248	TI-119.248	5
50	119.250	TI-119.250	5
55	119.255	TI-119.255	5
60	119.260	TI-119.260	5
65	119.265	TI-119.265	5
70	119.270	TI-119.270	6
75	119.275	TI-119.275	6
80	119.280	TI-119.280	6
85	119.285	TI-119.285	6
90	119.290	TI-119.290	6

7-028-01 L Wrench, Ø4.5mm for LISS System



7-028-02 T-Handle With Quick Coupling for LISS System




7-028-03 Tap with Quick Coupling End, Ø5.0mm for LISS System




7-028-04 Torque-Limiting Screwdriver Handle, 4.0Nm, for LISS System



7-028-05 Guide Sleeve, Ø7.3mm, for LISS System



- 7-028-06** Threaded Guide Pin, Ø2.0mm x Thread Length 11mm x Total Length 280mm, for LISS System



- 7-028-07** Trocar Ø7.2mm for LISS System



- 7-028-08** Lifting Drill Ø4.0mm for LISS System



- 7-028-09** Screwdriver, Hex 3.5mm, for LISS System



- 7-028-10** Drill sleeve for Locking Screw, Ø2.2mm, for LISS System



7-028-11 Drill Sleeve for Pin, Ø2.2mm, for LISS System



7-028-12 Drill Sleeve, Ø4.2mm, for LISS System



7-028-13 Depth Gauge, Ø2.2mm, for LISS System



7-028-14 Drill Bit Ø4.1 × 280mm for LISS System



7-028-15 Stopper for LISS System



7-028-16 Screw Caddy for 5.0mm Wise-Lock LISS System



7-029-01 5.0mm Wise-Lock Distal Femur (LISS) Guide, Left



7-029-02 5.0mm Wise-Lock Distal Femur (LISS) Guide, Right



7-029-03 Instrument Trays for 5.0mm Wise-Lock Distal Femur Plate(LISS) System



7-029-04 Container for 5.0mm Wise-Lock Distal Femur Plate(LISS) System



7-029 5.0mm Wise-Lock Distal Femur Plate (LISS) Instrument set



7-029 5.0mm Wise-Lock Distal Femur Plate (LISS) Instrument set

Codes	Set Consisting of:	Units
7-028-01	L Wrench, Ø4.5mm for LISS System	1
7-028-02	T-Handle With Quick Coupling for LISS System	1
7-028-03	Tap with Quick Coupling End, Ø5.0mm for LISS System	2
7-028-04	Torque-Limiting Screwdriver Handle, 4.0Nm, for LISS System	1
7-028-05	Guide Sleeve, Ø7.3mm, for LISS System	4
7-028-06	Threaded Guide Pin, Ø2.0mm x Thread Length 11mm x Total Length 280mm, for LISS System	4
7-028-07	Trocar Ø7.2mm for LISS System	1
7-028-08	Lifting Drill Ø4.0mm for LISS System	2
7-028-09	Screwdriver, Hex 3.5mm, for LISS System	2
7-028-10	Drill sleeve for Locking Screw, Ø2.2mm, for LISS System	1
7-028-11	Drill Sleeve for Pin, Ø2.2mm, for LISS System	2
7-028-12	Drill Sleeve, Ø4.2mm, for LISS System	2
7-028-13	Depth Gauge, Ø2.2mm, for LISS System	1
7-028-14	Drill Bit Ø4.1 x 280mm for LISS System	2
7-028-15	Stopper for LISS System	20
7-028-16	Screw Caddy for 5.0mm Wise-Lock LISS System	1
7-029-01	5.0mm Wise-Lock Distal Femur (LISS) Guide, Left	1
7-029-02	5.0mm Wise-Lock Distal Femur (LISS) Guide, Right	1
7-029-03	Instrument Trays for 5.0mm Wise-Lock Distal Femur Plate(LISS) System	2
7-029-04	Container for 5.0mm Wise-Lock Distal Femur Plate(LISS) System	1




5.0mm Wise-Lock Proximal Lateral Tibia Plate(LISS) System

5.0mm Wise-Lock Proximal Lateral Tibia Plate(LISS)


Compression and
Locking Hole



Right



Left



Holes	Left Direction	Right Direction
5	10-017-05LTI	10-017-05RTI
7	10-017-07LTI	10-017-07RTI
9	10-017-09LTI	10-017-09RTI
11	10-017-11LTI	10-017-11RTI
13	10-017-13LTI	10-017-13RTI

5.0mm Wise-Lock Screw, Self-Tapping, (Hex Head)

Length (mm)	Stainless Steel	Titanium
12	119.012	TI-119.012
14	119.014	TI-119.014
16	119.016	TI-119.016
18	119.018	TI-119.018
20	119.020	TI-119.020
22	119.022	TI-119.022
24	119.024	TI-119.024
26	119.026	TI-119.026
28	119.028	TI-119.028
30	119.030	TI-119.030
32	119.032	TI-119.032
34	119.034	TI-119.034
36	119.036	TI-119.036
38	119.038	TI-119.038
40	119.040	TI-119.040
42	119.042	TI-119.042
44	119.044	TI-119.044
46	119.046	TI-119.046
48	119.048	TI-119.048
50	119.050	TI-119.050
52	119.052	TI-119.052
54	119.054	TI-119.054
56	119.056	TI-119.056
58	119.058	TI-119.058
60	119.060	TI-119.060
65	119.065	TI-119.065
70	119.070	TI-119.070
75	119.075	TI-119.075
80	119.080	TI-119.080
85	119.085	TI-119.085
90	119.090	TI-119.090



7-028-16 Screw Caddy for 5.0mm Wise-Lock LISS System



**5.0mm Wise-Lock Screw, Self-Tapping,
(Hex Head)**

Length (mm)	Stainless Steel	Titanium	Unit
26	119.026	TI-119.026	4
28	119.028	TI-119.028	5
30	119.030	TI-119.030	5
32	119.032	TI-119.032	5
34	119.034	TI-119.034	5
36	119.036	TI-119.036	5
38	119.038	TI-119.038	5
40	119.040	TI-119.040	5
42	119.042	TI-119.042	5
44	119.044	TI-119.044	5
46	119.046	TI-119.046	5
48	119.048	TI-119.048	5
50	119.050	TI-119.050	5
52	119.052	TI-119.052	5
54	119.054	TI-119.054	5
56	119.056	TI-119.056	5
58	119.058	TI-119.058	5
60	119.060	TI-119.060	5
65	119.065	TI-119.065	5
70	119.070	TI-119.070	5
75	119.075	TI-119.075	5
80	119.080	TI-119.080	6
85	119.085	TI-119.085	6
90	119.090	TI-119.090	6

**5.0mm Wise-Lock Screw,
Self-Tapping & Self-Drilling,(Hex Head)**

Length (mm)	Stainless Steel	Titanium	Unit
26	119.226	TI-119.226	5
28	119.228	TI-119.228	5
30	119.230	TI-119.230	5
32	119.232	TI-119.232	5
34	119.234	TI-119.234	5
36	119.236	TI-119.236	5
38	119.238	TI-119.238	5
40	119.240	TI-119.240	5
42	119.242	TI-119.242	6
44	119.244	TI-119.244	5
46	119.246	TI-119.246	5
48	119.248	TI-119.248	5
50	119.250	TI-119.250	5
55	119.255	TI-119.255	5
60	119.260	TI-119.260	5
65	119.265	TI-119.265	5
70	119.270	TI-119.270	6
75	119.275	TI-119.275	6
80	119.280	TI-119.280	6
85	119.285	TI-119.285	6
90	119.290	TI-119.290	6

7-028-01 L Wrench, Ø4.5mm for LISS System



7-028-02 T-Handle With Quick Coupling for LISS System




7-028-03 Tap with Quick Coupling End, Ø5.0mm for LISS System



7-028-04 Torque-Limiting Screwdriver Handle, 4.0Nm, for LISS System



7-028-05 Guide Sleeve, Ø7.3mm, for LISS System



7-028-06 Threaded Guide Pin, Ø2.0mm x Thread Length 11mm x Total Length 280mm, for LISS System



7-028-07 Trocar Ø7.2mm for LISS System



7-028-08 Lifting Drill Ø4.0mm for LISS System




7-028-09 Screwdriver, Hex 3.5mm, for LISS System



7-028-10 Drill sleeve for Locking Screw, Ø2.2mm, for LISS System



7-028-11 Drill Sleeve for Pin, Ø2.2mm, for LISS System



7-028-12 Drill Sleeve, Ø4.2mm, for LISS System



7-028-13 Depth Gauge, Ø2.2mm, for LISS System



7-028-14 Drill Bit Ø4.1 × 280mm for LISS System



7-028-15 Stopper for LISS System



7-028-16 Screw Caddy for 5.0mm Wise-Lock LISS System



7-030-01 5.0mm Wise-Lock Proximal Lateral Tibia(LISS) Guide, Left



7-030-02 5.0mm Wise-Lock Proximal Lateral Tibia(LISS) Guide, Right



7-030-03 Trays for 5.0mm Wise-Lock Proximal Lateral Tibia Plate(LISS) Instrument set



7-030-04 Container for 5.0mm Wise-Lock Proximal Lateral Tibia Plate(LISS) Instrument set



7-030 5.0mm Wise-Lock Proximal Lateral Tibia Plate (LISS) Instrument set



7-030 5.0mm Wise-Lock Proximal Lateral Tibia Plate (LISS) Instrument set

Codes	Set Consisting of:	Units
7-028-01	L Wrench, Ø4.5mm for LISS System	1
7-028-02	T-Handle With Quick Coupling for LISS System	1
7-028-03	Tap with Quick Coupling End, Ø5.0mm for LISS System	2
7-028-04	Torque-Limiting Screwdriver Handle, 4.0Nm, for LISS System	1
7-028-05	Guide Sleeve, Ø7.3mm, for LISS System	4
7-028-06	Threaded Guide Pin, Ø2.0mm x Thread Length 11mm x Total Length 280mm, for LISS System	4
7-028-07	Trocar Ø7.2mm for LISS System	1
7-028-08	Lifting Drill Ø4.0mm for LISS System	2
7-028-09	Screwdriver, Hex 3.5mm, for LISS System	2
7-028-10	Drill sleeve for Locking Screw, Ø2.2mm, for LISS System	1
7-028-11	Drill Sleeve for Pin, Ø2.2mm, for LISS System	2
7-028-12	Drill Sleeve, Ø4.2mm, for LISS System	2
7-028-13	Depth Gauge, Ø2.2mm, for LISS System	1
7-028-14	Drill Bit Ø4.1 × 280mm for LISS System	2
7-028-15	Stopper for LISS System	20
7-028-16	Screw Caddy for 5.0mm Wise-Lock LISS System	1
7-030-01	5.0mm Wise-Lock Proximal Lateral Tibia(LISS) Guide, Left	1
7-030-02	5.0mm Wise-Lock Proximal Lateral Tibia(LISS) Guide, Right	1
7-030-03	Trays for 5.0mm Wise-Lock Proximal Lateral Tibia Plate(LISS) Instrument set	2
7-030-04	Container for 5.0mm Wise-Lock Proximal Lateral Tibia Plate(LISS) Instrument set	1