

7.3mm Cannulated Screws


7.3mm Cannulated Cancellous Screw, 16mm Thread, Self-Drilling


Length (mm)	Stainless Steel	Titanium
30	134.030	TI-134.030
35	134.035	TI-134.035
40	134.040	TI-134.040
45	134.045	TI-134.045
50	134.050	TI-134.050
55	134.055	TI-134.055
60	134.060	TI-134.060
65	134.065	TI-134.065
70	134.070	TI-134.070
75	134.075	TI-134.075
80	134.080	TI-134.080
85	134.085	TI-134.085
90	134.090	TI-134.090
95	134.095	TI-134.095
100	134.100	TI-134.100
105	134.105	TI-134.105
110	134.110	TI-134.110
* 115	134.115	TI-134.115
* 120	134.120	TI-134.120
* 125	134.125	TI-134.125
* 130	134.130	TI-134.130
* 135	134.135	TI-134.135
* 140	134.140	TI-134.140
* 145	134.145	TI-134.145
* 150	134.150	TI-134.150

* Sizes not available in Screw Caddy


7.3mm Cannulated Cancellous Screw, 32mm Thread, Self-Drilling

Length (mm)	Stainless Steel	Titanium
35	134.235	TI-134.235
40	134.240	TI-134.240
45	134.245	TI-134.245
50	134.250	TI-134.250
55	134.255	TI-134.255
60	134.260	TI-134.260
65	134.265	TI-134.265
70	134.270	TI-134.270
75	134.275	TI-134.275
80	134.280	TI-134.280
85	134.285	TI-134.285
90	134.290	TI-134.290
95	134.295	TI-134.295
100	134.300	TI-134.300
105	134.305	TI-134.305
110	134.310	TI-134.310
* 115	134.315	TI-134.315
* 120	134.320	TI-134.320
* 125	134.325	TI-134.325
* 130	134.330	TI-134.330
* 135	134.335	TI-134.335
* 140	134.340	TI-134.340
* 145	134.345	TI-134.345
* 150	134.350	TI-134.350

* Sizes not available in Screw Caddy



STERILE
7.3mm Cannulated Cancellous Screw, 16mm Thread, Self-Drilling


Length (mm)	Stainless Steel	Titanium
30	134.030-S	TI-134.030-S
35	134.035-S	TI-134.035-S
40	134.040-S	TI-134.040-S
45	134.045-S	TI-134.045-S
50	134.050-S	TI-134.050-S
55	134.055-S	TI-134.055-S
60	134.060-S	TI-134.060-S
65	134.065-S	TI-134.065-S
70	134.070-S	TI-134.070-S
75	134.075-S	TI-134.075-S
80	134.080-S	TI-134.080-S
85	134.085-S	TI-134.085-S
90	134.090-S	TI-134.090-S
95	134.095-S	TI-134.095-S
100	134.100-S	TI-134.100-S
105	134.105-S	TI-134.105-S
110	134.110-S	TI-134.110-S
* 115	134.115-S	TI-134.115-S
* 120	134.120-S	TI-134.120-S
* 125	134.125-S	TI-134.125-S
* 130	134.130-S	TI-134.130-S
* 135	134.135-S	TI-134.135-S
* 140	134.140-S	TI-134.140-S
* 145	134.145-S	TI-134.145-S
* 150	134.150-S	TI-134.150-S

* Sizes not available in Screw Caddy


7.3mm Cannulated Cancellous Screw, 32mm Thread, Self-Drilling

Length (mm)	Stainless Steel	Titanium
35	134.235-S	TI-134.235-S
40	134.240-S	TI-134.240-S
45	134.245-S	TI-134.245-S
50	134.250-S	TI-134.250-S
55	134.255-S	TI-134.255-S
60	134.260-S	TI-134.260-S
65	134.265-S	TI-134.265-S
70	134.270-S	TI-134.270-S
75	134.275-S	TI-134.275-S
80	134.280-S	TI-134.280-S
85	134.285-S	TI-134.285-S
90	134.290-S	TI-134.290-S
95	134.295-S	TI-134.295-S
100	134.300-S	TI-134.300-S
105	134.305-S	TI-134.305-S
110	134.310-S	TI-134.310-S
* 115	134.315-S	TI-134.315-S
* 120	134.320-S	TI-134.320-S
* 125	134.325-S	TI-134.325-S
* 130	134.330-S	TI-134.330-S
* 135	134.335-S	TI-134.335-S
* 140	134.340-S	TI-134.340-S
* 145	134.345-S	TI-134.345-S
* 150	134.350-S	TI-134.350-S

* Sizes not available in Screw Caddy




7.3mm Cannulated Cancellous Screw, Full Thread, Self-Drilling

Length (mm)	Stainless Steel	Titanium
20	148.020	TI-148.020
25	148.025	TI-148.025
30	148.030	TI-148.030
35	148.035	TI-148.035
40	148.040	TI-148.040
45	148.045	TI-148.045
50	148.050	TI-148.050
55	148.055	TI-148.055
60	148.060	TI-148.060
65	148.065	TI-148.065
70	148.070	TI-148.070
75	148.075	TI-148.075
80	148.080	TI-148.080
85	148.085	TI-148.085
90	148.090	TI-148.090
95	148.095	TI-148.095
100	148.100	TI-148.100
105	148.105	TI-148.105
110	148.110	TI-148.110
* 115	148.115	TI-148.115
* 120	148.120	TI-148.120
* 125	148.125	TI-148.125
* 130	148.130	TI-148.130



* Sizes not available in Screw Caddy

Washer For Large Fragment Screw

	Stainless Steel	Titanium
Code	100.120	TI-100.120



STERILE
7.3mm Cannulated Cancellous Screw, Full Thread, Self-Drilling

Length (mm)	Stainless Steel	Titanium
20	148.020-S	TI-148.020-S
25	148.025-S	TI-148.025-S
30	148.030-S	TI-148.030-S
35	148.035-S	TI-148.035-S
40	148.040-S	TI-148.040-S
45	148.045-S	TI-148.045-S
50	148.050-S	TI-148.050-S
55	148.055-S	TI-148.055-S
60	148.060-S	TI-148.060-S
65	148.065-S	TI-148.065-S
70	148.070-S	TI-148.070-S
75	148.075-S	TI-148.075-S
80	148.080-S	TI-148.080-S
85	148.085-S	TI-148.085-S
90	148.090-S	TI-148.090-S
95	148.095-S	TI-148.095-S
100	148.100-S	TI-148.100-S
105	148.105-S	TI-148.105-S
110	148.110-S	TI-148.110-S
* 115	148.115-S	TI-148.115-S
* 120	148.120-S	TI-148.120-S
* 125	148.125-S	TI-148.125-S
* 130	148.130-S	TI-148.130-S




* Sizes not available in Screw Caddy

Washer For Large Fragment Screw

	Stainless Steel	Titanium
Code	100.120-S	TI-100.120-S



ST-005-7.3 Screw Caddy for 7.3mm Cannulated Cancellous Screw System
(Tray Without Screw)




**7.3mm Cannulated Cancellous Screw,
16mm Thread, Self-Drilling**

Stainless Steel	Titanium	Length	Unit
134.030	TI-134.030	30mm	2
134.035	TI-134.035	35mm	2
134.040	TI-134.040	40mm	2
134.045	TI-134.045	45mm	2
134.050	TI-134.050	50mm	2
134.055	TI-134.055	55mm	2
134.060	TI-134.060	60mm	2
134.065	TI-134.065	65mm	2
134.070	TI-134.070	70mm	2
134.075	TI-134.075	75mm	4
134.080	TI-134.080	80mm	4
134.085	TI-134.085	85mm	4
134.090	TI-134.090	90mm	4
134.095	TI-134.095	95mm	2
134.100	TI-134.100	100mm	2
134.105	TI-134.105	105mm	2
134.110	TI-134.110	110mm	2

**7.3mm Cannulated Cancellous Screw,
Full Thread, Self-Drilling**

Stainless Steel	Titanium	Length	Unit
148.030	TI-148.030	30mm	2
148.035	TI-148.035	35mm	2
148.040	TI-148.040	40mm	2
148.045	TI-148.045	45mm	2
148.050	TI-148.050	50mm	2
148.055	TI-148.055	55mm	2
148.060	TI-148.060	60mm	2
148.065	TI-148.065	65mm	2
148.070	TI-148.070	70mm	2
148.075	TI-148.075	75mm	4
148.080	TI-148.080	80mm	4
148.085	TI-148.085	85mm	4
148.090	TI-148.090	90mm	4
148.095	TI-148.095	95mm	2
148.100	TI-148.100	100mm	2
148.105	TI-148.105	105mm	2
148.110	TI-148.110	110mm	2

**7.3mm Cannulated Cancellous Screw,
32mm Thread, Self-Drilling**

Stainless Steel	Titanium	Length	Unit
134.235	TI-134.235	35mm	2
134.240	TI-134.240	40mm	2
134.245	TI-134.245	45mm	2
134.250	TI-134.250	50mm	2
134.255	TI-134.255	55mm	2
134.260	TI-134.260	60mm	2
134.265	TI-134.265	65mm	2
134.270	TI-134.270	70mm	2
134.275	TI-134.275	75mm	4
134.280	TI-134.280	80mm	4
134.285	TI-134.285	85mm	4
134.290	TI-134.290	90mm	4
134.295	TI-134.295	95mm	2
134.300	TI-134.300	100mm	2
134.305	TI-134.305	105mm	2
134.310	TI-134.310	110mm	2

GW-2.8-250T Guide Wire with Threaded Tip, $\text{\O}2.8\text{mm}$ x Thread Length 10mm x Length 250mm



GW-2.8-250 Guide Wire, $\text{\O}2.8\text{mm}$ x Length 250mm



1532-01 Drill Bit Cannulated Plain Shank/Jacob Chuck End, $\text{\O}5.0\text{mm}$ x Length 250mm




1532-05 Stopper for Drill Bit $\text{\O}5.0\text{mm}$



1539-01 Allen Key, Hex 3mm, for $\text{\O}7.3\text{mm}$ Screw



3436-014 Direct Measuring Device for Cannulated Cancellous Screw




1538-01 Cannulated Countersink with T-Handle for Ø7.3mm Screw



1543-01 Wire Sleeve Ø10.0/2.9mm



1542-01 Drill Guide for Drill Bit Ø5.0mm



1541-01 Guide Sleeve Ø5.0/2.9mm



1537-01 Protection Sleeve for Ø10.0/2.9mm



1482-01 Wire Sleeve Bolt Key



1544-01 Wire Sleeve, Adjustable




1534-01 Cannulated Tap with T-Handle, for Ø7.3mm Screw



1533-01 Parallel Guide for Guide Wire Ø2.8mm



7-058-01 Screw Holding Forcep for Ø7.3mm Screw



1478-2.5 Stylet, Ø2.5mm



3433-03 Hexagonal Screwdriver, Cannulation Ø2.9mm, Hex 4.0mm



1540-01 Hexagonal Screwdriver, Hex 4.0mm for Ø7.3mm Screw



1540-05 Holding Sleeve for Cannulated Hexagonal Screwdriver, Hex 4.0mm




1545-01 Extraction Screw with T-Handle, for Ø7.3mm Screw



ST-005-7.3 Screw Caddy for 7.3mm Cannulated Cancellous Screw



IS006-7.3T Instrument Tray for 7.3mm Cannulated Screw Instrument Set



IS006-7.3B Container For 7.3mm Cannulated Cancellous Screw Instrument Set



IS006-7.3 7.3mm Cannulated Cancellous Screw Instrument Set



IS006-7.3 7.3mm Cannulated Cancellous Screw Instrument Set

Item Code	Item Name	Unit
GW-2.8-250T	Guide Wire with Threaded Tip, Ø2.8mm x Thread Length 10mm x Length 250mm	3
GW-2.8-250	Guide Wire, Ø2.8mm x Length 250mm	3
1532-01	Drill Bit Cannulated Plain Shank/Jacob Chuck End, Ø5.0mm x Length 250mm	2
1532-05	Stopper for Drill Bit Ø5.0mm	1
1539-01	Allen Key, Hex 3mm, for Ø7.3mm Screw	1
3436-014	Direct Measuring Device for Cannulated Cancellous Screw	1
1538-01	Cannulated Countersink with T-Handle for Ø7.3mm Screw	1
1543-01	Wire Sleeve Ø10.0/2.9mm	1
1542-01	Drill Guide for Drill Bit Ø5.0mm	1
1541-01	Guide Sleeve Ø5.0/2.9mm	1
1537-01	Protection Sleeve for Ø10.0/2.9mm	1
1482-01	Wire Sleeve Bolt Key	1
1544-01	Wire Sleeve, Adjustable	1
1534-01	Cannulated Tap with T-Handle, for Ø7.3mm Screw	1
1533-01	Parallel Guide for Guide Wire Ø2.8mm	1
7-058-01	Screw Holding Forcep for Ø7.3mm Screw	1
1478-2.5	Stylet, Ø2.5mm	1
3433-03	Hexagonal Screwdriver, Cannulation Ø2.9mm, Hex 4.0mm	1
1540-01	Hexagonal Screwdriver, Hex 4.0mm for Ø7.3mm Screw	1
1540-05	Holding Sleeve for Cannulated Hexagonal Screwdriver, Hex 4.0mm	1
1545-01	Extraction Screw with T-Handle, for Ø7.3mm Screw	1
ST-005-7.3	Screw Caddy for 7.3mm Cannulated Cancellous Screw	1
IS006-7.3T	Instrument Tray for 7.3mm Cannulated Screw Instrument Set	2
IS006-7.3B	Container For 7.3mm Cannulated Cancellous Screw Instrument Set	1