

AV-WISELOCK

AV-WISELOCK™

2.7/3.5mm Small Fragment



edition
2023

 **Auxein®**

about us

Auxein Medical is an integrated, research based, orthopaedic Implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in india and serves customers in over 75 countries worldwide.

Our Achievements



AV-WISELOCK

AV-WISELOCK™

AV-Wiselock 2.7/3.5mm Small Fragment

AV-Wiselock 2.7/3.5mm Ankle Plates

AV-Wiselock 2.7/3.5mm Elbow Plates

AV-Wiselock 2.7/3.5mm Proximal Tibia Plates

AV-Wiselock 2.7/3.5mm Foot Plates

PHELOS - 3.5mm AV-Wiselock Proximal Humerus Plate With Suture Cuts

3.5mm AV-Wiselock Medial Proximal Tibia Plate

3.5mm AV-Wiselock Extra-Articular Distal Humerus Plate

2.7/3.5mm AV Wiselock Superior Anterior Clavicle Plate with Lateral Extension

2.7/3.5mm

AV-WiseLock System

AV-WISELOCK

2.7/3.5mm ANKLE SYSTEM


INTRODUCTION

The Plates are part of Auxein 2.7/3.5mm AV-Wiselock Ankle System merging the AV-Wiselock locking technology with the conventional plating. The plates are anatomically pre-contoured to fit the bone anatomy.

The AV-Wiselock Ankle System is intended for fixation of the ankle in adults and adolescents (12-21) in which the growth plates have fused.


Specifically, Medial and Anteromedial Distal Tibia Plates are intended for fixation of osteotomies, fractures, nonunions, malunions, and replantations of bones and bone fragments of the diaphyseal and metaphyseal regions of the distal tibia, Distal Tibia T Plates and Distal Tibia L Plates are intended to buttress partial articular fractures and bone fragments of the distal tibia, and Lateral Distal Fibula Plates are intended for fixation of osteotomies, fractures, nonunions, malunions, and replantations of bones and bone fragments of the diaphyseal and metaphyseal regions of the distal fibula.

- The combi holes in the plate shaft combine dynamic compression unit with the AV-Wiselock locking screw holes.
- AV-Wiselock holes in the plate head is designed to provide fixed-angle construct at any desired angle in a range of 0° to 15° from the central axis.
- The K-wire holes accept wires up to 1.6mm in diameter to provisionally maintain the reduction.
- Enables choice of screw functionality in every combi-hole
- Choice of syndesmotic fixation placement with dedicated slots




Plates Available in

Titanium as per ISO 5832-3 and Stainless Steel as per ISO 5832-1



2.7/3.5mm AV-Wiselock Medial Distal Tibia Plate


Titanium				Stainless Steel			
Holes	Direction	Sterile	Non-Sterile	Holes	Direction	Sterile	Non-Sterile
4	Left	11-019-04LTI-S	11-019-04LTI	4	Left	11-019-04LSS-S	11-019-04LSS
4	Right	11-019-04RTI-S	11-019-04RTI	4	Right	11-019-04RSS-S	11-019-04RSS
6	Left	11-019-06LTI-S	11-019-06LTI	6	Left	11-019-06LSS-S	11-019-06LSS
6	Right	11-019-06RTI-S	11-019-06RTI	6	Right	11-019-06RSS-S	11-019-06RSS
8	Left	11-019-08LTI-S	11-019-08LTI	8	Left	11-019-08LSS-S	11-019-08LSS
8	Right	11-019-08RTI-S	11-019-08RTI	8	Right	11-019-08RSS-S	11-019-08RSS
10	Left	11-019-10LTI-S	11-019-10LTI	10	Left	11-019-10LSS-S	11-019-10LSS
10	Right	11-019-10RTI-S	11-019-10RTI	10	Right	11-019-10RSS-S	11-019-10RSS
12	Left	11-019-12LTI-S	11-019-12LTI	12	Left	11-019-12LSS-S	11-019-12LSS
12	Right	11-019-12RTI-S	11-019-12RTI	12	Right	11-019-12RSS-S	11-019-12RSS
14	Left	11-019-14LTI-S	11-019-14LTI	14	Left	11-019-14LSS-S	11-019-14LSS
14	Right	11-019-14RTI-S	11-019-14RTI	14	Right	11-019-14RSS-S	11-019-14RSS
16	Left	11-019-16LTI-S	11-019-16LTI	16	Left	11-019-16LSS-S	11-019-16LSS
16	Right	11-019-16RTI-S	11-019-16RTI	16	Right	11-019-16RSS-S	11-019-16RSS

2.7/3.5mm AV-Wiselock Anteromedial Distal Tibia Plate


Titanium			
Holes	Direction	Sterile	Non-Sterile
4	Left	11-020-04LTI-S	11-020-04LTI
4	Right	11-020-04RTI-S	11-020-04RTI
6	Left	11-020-06LTI-S	11-020-06LTI
6	Right	11-020-06RTI-S	11-020-06RTI
8	Left	11-020-08LTI-S	11-020-08LTI
8	Right	11-020-08RTI-S	11-020-08RTI
10	Left	11-020-10LTI-S	11-020-10LTI
10	Right	11-020-10RTI-S	11-020-10RTI
12	Left	11-020-12LTI-S	11-020-12LTI
12	Right	11-020-12RTI-S	11-020-12RTI
14	Left	11-020-14LTI-S	11-020-14LTI
14	Right	11-020-14RTI-S	11-020-14RTI
16	Left	11-020-16LTI-S	11-020-16LTI
16	Right	11-020-16RTI-S	11-020-16RTI

Stainless Steel			
Holes	Direction	Sterile	Non-Sterile
4	Left	11-020-04LSS-S	11-020-04LSS
4	Right	11-020-04RSS-S	11-020-04RSS
6	Left	11-020-06LSS-S	11-020-06LSS
6	Right	11-020-06RSS-S	11-020-06RSS
8	Left	11-020-08LSS-S	11-020-08LSS
8	Right	11-020-08RSS-S	11-020-08RSS
10	Left	11-020-10LSS-S	11-020-10LSS
10	Right	11-020-10RSS-S	11-020-10RSS
12	Left	11-020-12LSS-S	11-020-12LSS
12	Right	11-020-12RSS-S	11-020-12RSS
14	Left	11-020-14LSS-S	11-020-14LSS
14	Right	11-020-14RSS-S	11-020-14RSS
16	Left	11-020-16LSS-S	11-020-16LSS
16	Right	11-020-16RSS-S	11-020-16RSS



2.7mm AV-Wiselock Lateral Distal Fibula Plate

Titanium

Holes	Direction	Sterile	Non-Sterile
3	Left	11-021-03LTI-S	11-021-03LTI
3	Right	11-021-03RTI-S	11-021-03RTI
4	Left	11-021-04LTI-S	11-021-04LTI
4	Right	11-021-04RTI-S	11-021-04RTI
5	Left	11-021-05LTI-S	11-021-05LTI
5	Right	11-021-05RTI-S	11-021-05RTI
6	Left	11-021-06LTI-S	11-021-06LTI
6	Right	11-021-06RTI-S	11-021-06RTI
7	Left	11-021-07LTI-S	11-021-07LTI
7	Right	11-021-07RTI-S	11-021-07RTI
9	Left	11-021-09LTI-S	11-021-09LTI
9	Right	11-021-09RTI-S	11-021-09RTI
11	Left	11-021-11LTI-S	11-021-11LTI
11	Right	11-021-11RTI-S	11-021-11RTI
13	Left	11-021-13LTI-S	11-021-13LTI
13	Right	11-021-13RTI-S	11-021-13RTI
15	Left	11-021-15LTI-S	11-021-15LTI
15	Right	11-021-15RTI-S	11-021-15RTI

Stainless Steel

Holes	Direction	Sterile	Non-Sterile
3	Left	11-021-03LSS-S	11-021-03LSS
3	Right	11-021-03RSS-S	11-021-03RSS
4	Left	11-021-04LSS-S	11-021-04LSS
4	Right	11-021-04RSS-S	11-021-04RSS
5	Left	11-021-05LSS-S	11-021-05LSS
5	Right	11-021-05RSS-S	11-021-05RSS
6	Left	11-021-06LSS-S	11-021-06LSS
6	Right	11-021-06RSS-S	11-021-06RSS
7	Left	11-021-07LSS-S	11-021-07LSS
7	Right	11-021-07RSS-S	11-021-07RSS
9	Left	11-021-09LSS-S	11-021-09LSS
9	Right	11-021-09RSS-S	11-021-09RSS
11	Left	11-021-11LSS-S	11-021-11LSS
11	Right	11-021-11RSS-S	11-021-11RSS
13	Left	11-021-13LSS-S	11-021-13LSS
13	Right	11-021-13RSS-S	11-021-13RSS
15	Left	11-021-15LSS-S	11-021-15LSS
15	Right	11-021-15RSS-S	11-021-15RSS


2.7mm AV-Wiselock Distal Tibia L-Plate



Titanium			
Holes	Direction	Sterile	Non-Sterile
4	Left	11-022-04LTI-S	11-022-04LTI
4	Right	11-022-04RTI-S	11-022-04RTI
6	Left	11-022-06LTI-S	11-022-06LTI
6	Right	11-022-06RTI-S	11-022-06RTI


Stainless Steel			
Holes	Direction	Sterile	Non-Sterile
4	Left	11-022-04LSS-S	11-022-04LSS
4	Right	11-022-04RSS-S	11-022-04RSS
6	Left	11-022-06LSS-S	11-022-06LSS
6	Right	11-022-06RSS-S	11-022-06RSS

2.7mm AV-Wiselock Distal Tibia T-Plate



Titanium		
Holes	Sterile	Non-Sterile
4	11-023-04TI-S	11-023-04TI
6	11-023-06TI-S	11-023-06TI

Stainless Steel		
Holes	Sterile	Non-Sterile
4	11-023-04SS-S	11-023-04SS
6	11-023-06SS-S	11-023-06SS


2.7/3.5mm AV-Wiselock Anterolateral Distal Tibia Plate


Titanium			
Holes	Direction	Sterile	Non-Sterile
4	Left	11-025-04LTI-S	11-025-04LTI
4	Right	11-025-04RTI-S	11-025-04RTI
6	Left	11-025-06LTI-S	11-025-06LTI
6	Right	11-025-06RTI-S	11-025-06RTI
8	Left	11-025-08LTI-S	11-025-08LTI
8	Right	11-025-08RTI-S	11-025-08RTI
10	Left	11-025-10LTI-S	11-025-10LTI
10	Right	11-025-10RTI-S	11-025-10RTI
12	Left	11-025-12LTI-S	11-025-12LTI
12	Right	11-025-12RTI-S	11-025-12RTI
14	Left	11-025-14LTI-S	11-025-14LTI
14	Right	11-025-14RTI-S	11-025-14RTI
16	Left	11-025-16LTI-S	11-025-16LTI
16	Right	11-025-16RTI-S	11-025-16RTI
18	Left	11-025-18LTI-S	11-025-18LTI
18	Right	11-025-18RTI-S	11-025-18RTI

Stainless Steel			
Holes	Direction	Sterile	Non-Sterile
4	Left	11-025-04LSS-S	11-025-04LSS
4	Right	11-025-04RSS-S	11-025-04RSS
6	Left	11-025-06LSS-S	11-025-06LSS
6	Right	11-025-06RSS-S	11-025-06RSS
8	Left	11-025-08LSS-S	11-025-08LSS
8	Right	11-025-08RSS-S	11-025-08RSS
10	Left	11-025-10LSS-S	11-025-10LSS
10	Right	11-025-10RSS-S	11-025-10RSS
12	Left	11-025-12LSS-S	11-025-12LSS
12	Right	11-025-12RSS-S	11-025-12RSS
14	Left	11-025-14LSS-S	11-025-14LSS
14	Right	11-025-14RSS-S	11-025-14RSS
16	Left	11-025-16LSS-S	11-025-16LSS
16	Right	11-025-16RSS-S	11-025-16RSS
18	Left	11-025-18LSS-S	11-025-18LSS
18	Right	11-025-18RSS-S	11-025-18RSS


2.7mm AV-Wiselock Screw, Self-Tapping, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø2.7	6mm	TI-1472-006-S	TI-1472-006	1472-006-S	1472-006
Ø2.7	8mm	TI-1472-008-S	TI-1472-008	1472-008-S	1472-008
Ø2.7	10mm	TI-1472-010-S	TI-1472-010	1472-010-S	1472-010
Ø2.7	12mm	TI-1472-012-S	TI-1472-012	1472-012-S	1472-012
Ø2.7	14mm	TI-1472-014-S	TI-1472-014	1472-014-S	1472-014
Ø2.7	16mm	TI-1472-016-S	TI-1472-016	1472-016-S	1472-016
Ø2.7	18mm	TI-1472-018-S	TI-1472-018	1472-018-S	1472-018
Ø2.7	20mm	TI-1472-020-S	TI-1472-020	1472-020-S	1472-020
Ø2.7	22mm	TI-1472-022-S	TI-1472-022	1472-022-S	1472-022
Ø2.7	24mm	TI-1472-024-S	TI-1472-024	1472-024-S	1472-024
Ø2.7	26mm	TI-1472-026-S	TI-1472-026	1472-026-S	1472-026
Ø2.7	28mm	TI-1472-028-S	TI-1472-028	1472-028-S	1472-028
Ø2.7	30mm	TI-1472-030-S	TI-1472-030	1472-030-S	1472-030
Ø2.7	32mm	11-017-032TI-S	11-017-032TI	11-017-032SS-S	11-017-032SS
Ø2.7	34mm	11-017-034TI-S	11-017-034TI	11-017-034SS-S	11-017-034SS
Ø2.7	36mm	11-017-036TI-S	11-017-036TI	11-017-036SS-S	11-017-036SS
Ø2.7	38mm	11-017-038TI-S	11-017-038TI	11-017-038SS-S	11-017-038SS
Ø2.7	40mm	11-017-040TI-S	11-017-040TI	11-017-040SS-S	11-017-040SS
Ø2.7	42mm	11-017-042TI-S	11-017-042TI	11-017-042SS-S	11-017-042SS
Ø2.7	44mm	11-017-044TI-S	11-017-044TI	11-017-044SS-S	11-017-044SS
Ø2.7	46mm	11-017-046TI-S	11-017-046TI	11-017-046SS-S	11-017-046SS
Ø2.7	48mm	11-017-048TI-S	11-017-048TI	11-017-048SS-S	11-017-048SS
Ø2.7	50mm	11-017-050TI-S	11-017-050TI	11-017-050SS-S	11-017-050SS
Ø2.7	52mm	11-017-052TI-S	11-017-052TI	11-017-052SS-S	11-017-052SS
Ø2.7	54mm	11-017-054TI-S	11-017-054TI	11-017-054SS-S	11-017-054SS
Ø2.7	56mm	11-017-056TI-S	11-017-056TI	11-017-056SS-S	11-017-056SS
Ø2.7	58mm	11-017-058TI-S	11-017-058TI	11-017-058SS-S	11-017-058SS
Ø2.7	60mm	11-017-060TI-S	11-017-060TI	11-017-060SS-S	11-017-060SS
Ø2.7	62mm	11-017-062TI-S	11-017-062TI	11-017-062SS-S	11-017-062SS
Ø2.7	64mm	11-017-064TI-S	11-017-064TI	11-017-064SS-S	11-017-064SS
Ø2.7	66mm	11-017-066TI-S	11-017-066TI	11-017-066SS-S	11-017-066SS
Ø2.7	68mm	11-017-068TI-S	11-017-068TI	11-017-068SS-S	11-017-068SS
Ø2.7	70mm	11-017-070TI-S	11-017-070TI	11-017-070SS-S	11-017-070SS



3.5mm AV-Wiselock Screw, Self-Tapping, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø3.5	10mm	11-024-010TI-S	11-024-010TI	11-024-010SS-S	11-024-010SS
Ø3.5	12mm	11-024-012TI-S	11-024-012TI	11-024-012SS-S	11-024-012SS
Ø3.5	14mm	11-024-014TI-S	11-024-014TI	11-024-014SS-S	11-024-014SS
Ø3.5	16mm	11-024-016TI-S	11-024-016TI	11-024-016SS-S	11-024-016SS
Ø3.5	18mm	11-024-018TI-S	11-024-018TI	11-024-018SS-S	11-024-018SS
Ø3.5	20mm	11-024-020TI-S	11-024-020TI	11-024-020SS-S	11-024-020SS
Ø3.5	22mm	11-024-022TI-S	11-024-022TI	11-024-022SS-S	11-024-022SS
Ø3.5	24mm	11-024-024TI-S	11-024-024TI	11-024-024SS-S	11-024-024SS
Ø3.5	26mm	11-024-026TI-S	11-024-026TI	11-024-026SS-S	11-024-026SS
Ø3.5	28mm	11-024-028TI-S	11-024-028TI	11-024-028SS-S	11-024-028SS
Ø3.5	30mm	11-024-030TI-S	11-024-030TI	11-024-030SS-S	11-024-030SS
Ø3.5	32mm	11-024-032TI-S	11-024-032TI	11-024-032SS-S	11-024-032SS
Ø3.5	34mm	11-024-034TI-S	11-024-034TI	11-024-034SS-S	11-024-034SS
Ø3.5	36mm	11-024-036TI-S	11-024-036TI	11-024-036SS-S	11-024-036SS
Ø3.5	38mm	11-024-038TI-S	11-024-038TI	11-024-038SS-S	11-024-038SS
Ø3.5	40mm	11-024-040TI-S	11-024-040TI	11-024-040SS-S	11-024-040SS
Ø3.5	42mm	11-024-042TI-S	11-024-042TI	11-024-042SS-S	11-024-042SS
Ø3.5	44mm	11-024-044TI-S	11-024-044TI	11-024-044SS-S	11-024-044SS
Ø3.5	46mm	11-024-046TI-S	11-024-046TI	11-024-046SS-S	11-024-046SS
Ø3.5	48mm	11-024-048TI-S	11-024-048TI	11-024-048SS-S	11-024-048SS
Ø3.5	50mm	11-024-050TI-S	11-024-050TI	11-024-050SS-S	11-024-050SS
Ø3.5	52mm	11-024-052TI-S	11-024-052TI	11-024-052SS-S	11-024-052SS
Ø3.5	54mm	11-024-054TI-S	11-024-054TI	11-024-054SS-S	11-024-054SS
Ø3.5	56mm	11-024-056TI-S	11-024-056TI	11-024-056SS-S	11-024-056SS
Ø3.5	58mm	11-024-058TI-S	11-024-058TI	11-024-058SS-S	11-024-058SS
Ø3.5	60mm	11-024-060TI-S	11-024-060TI	11-024-060SS-S	11-024-060SS
Ø3.5	65mm	11-024-065TI-S	11-024-065TI	11-024-065SS-S	11-024-065SS
Ø3.5	70mm	11-024-070TI-S	11-024-070TI	11-024-070SS-S	11-024-070SS
Ø3.5	75mm	11-024-075TI-S	11-024-075TI	11-024-075SS-S	11-024-075SS
Ø3.5	80mm	11-024-080TI-S	11-024-080TI	11-024-080SS-S	11-024-080SS
Ø3.5	85mm	11-024-085TI-S	11-024-085TI	11-024-085SS-S	11-024-085SS
Ø3.5	90mm	11-024-090TI-S	11-024-090TI	11-024-090SS-S	11-024-090SS
Ø3.5	95mm	11-024-095TI-S	11-024-095TI	11-024-095SS-S	11-024-095SS




3.5mm Cortical Screw, Self-Tapping, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø3.5	10mm	11-002-010TI-S	11-002-010TI	11-002-010SS-S	11-002-010SS
Ø3.5	12mm	11-002-012TI-S	11-002-012TI	11-002-012SS-S	11-002-012SS
Ø3.5	14mm	11-002-014TI-S	11-002-014TI	11-002-014SS-S	11-002-014SS
Ø3.5	16mm	11-002-016TI-S	11-002-016TI	11-002-016SS-S	11-002-016SS
Ø3.5	18mm	11-002-018TI-S	11-002-018TI	11-002-018SS-S	11-002-018SS
Ø3.5	20mm	11-002-020TI-S	11-002-020TI	11-002-020SS-S	11-002-020SS
Ø3.5	22mm	11-002-022TI-S	11-002-022TI	11-002-022SS-S	11-002-022SS
Ø3.5	24mm	11-002-024TI-S	11-002-024TI	11-002-024SS-S	11-002-024SS
Ø3.5	26mm	11-002-026TI-S	11-002-026TI	11-002-026SS-S	11-002-026SS
Ø3.5	28mm	11-002-028TI-S	11-002-028TI	11-002-028SS-S	11-002-028SS
Ø3.5	30mm	11-002-030TI-S	11-002-030TI	11-002-030SS-S	11-002-030SS
Ø3.5	32mm	11-002-032TI-S	11-002-032TI	11-002-032SS-S	11-002-032SS
Ø3.5	34mm	11-002-034TI-S	11-002-034TI	11-002-034SS-S	11-002-034SS
Ø3.5	36mm	11-002-036TI-S	11-002-036TI	11-002-036SS-S	11-002-036SS
Ø3.5	38mm	11-002-038TI-S	11-002-038TI	11-002-038SS-S	11-002-038SS
Ø3.5	40mm	11-002-040TI-S	11-002-040TI	11-002-040SS-S	11-002-040SS
Ø3.5	42mm	11-002-042TI-S	11-002-042TI	11-002-042SS-S	11-002-042SS
Ø3.5	44mm	11-002-044TI-S	11-002-044TI	11-002-044SS-S	11-002-044SS
Ø3.5	46mm	11-002-046TI-S	11-002-046TI	11-002-046SS-S	11-002-046SS
Ø3.5	48mm	11-002-048TI-S	11-002-048TI	11-002-048SS-S	11-002-048SS
Ø3.5	50mm	11-002-050TI-S	11-002-050TI	11-002-050SS-S	11-002-050SS
Ø3.5	55mm	11-002-055TI-S	11-002-055TI	11-002-055SS-S	11-002-055SS
Ø3.5	60mm	11-002-060TI-S	11-002-060TI	11-002-060SS-S	11-002-060SS
Ø3.5	65mm	11-002-065TI-S	11-002-065TI	11-002-065SS-S	11-002-065SS
Ø3.5	70mm	11-002-070TI-S	11-002-070TI	11-002-070SS-S	11-002-070SS
Ø3.5	75mm	11-002-075TI-S	11-002-075TI	11-002-075SS-S	11-002-075SS
Ø3.5	80mm	11-002-080TI-S	11-002-080TI	11-002-080SS-S	11-002-080SS
Ø3.5	85mm	11-002-085TI-S	11-002-085TI	11-002-085SS-S	11-002-085SS
Ø3.5	90mm	11-002-090TI-S	11-002-090TI	11-002-090SS-S	11-002-090SS
Ø3.5	95mm	11-002-095TI-S	11-002-095TI	11-002-095SS-S	11-002-095SS




2.7mm Cortical Screw, Self Tapping, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø2.7	6mm	TI-1242.06-S	TI-1242.06	1242.06-S	1242.06
Ø2.7	8mm	TI-1242.08-S	TI-1242.08	1242.08-S	1242.08
Ø2.7	10mm	TI-1242.10-S	TI-1242.10	1242.10-S	1242.10
Ø2.7	12mm	TI-1242.12-S	TI-1242.12	1242.12-S	1242.12
Ø2.7	14mm	TI-1242.14-S	TI-1242.14	1242.14-S	1242.14
Ø2.7	16mm	TI-1242.16-S	TI-1242.16	1242.16-S	1242.16
Ø2.7	18mm	TI-1242.18-S	TI-1242.18	1242.18-S	1242.18
Ø2.7	20mm	TI-1242.20-S	TI-1242.20	1242.20-S	1242.20
Ø2.7	22mm	TI-1242.22-S	TI-1242.22	1242.22-S	1242.22
Ø2.7	24mm	TI-1242.24-S	TI-1242.24	1242.24-S	1242.24
Ø2.7	26mm	TI-1242.26-S	TI-1242.26	1242.26-S	1242.26
Ø2.7	28mm	TI-1242.28-S	TI-1242.28	1242.28-S	1242.28
Ø2.7	30mm	TI-1242.30-S	TI-1242.30	1242.30-S	1242.30
Ø2.7	32mm	TI-1242.32-S	TI-1242.32	1242.32-S	1242.32
Ø2.7	34mm	TI-1242.34-S	TI-1242.34	1242.34-S	1242.34
Ø2.7	36mm	TI-1242.36-S	TI-1242.36	1242.36-S	1242.36
Ø2.7	38mm	TI-1242.38-S	TI-1242.38	1242.38-S	1242.38
Ø2.7	40mm	TI-1242.40-S	TI-1242.40	1242.40-S	1242.40
Ø2.7	42mm	TI-1242.42-S	TI-1242.42	1242.42-S	1242.42
Ø2.7	44mm	TI-1242.44-S	TI-1242.44	1242.44-S	1242.44
Ø2.7	46mm	TI-1242.46-S	TI-1242.46	1242.46-S	1242.46
Ø2.7	48mm	TI-1242.48-S	TI-1242.48	1242.48-S	1242.48
Ø2.7	50mm	TI-1242.50-S	TI-1242.50	1242.50-S	1242.50
Ø2.7	52mm	TI-1242.52-S	TI-1242.52	1242.52-S	1242.52
Ø2.7	54mm	TI-1242.54-S	TI-1242.54	1242.54-S	1242.54
Ø2.7	56mm	TI-1242.56-S	TI-1242.56	1242.56-S	1242.56
Ø2.7	58mm	TI-1242.58-S	TI-1242.58	1242.58-S	1242.58
Ø2.7	60mm	TI-1242.60-S	TI-1242.60	1242.60-S	1242.60
Ø2.7	62mm	TI-1242.62-S	TI-1242.62	1242.62-S	1242.62
Ø2.7	64mm	TI-1242.64-S	TI-1242.64	1242.64-S	1242.64
Ø2.7	66mm	TI-1242.66-S	TI-1242.66	1242.66-S	1242.66
Ø2.7	68mm	TI-1242.68-S	TI-1242.68	1242.68-S	1242.68
Ø2.7	70mm	TI-1242.70-S	TI-1242.70	1242.70-S	1242.70




2.7mm Wise-Lock Screw, Self Tapping, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø2.7	6mm	TI-1241.06-S	TI-1241.06	11-012-006SS-S	11-012-006SS
Ø2.7	7mm	TI-1241.07-S	TI-1241.07	11-012-007SS-S	11-012-007SS
Ø2.7	8mm	TI-1241.08-S	TI-1241.08	11-012-008SS-S	11-012-008SS
Ø2.7	9mm	TI-1241.09-S	TI-1241.09	11-012-009SS-S	11-012-009SS
Ø2.7	10mm	TI-1241.10-S	TI-1241.10	11-012-010SS-S	11-012-010SS
Ø2.7	12mm	TI-1241.12-S	TI-1241.12	11-012-012SS-S	11-012-012SS
Ø2.7	14mm	TI-1241.14-S	TI-1241.14	11-012-014SS-S	11-012-014SS
Ø2.7	16mm	TI-1241.16-S	TI-1241.16	11-012-016SS-S	11-012-016SS
Ø2.7	18mm	TI-1241.18-S	TI-1241.18	11-012-018SS-S	11-012-018SS
Ø2.7	20mm	TI-1241.20-S	TI-1241.20	11-012-020SS-S	11-012-020SS
Ø2.7	22mm	TI-1241.22-S	TI-1241.22	11-012-022SS-S	11-012-022SS
Ø2.7	24mm	TI-1241.24-S	TI-1241.24	11-012-024SS-S	11-012-024SS
Ø2.7	26mm	TI-1241.26-S	TI-1241.26	11-012-026SS-S	11-012-026SS
Ø2.7	28mm	TI-1241.28-S	TI-1241.28	11-012-028SS-S	11-012-028SS
Ø2.7	30mm	TI-1241.30-S	TI-1241.30	11-012-030SS-S	11-012-030SS
Ø2.7	32mm	TI-1241.32-S	TI-1241.32	11-012-032SS-S	11-012-032SS
Ø2.7	34mm	TI-1241.34-S	TI-1241.34	11-012-034SS-S	11-012-034SS
Ø2.7	36mm	TI-1241.36-S	TI-1241.36	11-012-036SS-S	11-012-036SS
Ø2.7	38mm	TI-1241.38-S	TI-1241.38	11-012-038SS-S	11-012-038SS
Ø2.7	40mm	TI-1241.40-S	TI-1241.40	11-012-040SS-S	11-012-040SS
Ø2.7	42mm	TI-1241.42-S	TI-1241.42	11-012-042SS-S	11-012-042SS
Ø2.7	44mm	TI-1241.44-S	TI-1241.44	11-012-044SS-S	11-012-044SS
Ø2.7	46mm	TI-1241.46-S	TI-1241.46	11-012-046SS-S	11-012-046SS
Ø2.7	48mm	TI-1241.48-S	TI-1241.48	11-012-048SS-S	11-012-048SS
Ø2.7	50mm	TI-1241.50-S	TI-1241.50	11-012-050SS-S	11-012-050SS
Ø2.7	52mm	TI-1241.52-S	TI-1241.52	11-012-052SS-S	11-012-052SS
Ø2.7	54mm	TI-1241.54-S	TI-1241.54	11-012-054SS-S	11-012-054SS
Ø2.7	56mm	TI-1241.56-S	TI-1241.56	11-012-056SS-S	11-012-056SS
Ø2.7	58mm	TI-1241.58-S	TI-1241.58	11-012-058SS-S	11-012-058SS
Ø2.7	60mm	TI-1241.60-S	TI-1241.60	11-012-060SS-S	11-012-060SS
Ø2.7	62mm	11-012-062TI-S	11-012-062TI	11-012-062SS-S	11-012-062SS
Ø2.7	64mm	11-012-064TI-S	11-012-064TI	11-012-064SS-S	11-012-064SS
Ø2.7	66mm	11-012-066TI-S	11-012-066TI	11-012-066SS-S	11-012-066SS
Ø2.7	68mm	11-012-068TI-S	11-012-068TI	11-012-068SS-S	11-012-068SS
Ø2.7	70mm	11-012-070TI-S	11-012-070TI	11-012-070SS-S	11-012-070SS



3.5mm Wise-Lock Screw, Self Tapping, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø3.5	10mm	TI-1473-010-S	TI-1473-010	1473-010-S	1473-010
Ø3.5	12mm	TI-1473-012-S	TI-1473-012	1473-012-S	1473-012
Ø3.5	14mm	TI-1473-014-S	TI-1473-014	1473-014-S	1473-014
Ø3.5	16mm	TI-1473-016-S	TI-1473-016	1473-016-S	1473-016
Ø3.5	18mm	TI-1473-018-S	TI-1473-018	1473-018-S	1473-018
Ø3.5	20mm	TI-1473-020-S	TI-1473-020	1473-020-S	1473-020
Ø3.5	22mm	TI-1473-022-S	TI-1473-022	1473-022-S	1473-022
Ø3.5	24mm	TI-1473-024-S	TI-1473-024	1473-024-S	1473-024
Ø3.5	26mm	TI-1473-026-S	TI-1473-026	1473-026-S	1473-026
Ø3.5	28mm	TI-1473-028-S	TI-1473-028	1473-028-S	1473-028
Ø3.5	30mm	TI-1473-030-S	TI-1473-030	1473-030-S	1473-030
Ø3.5	32mm	TI-1473-032-S	TI-1473-032	1473-032-S	1473-032
Ø3.5	34mm	TI-1473-034-S	TI-1473-034	1473-034-S	1473-034
Ø3.5	36mm	TI-1473-036-S	TI-1473-036	1473-036-S	1473-036
Ø3.5	38mm	TI-1473-038-S	TI-1473-038	1473-038-S	1473-038
Ø3.5	40mm	TI-1473-040-S	TI-1473-040	1473-040-S	1473-040
Ø3.5	42mm	TI-1473-042-S	TI-1473-042	1473-042-S	1473-042
Ø3.5	44mm	TI-1473-044-S	TI-1473-044	1473-044-S	1473-044
Ø3.5	46mm	TI-1473-046-S	TI-1473-046	1473-046-S	1473-046
Ø3.5	48mm	TI-1473-048-S	TI-1473-048	1473-048-S	1473-048
Ø3.5	50mm	TI-1473-050-S	TI-1473-050	1473-050-S	1473-050
Ø3.5	52mm	TI-1473-052-S	TI-1473-052	1473-052-S	1473-052
Ø3.5	54mm	TI-1473-054-S	TI-1473-054	1473-054-S	1473-054
Ø3.5	56mm	TI-1473-056-S	TI-1473-056	1473-056-S	1473-056
Ø3.5	58mm	TI-1473-058-S	TI-1473-058	1473-058-S	1473-058
Ø3.5	60mm	TI-1473-060-S	TI-1473-060	1473-060-S	1473-060
Ø3.5	65mm	TI-1473-065-S	TI-1473-065	1473-065-S	1473-065
Ø3.5	70mm	TI-1473-070-S	TI-1473-070	1473-070-S	1473-070
Ø3.5	75mm	TI-1473-075-S	TI-1473-075	1473-075-S	1473-075
Ø3.5	80mm	TI-1473-080-S	TI-1473-080	1473-080-S	1473-080
Ø3.5	85mm	TI-1473-085-S	TI-1473-085	1473-085-S	1473-085
Ø3.5	90mm	TI-1473-090-S	TI-1473-090	1473-090-S	1473-090



AV-WISELOCK

2.7/3.5mm Elbow System-Locking

Designed to reduce the risk of soft tissue irritation

by taking into account the poor soft tissue coverage of the elbow, thereby improving the plate profile with an anatomical plate contour and enabling the screw heads to sit flush in the plate hole.

Offers simplicity and adaptability:



AV-Wiselock technology combines the simplicity of the established Wise-Lock screw insertion technique with the possibility of adapting screw angulations. Adaptability is enhanced during screw insertion by offering predefined screw angles while allowing $\pm 15^\circ$ off-axis angulation, if needed.

The Auxein's AV-Wiselock 2.7/3.5mm Elbow Plates is intended for fixation of fractures of the distal humerus, olecranon and ulna in adults and adolescents (12-21) in which the growth plates have fused.

Specifically:



Distal humerus plates are indicated for intra-articular fractures, comminuted supracondylar fractures, osteotomies, malunions and non-unions of the distal humerus.

Olecranon and Proximal ulna plates are indicated for fractures, osteotomies, malunions and non-unions of the olecranon and proximal ulna.



Plates Available in

Titanium as per ISO 5832-3 and Stainless Steel as per ISO 5832-1

2.7/3.5mm AV-Wiselock Proximal Olecranon Plate


Titanium

Holes	Direction	Sterile	Non-Sterile
2	Left	11-003-02TIL-S	11-003-02TIL
2	Right	11-003-02TIR-S	11-003-02TIR

Stainless Steel

Holes	Direction	Sterile	Non-Sterile
2	Left	11-003-02SSL-S	11-003-02SSL
2	Right	11-003-02SSR-S	11-003-02SSR


2.7/3.5mm AV-Wiselock Olecranon Plate


Titanium

Holes	Direction	Sterile	Non-Sterile
2	Left	11-004-02TIL-S	11-004-02TIL
2	Right	11-004-02TIR-S	11-004-02TIR
4	Left	11-004-04TIL-S	11-004-04TIL
4	Right	11-004-04TIR-S	11-004-04TIR
6	Left	11-004-06TIL-S	11-004-06TIL
6	Right	11-004-06TIR-S	11-004-06TIR
8	Left	11-004-08TIL-S	11-004-08TIL
8	Right	11-004-08TIR-S	11-004-08TIR

Stainless Steel

Holes	Direction	Sterile	Non-Sterile
2	Left	11-004-02SSL-S	11-004-02SSL
2	Right	11-004-02SSR-S	11-004-02SSR
4	Left	11-004-04SSL-S	11-004-04SSL
4	Right	11-004-04SSR-S	11-004-04SSR
6	Left	11-004-06SSL-S	11-004-06SSL
6	Right	11-004-06SSR-S	11-004-06SSR
8	Left	11-004-08SSL-S	11-004-08SSL
8	Right	11-004-08SSR-S	11-004-08SSR


2.7/3.5mm AV-Wiselock Extra Articular Proximal Ulna Plate


Titanium

Holes	Direction	Sterile	Non-Sterile
6	Left	11-005-06TIL-S	11-005-06TIL
6	Right	11-005-06TIR-S	11-005-06TIR
8	Left	11-005-08TIL-S	11-005-08TIL
8	Right	11-005-08TIR-S	11-005-08TIR
10	Left	11-005-10TIL-S	11-005-10TIL
10	Right	11-005-10TIR-S	11-005-10TIR
12	Left	11-005-12TIL-S	11-005-12TIL
12	Right	11-005-12TIR-S	11-005-12TIR

Stainless Steel

Holes	Direction	Sterile	Non-Sterile
6	Left	11-005-06SSL-S	11-005-06SSL
6	Right	11-005-06SSR-S	11-005-06SSR
8	Left	11-005-08SSL-S	11-005-08SSL
8	Right	11-005-08SSR-S	11-005-08SSR
10	Left	11-005-10SSL-S	11-005-10SSL
10	Right	11-005-10SSR-S	11-005-10SSR
12	Left	11-005-12SSL-S	11-005-12SSL
12	Right	11-005-12SSR-S	11-005-12SSR


2.7/3.5mm AV-Wiselock Dorsolateral Distal Humerus Plate with Lateral Support


Titanium

Holes	Direction	Sterile	Non-Sterile
3	Left	11-006-03TIL-S	11-006-03TIL
3	Right	11-006-03TIR-S	11-006-03TIR
4	Left	11-006-04TIL-S	11-006-04TIL
4	Right	11-006-04TIR-S	11-006-04TIR
7	Left	11-006-07TIL-S	11-006-07TIL
7	Right	11-006-07TIR-S	11-006-07TIR
9	Left	11-006-09TIL-S	11-006-09TIL
9	Right	11-006-09TIR-S	11-006-09TIR
11	Left	11-006-11TIL-S	11-006-11TIL
11	Right	11-006-11TIR-S	11-006-11TIR
13	Left	11-006-13TIL-S	11-006-13TIL
13	Right	11-006-13TIL-S	11-006-13TIL

Stainless Steel

Holes	Direction	Sterile	Non-Sterile
3	Left	11-006-03SSL-S	11-006-03SSL
3	Right	11-006-03SSR-S	11-006-03SSR
4	Left	11-006-04SSL-S	11-006-04SSL
4	Right	11-006-04SSR-S	11-006-04SSR
7	Left	11-006-07SSL-S	11-006-07SSL
7	Right	11-006-07SSR-S	11-006-07SSR
9	Left	11-006-09SSL-S	11-006-09SSL
9	Right	11-006-09SSR-S	11-006-09SSR
11	Left	11-006-11SSL-S	11-006-11SSL
11	Right	11-006-11SSR-S	11-006-11SSR
13	Left	11-006-13SSL-S	11-006-13SSL
13	Right	11-006-13SSL-S	11-006-13SSL


2.7/3.5mm AV-Wiselock Dorsolateral Distal Humerus Plate


Titanium

Holes	Direction	Sterile	Non-Sterile
3	Left	11-007-03TIL-S	11-007-03TIL
3	Right	11-007-03TIR-S	11-007-03TIR
4	Left	11-007-04TIL-S	11-007-04TIL
4	Right	11-007-04TIR-S	11-007-04TIR
7	Left	11-007-07TIL-S	11-007-07TIL
7	Right	11-007-07TIR-S	11-007-07TIR
9	Left	11-007-09TIL-S	11-007-09TIL
9	Right	11-007-09TIR-S	11-007-09TIR
11	Left	11-007-11TIL-S	11-007-11TIL
11	Right	11-007-11TIR-S	11-007-11TIR
13	Left	11-007-13TIL-S	11-007-13TIL
13	Right	11-007-13TIR-S	11-007-13TIR

Stainless Steel

Holes	Direction	Sterile	Non-Sterile
3	Left	11-007-03SSL-S	11-007-03SSL
3	Right	11-007-03SSR-S	11-007-03SSR
4	Left	11-007-04SSL-S	11-007-04SSL
4	Right	11-007-04SSR-S	11-007-04SSR
7	Left	11-007-07SSL-S	11-007-07SSL
7	Right	11-007-07SSR-S	11-007-07SSR
9	Left	11-007-09SSL-S	11-007-09SSL
9	Right	11-007-09SSR-S	11-007-09SSR
11	Left	11-007-11SSL-S	11-007-11SSL
11	Right	11-007-11SSR-S	11-007-11SSR
13	Left	11-007-13SSL-S	11-007-13SSL
13	Right	11-007-13SSR-S	11-007-13SSR


2.7/3.5mm AV-Wiselock Medial Distal Humerus Plate


Titanium

Holes	Direction	Sterile	Non-Sterile
1	Left	11-008-01TIL-S	11-008-01TIL
1	Right	11-008-01TIR-S	11-008-01TIR
2	Left	11-008-02TIL-S	11-008-02TIL
2	Right	11-008-02TIR-S	11-008-02TIR
4	Left	11-008-04TIL-S	11-008-04TIL
4	Right	11-008-04TIR-S	11-008-04TIR
6	Left	11-008-06TIL-S	11-008-06TIL
6	Right	11-008-06TIR-S	11-008-06TIR
8	Left	11-008-08TIL-S	11-008-08TIL
8	Right	11-008-08TIR-S	11-008-08TIR
10	Left	11-008-10TIL-S	11-008-10TIL
10	Right	11-008-10TIR-S	11-008-10TIR


Stainless Steel

Holes	Direction	Sterile	Non-Sterile
1	Left	11-008-01SSL-S	11-008-01SSL
1	Right	11-008-01SSR-S	11-008-01SSR
2	Left	11-008-02SSL-S	11-008-02SSL
2	Right	11-008-02SSR-S	11-008-02SSR
4	Left	11-008-04SSL-S	11-008-04SSL
4	Right	11-008-04SSR-S	11-008-04SSR
6	Left	11-008-06SSL-S	11-008-06SSL
6	Right	11-008-06SSR-S	11-008-06SSR
8	Left	11-008-08SSL-S	11-008-08SSL
8	Right	11-008-08SSR-S	11-008-08SSR
10	Left	11-008-10SSL-S	11-008-10SSL
10	Right	11-008-10SSR-S	11-008-10SSR

2.7/3.5mm AV-Wiselock Extended Medial Distal Humerus Plate


Titanium			
Holes	Direction	Sterile	Non-Sterile
1	Left	11-009-01TIL-S	11-009-01TIL
1	Right	11-009-01TIR-S	11-009-01TIR
2	Left	11-009-02TIL-S	11-009-02TIL
2	Right	11-009-02TIR-S	11-009-02TIR
4	Left	11-009-04TIL-S	11-009-04TIL
4	Right	11-009-04TIR-S	11-009-04TIR
6	Left	11-009-06TIL-S	11-009-06TIL
6	Right	11-009-06TIR-S	11-009-06TIR
8	Left	11-009-08TIL-S	11-009-08TIL
8	Right	11-009-08TIR-S	11-009-08TIR
10	Left	11-009-10TIL-S	11-009-10TIL
10	Right	11-009-10TIR-S	11-009-10TIR

Stainless Steel			
Holes	Direction	Sterile	Non-Sterile
1	Left	11-009-01SSL-S	11-009-01SSL
1	Right	11-009-01SSR-S	11-009-01SSR
2	Left	11-009-02SSL-S	11-009-02SSL
2	Right	11-009-02SSR-S	11-009-02SSR
4	Left	11-009-04SSL-S	11-009-04SSL
4	Right	11-009-04SSR-S	11-009-04SSR
6	Left	11-009-06SSL-S	11-009-06SSL
6	Right	11-009-06SSR-S	11-009-06SSR
8	Left	11-009-08SSL-S	11-009-08SSL
8	Right	11-009-08SSR-S	11-009-08SSR
10	Left	11-009-10SSL-S	11-009-10SSL
10	Right	11-009-10SSR-S	11-009-10SSR

2.7/3.5mm AV-Wiselock Lateral Distal Humerus Plate


Titanium			
Holes	Direction	Sterile	Non-Sterile
1	Left	11-010-01TIL-S	11-010-01TIL
1	Right	11-010-01TIR-S	11-010-01TIR
2	Left	11-010-02TIL-S	11-010-02TIL
2	Right	11-010-02TIR-S	11-010-02TIR
5	Left	11-010-05TIL-S	11-010-05TIL
5	Right	11-010-05TIR-S	11-010-05TIR
7	Left	11-010-07TIL-S	11-010-07TIL
7	Right	11-010-07TIR-S	11-010-07TIR
9	Left	11-010-09TIL-S	11-010-09TIL
9	Right	11-010-09TIR-S	11-010-09TIR
11	Left	11-010-11TIL-S	11-010-11TIL
11	Right	11-010-11TIR-S	11-010-11TIR


Stainless Steel			
Holes	Direction	Sterile	Non-Sterile
1	Left	11-010-01SSL-S	11-010-01SSL
1	Right	11-010-01SSR-S	11-010-01SSR
2	Left	11-010-02SSL-S	11-010-02SSL
2	Right	11-010-02SSR-S	11-010-02SSR
5	Left	11-010-05SSL-S	11-010-05SSL
5	Right	11-010-05SSR-S	11-010-05SSR
7	Left	11-010-07SSL-S	11-010-07SSL
7	Right	11-010-07SSR-S	11-010-07SSR
9	Left	11-010-09SSL-S	11-010-09SSL
9	Right	11-010-09SSR-S	11-010-09SSR
11	Left	11-010-11SSL-S	11-010-11SSL
11	Right	11-010-11SSR-S	11-010-11SSR

AV-WISELOCK

2.7mm Elbow System


All elbow plates have been developed through Anatomical Comprehensive analyses with innovative approaches

- Anatomic plate designs support fracture reduction
- Plate design derived from comprehensive analyses of human Anatomical
- Low profile plates to reduce soft tissue irritation
- Rounded plates with no sharp edges.
- User-friendly locking technology for stable treatment also of complex fractures with poor bone quality
- AV-Wiselock Multidirectional locking technology
- AV-Wiselock Variable angle of $\pm 15^\circ$ in all directions in Locking screw hole
- All plate holes (apart from oblong holes and compression holes) are compatible with AV-Wiselock or cortical screws
- Plates can be contoured to the individual anatomy of the patient
- K-wire holes for temporary fixation of the plate to the bone




Plates Available in

Titanium as per ISO 5832-3 and Stainless Steel as per ISO 5832-1

2.7mm AV-Wiselock Olecranon Tension M Plate



Titanium			
Holes		Sterile	Non-Sterile
6		11-056-01TI-S	11-056-01TI

Stainless Steel			
Holes		Sterile	Non-Sterile
6		11-056-01SS-S	11-056-01SS

2.7mm AV-Wiselock Medial Distal Humerus M Plate


Titanium				
Long	Holes	Direction	Sterile	Non-Sterile
	14	Left	11-056-06TIL-S	11-056-06TIL
	14	Right	11-056-06TIR-S	11-056-06TIR
Medium	Holes	Direction	Sterile	Non-Sterile
	10	Left	11-056-07TIL-S	11-056-07TIL
	10	Right	11-056-07TIR-S	11-056-07TIR
Short	Holes	Direction	Sterile	Non-Sterile
	9	Left	11-056-08TIL-S	11-056-08TIL
	9	Right	11-056-08TIR-S	11-056-08TIR

Stainless Steel				
Long	Holes	Direction	Sterile	Non-Sterile
	14	Left	11-056-06SSL-S	11-056-06SSL
	14	Right	11-056-06SSR-S	11-056-06SSR
Medium	Holes	Direction	Sterile	Non-Sterile
	10	Left	11-056-07SSL-S	11-056-07SSL
	10	Right	11-056-07SSR-S	11-056-07SSR
Short	Holes	Direction	Sterile	Non-Sterile
	9	Left	11-056-08SSL-S	11-056-08SSL
	9	Right	11-056-08SSR-S	11-056-08SSR

2.7mm AV-Wiselock Lateral Distal Humerus M Plate

Titanium
Stainless Steel
Long

Holes	Direction	Sterile	Non-Sterile
15	Left	11-056-09TIL-S	11-056-09TIL
15	Right	11-056-09TIR-S	11-056-09TIR

Holes	Direction	Sterile	Non-Sterile
15	Left	11-056-09SSL-S	11-056-09SSL
15	Right	11-056-09SSR-S	11-056-09SSR

Medium

Holes	Direction	Sterile	Non-Sterile
10	Left	11-056-10TIL-S	11-056-10TIL
10	Right	11-056-10TIR-S	11-056-10TIR

Holes	Direction	Sterile	Non-Sterile
10	Left	11-056-10SSL-S	11-056-10SSL
10	Right	11-056-10SSR-S	11-056-10SSR

Short

Holes	Direction	Sterile	Non-Sterile
9.	Left	11-056-11TIL-S	11-056-11TIL
9	Right	11-056-11TIR-S	11-056-11TIR

Holes	Direction	Sterile	Non-Sterile
9.	Left	11-056-11SSL-S	11-056-11SSL
9	Right	11-056-11SSR-S	11-056-11SSR

Washer For Small Fragment Screw
Titanium

Sterile	Non-Sterile
100.110TI-S	100.110TI

Stainless Steel

Sterile	Non-Sterile
100.110-S	100.110




2.4mm Cortical Screws, Self Tapping, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø2.4	6mm	TI-1144.06-S	TI-1144.06	1144.06-S	1144.06
Ø2.4	8mm	TI-1144.08-S	TI-1144.08	1144.08-S	1144.08
Ø2.4	10mm	TI-1144.10-S	TI-1144.10	1144.10-S	1144.10
Ø2.4	12mm	TI-1144.12-S	TI-1144.12	1144.12-S	1144.12
Ø2.4	14mm	TI-1144.14-S	TI-1144.14	1144.14-S	1144.14
Ø2.4	16mm	TI-1144.16-S	TI-1144.16	1144.16-S	1144.16
Ø2.4	18mm	TI-1144.18-S	TI-1144.18	1144.18-S	1144.18
Ø2.4	20mm	TI-1144.20-S	TI-1144.20	1144.20-S	1144.20
Ø2.4	22mm	TI-1144.22-S	TI-1144.22	1144.22-S	1144.22
Ø2.4	24mm	TI-1144.24-S	TI-1144.24	1144.24-S	1144.24
Ø2.4	26mm	TI-1144.26-S	TI-1144.26	1144.26-S	1144.26
Ø2.4	28mm	TI-1144.28-S	TI-1144.28	1144.28-S	1144.28
Ø2.4	30mm	TI-1144.30-S	TI-1144.30	1144.30-S	1144.30
Ø2.4	32mm	TI-1144.32-S	TI-1144.32	1144.32-S	1144.32
Ø2.4	34mm	TI-1144.34-S	TI-1144.34	1144.34-S	1144.34
Ø2.4	36mm	TI-1144.36-S	TI-1144.36	1144.36-S	1144.36
Ø2.4	38mm	TI-1144.38-S	TI-1144.38	1144.38-S	1144.38
Ø2.4	40mm	TI-1144.40-S	TI-1144.40	1144.40-S	1144.40




2.7mm AV-Wiselock Screw, Self-Tapping, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø2.7	6mm	TI-1472-006-S	TI-1472-006	1472-006-S	1472-006
Ø2.7	8mm	TI-1472-008-S	TI-1472-008	1472-008-S	1472-008
Ø2.7	10mm	TI-1472-010-S	TI-1472-010	1472-010-S	1472-010
Ø2.7	12mm	TI-1472-012-S	TI-1472-012	1472-012-S	1472-012
Ø2.7	14mm	TI-1472-014-S	TI-1472-014	1472-014-S	1472-014
Ø2.7	16mm	TI-1472-016-S	TI-1472-016	1472-016-S	1472-016
Ø2.7	18mm	TI-1472-018-S	TI-1472-018	1472-018-S	1472-018
Ø2.7	20mm	TI-1472-020-S	TI-1472-020	1472-020-S	1472-020
Ø2.7	22mm	TI-1472-022-S	TI-1472-022	1472-022-S	1472-022
Ø2.7	24mm	TI-1472-024-S	TI-1472-024	1472-024-S	1472-024
Ø2.7	26mm	TI-1472-026-S	TI-1472-026	1472-026-S	1472-026
Ø2.7	28mm	TI-1472-028-S	TI-1472-028	1472-028-S	1472-028
Ø2.7	30mm	TI-1472-030-S	TI-1472-030	1472-030-S	1472-030
Ø2.7	32mm	11-017-032TI-S	11-017-032TI	11-017-032SS-S	11-017-032SS
Ø2.7	34mm	11-017-034TI-S	11-017-034TI	11-017-034SS-S	11-017-034SS
Ø2.7	36mm	11-017-036TI-S	11-017-036TI	11-017-036SS-S	11-017-036SS
Ø2.7	38mm	11-017-038TI-S	11-017-038TI	11-017-038SS-S	11-017-038SS
Ø2.7	40mm	11-017-040TI-S	11-017-040TI	11-017-040SS-S	11-017-040SS
Ø2.7	42mm	11-017-042TI-S	11-017-042TI	11-017-042SS-S	11-017-042SS
Ø2.7	44mm	11-017-044TI-S	11-017-044TI	11-017-044SS-S	11-017-044SS
Ø2.7	46mm	11-017-046TI-S	11-017-046TI	11-017-046SS-S	11-017-046SS
Ø2.7	48mm	11-017-048TI-S	11-017-048TI	11-017-048SS-S	11-017-048SS
Ø2.7	50mm	11-017-050TI-S	11-017-050TI	11-017-050SS-S	11-017-050SS
Ø2.7	52mm	11-017-052TI-S	11-017-052TI	11-017-052SS-S	11-017-052SS
Ø2.7	54mm	11-017-054TI-S	11-017-054TI	11-017-054SS-S	11-017-054SS
Ø2.7	56mm	11-017-056TI-S	11-017-056TI	11-017-056SS-S	11-017-056SS
Ø2.7	58mm	11-017-058TI-S	11-017-058TI	11-017-058SS-S	11-017-058SS
Ø2.7	60mm	11-017-060TI-S	11-017-060TI	11-017-060SS-S	11-017-060SS
Ø2.7	62mm	11-017-062TI-S	11-017-062TI	11-017-062SS-S	11-017-062SS
Ø2.7	64mm	11-017-064TI-S	11-017-064TI	11-017-064SS-S	11-017-064SS
Ø2.7	66mm	11-017-066TI-S	11-017-066TI	11-017-066SS-S	11-017-066SS
Ø2.7	68mm	11-017-068TI-S	11-017-068TI	11-017-068SS-S	11-017-068SS
Ø2.7	70mm	11-017-070TI-S	11-017-070TI	11-017-070SS-S	11-017-070SS




3.5mm Wise-Lock Screw, Self Tapping, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø3.5	10mm	TI-1473-010-S	TI-1473-010	1473-010-S	1473-010
Ø3.5	12mm	TI-1473-012-S	TI-1473-012	1473-012-S	1473-012
Ø3.5	14mm	TI-1473-014-S	TI-1473-014	1473-014-S	1473-014
Ø3.5	16mm	TI-1473-016-S	TI-1473-016	1473-016-S	1473-016
Ø3.5	18mm	TI-1473-018-S	TI-1473-018	1473-018-S	1473-018
Ø3.5	20mm	TI-1473-020-S	TI-1473-020	1473-020-S	1473-020
Ø3.5	22mm	TI-1473-022-S	TI-1473-022	1473-022-S	1473-022
Ø3.5	24mm	TI-1473-024-S	TI-1473-024	1473-024-S	1473-024
Ø3.5	26mm	TI-1473-026-S	TI-1473-026	1473-026-S	1473-026
Ø3.5	28mm	TI-1473-028-S	TI-1473-028	1473-028-S	1473-028
Ø3.5	30mm	TI-1473-030-S	TI-1473-030	1473-030-S	1473-030
Ø3.5	32mm	TI-1473-032-S	TI-1473-032	1473-032-S	1473-032
Ø3.5	34mm	TI-1473-034-S	TI-1473-034	1473-034-S	1473-034
Ø3.5	36mm	TI-1473-036-S	TI-1473-036	1473-036-S	1473-036
Ø3.5	38mm	TI-1473-038-S	TI-1473-038	1473-038-S	1473-038
Ø3.5	40mm	TI-1473-040-S	TI-1473-040	1473-040-S	1473-040
Ø3.5	42mm	TI-1473-042-S	TI-1473-042	1473-042-S	1473-042
Ø3.5	44mm	TI-1473-044-S	TI-1473-044	1473-044-S	1473-044
Ø3.5	46mm	TI-1473-046-S	TI-1473-046	1473-046-S	1473-046
Ø3.5	48mm	TI-1473-048-S	TI-1473-048	1473-048-S	1473-048
Ø3.5	50mm	TI-1473-050-S	TI-1473-050	1473-050-S	1473-050
Ø3.5	52mm	TI-1473-052-S	TI-1473-052	1473-052-S	1473-052
Ø3.5	54mm	TI-1473-054-S	TI-1473-054	1473-054-S	1473-054
Ø3.5	56mm	TI-1473-056-S	TI-1473-056	1473-056-S	1473-056
Ø3.5	58mm	TI-1473-058-S	TI-1473-058	1473-058-S	1473-058
Ø3.5	60mm	TI-1473-060-S	TI-1473-060	1473-060-S	1473-060
Ø3.5	65mm	TI-1473-065-S	TI-1473-065	1473-065-S	1473-065
Ø3.5	70mm	TI-1473-070-S	TI-1473-070	1473-070-S	1473-070
Ø3.5	75mm	TI-1473-075-S	TI-1473-075	1473-075-S	1473-075
Ø3.5	80mm	TI-1473-080-S	TI-1473-080	1473-080-S	1473-080
Ø3.5	85mm	TI-1473-085-S	TI-1473-085	1473-085-S	1473-085
Ø3.5	90mm	TI-1473-090-S	TI-1473-090	1473-090-S	1473-090




3.5mm Cortical Screw, Self Tapping (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø3.5	10mm	11-002-010TI-S	11-002-010TI	11-002-010SS-S	11-002-010SS
Ø3.5	12mm	11-002-012TI-S	11-002-012TI	11-002-012SS-S	11-002-012SS
Ø3.5	14mm	11-002-014TI-S	11-002-014TI	11-002-014SS-S	11-002-014SS
Ø3.5	16mm	11-002-016TI-S	11-002-016TI	11-002-016SS-S	11-002-016SS
Ø3.5	18mm	11-002-018TI-S	11-002-018TI	11-002-018SS-S	11-002-018SS
Ø3.5	20mm	11-002-020TI-S	11-002-020TI	11-002-020SS-S	11-002-020SS
Ø3.5	22mm	11-002-022TI-S	11-002-022TI	11-002-022SS-S	11-002-022SS
Ø3.5	24mm	11-002-024TI-S	11-002-024TI	11-002-024SS-S	11-002-024SS
Ø3.5	26mm	11-002-026TI-S	11-002-026TI	11-002-026SS-S	11-002-026SS
Ø3.5	28mm	11-002-028TI-S	11-002-028TI	11-002-028SS-S	11-002-028SS
Ø3.5	30mm	11-002-030TI-S	11-002-030TI	11-002-030SS-S	11-002-030SS
Ø3.5	32mm	11-002-032TI-S	11-002-032TI	11-002-032SS-S	11-002-032SS
Ø3.5	34mm	11-002-034TI-S	11-002-034TI	11-002-034SS-S	11-002-034SS
Ø3.5	36mm	11-002-036TI-S	11-002-036TI	11-002-036SS-S	11-002-036SS
Ø3.5	38mm	11-002-038TI-S	11-002-038TI	11-002-038SS-S	11-002-038SS
Ø3.5	40mm	11-002-040TI-S	11-002-040TI	11-002-040SS-S	11-002-040SS
Ø3.5	42mm	11-002-042TI-S	11-002-042TI	11-002-042SS-S	11-002-042SS
Ø3.5	44mm	11-002-044TI-S	11-002-044TI	11-002-044SS-S	11-002-044SS
Ø3.5	46mm	11-002-046TI-S	11-002-046TI	11-002-046SS-S	11-002-046SS
Ø3.5	48mm	11-002-048TI-S	11-002-048TI	11-002-048SS-S	11-002-048SS
Ø3.5	50mm	11-002-050TI-S	11-002-050TI	11-002-050SS-S	11-002-050SS
Ø3.5	55mm	11-002-055TI-S	11-002-055TI	11-002-055SS-S	11-002-055SS
Ø3.5	60mm	11-002-060TI-S	11-002-060TI	11-002-060SS-S	11-002-060SS
Ø3.5	65mm	11-002-065TI-S	11-002-065TI	11-002-065SS-S	11-002-065SS
Ø3.5	70mm	11-002-070TI-S	11-002-070TI	11-002-070SS-S	11-002-070SS
Ø3.5	75mm	11-002-075TI-S	11-002-075TI	11-002-075SS-S	11-002-075SS
Ø3.5	80mm	11-002-080TI-S	11-002-080TI	11-002-080SS-S	11-002-080SS
Ø3.5	85mm	11-002-085TI-S	11-002-085TI	11-002-085SS-S	11-002-085SS
Ø3.5	90mm	11-002-090TI-S	11-002-090TI	11-002-090SS-S	11-002-090SS




2.7mm Cortical Screw, Self-Tapping, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø2.7	10mm	TI-1242.10-S	TI-1242.10	1242.10-S	1242.10
Ø2.7	12mm	TI-1242.12-S	TI-1242.12	1242.12-S	1242.12
Ø2.7	14mm	TI-1242.14-S	TI-1242.14	1242.14-S	1242.14
Ø2.7	16mm	TI-1242.16-S	TI-1242.16	1242.16-S	1242.16
Ø2.7	18mm	TI-1242.18-S	TI-1242.18	1242.18-S	1242.18
Ø2.7	20mm	TI-1242.20-S	TI-1242.20	1242.20-S	1242.20
Ø2.7	22mm	TI-1242.22-S	TI-1242.22	1242.22-S	1242.22
Ø2.7	24mm	TI-1242.24-S	TI-1242.24	1242.24-S	1242.24
Ø2.7	26mm	TI-1242.26-S	TI-1242.26	1242.26-S	1242.26
Ø2.7	28mm	TI-1242.28-S	TI-1242.28	1242.28-S	1242.28
Ø2.7	30mm	TI-1242.30-S	TI-1242.30	1242.30-S	1242.30
Ø2.7	32mm	TI-1242.32-S	TI-1242.32	1242.32-S	1242.32
Ø2.7	34mm	TI-1242.34-S	TI-1242.34	1242.34-S	1242.34
Ø2.7	36mm	TI-1242.36-S	TI-1242.36	1242.36-S	1242.36
Ø2.7	38mm	TI-1242.38-S	TI-1242.38	1242.38-S	1242.38
Ø2.7	40mm	TI-1242.40-S	TI-1242.40	1242.40-S	1242.40
Ø2.7	42mm	TI-1242.42-S	TI-1242.42	1242.42-S	1242.42
Ø2.7	44mm	TI-1242.44-S	TI-1242.44	1242.44-S	1242.44
Ø2.7	46mm	TI-1242.46-S	TI-1242.46	1242.46-S	1242.46
Ø2.7	48mm	TI-1242.48-S	TI-1242.48	1242.48-S	1242.48
Ø2.7	50mm	TI-1242.50-S	TI-1242.50	1242.50-S	1242.50
Ø2.7	52mm	TI-1242.52-S	TI-1242.52	1242.52-S	1242.52
Ø2.7	54mm	TI-1242.54-S	TI-1242.54	1242.54-S	1242.54
Ø2.7	56mm	TI-1242.56-S	TI-1242.56	1242.56-S	1242.56
Ø2.7	58mm	TI-1242.58-S	TI-1242.58	1242.58-S	1242.58
Ø2.7	60mm	TI-1242.60-S	TI-1242.60	1242.60-S	1242.60
Ø2.7	62mm	TI-1242.62-S	TI-1242.62	1242.62-S	1242.62
Ø2.7	64mm	TI-1242.64-S	TI-1242.64	1242.64-S	1242.64
Ø2.7	66mm	TI-1242.66-S	TI-1242.66	1242.66-S	1242.66
Ø2.7	68mm	TI-1242.68-S	TI-1242.68	1242.68-S	1242.68
Ø2.7	70mm	TI-1242.70-S	TI-1242.70	1242.70-S	1242.70




2.7mm Cortical Screw, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø2.7	6mm	11-013-006TI-S	11-013-006TI	11-013-006SS-S	11-013-006SS
Ø2.7	8mm	11-013-008TI-S	11-013-008TI	11-013-008SS-S	11-013-008SS
Ø2.7	10mm	11-013-010TI-S	11-013-010TI	11-013-010SS-S	11-013-010SS
Ø2.7	12mm	11-013-012TI-S	11-013-012TI	11-013-012SS-S	11-013-012SS
Ø2.7	14mm	11-013-014TI-S	11-013-014TI	11-013-014SS-S	11-013-014SS
Ø2.7	16mm	11-013-016TI-S	11-013-016TI	11-013-016SS-S	11-013-016SS
Ø2.7	18mm	11-013-018TI-S	11-013-018TI	11-013-018SS-S	11-013-018SS
Ø2.7	20mm	11-013-020TI-S	11-013-020TI	11-013-020SS-S	11-013-020SS
Ø2.7	22mm	11-013-022TI-S	11-013-022TI	11-013-022SS-S	11-013-022SS
Ø2.7	24mm	11-013-024TI-S	11-013-024TI	11-013-024SS-S	11-013-024SS
Ø2.7	26mm	11-013-026TI-S	11-013-026TI	11-013-026SS-S	11-013-026SS
Ø2.7	28mm	11-013-028TI-S	11-013-028TI	11-013-028SS-S	11-013-028SS
Ø2.7	30mm	11-013-030TI-S	11-013-030TI	11-013-030SS-S	11-013-030SS
Ø2.7	32mm	11-013-032TI-S	11-013-032TI	11-013-032SS-S	11-013-032SS
Ø2.7	34mm	11-013-034TI-S	11-013-034TI	11-013-034SS-S	11-013-034SS
Ø2.7	36mm	11-013-036TI-S	11-013-036TI	11-013-036SS-S	11-013-036SS
Ø2.7	38mm	11-013-038TI-S	11-013-038TI	11-013-038SS-S	11-013-038SS
Ø2.7	40mm	11-013-040TI-S	11-013-040TI	11-013-040SS-S	11-013-040SS
Ø2.7	42mm	11-013-042TI-S	11-013-042TI	11-013-042SS-S	11-013-042SS
Ø2.7	44mm	11-013-044TI-S	11-013-044TI	11-013-044SS-S	11-013-044SS
Ø2.7	46mm	11-013-046TI-S	11-013-046TI	11-013-046SS-S	11-013-046SS
Ø2.7	48mm	11-013-048TI-S	11-013-048TI	11-013-048SS-S	11-013-048SS
Ø2.7	50mm	11-013-050TI-S	11-013-050TI	11-013-050SS-S	11-013-050SS
Ø2.7	52mm	11-013-052TI-S	11-013-052TI	11-013-052SS-S	11-013-052SS
Ø2.7	54mm	11-013-054TI-S	11-013-054TI	11-013-054SS-S	11-013-054SS
Ø2.7	56mm	11-013-056TI-S	11-013-056TI	11-013-056SS-S	11-013-056SS
Ø2.7	58mm	11-013-058TI-S	11-013-058TI	11-013-058SS-S	11-013-058SS
Ø2.7	60mm	11-013-060TI-S	11-013-060TI	11-013-060SS-S	11-013-060SS
Ø2.7	62mm	11-013-062TI-S	11-013-062TI	11-013-062SS-S	11-013-062SS
Ø2.7	64mm	11-013-064TI-S	11-013-064TI	11-013-064SS-S	11-013-064SS
Ø2.7	66mm	11-013-066TI-S	11-013-066TI	11-013-066SS-S	11-013-066SS
Ø2.7	68mm	11-013-068TI-S	11-013-068TI	11-013-068SS-S	11-013-068SS
Ø2.7	70mm	11-013-070TI-S	11-013-070TI	11-013-070SS-S	11-013-070SS




2.4mm Cortical Screw, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø2.4	6mm	11-015-006TI-S	11-015-006TI	11-015-006SS-S	11-015-006SS
Ø2.4	8mm	11-015-008TI-S	11-015-008TI	11-015-008SS-S	11-015-008SS
Ø2.4	10mm	11-015-010TI-S	11-015-010TI	11-015-010SS-S	11-015-010SS
Ø2.4	12mm	11-015-012TI-S	11-015-012TI	11-015-012SS-S	11-015-012SS
Ø2.4	14mm	11-015-014TI-S	11-015-014TI	11-015-014SS-S	11-015-014SS
Ø2.4	16mm	11-015-016TI-S	11-015-016TI	11-015-016SS-S	11-015-016SS
Ø2.4	18mm	11-015-018TI-S	11-015-018TI	11-015-018SS-S	11-015-018SS
Ø2.4	20mm	11-015-020TI-S	11-015-020TI	11-015-020SS-S	11-015-020SS
Ø2.4	22mm	11-015-022TI-S	11-015-022TI	11-015-022SS-S	11-015-022SS
Ø2.4	24mm	11-015-024TI-S	11-015-024TI	11-015-024SS-S	11-015-024SS
Ø2.4	26mm	11-015-026TI-S	11-015-026TI	11-015-026SS-S	11-015-026SS
Ø2.4	28mm	11-015-028TI-S	11-015-028TI	11-015-028SS-S	11-015-028SS
Ø2.4	30mm	11-015-030TI-S	11-015-030TI	11-015-030SS-S	11-015-030SS
Ø2.4	32mm	11-015-032TI-S	11-015-032TI	11-015-032SS-S	11-015-032SS
Ø2.4	34mm	11-015-034TI-S	11-015-034TI	11-015-034SS-S	11-015-034SS
Ø2.4	36mm	11-015-036TI-S	11-015-036TI	11-015-036SS-S	11-015-036SS
Ø2.4	38mm	11-015-038TI-S	11-015-038TI	11-015-038SS-S	11-015-038SS
Ø2.4	40mm	11-015-040TI-S	11-015-040TI	11-015-040SS-S	11-015-040SS




2.7mm Wise-Lock Screw, Self Tapping, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø2.7	6mm	TI-1241.06-S	TI-1241.06	11-012-006SS-S	11-012-006SS
Ø2.7	7mm	TI-1241.07-S	TI-1241.07	11-012-007SS-S	11-012-007SS
Ø2.7	8mm	TI-1241.08-S	TI-1241.08	11-012-008SS-S	11-012-008SS
Ø2.7	9mm	TI-1241.09-S	TI-1241.09	11-012-009SS-S	11-012-009SS
Ø2.7	10mm	TI-1241.10-S	TI-1241.10	11-012-010SS-S	11-012-010SS
Ø2.7	12mm	TI-1241.12-S	TI-1241.12	11-012-012SS-S	11-012-012SS
Ø2.7	14mm	TI-1241.14-S	TI-1241.14	11-012-014SS-S	11-012-014SS
Ø2.7	16mm	TI-1241.16-S	TI-1241.16	11-012-016SS-S	11-012-016SS
Ø2.7	18mm	TI-1241.18-S	TI-1241.18	11-012-018SS-S	11-012-018SS
Ø2.7	20mm	TI-1241.20-S	TI-1241.20	11-012-020SS-S	11-012-020SS
Ø2.7	22mm	TI-1241.22-S	TI-1241.22	11-012-022SS-S	11-012-022SS
Ø2.7	24mm	TI-1241.24-S	TI-1241.24	11-012-024SS-S	11-012-024SS
Ø2.7	26mm	TI-1241.26-S	TI-1241.26	11-012-026SS-S	11-012-026SS
Ø2.7	28mm	TI-1241.28-S	TI-1241.28	11-012-028SS-S	11-012-028SS
Ø2.7	30mm	TI-1241.30-S	TI-1241.30	11-012-030SS-S	11-012-030SS
Ø2.7	32mm	TI-1241.32-S	TI-1241.32	11-012-032SS-S	11-012-032SS
Ø2.7	34mm	TI-1241.34-S	TI-1241.34	11-012-034SS-S	11-012-034SS
Ø2.7	36mm	TI-1241.36-S	TI-1241.36	11-012-036SS-S	11-012-036SS
Ø2.7	38mm	TI-1241.38-S	TI-1241.38	11-012-038SS-S	11-012-038SS
Ø2.7	40mm	TI-1241.40-S	TI-1241.40	11-012-040SS-S	11-012-040SS
Ø2.7	42mm	TI-1241.42-S	TI-1241.42	11-012-042SS-S	11-012-042SS
Ø2.7	44mm	TI-1241.44-S	TI-1241.44	11-012-044SS-S	11-012-044SS
Ø2.7	46mm	TI-1241.46-S	TI-1241.46	11-012-046SS-S	11-012-046SS
Ø2.7	48mm	TI-1241.48-S	TI-1241.48	11-012-048SS-S	11-012-048SS
Ø2.7	50mm	TI-1241.50-S	TI-1241.50	11-012-050SS-S	11-012-050SS
Ø2.7	52mm	TI-1241.52-S	TI-1241.52	11-012-052SS-S	11-012-052SS
Ø2.7	54mm	TI-1241.54-S	TI-1241.54	11-012-054SS-S	11-012-054SS
Ø2.7	56mm	TI-1241.56-S	TI-1241.56	11-012-056SS-S	11-012-056SS
Ø2.7	58mm	TI-1241.58-S	TI-1241.58	11-012-058SS-S	11-012-058SS
Ø2.7	60mm	TI-1241.60-S	TI-1241.60	11-012-060SS-S	11-012-060SS
Ø2.7	62mm	11-012-062TI-S	11-012-062TI	11-012-062SS-S	11-012-062SS
Ø2.7	64mm	11-012-064TI-S	11-012-064TI	11-012-064SS-S	11-012-064SS
Ø2.7	66mm	11-012-066TI-S	11-012-066TI	11-012-066SS-S	11-012-066SS
Ø2.7	68mm	11-012-068TI-S	11-012-068TI	11-012-068SS-S	11-012-068SS
Ø2.7	70mm	11-012-070TI-S	11-012-070TI	11-012-070SS-S	11-012-070SS



3.5mm Cortical Screw, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø3.5	10mm	11-016-010TI-S	11-016-010TI	11-016-010SS-S	11-016-010SS
Ø3.5	12mm	11-016-012TI-S	11-016-012TI	11-016-012SS-S	11-016-012SS
Ø3.5	14mm	11-016-014TI-S	11-016-014TI	11-016-014SS-S	11-016-014SS
Ø3.5	16mm	11-016-016TI-S	11-016-016TI	11-016-016SS-S	11-016-016SS
Ø3.5	18mm	11-016-018TI-S	11-016-018TI	11-016-018SS-S	11-016-018SS
Ø3.5	20mm	11-016-020TI-S	11-016-020TI	11-016-020SS-S	11-016-020SS
Ø3.5	22mm	11-016-022TI-S	11-016-022TI	11-016-022SS-S	11-016-022SS
Ø3.5	24mm	11-016-024TI-S	11-016-024TI	11-016-024SS-S	11-016-024SS
Ø3.5	26mm	11-016-026TI-S	11-016-026TI	11-016-026SS-S	11-016-026SS
Ø3.5	28mm	11-016-028TI-S	11-016-028TI	11-016-028SS-S	11-016-028SS
Ø3.5	30mm	11-016-030TI-S	11-016-030TI	11-016-030SS-S	11-016-030SS
Ø3.5	32mm	11-016-032TI-S	11-016-032TI	11-016-032SS-S	11-016-032SS
Ø3.5	34mm	11-016-034TI-S	11-016-034TI	11-016-034SS-S	11-016-034SS
Ø3.5	36mm	11-016-036TI-S	11-016-036TI	11-016-036SS-S	11-016-036SS
Ø3.5	38mm	11-016-038TI-S	11-016-038TI	11-016-038SS-S	11-016-038SS
Ø3.5	40mm	11-016-040TI-S	11-016-040TI	11-016-040SS-S	11-016-040SS
Ø3.5	42mm	11-016-042TI-S	11-016-042TI	11-016-042SS-S	11-016-042SS
Ø3.5	44mm	11-016-044TI-S	11-016-044TI	11-016-044SS-S	11-016-044SS
Ø3.5	46mm	11-016-046TI-S	11-016-046TI	11-016-046SS-S	11-016-046SS
Ø3.5	48mm	11-016-048TI-S	11-016-048TI	11-016-048SS-S	11-016-048SS
Ø3.5	50mm	11-016-050TI-S	11-016-050TI	11-016-050SS-S	11-016-050SS
Ø3.5	55mm	11-016-055TI-S	11-016-055TI	11-016-055SS-S	11-016-055SS
Ø3.5	60mm	11-016-060TI-S	11-016-060TI	11-016-060SS-S	11-016-060SS
Ø3.5	65mm	11-016-065TI-S	11-016-065TI	11-016-065SS-S	11-016-065SS
Ø3.5	70mm	11-016-070TI-S	11-016-070TI	11-016-070SS-S	11-016-070SS
Ø3.5	75mm	11-016-075TI-S	11-016-075TI	11-016-075SS-S	11-016-075SS
Ø3.5	80mm	11-016-080TI-S	11-016-080TI	11-016-080SS-S	11-016-080SS
Ø3.5	85mm	11-016-085TI-S	11-016-085TI	11-016-085SS-S	11-016-085SS
Ø3.5	90mm	11-016-090TI-S	11-016-090TI	11-016-090SS-S	11-016-090SS



AV-WISELOCK

3.5MM PROXIMAL TIBIA SYSTEM

Small Bend & Large Bend

The 3.5mm AV-Wiselock Proximal Tibia Plating System combines AV-Wiselock locking screw technology with traditional locking system.

The plates are designed in two bends (small bend and large bend) to accommodate varying anatomies of tibia with shaft hole varying from 4 to 14 shaft holes.


The 3.5mm AV-Wiselock Proximal Tibia Plates are intended to treat fractures of the proximal tibia in adults and adolescents in which the growth plates have fused. Plates can also be used for treatment of nonunions, malunions, tibial osteotomies and osteopenic bone.

The Av-Wiselock Proximal Tibia Plates Are Intended To Treat Fractures Of The Proximal Tibia In Adults And Adolescents In Which The Growth Plates Have Fused Including:

Simple, Comminuted, Lateral Wedge, Depression, Medial Wedge, Bicondylar Combination Of Lateral Wedge And Depression, Periprosthetic, And Fractures With Associated Shaft Fractures. Plates Can Also Be Used For Treatment Of Nonunions, Malunions, Tibial Osteotomies And Osteopenic Bone.

Plates Available in


Titanium as per ISO 5832-3 and Stainless Steel as per ISO 5832-1

3.5mm AV-Wiselock Proximal Tibia Plate, Small Bend

Titanium

Holes	Direction	Sterile	Non-Sterile
4	Left	11-026-04LTI-S	11-026-04LTI
4	Right	11-026-04RTI-S	11-026-04RTI
6	Left	11-026-06LTI-S	11-026-06LTI
6	Right	11-026-06RTI-S	11-026-06RTI
8	Left	11-026-08LTI-S	11-026-08LTI
8	Right	11-026-08RTI-S	11-026-08RTI
10	Left	11-026-10LTI-S	11-026-10LTI
10	Right	11-026-10RTI-S	11-026-10RTI
12	Left	11-026-12LTI-S	11-026-12LTI
12	Right	11-026-12RTI-S	11-026-12RTI
14	Left	11-026-14LTI-S	11-026-14LTI
14	Right	11-026-14RTI-S	11-026-14RTI

Stainless Steel

Holes	Direction	Sterile	Non-Sterile
4	Left	11-026-04LSS-S	11-026-04LSS
4	Right	11-026-04RSS-S	11-026-04RSS
6	Left	11-026-06LSS-S	11-026-06LSS
6	Right	11-026-06RSS-S	11-026-06RSS
8	Left	11-026-08LSS-S	11-026-08LSS
8	Right	11-026-08RSS-S	11-026-08RSS
10	Left	11-026-10LSS-S	11-026-10LSS
10	Right	11-026-10RSS-S	11-026-10RSS
12	Left	11-026-12LSS-S	11-026-12LSS
12	Right	11-026-12RSS-S	11-026-12RSS
14	Left	11-026-14LSS-S	11-026-14LSS
14	Right	11-026-14RSS-S	11-026-14RSS


3.5mm AV-Wiselock Proximal Tibia Plate, Large Bend


Titanium			
Holes	Direction	Sterile	Non-Sterile
4	Left	11-027-04LTI-S	11-027-04LTI
4	Right	11-027-04RTI-S	11-027-04RTI
6	Left	11-027-06LTI-S	11-027-06LTI
6	Right	11-027-06RTI-S	11-027-06RTI
8	Left	11-027-08LTI-S	11-027-08LTI
8	Right	11-027-08RTI-S	11-027-08RTI
10	Left	11-027-10LTI-S	11-027-10LTI
10	Right	11-027-10RTI-S	11-027-10RTI
12	Left	11-027-12LTI-S	11-027-12LTI
12	Right	11-027-12RTI-S	11-027-12RTI
14	Left	11-027-14LTI-S	11-027-14LTI
14	Right	11-027-14RTI-S	11-027-14RTI

Stainless Steel			
Holes	Direction	Sterile	Non-Sterile
4	Left	11-027-04LSS-S	11-027-04LSS
4	Right	11-027-04RSS-S	11-027-04RSS
6	Left	11-027-06LSS-S	11-027-06LSS
6	Right	11-027-06RSS-S	11-027-06RSS
8	Left	11-027-08LSS-S	11-027-08LSS
8	Right	11-027-08RSS-S	11-027-08RSS
10	Left	11-027-10LSS-S	11-027-10LSS
10	Right	11-027-10RSS-S	11-027-10RSS
12	Left	11-027-12LSS-S	11-027-12LSS
12	Right	11-027-12RSS-S	11-027-12RSS
14	Left	11-027-14LSS-S	11-027-14LSS
14	Right	11-027-14RSS-S	11-027-14RSS


3.5mm AV-Wiselock Screw, Self-Tapping

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø3.5	10mm	11-024-010TI-S	11-024-010TI	11-024-010SS-S	11-024-010SS
Ø3.5	12mm	11-024-012TI-S	11-024-012TI	11-024-012SS-S	11-024-012SS
Ø3.5	14mm	11-024-014TI-S	11-024-014TI	11-024-014SS-S	11-024-014SS
Ø3.5	16mm	11-024-016TI-S	11-024-016TI	11-024-016SS-S	11-024-016SS
Ø3.5	18mm	11-024-018TI-S	11-024-018TI	11-024-018SS-S	11-024-018SS
Ø3.5	20mm	11-024-020TI-S	11-024-020TI	11-024-020SS-S	11-024-020SS
Ø3.5	22mm	11-024-022TI-S	11-024-022TI	11-024-022SS-S	11-024-022SS
Ø3.5	24mm	11-024-024TI-S	11-024-024TI	11-024-024SS-S	11-024-024SS
Ø3.5	26mm	11-024-026TI-S	11-024-026TI	11-024-026SS-S	11-024-026SS
Ø3.5	28mm	11-024-028TI-S	11-024-028TI	11-024-028SS-S	11-024-028SS
Ø3.5	30mm	11-024-030TI-S	11-024-030TI	11-024-030SS-S	11-024-030SS
Ø3.5	32mm	11-024-032TI-S	11-024-032TI	11-024-032SS-S	11-024-032SS
Ø3.5	34mm	11-024-034TI-S	11-024-034TI	11-024-034SS-S	11-024-034SS
Ø3.5	36mm	11-024-036TI-S	11-024-036TI	11-024-036SS-S	11-024-036SS
Ø3.5	38mm	11-024-038TI-S	11-024-038TI	11-024-038SS-S	11-024-038SS
Ø3.5	40mm	11-024-040TI-S	11-024-040TI	11-024-040SS-S	11-024-040SS
Ø3.5	42mm	11-024-042TI-S	11-024-042TI	11-024-042SS-S	11-024-042SS
Ø3.5	44mm	11-024-044TI-S	11-024-044TI	11-024-044SS-S	11-024-044SS
Ø3.5	46mm	11-024-046TI-S	11-024-046TI	11-024-046SS-S	11-024-046SS
Ø3.5	48mm	11-024-048TI-S	11-024-048TI	11-024-048SS-S	11-024-048SS
Ø3.5	50mm	11-024-050TI-S	11-024-050TI	11-024-050SS-S	11-024-050SS
Ø3.5	52mm	11-024-052TI-S	11-024-052TI	11-024-052SS-S	11-024-052SS
Ø3.5	54mm	11-024-054TI-S	11-024-054TI	11-024-054SS-S	11-024-054SS
Ø3.5	56mm	11-024-056TI-S	11-024-056TI	11-024-056SS-S	11-024-056SS
Ø3.5	58mm	11-024-058TI-S	11-024-058TI	11-024-058SS-S	11-024-058SS
Ø3.5	60mm	11-024-060TI-S	11-024-060TI	11-024-060SS-S	11-024-060SS
Ø3.5	65mm	11-024-065TI-S	11-024-065TI	11-024-065SS-S	11-024-065SS
Ø3.5	70mm	11-024-070TI-S	11-024-070TI	11-024-070SS-S	11-024-070SS
Ø3.5	75mm	11-024-075TI-S	11-024-075TI	11-024-075SS-S	11-024-075SS
Ø3.5	80mm	11-024-080TI-S	11-024-080TI	11-024-080SS-S	11-024-080SS
Ø3.5	85mm	11-024-085TI-S	11-024-085TI	11-024-085SS-S	11-024-085SS
Ø3.5	90mm	11-024-090TI-S	11-024-090TI	11-024-090SS-S	11-024-090SS
Ø3.5	95mm	11-024-095TI-S	11-024-095TI	11-024-095SS-S	11-024-095SS



3.5mm Wise-Lock Screw, Self Tapping, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø3.5	10mm	TI-1473-010-S	TI-1473-010	1473-010-S	1473-010
Ø3.5	12mm	TI-1473-012-S	TI-1473-012	1473-012-S	1473-012
Ø3.5	14mm	TI-1473-014-S	TI-1473-014	1473-014-S	1473-014
Ø3.5	16mm	TI-1473-016-S	TI-1473-016	1473-016-S	1473-016
Ø3.5	18mm	TI-1473-018-S	TI-1473-018	1473-018-S	1473-018
Ø3.5	20mm	TI-1473-020-S	TI-1473-020	1473-020-S	1473-020
Ø3.5	22mm	TI-1473-022-S	TI-1473-022	1473-022-S	1473-022
Ø3.5	24mm	TI-1473-024-S	TI-1473-024	1473-024-S	1473-024
Ø3.5	26mm	TI-1473-026-S	TI-1473-026	1473-026-S	1473-026
Ø3.5	28mm	TI-1473-028-S	TI-1473-028	1473-028-S	1473-028
Ø3.5	30mm	TI-1473-030-S	TI-1473-030	1473-030-S	1473-030
Ø3.5	32mm	TI-1473-032-S	TI-1473-032	1473-032-S	1473-032
Ø3.5	34mm	TI-1473-034-S	TI-1473-034	1473-034-S	1473-034
Ø3.5	36mm	TI-1473-036-S	TI-1473-036	1473-036-S	1473-036
Ø3.5	38mm	TI-1473-038-S	TI-1473-038	1473-038-S	1473-038
Ø3.5	40mm	TI-1473-040-S	TI-1473-040	1473-040-S	1473-040
Ø3.5	42mm	TI-1473-042-S	TI-1473-042	1473-042-S	1473-042
Ø3.5	44mm	TI-1473-044-S	TI-1473-044	1473-044-S	1473-044
Ø3.5	46mm	TI-1473-046-S	TI-1473-046	1473-046-S	1473-046
Ø3.5	48mm	TI-1473-048-S	TI-1473-048	1473-048-S	1473-048
Ø3.5	50mm	TI-1473-050-S	TI-1473-050	1473-050-S	1473-050
Ø3.5	52mm	TI-1473-052-S	TI-1473-052	1473-052-S	1473-052
Ø3.5	54mm	TI-1473-054-S	TI-1473-054	1473-054-S	1473-054
Ø3.5	56mm	TI-1473-056-S	TI-1473-056	1473-056-S	1473-056
Ø3.5	58mm	TI-1473-058-S	TI-1473-058	1473-058-S	1473-058
Ø3.5	60mm	TI-1473-060-S	TI-1473-060	1473-060-S	1473-060
Ø3.5	65mm	TI-1473-065-S	TI-1473-065	1473-065-S	1473-065
Ø3.5	70mm	TI-1473-070-S	TI-1473-070	1473-070-S	1473-070
Ø3.5	75mm	TI-1473-075-S	TI-1473-075	1473-075-S	1473-075
Ø3.5	80mm	TI-1473-080-S	TI-1473-080	1473-080-S	1473-080
Ø3.5	85mm	TI-1473-085-S	TI-1473-085	1473-085-S	1473-085
Ø3.5	90mm	TI-1473-090-S	TI-1473-090	1473-090-S	1473-090
Ø3.5	95mm	TI-1473-095-S	TI-1473-095	1473-095-S	1473-095



3.5mm Cortical Screw, Self Tapping (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø3.5	10mm	11-002-010TI-S	11-002-010TI	11-002-010SS-S	11-002-010SS
Ø3.5	12mm	11-002-012TI-S	11-002-012TI	11-002-012SS-S	11-002-012SS
Ø3.5	14mm	11-002-014TI-S	11-002-014TI	11-002-014SS-S	11-002-014SS
Ø3.5	16mm	11-002-016TI-S	11-002-016TI	11-002-016SS-S	11-002-016SS
Ø3.5	18mm	11-002-018TI-S	11-002-018TI	11-002-018SS-S	11-002-018SS
Ø3.5	20mm	11-002-020TI-S	11-002-020TI	11-002-020SS-S	11-002-020SS
Ø3.5	22mm	11-002-022TI-S	11-002-022TI	11-002-022SS-S	11-002-022SS
Ø3.5	24mm	11-002-024TI-S	11-002-024TI	11-002-024SS-S	11-002-024SS
Ø3.5	26mm	11-002-026TI-S	11-002-026TI	11-002-026SS-S	11-002-026SS
Ø3.5	28mm	11-002-028TI-S	11-002-028TI	11-002-028SS-S	11-002-028SS
Ø3.5	30mm	11-002-030TI-S	11-002-030TI	11-002-030SS-S	11-002-030SS
Ø3.5	32mm	11-002-032TI-S	11-002-032TI	11-002-032SS-S	11-002-032SS
Ø3.5	34mm	11-002-034TI-S	11-002-034TI	11-002-034SS-S	11-002-034SS
Ø3.5	36mm	11-002-036TI-S	11-002-036TI	11-002-036SS-S	11-002-036SS
Ø3.5	38mm	11-002-038TI-S	11-002-038TI	11-002-038SS-S	11-002-038SS
Ø3.5	40mm	11-002-040TI-S	11-002-040TI	11-002-040SS-S	11-002-040SS
Ø3.5	42mm	11-002-042TI-S	11-002-042TI	11-002-042SS-S	11-002-042SS
Ø3.5	44mm	11-002-044TI-S	11-002-044TI	11-002-044SS-S	11-002-044SS
Ø3.5	46mm	11-002-046TI-S	11-002-046TI	11-002-046SS-S	11-002-046SS
Ø3.5	48mm	11-002-048TI-S	11-002-048TI	11-002-048SS-S	11-002-048SS
Ø3.5	50mm	11-002-050TI-S	11-002-050TI	11-002-050SS-S	11-002-050SS
Ø3.5	55mm	11-002-055TI-S	11-002-055TI	11-002-055SS-S	11-002-055SS
Ø3.5	60mm	11-002-060TI-S	11-002-060TI	11-002-060SS-S	11-002-060SS
Ø3.5	65mm	11-002-065TI-S	11-002-065TI	11-002-065SS-S	11-002-065SS
Ø3.5	70mm	11-002-070TI-S	11-002-070TI	11-002-070SS-S	11-002-070SS
Ø3.5	75mm	11-002-075TI-S	11-002-075TI	11-002-075SS-S	11-002-075SS
Ø3.5	80mm	11-002-080TI-S	11-002-080TI	11-002-080SS-S	11-002-080SS
Ø3.5	85mm	11-002-085TI-S	11-002-085TI	11-002-085SS-S	11-002-085SS
Ø3.5	90mm	11-002-090TI-S	11-002-090TI	11-002-090SS-S	11-002-090SS
Ø3.5	95mm	11-002-095TI-S	11-002-095TI	11-002-095SS-S	11-002-095SS



AV-WISELOCK

2.7mm Foot Hallux System

2.7mm AV-Wiselock Foot systems offers the different multidirectional and angular stable plates with different designs, lengths and widths for a range of applications in the foot.

Indication

AV-Wiselock foot System is indicated for fixation of osteotomies, fusions, fractures, nonunions, malunions, and replantations of small bones and small bone fragments in adult and adolescent (12–21 years) patients, including the Foot and Ankle.


Plate Features

- Low-Profile plates with better anatomical Fit
- Oblong hole for variable positioning of the plate
- AV-Wiselock– variable angle of $\pm 15^\circ$ in all directions in each screw hole
- Pre-angled Variable Angle holes for oriented screw placement
- Rounded edges and a smooth surface for soft tissue protection
- Universal Variable angle drill guides available for easy angulation of screws
- K-wire holes and slot for temporary plate fixation and compression

Plates Available in

Titanium as per ISO 5832-3 and Stainless Steel as per ISO 5832-1

2.7mm AV-Wiselock MTP Fusion M Plate




Titanium

Holes	Degree	Direction	Sterile	Non-Sterile
7	0	Left	11-057-01TIL-S	11-057-01TIL
7	0	Right	11-057-01TIR-S	11-057-01TIR
7	5	Left	11-057-02TIL-S	11-057-02TIL
7	5	Right	11-057-02TIR-S	11-057-02TIR
7	10	Left	11-057-03TIL-S	11-057-03TIL
7	10	Right	11-057-03TIR-S	11-057-03TIR

Stainless Steel

Holes	Degree	Direction	Sterile	Non-Sterile
7	0	Left	11-057-01SSL-S	11-057-01SSL
7	0	Right	11-057-01SSR-S	11-057-01SSR
7	5	Left	11-057-02SSL-S	11-057-02SSL
7	5	Right	11-057-02SSR-S	11-057-02SSR
7	10	Left	11-057-03SSL-S	11-057-03SSL
7	10	Right	11-057-03SSR-S	11-057-03SSR

2.7mm AV-Wiselock MTP Revision M Plate




Titanium

Holes	Degree	Direction	Sterile	Non-Sterile
9	5	Left	11-057-04TIL-S	11-057-04TIL
9	5	Right	11-057-04TIR-S	11-057-04TIR
9	10	Left	11-057-05TIL-S	11-057-05TIL
9	10	Right	11-057-05TIR-S	11-057-05TIR



Stainless Steel

Holes	Degree	Direction	Sterile	Non-Sterile
9	5	Left	11-057-04SSL-S	11-057-04SSL
9	5	Right	11-057-04SSR-S	11-057-04SSR
9	10	Left	11-057-05SSL-S	11-057-05SSL
9	10	Right	11-057-05SSR-S	11-057-05SSR

2.7mm AV-Wiselock TMT-I Medial Fusion M Plate


Titanium			
Hole	Direction	Sterile	Non-Sterile
7	Left	11-057-06TIL-S	11-057-06TIL
7	Right	11-057-06TIR-S	11-057-06TIR

Stainless Steel			
Hole	Direction	Sterile	Non-Sterile
7	Left	11-057-06SSL-S	11-057-06SSL
7	Right	11-057-06SSR-S	11-057-06SSR


2.7mm AV-Wiselock TMT-I Plantar Fusion M Plate


Titanium			
Hole	Direction	Sterile	Non-Sterile
6	Left	11-057-07TIL-S	11-057-07TIL
6	Right	11-057-07TIR-S	11-057-07TIR

Stainless Steel			
Hole	Direction	Sterile	Non-Sterile
6	Left	11-057-07SSL-S	11-057-07SSL
6	Right	11-057-07SSR-S	11-057-07SSR


2.7mm Cortical Screw, Self-Tapping, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø2.7	8mm	TI-1242.08-S	TI-1242.08	1242.08-S	1242.08
Ø2.7	10mm	TI-1242.10-S	TI-1242.10	1242.10-S	1242.10
Ø2.7	12mm	TI-1242.12-S	TI-1242.12	1242.12-S	1242.12
Ø2.7	14mm	TI-1242.14-S	TI-1242.14	1242.14-S	1242.14
Ø2.7	16mm	TI-1242.16-S	TI-1242.16	1242.16-S	1242.16
Ø2.7	18mm	TI-1242.18-S	TI-1242.18	1242.18-S	1242.18
Ø2.7	20mm	TI-1242.20-S	TI-1242.20	1242.20-S	1242.20
Ø2.7	22mm	TI-1242.22-S	TI-1242.22	1242.22-S	1242.22
Ø2.7	24mm	TI-1242.24-S	TI-1242.24	1242.24-S	1242.24
Ø2.7	26mm	TI-1242.26-S	TI-1242.26	1242.26-S	1242.26
Ø2.7	28mm	TI-1242.28-S	TI-1242.28	1242.28-S	1242.28
Ø2.7	30mm	TI-1242.30-S	TI-1242.30	1242.30-S	1242.30
Ø2.7	32mm	TI-1242.32-S	TI-1242.32	1242.32-S	1242.32
Ø2.7	34mm	TI-1242.34-S	TI-1242.34	1242.34-S	1242.34
Ø2.7	36mm	TI-1242.36-S	TI-1242.36	1242.36-S	1242.36
Ø2.7	38mm	TI-1242.38-S	TI-1242.38	1242.38-S	1242.38
Ø2.7	40mm	TI-1242.40-S	TI-1242.40	1242.40-S	1242.40
Ø2.7	42mm	TI-1242.42-S	TI-1242.42	1242.42-S	1242.42
Ø2.7	44mm	TI-1242.44-S	TI-1242.44	1242.44-S	1242.44
Ø2.7	46mm	TI-1242.46-S	TI-1242.46	1242.46-S	1242.46




2.7mm AV-Wiselock Screw, Self Tapping, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø2.7	8mm	TI-1472-008-S	TI-1472-008	1472-008-S	1472-008
Ø2.7	10mm	TI-1472-010-S	TI-1472-010	1472-010-S	1472-010
Ø2.7	12mm	TI-1472-012-S	TI-1472-012	1472-012-S	1472-012
Ø2.7	14mm	TI-1472-014-S	TI-1472-014	1472-014-S	1472-014
Ø2.7	16mm	TI-1472-016-S	TI-1472-016	1472-016-S	1472-016
Ø2.7	18mm	TI-1472-018-S	TI-1472-018	1472-018-S	1472-018
Ø2.7	20mm	TI-1472-020-S	TI-1472-020	1472-020-S	1472-020
Ø2.7	22mm	TI-1472-022-S	TI-1472-022	1472-022-S	1472-022
Ø2.7	24mm	TI-1472-024-S	TI-1472-024	1472-024-S	1472-024
Ø2.7	26mm	TI-1472-026-S	TI-1472-026	1472-026-S	1472-026
Ø2.7	28mm	TI-1472-028-S	TI-1472-028	1472-028-S	1472-028
Ø2.7	30mm	TI-1472-030-S	TI-1472-030	1472-030-S	1472-030
Ø2.7	32mm	11-017-032TI-S	11-017-032TI	11-017-032SS-S	11-017-032SS
Ø2.7	34mm	11-017-034TI-S	11-017-034TI	11-017-034SS-S	11-017-034SS
Ø2.7	36mm	11-017-036TI-S	11-017-036TI	11-017-036SS-S	11-017-036SS
Ø2.7	38mm	11-017-038TI-S	11-017-038TI	11-017-038SS-S	11-017-038SS
Ø2.7	40mm	11-017-040TI-S	11-017-040TI	11-017-040SS-S	11-017-040SS
Ø2.7	42mm	11-017-042TI-S	11-017-042TI	11-017-042SS-S	11-017-042SS
Ø2.7	44mm	11-017-044TI-S	11-017-044TI	11-017-044SS-S	11-017-044SS
Ø2.7	46mm	11-017-046TI-S	11-017-046TI	11-017-046SS-S	11-017-046SS



2.7mm Wise-Lock Screw, Self-Tapping, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø2.7	8mm	TI-1241.08-S	TI-1241.08	11-012-008SS-S	11-012-008SS
Ø2.7	9mm	TI-1241.09-S	TI-1241.09	11-012-009SS-S	11-012-009SS
Ø2.7	10mm	TI-1241.10-S	TI-1241.10	11-012-010SS-S	11-012-010SS
Ø2.7	12mm	TI-1241.12-S	TI-1241.12	11-012-012SS-S	11-012-012SS
Ø2.7	14mm	TI-1241.14-S	TI-1241.14	11-012-014SS-S	11-012-014SS
Ø2.7	16mm	TI-1241.16-S	TI-1241.16	11-012-016SS-S	11-012-016SS
Ø2.7	18mm	TI-1241.18-S	TI-1241.18	11-012-018SS-S	11-012-018SS
Ø2.7	20mm	TI-1241.20-S	TI-1241.20	11-012-020SS-S	11-012-020SS
Ø2.7	22mm	TI-1241.22-S	TI-1241.22	11-012-022SS-S	11-012-022SS
Ø2.7	24mm	TI-1241.24-S	TI-1241.24	11-012-024SS-S	11-012-024SS
Ø2.7	26mm	TI-1241.26-S	TI-1241.26	11-012-026SS-S	11-012-026SS
Ø2.7	28mm	TI-1241.28-S	TI-1241.28	11-012-028SS-S	11-012-028SS
Ø2.7	30mm	TI-1241.30-S	TI-1241.30	11-012-030SS-S	11-012-030SS



4.0mm Cancellous Screw, Short Thread

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø4.0	10mm	TI-109.010-S	TI-109.010	109.010-S	109.010
Ø4.0	12mm	TI-109.012-S	TI-109.012	109.012-S	109.012
Ø4.0	14mm	TI-109.014-S	TI-109.014	109.014-S	109.014
Ø4.0	16mm	TI-109.016-S	TI-109.016	109.016-S	109.016
Ø4.0	18mm	TI-109.018-S	TI-109.018	109.018-S	109.018
Ø4.0	20mm	TI-109.020-S	TI-109.020	109.020-S	109.020
Ø4.0	22mm	TI-109.022-S	TI-109.022	109.022-S	109.022
Ø4.0	24mm	TI-109.024-S	TI-109.024	109.024-S	109.024
Ø4.0	26mm	TI-109.026-S	TI-109.026	109.026-S	109.026
Ø4.0	28mm	TI-109.028-S	TI-109.028	109.028-S	109.028
Ø4.0	30mm	TI-109.030-S	TI-109.030	109.030-S	109.030
Ø4.0	32mm	TI-109.032-S	TI-109.032	109.032-S	109.032
Ø4.0	34mm	TI-109.034-S	TI-109.034	109.034-S	109.034
Ø4.0	35mm	TI-109.035-S	TI-109.035	109.035-S	109.035
Ø4.0	36mm	TI-109.036-S	TI-109.036	109.036-S	109.036
Ø4.0	38mm	TI-109.038-S	TI-109.038	109.038-S	109.038
Ø4.0	40mm	TI-109.040-S	TI-109.040	109.040-S	109.040
Ø4.0	42mm	TI-109.042-S	TI-109.042	109.042-S	109.042
Ø4.0	44mm	TI-109.044-S	TI-109.044	109.044-S	109.044
Ø4.0	45mm	TI-109.045-S	TI-109.045	109.045-S	109.045



AV-WISELOCK

3.5mm Foot System

Plate Features


- Low-Profile plates with better anatomical Fit
- AV-Wiselock– variable angle of $\pm 15^\circ$ in all directions in each screw hole
- Pre-angled Variable Angle holes for oriented screw placement
- Rounded edges and a smooth surface for soft tissue protection
- Universal Variable angle drill guides available for easy angulation of screws
- K-wire holes to assist with temporary plate fixation



Plates Available in

Titanium as per ISO 5832-3 and Stainless Steel as per ISO 5832-1

3.5mm AV-Wiselock Medial Column Fusion Plate





Titanium

Hole	Direction	Sterile	Non-Sterile
8	Left	11-093-08TIL-S	11-093-08TIL
8	Right	11-093-08TIR-S	11-093-08TIR

Stainless Steel

Hole	Direction	Sterile	Non-Sterile
8	Left	11-093-08SSL-S	11-093-08SSL
8	Right	11-093-08SSR-S	11-093-08SSR

3.5mm AV-Wiselock Medial Column Fusion Plate, Talus Extension




Titanium

Hole	Direction	Sterile	Non-Sterile
9	Left	11-094-09TIL-S	11-094-09TIL
9	Right	11-094-09TIR-S	11-094-09TIR

Stainless Steel



Hole	Direction	Sterile	Non-Sterile
9	Left	11-094-09SSL-S	11-094-09SSL
9	Right	11-094-09SSR-S	11-094-09SSR

3.5mm AV-Wiselock Medial Column Fusion Plantar Plate



Titanium			
Hole	Direction	Sterile	Non-Sterile
7	Left	11-095-07TIL-S	11-095-07TIL
7	Right	11-095-07TIR-S	11-095-07TIR

Stainless Steel			
Hole	Direction	Sterile	Non-Sterile
7	Left	11-095-07SSL-S	11-095-07SSL
7	Right	11-095-07SSR-S	11-095-07SSR

3.5mm AV-Wiselock Calcaneus M Plate


Small	Titanium			
	Hole	Direction	Sterile	Non-Sterile
	12	Left	11-060-01TIL-S	11-060-01TIL
	12	Right	11-060-01TIR-S	11-060-01TIR

Stainless Steel			
Hole	Direction	Sterile	Non-Sterile
12	Left	11-060-01SSL-S	11-060-01SSL
12	Right	11-060-01SSR-S	11-060-01SSR

Medium	Titanium			
	Hole	Direction	Sterile	Non-Sterile
	13	Left	11-060-02TIL-S	11-060-02TIL
	13	Right	11-060-02TIR-S	11-060-02TIR

Stainless Steel			
Hole	Direction	Sterile	Non-Sterile
13	Left	11-060-02SSL-S	11-060-02SSL
13	Right	11-060-02SSR-S	11-060-02SSR

Large	Titanium			
	Hole	Direction	Sterile	Non-Sterile
	13	Left	11-060-03TIL-S	11-060-03TIL
	13	Right	11-060-03TIR-S	11-060-03TIR

Stainless Steel			
Hole	Direction	Sterile	Non-Sterile
13	Left	11-060-03SSL-S	11-060-03SSL
13	Right	11-060-03SSR-S	11-060-03SSR

AV-WISELOCK

PHEELOS - 3.5mm Small Fragment System

FEATURES

- Proximal screw holes for AV-Wiselock screws Dia 3.5 mm enable an angular stable construct to enhance the grip in osteoporotic bone and multi-fragment fractures
- Carefully apply for osteoporotic bone
- Optimal screw placement
- Proximal holes for suturing to help maintain fracture reduction
- Distal Locking long holes for maximum adaptability

INDICATIONS


- Dislocated two, three and four fragment fractures of the proximal humerus, including fractures involving osteopenic bone
- Pseudarthroses in the proximal humerus
- Osteotomies in the proximal humerus



Plates Available in

Titanium as per ISO 5832-3 and Stainless Steel as per ISO 5832-1



PHEELOS - 3.5mm AV-Wiselock Proximal Humerus Plate With Suture Cuts

Short

Titanium

Hole	Sterile	Non-Sterile
2	TI-11-075-02-S	TI-11-075-02
3	TI-11-075-03-S	TI-11-075-03
4	TI-11-075-04-S	TI-11-075-04
5	TI-11-075-05-S	TI-11-075-05

Stainless Steel

Hole	Sterile	Non-Sterile
2	11-075-02SS-S	11-075-02SS
3	11-075-03SS-S	11-075-03SS
4	11-075-04SS-S	11-075-04SS
5	11-075-05SS-S	11-075-05SS


Long

Titanium

Hole	Sterile	Non-Sterile
5	TI-11-076-05-S	TI-11-076-05
6	TI-11-076-06-S	TI-11-076-06
8	TI-11-076-08-S	TI-11-076-08
10	TI-11-076-10-S	TI-11-076-10
12	TI-11-076-12-S	TI-11-076-12

Stainless Steel

Hole	Sterile	Non-Sterile
5	11-076-05SS-S	11-076-05SS
6	11-076-06SS-S	11-076-06SS
8	11-076-08SS-S	11-076-08SS
10	11-076-10SS-S	11-076-10SS
12	11-076-12SS-S	11-076-12SS


3.5mm AV-Wiselock Medial Proximal Tibia Plate


Titanium

Hole	Direction	Sterile	Non-Sterile
4	Left	11-080-04TIL-S	11-080-04TIL
4	Right	11-080-04TIR-S	11-080-04TIR
6	Left	11-080-06TIL-S	11-080-06TIL
6	Right	11-080-06TIR-S	11-080-06TIR
8	Left	11-080-08TIL-S	11-080-08TIL
8	Right	11-080-08TIR-S	11-080-08TIR
10	Left	11-080-10TIL-S	11-080-10TIL
10	Right	11-080-10TIR-S	11-080-10TIR
12	Left	11-080-12TIL-S	11-080-12TIL
12	Right	11-080-12TIR-S	11-080-12TIR
14	Left	11-080-14TIL-S	11-080-14TIL
14	Right	11-080-14TIR-S	11-080-14TIR
16	Left	11-080-16TIL-S	11-080-16TIL
16	Right	11-080-16TIR-S	11-080-16TIR
18	Left	11-080-18TIL-S	11-080-18TIL
18	Right	11-080-18TIR-S	11-080-18TIR
20	Left	11-080-20TIL-S	11-080-20TIL
20	Right	11-080-20TIR-S	11-080-20TIR

Stainless Steel

Hole	Direction	Sterile	Non-Sterile
4	Left	11-080-04SSL-S	11-080-04SSL
4	Right	11-080-04SSR-S	11-080-04SSR
6	Left	11-080-06SSL-S	11-080-06SSL
6	Right	11-080-06SSR-S	11-080-06SSR
8	Left	11-080-08SSL-S	11-080-08SSL
8	Right	11-080-08SSR-S	11-080-08SSR
10	Left	11-080-10SSL-S	11-080-10SSL
10	Right	11-080-10SSR-S	11-080-10SSR
12	Left	11-080-12SSL-S	11-080-12SSL
12	Right	11-080-12SSR-S	11-080-12SSR
14	Left	11-080-14SSL-S	11-080-14SSL
14	Right	11-080-14SSR-S	11-080-14SSR
16	Left	11-080-16SSL-S	11-080-16SSL
16	Right	11-080-16SSR-S	11-080-16SSR
18	Left	11-080-18SSL-S	11-080-18SSL
18	Right	11-080-18SSR-S	11-080-18SSR
20	Left	11-080-20SSL-S	11-080-20SSL
20	Right	11-080-20SSR-S	11-080-20SSR


3.5mm AV-Wiselock Extra-Articular Distal Humerus Plate


Titanium

Hole	Direction	Sterile	Non-Sterile
4	Left	11-096-04TIL-S	11-096-04TIL
4	Right	11-096-04TIR-S	11-096-04TIR
6	Left	11-096-06TIL-S	11-096-06TIL
6	Right	11-096-06TIR-S	11-096-06TIR
8	Left	11-096-08TIL-S	11-096-08TIL
8	Right	11-096-08TIR-S	11-096-08TIR
10	Left	11-096-10TIL-S	11-096-10TIL
10	Right	11-096-10TIR-S	11-096-10TIR
12	Left	11-096-12TIL-S	11-096-12TIL
12	Right	11-096-12TIR-S	11-096-12TIR
14	Left	11-096-14TIL-S	11-096-14TIL
14	Right	11-096-14TIR-S	11-096-14TIR

Stainless Steel

Hole	Direction	Sterile	Non-Sterile
4	Left	11-096-04SSL-S	11-096-04SSL
4	Right	11-096-04SSR-S	11-096-04SSR
6	Left	11-096-06SSL-S	11-096-06SSL
6	Right	11-096-06SSR-S	11-096-06SSR
8	Left	11-096-08SSL-S	11-096-08SSL
8	Right	11-096-08SSR-S	11-096-08SSR
10	Left	11-096-10SSL-S	11-096-10SSL
10	Right	11-096-10SSR-S	11-096-10SSR
12	Left	11-096-12SSL-S	11-096-12SSL
12	Right	11-096-12SSR-S	11-096-12SSR
14	Left	11-096-14SSL-S	11-096-14SSL
14	Right	11-096-14SSR-S	11-096-14SSR


2.7/3.5mm AV Wiselock Superior Anterior Clavicle Plate with Lateral Extension


Titanium			
Hole	Direction	Sterile	Non-Sterile
3	Left	11-087-03LTI-S	11-087-03LTI
3	Right	11-087-03RTI-S	11-087-03RTI
4	Left	11-087-04LTI-S	11-087-04LTI
4	Right	11-087-04RTI-S	11-087-04RTI
5	Left	11-087-05LTI-S	11-087-05LTI
5	Right	11-087-05RTI-S	11-087-05RTI
6	Left	11-087-06LTI-S	11-087-06LTI
6	Right	11-087-06RTI-S	11-087-06RTI
7	Left	11-087-07LTI-S	11-087-07LTI
7	Right	11-087-07RTI-S	11-087-07RTI
8	Left	11-087-08LTI-S	11-087-08LTI
8	Right	11-087-08RTI-S	11-087-08RTI

Stainless Steel			
Hole	Direction	Sterile	Non-Sterile
3	Left	11-087-03LSS-S	11-087-03LSS
3	Right	11-087-03RSS-S	11-087-03RSS
4	Left	11-087-04LSS-S	11-087-04LSS
4	Right	11-087-04RSS-S	11-087-04RSS
5	Left	11-087-05LSS-S	11-087-05LSS
5	Right	11-087-05RSS-S	11-087-05RSS
6	Left	11-087-06LSS-S	11-087-06LSS
6	Right	11-087-06RSS-S	11-087-06RSS
7	Left	11-087-07LSS-S	11-087-07LSS
7	Right	11-087-07RSS-S	11-087-07RSS
8	Left	11-087-08LSS-S	11-087-08LSS
8	Right	11-087-08RSS-S	11-087-08RSS


3.5mm AV-Wiselock Screw, Self-Tapping

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø3.5	10mm	11-024-010TI-S	11-024-010TI	11-024-010SS-S	11-024-010SS
Ø3.5	12mm	11-024-012TI-S	11-024-012TI	11-024-012SS-S	11-024-012SS
Ø3.5	14mm	11-024-014TI-S	11-024-014TI	11-024-014SS-S	11-024-014SS
Ø3.5	16mm	11-024-016TI-S	11-024-016TI	11-024-016SS-S	11-024-016SS
Ø3.5	18mm	11-024-018TI-S	11-024-018TI	11-024-018SS-S	11-024-018SS
Ø3.5	20mm	11-024-020TI-S	11-024-020TI	11-024-020SS-S	11-024-020SS
Ø3.5	22mm	11-024-022TI-S	11-024-022TI	11-024-022SS-S	11-024-022SS
Ø3.5	24mm	11-024-024TI-S	11-024-024TI	11-024-024SS-S	11-024-024SS
Ø3.5	26mm	11-024-026TI-S	11-024-026TI	11-024-026SS-S	11-024-026SS
Ø3.5	28mm	11-024-028TI-S	11-024-028TI	11-024-028SS-S	11-024-028SS
Ø3.5	30mm	11-024-030TI-S	11-024-030TI	11-024-030SS-S	11-024-030SS
Ø3.5	32mm	11-024-032TI-S	11-024-032TI	11-024-032SS-S	11-024-032SS
Ø3.5	34mm	11-024-034TI-S	11-024-034TI	11-024-034SS-S	11-024-034SS
Ø3.5	36mm	11-024-036TI-S	11-024-036TI	11-024-036SS-S	11-024-036SS
Ø3.5	38mm	11-024-038TI-S	11-024-038TI	11-024-038SS-S	11-024-038SS
Ø3.5	40mm	11-024-040TI-S	11-024-040TI	11-024-040SS-S	11-024-040SS
Ø3.5	42mm	11-024-042TI-S	11-024-042TI	11-024-042SS-S	11-024-042SS
Ø3.5	44mm	11-024-044TI-S	11-024-044TI	11-024-044SS-S	11-024-044SS
Ø3.5	46mm	11-024-046TI-S	11-024-046TI	11-024-046SS-S	11-024-046SS
Ø3.5	48mm	11-024-048TI-S	11-024-048TI	11-024-048SS-S	11-024-048SS
Ø3.5	50mm	11-024-050TI-S	11-024-050TI	11-024-050SS-S	11-024-050SS
Ø3.5	52mm	11-024-052TI-S	11-024-052TI	11-024-052SS-S	11-024-052SS
Ø3.5	54mm	11-024-054TI-S	11-024-054TI	11-024-054SS-S	11-024-054SS
Ø3.5	56mm	11-024-056TI-S	11-024-056TI	11-024-056SS-S	11-024-056SS
Ø3.5	58mm	11-024-058TI-S	11-024-058TI	11-024-058SS-S	11-024-058SS
Ø3.5	60mm	11-024-060TI-S	11-024-060TI	11-024-060SS-S	11-024-060SS
Ø3.5	65mm	11-024-065TI-S	11-024-065TI	11-024-065SS-S	11-024-065SS
Ø3.5	70mm	11-024-070TI-S	11-024-070TI	11-024-070SS-S	11-024-070SS
Ø3.5	75mm	11-024-075TI-S	11-024-075TI	11-024-075SS-S	11-024-075SS
Ø3.5	80mm	11-024-080TI-S	11-024-080TI	11-024-080SS-S	11-024-080SS
Ø3.5	85mm	11-024-085TI-S	11-024-085TI	11-024-085SS-S	11-024-085SS
Ø3.5	90mm	11-024-090TI-S	11-024-090TI	11-024-090SS-S	11-024-090SS
Ø3.5	95mm	11-024-095TI-S	11-024-095TI	11-024-095SS-S	11-024-095SS



3.5mm Wise-Lock Screw, Self Tapping, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø3.5	10mm	TI-1473-010-S	TI-1473-010	1473-010-S	1473-010
Ø3.5	12mm	TI-1473-012-S	TI-1473-012	1473-012-S	1473-012
Ø3.5	14mm	TI-1473-014-S	TI-1473-014	1473-014-S	1473-014
Ø3.5	16mm	TI-1473-016-S	TI-1473-016	1473-016-S	1473-016
Ø3.5	18mm	TI-1473-018-S	TI-1473-018	1473-018-S	1473-018
Ø3.5	20mm	TI-1473-020-S	TI-1473-020	1473-020-S	1473-020
Ø3.5	22mm	TI-1473-022-S	TI-1473-022	1473-022-S	1473-022
Ø3.5	24mm	TI-1473-024-S	TI-1473-024	1473-024-S	1473-024
Ø3.5	26mm	TI-1473-026-S	TI-1473-026	1473-026-S	1473-026
Ø3.5	28mm	TI-1473-028-S	TI-1473-028	1473-028-S	1473-028
Ø3.5	30mm	TI-1473-030-S	TI-1473-030	1473-030-S	1473-030
Ø3.5	32mm	TI-1473-032-S	TI-1473-032	1473-032-S	1473-032
Ø3.5	34mm	TI-1473-034-S	TI-1473-034	1473-034-S	1473-034
Ø3.5	36mm	TI-1473-036-S	TI-1473-036	1473-036-S	1473-036
Ø3.5	38mm	TI-1473-038-S	TI-1473-038	1473-038-S	1473-038
Ø3.5	40mm	TI-1473-040-S	TI-1473-040	1473-040-S	1473-040
Ø3.5	42mm	TI-1473-042-S	TI-1473-042	1473-042-S	1473-042
Ø3.5	44mm	TI-1473-044-S	TI-1473-044	1473-044-S	1473-044
Ø3.5	46mm	TI-1473-046-S	TI-1473-046	1473-046-S	1473-046
Ø3.5	48mm	TI-1473-048-S	TI-1473-048	1473-048-S	1473-048
Ø3.5	50mm	TI-1473-050-S	TI-1473-050	1473-050-S	1473-050
Ø3.5	52mm	TI-1473-052-S	TI-1473-052	1473-052-S	1473-052
Ø3.5	54mm	TI-1473-054-S	TI-1473-054	1473-054-S	1473-054
Ø3.5	56mm	TI-1473-056-S	TI-1473-056	1473-056-S	1473-056
Ø3.5	58mm	TI-1473-058-S	TI-1473-058	1473-058-S	1473-058
Ø3.5	60mm	TI-1473-060-S	TI-1473-060	1473-060-S	1473-060
Ø3.5	65mm	TI-1473-065-S	TI-1473-065	1473-065-S	1473-065
Ø3.5	70mm	TI-1473-070-S	TI-1473-070	1473-070-S	1473-070
Ø3.5	75mm	TI-1473-075-S	TI-1473-075	1473-075-S	1473-075
Ø3.5	80mm	TI-1473-080-S	TI-1473-080	1473-080-S	1473-080
Ø3.5	85mm	TI-1473-085-S	TI-1473-085	1473-085-S	1473-085
Ø3.5	90mm	TI-1473-090-S	TI-1473-090	1473-090-S	1473-090
Ø3.5	95mm	TI-1473-095-S	TI-1473-095	1473-095-S	1473-095



3.5mm Cortical Screw, Self Tapping, (Star Head)


Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø3.5	10mm	11-002-010TI-S	11-002-010TI	11-002-010SS-S	11-002-010SS
Ø3.5	12mm	11-002-012TI-S	11-002-012TI	11-002-012SS-S	11-002-012SS
Ø3.5	14mm	11-002-014TI-S	11-002-014TI	11-002-014SS-S	11-002-014SS
Ø3.5	16mm	11-002-016TI-S	11-002-016TI	11-002-016SS-S	11-002-016SS
Ø3.5	18mm	11-002-018TI-S	11-002-018TI	11-002-018SS-S	11-002-018SS
Ø3.5	20mm	11-002-020TI-S	11-002-020TI	11-002-020SS-S	11-002-020SS
Ø3.5	22mm	11-002-022TI-S	11-002-022TI	11-002-022SS-S	11-002-022SS
Ø3.5	24mm	11-002-024TI-S	11-002-024TI	11-002-024SS-S	11-002-024SS
Ø3.5	26mm	11-002-026TI-S	11-002-026TI	11-002-026SS-S	11-002-026SS
Ø3.5	28mm	11-002-028TI-S	11-002-028TI	11-002-028SS-S	11-002-028SS
Ø3.5	30mm	11-002-030TI-S	11-002-030TI	11-002-030SS-S	11-002-030SS
Ø3.5	32mm	11-002-032TI-S	11-002-032TI	11-002-032SS-S	11-002-032SS
Ø3.5	34mm	11-002-034TI-S	11-002-034TI	11-002-034SS-S	11-002-034SS
Ø3.5	36mm	11-002-036TI-S	11-002-036TI	11-002-036SS-S	11-002-036SS
Ø3.5	38mm	11-002-038TI-S	11-002-038TI	11-002-038SS-S	11-002-038SS
Ø3.5	40mm	11-002-040TI-S	11-002-040TI	11-002-040SS-S	11-002-040SS
Ø3.5	42mm	11-002-042TI-S	11-002-042TI	11-002-042SS-S	11-002-042SS
Ø3.5	44mm	11-002-044TI-S	11-002-044TI	11-002-044SS-S	11-002-044SS
Ø3.5	46mm	11-002-046TI-S	11-002-046TI	11-002-046SS-S	11-002-046SS
Ø3.5	48mm	11-002-048TI-S	11-002-048TI	11-002-048SS-S	11-002-048SS
Ø3.5	50mm	11-002-050TI-S	11-002-050TI	11-002-050SS-S	11-002-050SS
Ø3.5	55mm	11-002-055TI-S	11-002-055TI	11-002-055SS-S	11-002-055SS
Ø3.5	60mm	11-002-060TI-S	11-002-060TI	11-002-060SS-S	11-002-060SS
Ø3.5	65mm	11-002-065TI-S	11-002-065TI	11-002-065SS-S	11-002-065SS
Ø3.5	70mm	11-002-070TI-S	11-002-070TI	11-002-070SS-S	11-002-070SS
Ø3.5	75mm	11-002-075TI-S	11-002-075TI	11-002-075SS-S	11-002-075SS
Ø3.5	80mm	11-002-080TI-S	11-002-080TI	11-002-080SS-S	11-002-080SS
Ø3.5	85mm	11-002-085TI-S	11-002-085TI	11-002-085SS-S	11-002-085SS
Ø3.5	90mm	11-002-090TI-S	11-002-090TI	11-002-090SS-S	11-002-090SS
Ø3.5	95mm	11-002-095TI-S	11-002-095TI	11-002-095SS-S	11-002-095SS



7-012-01 Conical Extraction Screw, Ø3.5mm, for 2.7/3.5mm AV-Wiselock System



7-012-02 Drill Bit with Quick Coupling with Stopper, Ø2.8mm x Length 165mm, for 2.7/3.5mm AV-Wiselock System



7-012-03 Star Screwdriver Shaft, T8, for 2.7/3.5mm AV-Wiselock System



7-012-04 Depth Gauge measuring upto 60mm for 2.7/3.5mm AV-Wiselock System



7-012-05 Drill Bit with Quick Coupling, Ø2.0mm x Length 140mm, for 2.7/3.5mm AV-Wiselock System



7-012-06 Bone Tap Quick Coupling for Cortical Screws, Ø3.5mm, for 2.7/3.5mm AV-Wiselock System



7-012-07 Screwdriver Shaft, T15, for 2.7/3.5mm AV-Wiselock System



7-012-08 Threaded Drill Sleeve for Ø2.8mm Drill Bit - 2.7/3.5mm AV-Wiselock System



7-012-09 Self-Centering Double Drill Guide, Ø2.5/3.5mm, for 2.7/3.5mm AV-Wiselock System



7-012-10 Torque Limiting Attachment, 1.2Nm, for 2.7/3.5mm AV-Wiselock System



7-012-11 Torque Limiting Attachment, 1.5Nm, for 2.7/3.5mm AV-Wiselock System



7-012-12 Drill Bit with Quick Coupling End, Ø2.5mm x Length 180mm - 2.7/3.5mm AV-Wiselock System



- 7-012-13** Drill Bit with Quick Coupling End, Ø3.5mm x Length 110mm - 2.7/3.5mm AV-Wiselock System



- 7-012-14** 2.0mm Variable Angle Guide for 2.7/3.5mm AV-Wiselock System



- 7-012-15** Self-Centering Double Drill Guide, Ø2.0/2.7mm, for 2.7/3.5mm AV-Wiselock System



- 7-012-16** 2.8mm Variable Angle Double Drill Guide



- 7-012-17** Handle with Quick Coupling for 2.7/3.5mm AV-Wiselock System



- 7-012-18** Handle for Torque Limiting Attachment, for 2.7/3.5mm AV-Wiselock System



7-012-19 2.5Nm Torque Limiting Handle for 2.7/3.5mm AV-Wiselock System



7-012-20 Depth Gauge measuring upto 100mm, for 2.7/3.5mm AV-Wiselock System



7-012-21 Percutaneous Depth Gauge Adaptor for 3.5mm Screws - 2.7/3.5mm AV-Wiselock System



7-012-22 Kirschner Wire with One Sided Trocar Tip Ø1.6mm x Length 150mm



7-012-23 Periarticular Reduction Forcep, Small, for 2.7/3.5mm AV-Wiselock System



7-012-24 Periarticular Reduction Forcep, Medium, for 2.7/3.5mm AV-Wiselock System




7-012-25 Large Compression Forcep with Speed Lock for 2.7/3.5mm AV-Wiselock System



7-012-26 Large Distraction Forcep with Speed Lock for 2.7/3.5mm AV-Wiselock System



7-012-27 3.5mm Post for AV-Wiselock Locking Hole



7-012-28 2.8mm Compression Wire, 10mm Thread for 2.7/3.5mm AV-Wiselock System



7-012-29 2.8mm Compression Wire, 15mm Thread for 2.7/3.5mm AV-Wiselock System



7-012-30 2.8mm Compression Wire, 20mm Thread for 2.7/3.5mm AV-Wiselock System



7-012-31 2.8mm Compression Wire, 25mm Thread for 2.7/3.5mm AV-Wiselock System



7-012-32 2.8mm Compression Wire, 30mm Thread for 2.7/3.5mm AV-Wiselock System



7-012-33 2.8mm Compression Wire, 35mm Thread for 2.7/3.5mm AV-Wiselock System



7-012-34 2.8mm Compression Wire, 40mm Thread for 2.7/3.5mm AV-Wiselock System



7-012-35 2.8mm Compression Wire, 45mm Thread for 2.7/3.5mm AV-Wiselock System



7-012-36 2.8mm Compression Wire, 50mm Thread for 2.7/3.5mm AV-Wiselock System



7-012-37 2.8mm Compression Wire, 55mm Thread for 2.7/3.5mm AV-Wiselock System



7-012-38 2.8mm Compression Wire, 60mm Thread for 2.7/3.5mm AV-Wiselock System




7-012-39 Guide Wire, Ø2.0mm x Length 150mm



7-012-40 Guide Wire with Threaded Tip, Ø2.0mm x Thread Length 10mm x Length 150mm



7-012-41 Threaded Drill Sleeve for Drill Bit 2.0mm - 2.7/3.5mm AV-Wiselock System



7-012-42 Tap, HA 2.7mm for 2.7/3.5mm AV-Wiselock System



7-012-43 Screwdriver Star Head, T8, for 2.7/3.5mm AV-Wiselock System



7-012-44 T-Handle with Quick Coupling for 2.7/3.5mm AV-Wiselock System



7-012-45 Self-Centering Bone Holding Forcep, Speed Lock, 190mm for 2.7/3.5mm AV-Wiselock System



7-012-46 Reduction Forcep, Serrated Jaws, Speed Lock, 160mm for 2.7/3.5mm AV-Wiselock System



7-012-47 Reduction Forcep, Pointed, Ratchet Lock, 180mm for 2.7/3.5mm AV-Wiselock System



7-012-48 Bending Iron 1 for 2.7/3.5mm AV-Wiselock System



7-012-49 Bending Iron 2 for 2.7/3.5mm AV-Wiselock System



7-012-50 Ø1.2mm Kirschner Wire for 2.7/3.5mm AV-Wiselock System



7-012-51 Hohmann Retractor, 6.5mm for 2.7/3.5mm AV-Wiselock System



7-012-52 Hohmann Retractor, 8.5mm for 2.7/3.5mm AV-Wiselock System



7-012-53 Hohmann Retractor, 15.5mm for 2.7/3.5mm AV-Wiselock System



7-012-54 Periosteal Elevator with Silicon Handle, Curved, 6mm for 2.7/3.5mm AV-Wiselock System



7-012-55 Screw Holding Forcep for 2.7/3.5mm AV-Wiselock System



7-012-56 Countersink for 2.7/3.5mm AV-Wiselock System



7-012-57 Screw Holding Sleeve for 2.7mm Cortical Screw



7-012-60 Screw Holding Sleeve for 2.7mm Locking Screw



7-012-61 Screw Holding Sleeve for 3.5mm Cortical Screw




7-012-62 Screw Holding Sleeve for 3.5mm Locking Screw




7-050-05


T- Handle for Threaded Drill Sleeve


7-012-58

Instrument Trays for 2.7/3.5mm AV-Wiselock System


7-012-59

Container for 2.7/3.5mm AV-Wiselock System Instrument Set



7-158-01 Pheelos Aiming Device with Nut for AV-Wiselock System



3300-02 Centering Sleeve, for Pheelos Aiming Device



7-158-02 Drill Sleeve, for Pheelos Aiming Device for AV-Wiselock System



3300-04 Centering Sleeve for Kirschner Wire Ø 1.6mm



3300-05 Pheelos Direct Measuring Device



451-1.6-150 Kirschner Wire with Trocar Tip Both End Ø1.6mm x Length 150mm



7-158 PHEELOS Instrument Set for AV-Wiselock System

Codes	Set consisting of:	Units
7-158-01	Pheelos Aiming Device with Nut for AV-Wiselock System	1
7-158-02	Drill Sleeve, for Pheelos Aiming Device for AV-Wiselock System	1
3300-02	Centering Sleeve, for Pheelos Aiming Device	1
3300-04	Centering Sleeve for Kirschner Wire Ø 1.6mm	1
3300-05	Pheelos Direct Measuring Device	1
451-1.6-150	Kirschner Wire with Trocar Tip Both End Ø1.6mm x Length 150mm	1

7-012 : 2.7/3.5mm AV-Wiselock System Instrument Set



7-012 : 2.7/3.5mm AV-Wiselock System Instrument Set

Code	Set Consisting of	Units
7-012-01	Conical Extraction Screw, Ø3.5mm, for 2.7/3.5mm AV-Wiselock System	1
7-012-02	Drill Bit with Quick Coupling with Stopper, Ø2.8mm x Length 165mm, for 2.7/3.5mm AV-Wiselock System	2
7-012-03	Star Screwdriver Shaft, T8, for 2.7/3.5mm AV-Wiselock System	2
7-012-04	Depth Gauge measuring upto 60mm for 2.7/3.5mm AV-Wiselock System	1
7-012-05	Drill Bit with Quick Coupling, Ø2.0mm x Length 140mm, for 2.7/3.5mm AV-Wiselock System	2
7-012-06	Bone Tap Quick Coupling for Cortical Screws, Ø3.5mm, for 2.7/3.5mm AV-Wiselock System	1
7-012-07	Screwdriver Shaft, T15, for 2.7/3.5mm AV-Wiselock System	2
7-012-08	Threaded Drill Sleeve for Ø2.8mm Drill Bit - 2.7/3.5mm AV-Wiselock System	3
7-012-09	Self-Centering Double Drill Guide, Ø2.5/3.5mm, for 2.7/3.5mm AV-Wiselock System	1
7-012-10	Torque Limiting Attachment, 1.2Nm, for 2.7/3.5mm AV-Wiselock System	1
7-012-11	Torque Limiting Attachment, 1.5Nm, for 2.7/3.5mm AV-Wiselock System	1
7-012-12	Drill Bit with Quick Coupling End, Ø2.5mm x Length 180mm - 2.7/3.5mm AV-Wiselock System	2
7-012-13	Drill Bit with Quick Coupling End, Ø3.5mm x Length 110mm - 2.7/3.5mm AV-Wiselock System	1
7-012-14	2.0mm Variable Angle Guide for 2.7/3.5mm AV-Wiselock System	1
7-012-15	Self-Centering Double Drill Guide, Ø2.0/2.7mm, for 2.7/3.5mm AV-Wiselock System	1
7-012-16	2.8mm Variable Angle Double Drill Guide	1
7-012-17	Handle with Quick Coupling for 2.7/3.5mm AV-Wiselock System	1
7-012-18	Handle for Torque Limiting Attachment, for 2.7/3.5mm AV-Wiselock System	1
7-012-19	2.5Nm Torque Limiting Handle for 2.7/3.5mm AV-Wiselock System	1
7-012-20	Depth Gauge measuring upto 100mm, for 2.7/3.5mm AV-Wiselock System	1
7-012-21	Percutaneous Depth Gauge Adaptor for 3.5mm Screws - 2.7/3.5mm AV-Wiselock System	1
7-012-22	Kirschner Wire with One Sided Trocar Tip Ø1.6mm x Length 150mm	5
7-012-23	Periarticular Reduction Forcep, Small, for 2.7/3.5mm AV-Wiselock System	1
7-012-24	Periarticular Reduction Forcep, Medium, for 2.7/3.5mm AV-Wiselock System	1
7-012-25	Large Compression Forcep with Speed Lock for 2.7/3.5mm AV-Wiselock System	1
7-012-26	Large Distraction Forcep with Speed Lock for 2.7/3.5mm AV-Wiselock System	1
7-012-27	3.5mm Post for AV-Wiselock Locking Hole	2
7-012-28	2.8mm Compression Wire, 10mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-29	2.8mm Compression Wire, 15mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-30	2.8mm Compression Wire, 20mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-31	2.8mm Compression Wire, 25mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-32	2.8mm Compression Wire, 30mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-33	2.8mm Compression Wire, 35mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-34	2.8mm Compression Wire, 40mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-35	2.8mm Compression Wire, 45mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-36	2.8mm Compression Wire, 50mm Thread for 2.7/3.5mm AV-Wiselock System	2

Code	Set Consisting of	Units
7-012-37	2.8mm Compression Wire, 55mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-38	2.8mm Compression Wire, 60mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-39	Guide Wire, Ø2.0mm x Length 150mm	5
7-012-40	Guide Wire with Threaded Tip, Ø2.0mm x Thread Length 10mm x Length 150mm	2
7-012-41	Threaded Drill Sleeve for Drill Bit 2.0mm - 2.7/3.5mm AV-Wiselock System	3
7-012-42	Tap, HA 2.7mm for 2.7/3.5mm AV-Wiselock System	1
7-012-43	Screwdriver Star Head, T8, for 2.7/3.5mm AV-Wiselock System	1
7-012-44	T-Handle with Quick Coupling for 2.7/3.5mm AV-Wiselock System	1
7-012-45	Self-Centering Bone Holding Forcep, Speed Lock, 190mm for 2.7/3.5mm AV-Wiselock System	2
7-012-46	Reduction Forcep, Serrated Jaws, Speed Lock, 160mm for 2.7/3.5mm AV-Wiselock System	1
7-012-47	Reduction Forcep, Pointed, Ratchet Lock, 180mm for 2.7/3.5mm AV-Wiselock System	1
7-012-48	Bending Iron 1 for 2.7/3.5mm AV-Wiselock System	1
7-012-49	Bending Iron 2 for 2.7/3.5mm AV-Wiselock System	1
7-012-50	Ø1.2mm Kirschner Wire for 2.7/3.5mm AV-Wiselock System	3
7-012-51	Hohmann Retractor, 6.5mm for 2.7/3.5mm AV-Wiselock System	1
7-012-52	Hohmann Retractor, 8.5mm for 2.7/3.5mm AV-Wiselock System	1
7-012-53	Hohmann Retractor, 15.5mm for 2.7/3.5mm AV-Wiselock System	2
7-012-54	Periosteal Elevator with Silicon Handle, Curved, 6mm for 2.7/3.5mm AV-Wiselock System	1
7-012-55	Screw Holding Forcep for 2.7/3.5mm AV-Wiselock System	1
7-012-56	Countersink for 2.7/3.5mm AV-Wiselock System	1
7-012-57	Screw Holding Sleeve for 2.7mm Cortical Screw	1
7-012-60	Screw Holding Sleeve for 2.7mm Locking Screw	1
7-012-61	Screw Holding Sleeve for 3.5mm Cortical Screw	1
7-012-62	Screw Holding Sleeve for 3.5mm Locking Screw	1
7-050-05	T- Handle for Threaded Drill Sleeve	1
7-012-58	Instrument Trays for 2.7/3.5mm AV-Wiselock System	3
7-012-59	Container for 2.7/3.5mm AV-Wiselock System Instrument Set	1

**USA**

Auxein Inc.
1500 Nw 89th Court, Suite 107-108
Doral, Florida 33172
Tel: +1 305 395 6062
E Fax: +1 305 395 6262
Email: USoffice@auxein.com

MEXICO

Auxein México, S.A. de C.V.
Tepic 139 int 801, Colonia Roma Sur,
Alcaldía Cuauhtémoc, CDMX,
México, C.P. 06760
Tel: +521 55 7261 0318
Email: info@auxein.mx

INDIA

Auxein Medical Pvt. Ltd.
Plot No. 168-169-170, Phase-4,
Kundli Industrial Area,
HSIIDC, Sector-57, Sonapat - 131028, Haryana
Tel: +91 99106 43638 | Fax: +91 86077 70197
Email: info@auxein.com