


AV-WISELOCK™

Distal Radius and Volar Rim
System - Locking



about us

Auxein Medical is an integrated, research based, orthopaedic Implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in india and serves customers in over 75 countries worldwide.

Our Achievements



AV-WISELOCK

2.4mm Distal Radius System

Plate Features

- 2.4mm Variable Angle Distal Radius System are used for displaced extra-articular and intra-articular distal radius fractures and corrective osteotomies of the distal radius.
- 2.4mm AV-Wise-Lock Distal Radius System are used for displaced extra-articular and intra-articular distal radius fractures and corrective osteotomies of the distal radius.

Plate Features

- AV-Wiselock Variable Angle of $\pm 15^\circ$ in all directions in each screw hole
- Pre-angled Variable Angle holes for oriented screw placement
- Rounded edges and a smooth surface for soft tissue protection
- Oblong hole for variable positioning of the plate
- Universal Variable angle drill guides available for rapid and easy angulation of screws
- K-wire holes to assist with temporary plate fixation



Anatomical Contoured
Buttrressing of distal fragments
due to anatomically contoured
plate.



Guiding Block
Allows guided drilling and screw
insertion in the predefined
nominal angle





**Oblong AV-Wiselock
Combi Hole**
Plate positioning on the bone

Plates Available in

Titanium as per ISO 5832-3 and Stainless Steel as per ISO 5832-1

2.4mm AV-Wiselock Distal Radius Volar Plate, Extra-articular



5 Head Holes

Titanium		
Holes	Sterile	Non-Sterile
3	11-043-03LTI-S	11-043-03LTI
3	11-043-03RTI-S	11-043-03RTI
5	11-043-05LTI-S	11-043-05LTI
5	11-043-05RTI-S	11-043-05RTI

Stainless Steel		
Holes	Sterile	Non-Sterile
3	11-043-03LSS-S	11-043-03LSS
3	11-043-03RSS-S	11-043-03RSS
5	11-043-05LSS-S	11-043-05LSS
5	11-043-05RSS-S	11-043-05RSS

4 Head Holes

Titanium		
Holes	Sterile	Non-Sterile
3	11-044-03LTI-S	11-044-03LTI
3	11-044-03RTI-S	11-044-03RTI
5	11-044-05LTI-S	11-044-05LTI
5	11-044-05RTI-S	11-044-05RTI

Stainless Steel		
Holes	Sterile	Non-Sterile
3	11-044-03LSS-S	11-044-03LSS
3	11-044-03RSS-S	11-044-03RSS
5	11-044-05LSS-S	11-044-05LSS
5	11-044-05RSS-S	11-044-05RSS


2.4mm AV-Wiselock Dorsal Distal Radius Plate, Radial Column



Titanium		
Holes	Sterile	Non-Sterile
5	11-045-05TI-S	11-045-05TI
6	11-045-06TI-S	11-045-06TI

Stainless Steel		
Holes	Sterile	Non-Sterile
5	11-045-05SS-S	11-045-05SS
6	11-045-06SS-S	11-045-06SS


2.4mm AV-Wiselock Dorsal Distal Radius Plate, Intermediate Column, (2 Head Holes)



Titanium			
Holes	Direction	Sterile	Non-Sterile
3	Left	11-046-03LTI-S	11-046-03LTI
3	Right	11-046-03RTI-S	11-046-03RTI
4	Left	11-046-04LTI-S	11-046-04LTI
4	Right	11-046-04RTI-S	11-046-04RTI

Stainless Steel			
Holes	Direction	Sterile	Non-Sterile
3	Left	11-046-03LSS-S	11-046-03LSS
3	Right	11-046-03RSS-S	11-046-03RSS
4	Left	11-046-04LSS-S	11-046-04LSS
4	Right	11-046-04RSS-S	11-046-04RSS

2.4mm AV-Wiselock "L" Dorsal Distal Radius Plate



2 Head Holes

Titanium			
Holes	Direction	Sterile	Non-Sterile
3	Left	11-047-03LTI-S	11-047-03LTI
3	Right	11-047-03RTI-S	11-047-03RTI
5	Left	11-047-05LTI-S	11-047-05LTI
5	Right	11-047-05RTI-S	11-047-05RTI

Stainless Steel			
Holes	Direction	Sterile	Non-Sterile
3	Left	11-047-03LSS-S	11-047-03LSS
3	Right	11-047-03RSS-S	11-047-03RSS
5	Left	11-047-05LSS-S	11-047-05LSS
5	Right	11-047-05RSS-S	11-047-05RSS

3 Head Holes

Titanium			
Holes	Direction	Sterile	Non-Sterile
3	Left	11-048-03LTI-S	11-048-03LTI
3	Right	11-048-03RTI-S	11-048-03RTI
5	Left	11-048-05LTI-S	11-048-05LTI
5	Right	11-048-05RTI-S	11-048-05RTI

Stainless Steel			
Holes	Direction	Sterile	Non-Sterile
3	Left	11-048-03LSS-S	11-048-03LSS
3	Right	11-048-03RSS-S	11-048-03RSS
5	Left	11-048-05LSS-S	11-048-05LSS
5	Right	11-048-05RSS-S	11-048-05RSS

Oblique (3 Head Holes)

Titanium			
Holes	Direction	Sterile	Non-Sterile
3	Left	11-049-03LTI-S	11-049-03LTI
3	Right	11-049-03RTI-S	11-049-03RTI
5	Left	11-049-05LTI-S	11-049-05LTI
5	Right	11-049-05RTI-S	11-049-05RTI

Stainless Steel			
Holes	Direction	Sterile	Non-Sterile
3	Left	11-049-03LSS-S	11-049-03LSS
3	Right	11-049-03RSS-S	11-049-03RSS
5	Left	11-049-05LSS-S	11-049-05LSS
5	Right	11-049-05RSS-S	11-049-05RSS

2.4mm AV-Wiselock "T" Dorsal Distal Radius Plate, (3 Head Holes)



Titanium		
Holes	Sterile	Non-Sterile
3	11-050-03TI-S	11-050-03TI
5	11-050-05TI-S	11-050-05TI

Stainless Steel		
Holes	Sterile	Non-Sterile
3	11-050-03SS-S	11-050-03SS
5	11-050-05SS-S	11-050-05SS

2.4mm Variable Angle Two-Column Volar Distal Radius Plate, **Titanium**



Width 19.5mm



Narrow, 6 Head Holes

Width 22mm



Standard, 6 Head Holes

Width 25.5mm



Wide, 7 Head Holes

Narrow

Holes	Direction	Sterile	Non-Sterile
2	Left	TI-546.002L-S	TI-546.002L
2	Right	TI-546.002R-S	TI-546.002R
3	Left	TI-546.003L-S	TI-546.003L
3	Right	TI-546.003R-S	TI-546.003R
4	Left	TI-546.004L-S	TI-546.004L
4	Right	TI-546.004R-S	TI-546.004R
5	Left	TI-546.005L-S	TI-546.005L
5	Right	TI-546.005R-S	TI-546.005R
6	Left	TI-546.006L-S	TI-546.006L
6	Right	TI-546.006R-S	TI-546.006R
7	Left	TI-546.007L-S	TI-546.007L
7	Right	TI-546.007R-S	TI-546.007R
8	Left	TI-546.008L-S	TI-546.008L
8	Right	TI-546.008R-S	TI-546.008R

Standard

Holes	Direction	Sterile	Non-Sterile
2	Left	TI-546.202L-S	TI-546.202L
2	Right	TI-546.202R-S	TI-546.202R
3	Left	TI-546.203L-S	TI-546.203L
3	Right	TI-546.203R-S	TI-546.203R
4	Left	TI-546.204L-S	TI-546.204L
4	Right	TI-546.204R-S	TI-546.204R
5	Left	TI-546.205L-S	TI-546.205L
5	Right	TI-546.205R-S	TI-546.205R
6	Left	TI-546.206L-S	TI-546.206L
6	Right	TI-546.206R-S	TI-546.206R
7	Left	TI-546.207L-S	TI-546.207L
7	Right	TI-546.207R-S	TI-546.207R
8	Left	TI-546.208L-S	TI-546.208L
8	Right	TI-546.208R-S	TI-546.208R

Wide

Holes	Direction	Sterile	Non-Sterile
2	Left	TI-548.002L-S	TI-548.002L
2	Right	TI-548.002R-S	TI-548.002R
3	Left	TI-548.003L-S	TI-548.003L
3	Right	TI-548.003R-S	TI-548.003R
4	Left	TI-548.004L-S	TI-548.004L
4	Right	TI-548.004R-S	TI-548.004R
5	Left	TI-548.005L-S	TI-548.005L
5	Right	TI-548.005R-S	TI-548.005R
6	Left	TI-548.006L-S	TI-548.006L
6	Right	TI-548.006R-S	TI-548.006R
7	Left	TI-548.007L-S	TI-548.007L
7	Right	TI-548.007R-S	TI-548.007R
8	Left	TI-548.008L-S	TI-548.008L
8	Right	TI-548.008R-S	TI-548.008R

2.4mm Variable Angle Two-Column Volar Distal Radius Plate, **Stainless Steel**



Width 19.5mm



Narrow, 6 Head Holes

Width 22mm



Standard, 6 Head Holes

Width 25.5mm



Wide, 7 Head Holes

Narrow

Holes	Direction	Sterile	Non-Sterile
2	Left	546.002L-S	546.002L
2	Right	546.002R-S	546.002R
3	Left	546.003L-S	546.003L
3	Right	546.003R-S	546.003R
4	Left	546.004L-S	546.004L
4	Right	546.004R-S	546.004R
5	Left	546.005L-S	546.005L
5	Right	546.005R-S	546.005R
6	Left	546.006L-S	546.006L
6	Right	546.006R-S	546.006R
7	Left	546.007L-S	546.007L
7	Right	546.007R-S	546.007R
8	Left	546.008L-S	546.008L
8	Right	546.008R-S	546.008R

Standard

Holes	Direction	Sterile	Non-Sterile
2	Left	546.202L-S	546.202L
2	Right	546.202R-S	546.202R
3	Left	546.203L-S	546.203L
3	Right	546.203R-S	546.203R
4	Left	546.204L-S	546.204L
4	Right	546.204R-S	546.204R
5	Left	546.205L-S	546.205L
5	Right	546.205R-S	546.205R
6	Left	546.206L-S	546.206L
6	Right	546.206R-S	546.206R
7	Left	546.207L-S	546.207L
7	Right	546.207R-S	546.207R
8	Left	546.208L-S	546.208L
8	Right	546.208R-S	546.208R

Wide

Holes	Direction	Sterile	Non-Sterile
2	Left	548.002L-S	548.002L
2	Right	548.002R-S	548.002R
3	Left	548.003L-S	548.003L
3	Right	548.003R-S	548.003R
4	Left	548.004L-S	548.004L
4	Right	548.004R-S	548.004R
5	Left	548.005L-S	548.005L
5	Right	548.005R-S	548.005R
6	Left	548.006L-S	548.006L
6	Right	548.006R-S	548.006R
7	Left	548.007L-S	548.007L
7	Right	548.007R-S	548.007R
8	Left	548.008L-S	548.008L
8	Right	548.008R-S	548.008R

Aiming Block For Variable Angle Two-Column Volar Distal Radius Plate, Narrow
Titanium

Holes	Direction	Width	Sterile	Non-Sterile
6	Left	19.5mm	TI-650-19.5L-S	TI-650-19.5L
6	Right	19.5mm	TI-650-19.5R-S	TI-650-19.5R

Stainless Steel

Holes	Direction	Width	Sterile	Non-Sterile
6	Left	19.5mm	650-19.5L-S	650-19.5L
6	Right	19.5mm	650-19.5R-S	650-19.5R


Aiming Block For Variable Angle Two-Column Volar Distal Radius Plate, Standard
Titanium

Holes	Direction	Width	Sterile	Non-Sterile
6	Left	22mm	TI-650-22L-S	TI-650-22L
6	Right	22mm	TI-650-22R-S	TI-650-22R

Stainless Steel

Holes	Direction	Width	Sterile	Non-Sterile
6	Left	22mm	650-22L-S	650-22L
6	Right	22mm	650-22R-S	650-22R


Aiming Block For Variable Angle Two-Column Volar Distal Radius Plate, Wide
Titanium

Holes	Direction	Width	Sterile	Non-Sterile
6	Left	25.5mm	TI-650-25.5L-S	TI-650-25.5L
6	Right	25.5mm	TI-650-25.5R-S	TI-650-25.5R

Stainless Steel

Holes	Direction	Width	Sterile	Non-Sterile
6	Left	25.5mm	650-25.5L-S	650-25.5L
6	Right	25.5mm	650-25.5R-S	650-25.5R



Aiming Block Drill Sleeve

Titanium

Sterile	Non-Sterile
TI-650-01-S	TI-650-01

Stainless Steel

Sterile	Non-Sterile
650-01-S	650-01



Aiming Block Screw

Titanium


Sterile	Non-Sterile
TI-650-02-S	TI-650-02

Stainless Steel

Sterile	Non-Sterile
650-02-S	650-02



2.4mm Variable Angle Volar Rim Distal Radius Plate



6 Head Holes

Titanium			
Holes	Direction	Sterile	Non-Sterile
3	Left	TI-560-06-03L-S	TI-560-06-03L
3	Right	TI-560-06-03R-S	TI-560-06-03R
4	Left	TI-560-06-04L-S	TI-560-06-04L
4	Right	TI-560-06-04R-S	TI-560-06-04R
5	Left	TI-560-06-05L-S	TI-560-06-05L
5	Right	TI-560-06-05R-S	TI-560-06-05R
6	Left	TI-560-06-06L-S	TI-560-06-06L
6	Right	TI-560-06-06R-S	TI-560-06-06R
7	Left	TI-560-06-07L-S	TI-560-06-07L
7	Right	TI-560-06-07R-S	TI-560-06-07R
8	Left	TI-560-06-08L-S	TI-560-06-08L
8	Right	TI-560-06-08R-S	TI-560-06-08R

Stainless Steel			
Holes	Direction	Sterile	Non-Sterile
3	Left	560-06-03L-S	560-06-03L
3	Right	560-06-03R-S	560-06-03R
4	Left	560-06-04L-S	560-06-04L
4	Right	560-06-04R-S	560-06-04R
5	Left	560-06-05L-S	560-06-05L
5	Right	560-06-05R-S	560-06-05R
6	Left	560-06-06L-S	560-06-06L
6	Right	560-06-06R-S	560-06-06R
7	Left	560-06-07L-S	560-06-07L
7	Right	560-06-07R-S	560-06-07R
8	Left	560-06-08L-S	560-06-08L
8	Right	560-06-08R-S	560-06-08R

7 Head Holes

Titanium			
Holes	Direction	Sterile	Non-Sterile
3	Left	TI-560-07-03L-S	TI-560-07-03L
3	Right	TI-560-07-03R-S	TI-560-07-03R
4	Left	TI-560-07-04L-S	TI-560-07-04L
4	Right	TI-560-07-04R-S	TI-560-07-04R
5	Left	TI-560-07-05L-S	TI-560-07-05L
5	Right	TI-560-07-05R-S	TI-560-07-05R
6	Left	TI-560-07-06L-S	TI-560-07-06L
6	Right	TI-560-07-06R-S	TI-560-07-06R
7	Left	TI-560-07-07L-S	TI-560-07-07L
7	Right	TI-560-07-07R-S	TI-560-07-07R
8	Left	TI-560-07-08L-S	TI-560-07-08L
8	Right	TI-560-07-08R-S	TI-560-07-08R

Stainless Steel			
Holes	Direction	Sterile	Non-Sterile
3	Left	560-07-03L-S	560-07-03L
3	Right	560-07-03R-S	560-07-03R
4	Left	560-07-04L-S	560-07-04L
4	Right	560-07-04R-S	560-07-04R
5	Left	560-07-05L-S	560-07-05L
5	Right	560-07-05R-S	560-07-05R
6	Left	560-07-06L-S	560-07-06L
6	Right	560-07-06R-S	560-07-06R
7	Left	560-07-07L-S	560-07-07L
7	Right	560-07-07R-S	560-07-07R
8	Left	560-07-08L-S	560-07-08L
8	Right	560-07-08R-S	560-07-08R

Aiming Block For Variable Angle Volar Rim Distal Radius Plate, 6 Head Holes

Titanium

Direction	Sterile	Non-Sterile
Left	TI-651-06L-S	TI-651-06L
Right	TI-651-06R-S	TI-651-06R

Stainless Steel

Direction	Sterile	Non-Sterile
Left	651-06L-S	651-06L
Right	651-06R-S	651-06R



Aiming Block For Variable Angle Volar Rim Distal Radius Plate, 7 Head Holes

Titanium

Direction	Sterile	Non-Sterile
Left	TI-651-07L-S	TI-651-07L
Right	TI-651-07R-S	TI-651-07R

Stainless Steel

Direction	Sterile	Non-Sterile
Left	651-07L-S	651-07L
Right	651-07R-S	651-07R



Aiming Block Drill Sleeve

Titanium

Sterile	Non-Sterile
TI-650-01-S	TI-650-01

Stainless Steel

Sterile	Non-Sterile
650-01-S	650-01



Aiming Block Screw

Titanium

Sterile	Non-Sterile
TI-650-02-S	TI-650-02


Stainless Steel

Sterile	Non-Sterile
650-02-S	650-02




2.4mm Variable Angle Screw, Self-Tapping, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø2.4	6mm	TI-549.006-S	TI-549.006	549.006-S	549.006
Ø2.4	8mm	TI-549.008-S	TI-549.008	549.008-S	549.008
Ø2.4	10mm	TI-549.010-S	TI-549.010	549.010-S	549.010
Ø2.4	12mm	TI-549.012-S	TI-549.012	549.012-S	549.012
Ø2.4	14mm	TI-549.014-S	TI-549.014	549.014-S	549.014
Ø2.4	16mm	TI-549.016-S	TI-549.016	549.016-S	549.016
Ø2.4	18mm	TI-549.018-S	TI-549.018	549.018-S	549.018
Ø2.4	20mm	TI-549.020-S	TI-549.020	549.020-S	549.020
Ø2.4	22mm	TI-549.022-S	TI-549.022	549.022-S	549.022
Ø2.4	24mm	TI-549.024-S	TI-549.024	549.024-S	549.024
Ø2.4	26mm	TI-549.026-S	TI-549.026	549.026-S	549.026
Ø2.4	28mm	TI-549.028-S	TI-549.028	549.028-S	549.028
Ø2.4	30mm	TI-549.030-S	TI-549.030	549.030-S	549.030




2.4mm Wise-Lock Screw, Self-Tapping, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø2.4	6mm	TI-1145.06-S	TI-1145.06	1145.06-S	1145.06
Ø2.4	8mm	TI-1145.08-S	TI-1145.08	1145.08-S	1145.08
Ø2.4	10mm	TI-1145.10-S	TI-1145.10	1145.10-S	1145.10
Ø2.4	12mm	TI-1145.12-S	TI-1145.12	1145.12-S	1145.12
Ø2.4	14mm	TI-1145.14-S	TI-1145.14	1145.14-S	1145.14
Ø2.4	16mm	TI-1145.16-S	TI-1145.16	1145.16-S	1145.16
Ø2.4	18mm	TI-1145.18-S	TI-1145.18	1145.18-S	1145.18
Ø2.4	20mm	TI-1145.20-S	TI-1145.20	1145.20-S	1145.20
Ø2.4	22mm	TI-1145.22-S	TI-1145.22	1145.22-S	1145.22
Ø2.4	24mm	TI-1145.24-S	TI-1145.24	1145.24-S	1145.24
Ø2.4	26mm	TI-1145.26-S	TI-1145.26	1145.26-S	1145.26
Ø2.4	28mm	TI-1145.28-S	TI-1145.28	1145.28-S	1145.28
Ø2.4	30mm	TI-1145.30-S	TI-1145.30	1145.30-S	1145.30



2.4mm Cortical Screw, Self-Tapping, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø2.4	6mm	TI-1144.06-S	TI-1144.06	1144.06-S	1144.06
Ø2.4	8mm	TI-1144.08-S	TI-1144.08	1144.08-S	1144.08
Ø2.4	10mm	TI-1144.10-S	TI-1144.10	1144.10-S	1144.10
Ø2.4	12mm	TI-1144.12-S	TI-1144.12	1144.12-S	1144.12
Ø2.4	14mm	TI-1144.14-S	TI-1144.14	1144.14-S	1144.14
Ø2.4	16mm	TI-1144.16-S	TI-1144.16	1144.16-S	1144.16
Ø2.4	18mm	TI-1144.18-S	TI-1144.18	1144.18-S	1144.18
Ø2.4	20mm	TI-1144.20-S	TI-1144.20	1144.20-S	1144.20
Ø2.4	22mm	TI-1144.22-S	TI-1144.22	1144.22-S	1144.22
Ø2.4	24mm	TI-1144.24-S	TI-1144.24	1144.24-S	1144.24
Ø2.4	26mm	TI-1144.26-S	TI-1144.26	1144.26-S	1144.26
Ø2.4	28mm	TI-1144.28-S	TI-1144.28	1144.28-S	1144.28
Ø2.4	30mm	TI-1144.30-S	TI-1144.30	1144.30-S	1144.30



2.7mm Cortical Screw, Self-Tapping, (Star Head)

Dia	Length	Titanium		Stainless Steel	
		Sterile	Non-Sterile	Sterile	Non-Sterile
Ø2.4	6mm	TI-1242.06-S	TI-1242.06	1242.06-S	1242.06
Ø2.4	8mm	TI-1242.08-S	TI-1242.08	1242.08-S	1242.08
Ø2.4	10mm	TI-1242.10-S	TI-1242.10	1242.10-S	1242.10
Ø2.4	12mm	TI-1242.12-S	TI-1242.12	1242.12-S	1242.12
Ø2.4	14mm	TI-1242.14-S	TI-1242.14	1242.14-S	1242.14
Ø2.4	16mm	TI-1242.16-S	TI-1242.16	1242.16-S	1242.16
Ø2.4	18mm	TI-1242.18-S	TI-1242.18	1242.18-S	1242.18
Ø2.4	20mm	TI-1242.20-S	TI-1242.20	1242.20-S	1242.20
Ø2.4	22mm	TI-1242.22-S	TI-1242.22	1242.22-S	1242.22
Ø2.4	24mm	TI-1242.24-S	TI-1242.24	1242.24-S	1242.24
Ø2.4	26mm	TI-1242.26-S	TI-1242.26	1242.26-S	1242.26
Ø2.4	28mm	TI-1242.28-S	TI-1242.28	1242.28-S	1242.28
Ø2.4	30mm	TI-1242.30-S	TI-1242.30	1242.30-S	1242.30



1312-01 Depth Gauge measuring upto 60mm




1312-02 T-Handle With Quick Coupling



1312-03-15 Retractor, 15mm



1312-04-1.8 Drill Bit with Quick Coupling End, Ø1.8mm x Length 140mm



1312-05-2.0 Drill Bit with Quick Coupling End, Ø2.0mm x Length 140mm



1312-06-2.4 Drill Bit with Quick Coupling End, Ø2.4mm x Length 140mm



1312-07-2.7

Drill Bit with Quick Coupling End, Ø2.7mm x Length 140mm


1312-08-1.2

Guide Wire, Ø1.2mm x Length 150mm


1312-09-1.2

Guide Wire With Threaded Tip, Ø1.2mm x Thread Length 5mm x Length 150mm


1312-10-2.4

Tap, HA 2.4mm


1312-11-2.7

Tap, HA 2.7mm


1312-13-1.8

Threaded Drill Sleeve, Ø1.8/2.4mm (Double Start Thread)



1312-15 Self-Centering Double Drill Guide $\Phi 1.8/2.4\text{mm}$



1312-16 Self-Centering Double Drill Guide $\Phi 2.0/2.7\text{mm}$



1312-17 Depth Gauge with Protector measuring upto 50mm



1312-18 Star Screwdriver Shaft, T8



1312-21 Torque Limiting Attachment, 0.8Nm



7-019-03 Handle for Torque Limiting Attachment



1312-22-6.5

Retractor, 6.5mm


1312-23

Double Periosteal Elevator, Round


1312-24

Double Periosteal Elevator, Flat


1312-25

Reduction Forcep with Points, Small


1312-26

Reduction Forcep with Points, Large


1312-27

Reduction Forcep with Points, Serrated



1312-28

Bending Forcep


1312-30

Universal Variable Angle Drill Guide, Ø1.8mm


1312-20-2.7

Screw Holding Sleeve for 2.7mm Cortical Screw


1312-20-2.4

Screw Holding Sleeve for 2.4mm Cortical Screw


1312-30-2.4

Screw Holding Sleeve for 2.4mm Wise-lock Screw


1312-22

Quick Coupling Handle



1312-29

Sharp Hook




1312-31

Screwdriver Star Head




1312-35

Instrument trays for 2.4mm AV-Wiselock Distal Radius Instrument Set



1312-34

Container for 2.4mm AV-Wiselock Distal Radius Instrument



1312-000 2.4mm AV-Wiselock Distal Radius Instrument Set



Code	Set Consisting of	Units
1312-01	Depth Gauge measuring upto 60mm	1
1312-02	T-Handle With Quick Coupling	1
1312-03-15	Retractor, 15mm	1
1312-04-1.8	Drill Bit with Quick Coupling End, Ø1.8mm x Length 140mm	3
1312-05-2.0	Drill Bit with Quick Coupling End, Ø2.0mm x Length 140mm	2
1312-06-2.4	Drill Bit with Quick Coupling End, Ø2.4mm x Length 140mm	1
1312-07-2.7	Drill Bit with Quick Coupling End, Ø2.7mm x Length 140mm	1
1312-08-1.2	Guide Wire, Ø1.2mm x Length 150mm	5
1312-09-1.2	Guide Wire With Threaded Tip, Ø1.2mm x Thread Length 5mm x Length 150mm	2
1312-10-2.4	Tap, HA 2.4mm	1
1312-11-2.7	Tap, HA 2.7mm	1
1312-13-1.8	Threaded Drill Sleeve, Ø1.8/2.4mm (Double Start Thread)	3
1312-15	Self-Centering Double Drill Guide Ø1.8/2.4mm	1
1312-16	Self-Centering Double Drill Guide Ø2.0/2.7mm	1
1312-17	Depth Gauge with Protector measuring upto 50mm	1
1312-18	Star Screwdriver Shaft, T8	2
1312-21	Torque Limiting Attachment, 0.8Nm	1
7-019-03	Handle for Torque Limiting Attachment	1
1312-22-6.5	Retractor, 6.5mm	1
1312-23	Double Periosteal Elevator, Round	1
1312-24	Double Periosteal Elevator, Flat	1
1312-25	Reduction Forcep with Points, Small	1
1312-26	Reduction Forcep with Points, Large	1
1312-27	Reduction Forcep with Points, Serrated	1
1312-28	Bending Forcep	2
1312-30	Universal Variable Angle Drill Guide, Ø1.8mm	1
1312-20-2.7	Screw Holding Sleeve for 2.7mm Cortical Screw	1
1312-20-2.4	Screw Holding Sleeve for 2.4mm Cortical Screw	1
1312-30-2.4	Screw Holding Sleeve for 2.4mm Wise-lock Screw	1
1312-22	Quick Coupling Handle	1
1312-29	Sharp Hook	1
1312-31	Screwdriver Star Head	1
1312-35	Instrument trays for 2.4mm AV-Wiselock Distal Radius Instrument Set	2
1312-34	Container for 2.4mm AV-Wiselock Distal Radius Instrument	1

Implant Box Optional

1312-33	Implant Box with cover for 2.4mm Variable Angle Distal Radius system	1
10-050	Implant Box with cover for 2.4mm Variable Angle Volar Rim Distal Radius System	1
11-051	Implant Box with cover for 2.4mm AV-Wiselock Distal Radius Dorsal & Extra-Articular System	1
11-077	Implant Box with cover for 2.4mm AV-Wiselock Hand and Wrist - Arthrodesis System	1

**USA**

Auxein Inc.
1500 Nw 89th Court, Suite 107-108
Doral, Florida 33172
Tel: +1 305 395 6062
E Fax: +1 305 395 6262
Email: USoffice@auxein.com

MEXICO

Auxein México, S.A. de C.V.
Tepic 139 int 801, Colonia Roma Sur,
Alcaldía Cuauhtémoc, CDMX,
México, C.P. 06760
Tel: +521 55 7261 0318
Email: info@auxein.mx

INDIA

Auxein Medical Pvt. Ltd.
Plot No. 168-169-170, Phase-4,
Kundli Industrial Area,
HSIIDC, Sector-57, Sonapat - 131028, Haryana
Tel: +91 99106 43638 | Fax: +91 86077 70197
Email: info@auxein.com