

Mini Fragment

1.5mm Cortical Screw(Hex Head)

Drill Bit: 1.1mm



Length (mm)	Stainless Steel	Titanium
6	101.006	TI-101.006
7	101.007	TI-101.007
8	101.008	TI-101.008
9	101.009	TI-101.009
10	101.010	TI-101.010
11	101.011	TI-101.011
12	101.012	TI-101.012
14	101.014	TI-101.014
16	101.016	TI-101.016
18	101.018	TI-101.018
20	101.020	TI-101.020

1.5mm Cortical Screw, Self -Tapping(Hex Head)

Drill Bit: 1.1mm



Length (mm)	Stainless Steel	Titanium
6	101.206	TI-101.206
7	101.207	TI-101.207
8	101.208	TI-101.208
9	101.209	TI-101.209
10	101.210	TI-101.210
11	101.211	TI-101.211
12	101.212	TI-101.212
14	101.214	TI-101.214
16	101.216	TI-101.216
18	101.218	TI-101.218
20	101.220	TI-101.220

2.0mm Cortical Screw(Hex Head)

Drill Bit : 1.5mm



Length (mm)	Stainless Steel	Titanium
6	102.006	TI-102.006
7	102.007	TI-102.007
8	102.008	TI-102.008
9	102.009	TI-102.009
10	102.010	TI-102.010
12	102.012	TI-102.012
14	102.014	TI-102.014
16	102.016	TI-102.016
18	102.018	TI-102.018
20	102.020	TI-102.020

2.0mm Cortical Screw, Self -Tapping(Hex Head)

Drill Bit : 1.5mm



Length (mm)	Stainless Steel	Titanium
6	102.306	TI-102.306
7	102.307	TI-102.307
8	102.308	TI-102.308
9	102.309	TI-102.309
10	102.310	TI-102.310
12	102.312	TI-102.312
14	102.314	TI-102.314
16	102.316	TI-102.316
18	102.318	TI-102.318
20	102.320	TI-102.320

2.7mm Cortical Screw(Hex Head)

Drill Bit : 2.0mm



Length (mm)	Stainless Steel	Titanium
6	103.006	TI-103.006
8	103.008	TI-103.008
10	103.010	TI-103.010
12	103.012	TI-103.012
14	103.014	TI-103.014
16	103.016	TI-103.016
18	103.018	TI-103.018
20	103.020	TI-103.020
22	103.022	TI-103.022
24	103.024	TI-103.024
26	103.026	TI-103.026
28	103.028	TI-103.028
30	103.030	TI-103.030

2.7mm Cortical Screw, Self -Tapping(Hex Head)

Drill Bit : 2.0mm



Length (mm)	Stainless Steel	Titanium
6	103.206	TI-103.206
8	103.208	TI-103.208
10	103.210	TI-103.210
12	103.212	TI-103.212
14	103.214	TI-103.214
16	103.216	TI-103.216
18	103.218	TI-103.218
20	103.220	TI-103.220
22	103.222	TI-103.222
24	103.224	TI-103.224
26	103.226	TI-103.226
28	103.228	TI-103.228
30	103.230	TI-103.230

2.0mm Mini "L" Plate



Holes	Length (mm)	Left Direction		Right Direction	
		Stainless Steel	Titanium	Stainless Steel	Titanium
2	18.60	906.100L	TI-906.100L	906.100R	TI-906.100R
4	30.20	906.100-04L	TI-906.100-04L	906.100-04R	TI-906.100-04R
6	41.80	906.100-06L	TI-906.100-06L	906.100-06R	TI-906.100-06R
8	53.40	906.100-08L	TI-906.100-08L	906.100-08R	TI-906.100-08R

2.0mm Mini Straight Plate



Holes	Length (mm)	Stainless Steel	Titanium
2	12.40	909.102	TI-909.102
3	18.40	909.103	TI-909.103
4	24.40	909.104	TI-909.104
5	30.40	909.105	TI-909.105
6	36.40	909.106	TI-909.106
7	42.40	909.107	TI-909.107
8	48.40	909.108	TI-909.108

2.0mm DCP Plate



Holes	Length (mm)	Stainless Steel	Titanium
4	22.00	908.104	TI-908.104
5	27.00	908.105	TI-908.105
6	32.00	908.106	TI-908.106
7	37.00	908.107	TI-908.107
8	42.00	908.108	TI-908.108

2.0mm Mini "T" Plate

Shaft with two holes & Head with two holes.



	Length (mm)	Stainless Steel	Titanium
Code	18.60	905.100	TI-905.100

2.7mm Mini "L" plate



Type	Length (mm)	Stainless Steel	Titanium
Left	32.80	705.100L	TI-705.100L
Right	32.80	705.100R	TI-705.100R

2.7mm "T" Plate



	Length(mm)	Stainless Steel	Titanium
Code	34.50	707.100	TI-707.100

2.7mm Quarter Tubular Plate



Holes	Length(mm)	Stainless Steel	Titanium
2	18.00	701.102	TI-701.102
3	26.00	701.103	TI-701.103
4	34.00	701.104	TI-701.104
5	42.00	701.105	TI-701.105
6	50.00	701.106	TI-701.106
7	58.00	701.107	TI-701.107
8	66.00	701.108	TI-701.108
9	74.00	701.109	TI-701.109
10	82.00	701.110	TI-701.110

2.7mm Dynamic Compression Plate



Holes	Length(mm)	Stainless Steel	Titanium
2	20.00	719.102	TI-719.102
3	28.00	719.103	TI-719.103
4	36.00	719.104	TI-719.104
5	44.00	719.105	TI-719.105
6	52.00	719.106	TI-719.106
7	60.00	719.107	TI-719.107
8	68.00	719.108	TI-719.108
9	76.00	719.109	TI-719.109
10	84.00	719.110	TI-719.110
12	100.0	719.112	TI-719.112

2.7mm Reconstruction Plate, Straight



Holes	Length(mm)	Stainless Steel	Titanium
2	16.20	726.102	TI-726.102
3	24.20	726.103	TI-726.103
4	32.20	726.104	TI-726.104
5	40.20	726.105	TI-726.105
6	48.20	726.106	TI-726.106
7	56.20	726.107	TI-726.107
8	64.20	726.108	TI-726.108
9	72.20	726.109	TI-726.109
10	80.20	726.110	TI-726.110
11	88.20	726.111	TI-726.111
12	96.20	726.112	TI-726.112
14	112.2	726.114	TI-726.114
16	128.2	726.116	TI-726.116
18	144.2	726.118	TI-726.118
20	160.2	726.120	TI-726.120
22	176.2	726.122	TI-726.122
24	192.2	726.124	TI-726.124

7-044-32

Implant Box for 1.5mm
Mini Fragment System


7-044-34

Implant Box for 2.0mm
Mini Fragment System


7-044-35

Implant Box for 2.7mm
Mini Fragment System



Instruments for 1.5mm System

7-044-01 Drill Bit, Ø1.1mm x Length 60mm



7-044-02 Double Drill Guide, Ø1.5/1.1mm



7-044-03 Tap Ø1.5mm, Quick Coupling


Instruments for 2.0mm System

7-044-05 Drill Bit, Ø1.5mm x Length 80mm



7-044-06 Double Drill Guide, Ø2.0/1.5mm




7-044-07 Tap Ø2.0mm, Quick Coupling



Instruments for 2.7mm System

7-044-08 Drill Bit, Ø2.0mm x Length 80mm



7-044-09 Double Drill Guide, Ø2.7/2.0mm



7-044-10 Tap Ø2.7mm, Quick Coupling



7-044-12 Countersink, for Screws 2.7mm


General Instruments

7-044-13 T-Handle with Quick Coupling for Mini Fragment



7-044-14 Mini Plate and Bone Holding Forcep



7-044-15

Mini Screw Holding Forcep

**7-044-16**

Mini Plate Bender (Pair)

**7-044-17**

Mini Periosteal Elevator, Curved, 3mm

**7-044-18**

Mini Periosteal Elevator, Curved, 5mm

**7-044-19**

Mini Reduction Forcep - Pointed

**7-044-20**

Mini Reduction Forcep - Serrated



7-044-21

Mini Hohmann Retractor, 6.0mm



7-044-22

Mini Hohmann Retractor, 8.5mm



7-044-23

Mini Hohmann Retractor, 15mm



7-044-24

Sharp Hook



7-044-25

Quick Coupling Handle, Straight



7-044-26

Depth Gauge measuring upto 50mm for Mini Fragment



7-044-27

Skin Hook, Single Pronged


7-044-28

Skin Hook, Double Pronged


7-044-29


Plate Cutter & Bending Forcep


7-044-30

Plate Holding Forcep


7-044-31

Plate Cutter


7-044-04

Screwdriver Shaft with Holding Sleeve, Hex 1.5mm, for 1.5mm Screw



- 7-044-11** Screwdriver Shaft with Holding Sleeve, Hex 1.5mm, for 2.0mm Screw



- 7-044-36** Screwdriver Shaft with Holding Sleeve, Hex 2.0mm, for 2.7mm Screw



- 7-044-32** Implant Box for 1.5mm Mini Fragment System



- 7-044-34** Implant Box for 2.0mm Mini Fragment System



- 7-044-35** Implant Box for 2.7mm Mini Fragment System



- 7-044-33** Container for Mini Fragment Instrument Set



7-044 : Mini Fragment Instrument Set



7-044 : Mini Fragment Instrument Set

Codes	Set of Consisting of	Units
	Instruments for 1.5mm System	
7-044-01	Drill Bit, Ø1.1mm x Length 60mm	3
7-044-02	Double Drill Guide, Ø1.5/1.1mm	1
7-044-03	Tap Ø1.5mm, Quick Coupling	1
	Instruments for 2.0mm System	
7-044-05	Drill Bit, Ø1.5mm x Length 80mm	3
7-044-06	Double Drill Guide, Ø2.0/1.5mm	1
7-044-07	Tap Ø2.0mm, Quick Coupling	1
	Instruments for 2.7mm System	
7-044-08	Drill Bit, Ø2.0mm x Length 80mm	3
7-044-09	Double Drill Guide, Ø2.7/2.0mm for Mini Fragment	1
7-044-10	Tap Ø2.7mm, Quick Coupling	1
7-044-12	Countersink, for Screws 2.7mm	1
	General Instrument	
7-044-13	T-Handle with Quick Coupling for Mini Fragment	1
7-044-14	Mini Plate and Bone Holding Forcep	1
7-044-15	Mini Screw Holding Forcep	1
7-044-16	Mini Plate Bender (Pair)	1
7-044-17	Mini Periosteal Elevator, Curved, 3mm	1
7-044-18	Mini Periosteal Elevator, Curved, 5mm	1
7-044-19	Mini Reduction Forcep - Pointed	1
7-044-20	Mini Reduction Forcep - Serrated	1
7-044-21	Mini Hohmann Retractor, 6.0mm	1
7-044-22	Mini Hohmann Retractor, 8.5mm	1
7-044-23	Mini Hohmann Retractor, 15mm	1
7-044-24	Sharp Hook	1
7-044-25	Quick Coupling Handle, Straight	1
7-044-26	Depth Gauge measuring upto 50mm for Mini Fragment	1
7-044-27	Skin Hook, Single Pronged	1
7-044-28	Skin Hook, Double Pronged	1
7-044-29	Plate Cutter & Bending Forcep	1
7-044-30	Plate Holding Forcep	1
7-044-31	Plate Cutter	1
7-044-04	Screwdriver Shaft with Holding Sleeve, Hex 1.5mm, for 1.5mm Screw	1
7-044-11	Screwdriver Shaft with Holding Sleeve, Hex 1.5mm, for 2.0mm Screw	1
7-044-36	Screwdriver Shaft with Holding Sleeve, Hex 2.0mm, for 2.7mm Screw	1
7-044-32	Implant Box for 1.5mm Mini Fragment System	1
7-044-34	Implant Box for 2.0mm Mini Fragment System	1
7-044-35	Implant Box for 2.7mm Mini Fragment System	1
7-044-33	Container for Mini Fragment Instrument Set	1