

PEDIATRIC HIP PLATE




about us

Auxein Medical is an integrated, research based, orthopaedic Implants & Instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and innovator products that exceed the highest expectations of our customers in term of quality and safety. The company has world-class manufacturing unit established in India and serves customers in over 60 countries worldwide.


Our Achievements

- US FDA 510K APPROVED •**
(Arthroscopy, Spine, Nailing, Locking Plates & Screws)
- EC Certificate - Class IIa & IIb •
- ISO 13485:2016 •
- ISO 9001:2015 •
- India FDA Approved •
- India MDR 2017 •

Pediatric Hip **Plate System**

3.5mm Wise-Lock Pediatric Hip Plate


Degree	Length	Holes	Width	Stainless Steel	Titanium
90°	75mm	3	18.5mm	1121.090	TI-1121.090
100°	75mm	3	18.5mm	1121.100	TI-1121.100
110°	75mm	3	18.5mm	1121.110	TI-1121.110
120°	75mm	4	18.5mm	1121.120	TI-1121.120
150°	60mm	3	18.5mm	1121.150	TI-1121.150

4.5/5.0mm Wise-Lock Paediatric Hip Plate


Degree	Length	Holes	Width	Stainless Steel	Titanium
90°	90mm	3	22.5mm	1122.090	TI-1122.090
100°	90mm	3	22.5mm	1122.100	TI-1122.100
110°	90mm	3	22.5mm	1122.110	TI-1122.110
120°	95mm	4	22.5mm	1122.120	TI-1122.120
150°	75mm	3	22.5mm	1122.150	TI-1122.150

3.5mm Wise-Lock Screw, Self-Tapping, Hex Head

Length (mm)	Stainless Steel	Titanium
10	117.010	TI-117.010
12	117.012	TI-117.012
14	117.014	TI-117.014
16	117.016	TI-117.016
18	117.018	TI-117.018
20	117.020	TI-117.020
22	117.022	TI-117.022
24	117.024	TI-117.024
26	117.026	TI-117.026
28	117.028	TI-117.028
30	117.030	TI-117.030
32	117.032	TI-117.032
34	117.034	TI-117.034
35	117.035	TI-117.035
36	117.036	TI-117.036
38	117.038	TI-117.038
40	117.040	TI-117.040
42	117.042	TI-117.042
44	117.044	TI-117.044
45	117.045	TI-117.045
46	117.046	TI-117.046
48	117.048	TI-117.048
50	117.050	TI-117.050
52	117.052	TI-117.052
54	117.054	TI-117.054
55	117.055	TI-117.055
56	117.056	TI-117.056
58	117.058	TI-117.058
60	117.060	TI-117.060
65	117.065	TI-117.065
70	117.070	TI-117.070
75	117.075	TI-117.075
80	117.080	TI-117.080



3.5mm Cortical Screw, Self-Tapping, Hex Head

Length (mm)	Stainless Steel	Titanium
10	104.210	TI-104.210
12	104.212	TI-104.212
14	104.214	TI-104.214
16	104.216	TI-104.216
18	104.218	TI-104.218
20	104.220	TI-104.220
22	104.222	TI-104.222
24	104.224	TI-104.224
26	104.226	TI-104.226
28	104.228	TI-104.228
30	104.230	TI-104.230
32	104.232	TI-104.232
34	104.234	TI-104.234
36	104.236	TI-104.236
38	104.238	TI-104.238
40	104.240	TI-104.240
42	104.242	TI-104.242
44	104.244	TI-104.244
46	104.246	TI-104.246
48	104.248	TI-104.248
50	104.250	TI-104.250
55	104.255	TI-104.255
60	104.260	TI-104.260
65	104.265	TI-104.265
70	104.270	TI-104.270
75	104.275	TI-104.275
80	104.280	TI-104.280
85	104.285	TI-104.285
90	104.290	TI-104.290



4.5mm Cortical Screws, Self-Tapping, Hex Head

Length (mm)	Stainless Steel	Titanium
20	105.220	TI-105.220
22	105.222	TI-105.222
24	105.224	TI-105.224
26	105.226	TI-105.226
28	105.228	TI-105.228
30	105.230	TI-105.230
32	105.232	TI-105.232
34	105.234	TI-105.234
36	105.236	TI-105.236
38	105.238	TI-105.238
40	105.240	TI-105.240
42	105.242	TI-105.242
44	105.244	TI-105.244
46	105.246	TI-105.246
48	105.248	TI-105.248
50	105.250	TI-105.250
52	105.252	TI-105.252
54	105.254	TI-105.254
55	105.255	TI-105.255
56	105.256	TI-105.256
58	105.258	TI-105.258
60	105.260	TI-105.260
65	105.265	TI-105.265
70	105.270	TI-105.270
80	105.280	TI-105.280
85	105.285	TI-105.285
90	105.290	TI-105.290



5.0mm Wise-Lock Screw, Self-Tapping, Hex Head

Length (mm)	Stainless Steel	Titanium
12	119.012	TI-119.012
14	119.014	TI-119.014
16	119.016	TI-119.016
18	119.018	TI-119.018
20	119.020	TI-119.020
22	119.022	TI-119.022
24	119.024	TI-119.024
26	119.026	TI-119.026
28	119.028	TI-119.028
30	119.030	TI-119.030
32	119.032	TI-119.032
34	119.034	TI-119.034
36	119.036	TI-119.036
38	119.038	TI-119.038
40	119.040	TI-119.040
42	119.042	TI-119.042
44	119.044	TI-119.044
46	119.046	TI-119.046
48	119.048	TI-119.048
50	119.050	TI-119.050
52	119.052	TI-119.052
54	119.054	TI-119.054
56	119.056	TI-119.056
58	119.058	TI-119.058
60	119.060	TI-119.060
65	119.065	TI-119.065
70	119.070	TI-119.070
75	119.075	TI-119.075
80	119.080	TI-119.080
85	119.085	TI-119.085
90	119.090	TI-119.090



1471-01 Guide Wire with Threaded Tip, Ø2.0mm x Thread Length 7.5mm x Length 150mm



1471-02 Positioner for Aiming Block



1471-03 Aiming Block, Ø3.5mm



1471-04 Aiming Block, Ø5.0mm



1471-05 Guide Wire, Ø2.8mm x Length 200mm



1471-06 Positioner for Osteotomy



1471-07 Reduction Forcep, Serrated Jaw, Curved



1471-08 Threaded Drill Sleeve for Drill Bit Ø2.8mm



1471-09 Threaded Drill Sleeve for Drill Bit Ø4.3mm



1471-10 Reduction Sleeve, Ø4.3/2.8mm



1471-11 Universal Depth Gauge for Pediatric Hip Plate



1471-13 Screwdriver, Hex 2.5mm for Pediatric Hip Plate



1471-13S Holding Sleeve for Screwdriver, Hex 2.5mm



1471-15 Screwdriver, Hex 3.5mm for Pediatric Hip Plate



1471-15S Holding Sleeve for Hexagonal Screwdriver, Hex 3.5mm




1471-16 Screwdriver Shaft, Hex 2.5mm



1471-17 Screwdriver Shaft, Hex 3.5mm



1471-18 T Handle with Quick Coupling for Pediatric Hip Plate



1471-21 Torque Limiting Screwdriver, 4.0Nm



1471-21S Holding Sleeve for Torque Limiting Screwdriver, 4.0Nm



1471-22 Torque Limiting Screwdriver, 1.5Nm



1471-22S Holding Sleeve for Torque Limiting Screwdriver, 1.5Nm



1471-23 Drill Bit with Quick Couling, Ø2.5mm x Length 112mm



1471-24 Drill Bit with Quick Couling, Ø2.8mm x Length 160mm



1471-25 Drill Bit with Quick Couling, Ø3.2mm x Length 147mm



1471-26 Drill Bit with Quick Couling, Ø4.3mm x Length 230mm



1471-27 Stopper for Drill Bit Ø2.8mm



1471-28 Stopper for Drill Bit Ø4.3mm



1471-29 Depth Gauge measuring upto 60mm for Pediatric Hip Plate



1471-30 Depth Gauge measuring upto 100mm for Pediatric Hip Plate



1471-31 Self-Centering Bone Holding Forcep for Pediatric Hip Plate**1471-32** Double Drill Guide, Ø2.5/3.5mm, for Pediatric Hip Plate**1471-33** Double Drill Guide, Ø3.2/4.5mm, for Pediatric Hip Plate**1471-34** Tap, HA3.5mm**1471-35** Tap, HA4.5mm**1471-36** Positioning Plate, Triangular, 90°/50°/40°

1471-37 Positioning Plate, Triangular, 80°/70°/30°




1471-38 Positioning Plate, Triangular, 100°/60°/20°



1471-39 Medialization Instrument for Pediatric Hip Plate






1471-41 Instrument Trays for Pediatric Hip Plate Instrument Set



1471-40 Container with PPSU Lid for Pediatric Hip Plate Instrument Set



INS-1471-000 Pediatric Hip Plate Instrument Set



INS-1471-000 Pediatric Hip Plate Instrument Set

Code	Set Consisting of	Units
1471-01	Guide Wire with Threaded Tip, Ø2.0mm x Thread Length 7.5mm x Length 150mm	4
1471-02	Positioner for Aiming Block	1
1471-03	Aiming Block, Ø3.5mm	1
1471-04	Aiming Block, Ø5.0mm	1
1471-05	Guide Wire, Ø2.8mm x Length 200mm	8
1471-06	Positioner for Osteotomy	1
1471-07	Reduction Forcep, Serrated Jaw, Curved	1
1471-08	Threaded Drill Sleeve for Drill Bit Ø2.8mm	3
1471-09	Threaded Drill Sleeve for Drill Bit Ø4.3mm	3
1471-10	Reduction Sleeve, Ø4.3/2.8mm	2
1471-11	Universal Depth Gauge for Pediatric Hip Plate	1
1471-13	Screwdriver, Hex 2.5mm for Pediatric Hip Plate	1
1471-13S	Holding Sleeve for Screwdriver, Hex 2.5mm	1
1471-15	Screwdriver, Hex 3.5mm for Pediatric Hip Plate	1
1471-15S	Holding Sleeve for Hexagonal Screwdriver, Hex 3.5mm	1
1471-16	Screwdriver Shaft, Hex 2.5mm	1
1471-17	Screwdriver Shaft, Hex 3.5mm	1
1471-18	T Handle with Quick Coupling for Pediatric Hip Plate	1
1471-21	Torque Limiting Screwdriver, 4.0Nm	1
1471-21S	Holding Sleeve for Torque Limiting Screwdriver, 4.0Nm	1
1471-22	Torque Limiting Screwdriver, 1.5Nm	1
1471-22S	Holding Sleeve for Torque Limiting Screwdriver, 1.5Nm	1
1471-23	Drill Bit with Quick Couling, Ø2.5mm x Length 112mm	2
1471-24	Drill Bit with Quick Couling, Ø2.8mm x Length 160mm	2
1471-25	Drill Bit with Quick Couling, Ø3.2mm x Length 147mm	2
1471-26	Drill Bit with Quick Couling, Ø4.3mm x Length 230mm	2
1471-27	Stopper for Drill Bit Ø2.8mm	2
1471-28	Stopper for Drill Bit Ø4.3mm	2
1471-29	Depth Gauge measuring upto 60mm for Pediatric Hip Plate	1
1471-30	Depth Gauge measuring upto 100mm for Pediatric Hip Plate	1
1471-31	Self-Centering Bone Holding Forcep for Pediatric Hip Plate	1
1471-32	Double Drill Guide, Ø2.5/3.5mm, for Pediatric Hip Plate	1
1471-33	Double Drill Guide, Ø3.2/4.5mm, for Pediatric Hip Plate	1
1471-34	Tap, HA3.5mm	1

Code	Set Consisting of	Units
1471-35	Tap, HA4.5mm	1
1471-36	Positioning Plate, Triangular, 90°/50°/40°	1
1471-37	Positioning Plate, Triangular, 80°/70°/30°	1
1471-38	Positioning Plate, Triangular, 100°/60°/20°	1
1471-39	Medialization Instrument for Pediatric Hip Plate	2
1471-41	Instrument Trays for Pediatric Hip Plate Instrument Set	3
1471-40	Container with PPSU Lid for Pediatric Hip Plate Instrument Set	1



USA

Auxein Inc.
1500 Nw 89th Court, Suite 107-108
Doral, Florida 33172
Tel: +1 305 395 6062
E Fax: +1 305 395 6262
Email: USoffice@auxein.com

MEXICO

Auxein México, S.A. de C.V.
Tepic 139 int 801, Colonia Roma Sur,
Alcaldía Cuauhtémoc, CDMX,
México, C.P. 06760
Tel: +521 55 7261 0318
Email: info@auxein.mx

INDIA

Auxein Medical Pvt. Ltd.
Plot No. 168-169-170, Phase-4,
Kundli Industrial Area,
HSIIDC, Sector-57, Sonapat - 131028, Haryana
Tel: +91 99106 43638 | Fax: +91 86077 70197
Email: info@auxein.com