


PELVIC C-CLAMP

PELVIC C-CLAMP



edition
2022

about us

Auxein Medical is an integrated, research based, orthopaedic Implants & Instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and innovator products that exceed the highest expectations of our customers in term of quality and safety. The company has world-class manufacturing unit established in India and serves customers in over 60 countries worldwide.

Our Achievements

US FDA 510K APPROVED •

(Arthroscopy, Spine, Nailing, Locking Plates & Screws)

EC Certificate - Class IIa & IIb •

ISO 13485:2016 •


ISO 9001:2015 •

India FDA Approved •

India MDR 2017 •

7-037-01


Pelvic C-Clamp Arm


7-037-02


Nail for Pelvic C-Clamp, Cannulated, Long


7-037-03

Nail for Pelvic C-Clamp, Cannulated, Short


7-037-04

Kirschner Wire with One Sided Trocar Tip Ø2.5mm x Length 280mm


7-037-05

Guide Handle for Ø2.5mm Kirschner Wire



7-037-06

Ratchet wrench for Nut, Hexagonal, 11mm


7-037-07

Hammer with Ø11mm Socket Wrench


7-037-08

Plier, Flat Nosed


7-037-09


Wire Cutter


7-037-11


Instrument Tray for Pelvic C-Clamp Instrument Set


7-037-10

Container for Pelvic C-Clamp Instrument Set



7-037 Pelvic C-Clamp Instrument set



7-037 Pelvic C-Clamp Instrument set

Code	Set Consisting of	Qty.
7-037-01	Pelvic C-Clamp Arm	1
7-037-02	Nail for Pelvic C-Clamp, Cannulated, Long	2
7-037-03	Nail for Pelvic C-Clamp, Cannulated, Short	2
7-037-04	Kirschner Wire with One Sided Trocar Tip Ø2.5mm x Length 280mm	10
7-037-05	Guide Handle for Ø2.5mm Kirschner Wire	1
7-037-06	Ratchet wrench for Nut, Hexagonal, 11mm	2
7-037-07	Hammer with Ø11mm Socket Wrench	1
7-037-08	Plier, Flat Nosed	1
7-037-09	Wire Cutter	1
7-037-11	Instrument Tray for Pelvic C-Clamp Instrument Set	1
7-037-10	Container for Pelvic C-Clamp Instrument Set	1



USA

Auxein Inc.
1500 Nw 89th Court, Suite 107-108
Doral, Florida 33172
Tel: +1 305 395 6062
E Fax: +1 305 395 6262
Email: USoffice@auxein.com

MEXICO

Auxein México, S.A. de C.V.
Tepic 139 int 801, Colonia Roma Sur,
Alcaldía Cuauhtémoc, CDMX,
México, C.P. 06760
Tel: +521 55 7261 0318
Email: info@auxein.mx

INDIA

Auxein Medical Pvt. Ltd.
Plot No. 168-169-170, Phase-4,
Kundli Industrial Area,
HSIIDC, Sector-57, Sonapat - 131028, Haryana
Tel: +91 99106 43638 | Fax: +91 86077 70197
Email: info@auxein.com