



PELVIC
PLATING SYSTEM

PELVIC

edition
2022

 **Auxein**[®]

about us

Auxein Medical is an integrated, research based, orthopaedic Implants & Instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and innovator products that exceed the highest expectations of our customers in term of quality and safety. The company has world-class manufacturing unit established in India and serves customers in over 60 countries worldwide.

Our Achievements

- US FDA 510K APPROVED •**
(Arthroscopy, Spine, Nailing, Locking Plates & Screws)
- EC Certificate - Class IIa & IIb •
- ISO 13485:2016 •
- ISO 9001:2015 •
- India FDA Approved •
- India MDR 2017 •

Orthopaedic Implants



Pelvic Plating System

Pelvic Plating System is a comprehensive set of plates, screws, and instrumentation for the treatment of pelvic ring and acetabular fractures. Designed to treat a wide variety of challenging pelvic fractures, the plates of the Pelvic Plating System are strategically precontoured where it may save time for the surgeon, and left noncontoured in some sections to allow for buttressing of fractures. Indication-specific plates are offered, as well as reconstruction-style plates to address a variety of fracture patterns. Enhancements to traditional pelvic instrumentation are designed to simplify surgical techniques.

Indications for Use:

- Fractures, fusions, and osteotomies of the acetabulum
- Fractures, fusions, and osteotomies of the sacrum
- Fractures, fusions, and osteotomies of the ilium
- Fractures, fusions, and osteotomies of the pelvic ring
- Sacroiliac joint dislocations
- Pubic symphysis disruptions



Golden for Right Direction



Blue for Left Direction

3.5mm Pubic Symphysis Plate



Holes	Titanium	Stainless Steel
4	13-005-04TI	13-005-04SS
6	13-005-06TI	13-005-06SS

3.5mm Posterior Wall Acetabular Fragment Plate



Holes	Titanium	Stainless Steel	Direction
9	13-006-09LTI	13-006-09LSS	Left
9	13-006-09RTI	13-006-09RSS	Right

3.5mm Curved Posterior Wall Acetabular Fragment Plate



Holes	Titanium	Stainless Steel	Direction
11	13-007-011LTI	13-007-011LSS	Left
11	13-007-011RTI	13-007-011RSS	Right

3.5mm Superior Sacro-iliac Plate



Holes	Titanium	Stainless Steel
4	13-008-04TI	13-008-04SS

3.5mm Posterior Wall Acetabular Plate



Holes	Titanium	Stainless Steel
8	13-009-08TI	13-009-08SS

3.5mm Acetabular Spring Plate



Holes	Titanium	Stainless Steel
2	13-010-02TI	13-010-02SS
3	13-010-03TI	13-010-03SS

3.5mm Reconstruction Pelvic Plate

Holes	Titanium	Stainless Steel
3	13-011-03TI	13-011-03SS
4	13-011-04TI	13-011-04SS
6	13-011-06TI	13-011-06SS
8	13-011-08TI	13-011-08SS
10	13-011-10TI	13-011-10SS
12	13-011-12TI	13-011-12SS
14	13-011-14TI	13-011-14SS
16	13-011-16TI	13-011-16SS



3.5mm Interlocking Reconstruction Plate

Holes	Titanium	Stainless Steel
11	13-012-11TI	13-012-11SS



3.5mm Quadrilateral Surface Plate

Holes	Titanium	Stainless Steel	Direction
2	13-013-02LTI	13-013-02LSS	Left
2	13-013-02RTI	13-013-02RSS	Right



3.5mm Intrapelvic Plate

Holes	Titanium	Stainless Steel	Direction
5	13-014-05LTI	13-014-05LSS	Left
5	13-014-05RTI	13-014-05RSS	Right
9	13-014-09LTI	13-014-09LSS	Left
9	13-014-09RTI	13-014-09RSS	Right



3.5mm Anterior Brim Plate

Holes	Titanium	Stainless Steel	Direction
12	13-015-12LTI	13-015-12LSS	Left
12	13-015-12RTI	13-015-12RSS	Right
14	13-015-14LTI	13-015-14LSS	Left
14	13-015-14RTI	13-015-14RSS	Right



3.5mm Low Profile Cortical Screw, (Hex Head)

Length	Titanium	Stainless Steel
10mm	13-017-010TI	13-017-010SS
12mm	13-017-012TI	13-017-012SS
14mm	13-017-014TI	13-017-014SS
16mm	13-017-016TI	13-017-016SS
18mm	13-017-018TI	13-017-018SS
20mm	13-017-020TI	13-017-020SS
22mm	13-017-022TI	13-017-022SS
24mm	13-017-024TI	13-017-024SS
26mm	13-017-026TI	13-017-026SS
28mm	13-017-028TI	13-017-028SS
30mm	13-017-030TI	13-017-030SS
32mm	13-017-032TI	13-017-032SS
34mm	13-017-034TI	13-017-034SS
36mm	13-017-036TI	13-017-036SS
38mm	13-017-038TI	13-017-038SS
40mm	13-017-040TI	13-017-040SS
45mm	13-017-045TI	13-017-045SS
50mm	13-017-050TI	13-017-050SS
55mm	13-017-055TI	13-017-055SS
60mm	13-017-060TI	13-017-060SS
65mm	13-017-065TI	13-017-065SS
70mm	13-017-070TI	13-017-070SS
75mm	13-017-075TI	13-017-075SS
80mm	13-017-080TI	13-017-080SS
85mm	13-017-085TI	13-017-085SS
90mm	13-017-090TI	13-017-090SS
95mm	13-017-095TI	13-017-095SS
100mm	13-017-100TI	13-017-100SS
105mm	13-017-105TI	13-017-105SS
110mm	13-017-110TI	13-017-110SS
115mm	13-017-115TI	13-017-115SS
120mm	13-017-120TI	13-017-120SS
125mm	13-017-125TI	13-017-125SS
130mm	13-017-130TI	13-017-130SS
135mm	13-017-135TI	13-017-135SS
140mm	13-017-140TI	13-017-140SS
145mm	13-017-145TI	13-017-145SS
150mm	13-017-150TI	13-017-150SS





4.0mm Low Profile Cancellous Screw, (Hex Head)


Length	Titanium	Stainless Steel
50mm	13-018-50TI	13-018-50SS
55mm	13-018-55TI	13-018-55SS
60mm	13-018-60TI	13-018-60SS
65mm	13-018-65TI	13-018-65SS
70mm	13-018-70TI	13-018-70SS
75mm	13-018-75TI	13-018-75SS
80mm	13-018-80TI	13-018-80SS
85mm	13-018-85TI	13-018-85SS
90mm	13-018-90TI	13-018-90SS
95mm	13-018-95TI	13-018-95SS
100mm	13-018-100TI	13-018-100SS
105mm	13-018-105TI	13-018-105SS
110mm	13-018-110TI	13-018-110SS
115mm	13-018-115TI	13-018-115SS
120mm	13-018-120TI	13-018-120SS
125mm	13-018-125TI	13-018-125SS
130mm	13-018-130TI	13-018-130SS
135mm	13-018-135TI	13-018-135SS
140mm	13-018-140TI	13-018-140SS
145mm	13-018-145TI	13-018-145SS
150mm	13-018-150TI	13-018-150SS



13-027 : Pelvic Implant Set



Screw Caddy



7-032-01

Template A, Small



7-032-02

Template B, Medium



7-032-03

Template C, Large



7-032-04

Template D, Large Pointed



7-032-05

Bone Retractor



7-032-06

Muscle Retractor



7-032-07 Periosteum Elevator, 14mm, Curved



7-032-08 Hammer for Pelvic Plate



7-032-09 Bending Forcep for Pelvic Plate



7-032-10 Reduction Forcep for Pelvic Plate



7-032-11 Multi-Angle Adjustment Forcep



7-032-12 Eccentric Multi-Angle Adjustment Forcep



7-032-13 Holding Forcep, Large



7-032-14 Holding Forcep, Small



7-032-15 Plate Cutting Forcep



7-032-16 Bone Retractor, Short



7-032-17 Bone Retractor, Long




7-032-18 Kirschner Wire, Ø1.2mm x Length 155mm



7-032-19 Kirschner Wire, Ø1.6mm x Length 215mm



7-032-20 Fixation Pin for Pelvic Plate



7-032-21 Drill Bit with Quick Coupling, Ø2.0mm x Length 145mm



7-032-22 Drill Bit with Quick Coupling, Ø2.8mm x Length 300mm



7-032-23 Quick Release Flexible Drill Bit, Ø3.5mm x Length 300mm



7-032-24 Drill Bit with Quick Coupling, Ø 3.5mm x Length 300mm



7-032-25 Drill Guide, Ø3.5mm



7-032-26 Drill Guide, Ø2.8mm



7-032-27 Soft Tissue Protector, Long, Ø7.8mm



7-032-28 Drill Guide, Ø3.5mm (for Quick Release Flexible Drill Bit)



7-032-29 Double Drill Guide, Ø2.0/2.8mm



7-032-30 Plate Bender, Large



7-032-31 Template for Pelvic Plate



7-032-32 Bone Graft Drill Assembly, Ø7.0mm



7-032-33 Depth Gauge measuring upto 150mm for Pelvic Plate



7-032-34 Reduction Clamp



7-032-35 Depth Gauge measuring upto 65mm for Pelvic Plate



7-032-36 Plate Reduction Clamp



7-032-37 Screwdriver Shaft Quick Coupling, Hex 2.5 x Length 150mm



7-032-38 Screwdriver Shaft Quick Coupling, Hex 2.5 x Length 228mm



7-032-39 Screwdriver Shaft Quick Coupling, Hex 2.5 x Length 250mm




7-032-40 Ratchet Handle for Pelvic Plate



7-032-41 2.8mm Offset Double Drill Guide







7-032-44 Instrument Trays for Container-I for Pelvic Plating System Instrument set



7-032-45 Instrument Trays for Container-II for Pelvic Plating System Instrument set



7-032 Pelvic Plating System Instrument set




7-032-44

Instrument Trays for Container-I
for Pelvic Plating System Instrument set

7-032-45

Instrument Trays for Container-II
for Pelvic Plating System Instrument set

7-032-42 Container-I for Pelvic Plating System Instrument set



7-032-43 Container-II for Pelvic Plating System Instrument set



7-032 Pelvic Instrument Set

Item Code	Item Name	Unit
7-032-01	Template A, Small	1
7-032-02	Template B, Medium	1
7-032-03	Template C, Large	1
7-032-04	Template D, Large Pointed	1
7-032-05	Bone Retractor	1
7-032-06	Muscle Retractor	1
7-032-07	Periosteum Elevator, 14mm, Curved	1
7-032-08	Hammer for Pelvic Plate	1
7-032-09	Bending Forcep for Pelvic Plate	1
7-032-10	Reduction Forcep for Pelvic Plate	1
7-032-11	Multi-Angle Adjustment Forcep	1
7-032-12	Eccentric Multi-Angle Adjustment Forcep	1
7-032-13	Holding Forcep, Large	1
7-032-14	Holding Forcep, Small	1
7-032-15	Plate Cutting Forcep	1
7-032-16	Bone Retractor, Short	1
7-032-17	Bone Retractor, Long	1
7-032-18	Kirschner Wire, Ø1.2mm x Length 155mm	5
7-032-19	Kirschner Wire, Ø1.6mm x Length 215mm	5
7-032-20	Fixation Pin for Pelvic Plate	5
7-032-21	Drill Bit with Quick Coupling, Ø2.0mm x Length 145mm	1
7-032-22	Drill Bit with Quick Coupling, Ø2.8mm x Length 300mm	2
7-032-23	Quick Release Flexible Drill Bit, Ø3.5mm x Length 300mm	2
7-032-24	Drill Bit with Quick Coupling, Ø 3.5mm x Length 300mm	2
7-032-25	Drill Guide, Ø3.5mm	1
7-032-26	Drill Guide, Ø2.8mm	1
7-032-27	Soft Tissue Protector, Long, Ø7.8mm	1
7-032-28	Drill Guide, Ø3.5mm (for Quick Release Flexible Drill Bit)	1
7-032-29	Double Drill Guide, Ø2.0/2.8mm	1
7-032-30	Plate Bender, Large	2
7-032-31	Template for Pelvic Plate	1
7-032-32	Bone Graft Drill Assembly, Ø7.0mm	1
7-032-33	Depth Gauge measuring upto 150mm for Pelvic Plate	1
7-032-34	Reduction Clamp	1
7-032-35	Depth Gauge measuring upto 65mm for Pelvic Plate	1

Item Code	Item Name	Unit
7-032-36	Plate Reduction Clamp	1
7-032-37	Screwdriver Shaft Quick Coupling, Hex 2.5 x Length 150mm	1
7-032-38	Screwdriver Shaft Quick Coupling, Hex 2.5 x Length 228mm	1
7-032-39	Screwdriver Shaft Quick Coupling, Hex 2.5 x Length 250mm	1
7-032-40	Ratchet Handle for Pelvic Plate	1
7-032-41	2.8mm Offset Double Drill Guide	1
7-032-44	Instrument Trays for Container-I for Pelvic Plating System Instrument set	3
7-032-45	Instrument Trays for Container-II for Pelvic Plating System Instrument set	2
7-032-42	Container-I for Pelvic Plating System Instrument set	1
7-032-43	Container-II for Pelvic Plating System Instrument set	1



USA

Auxein Inc.
1500 Nw 89th Court, Suite 107-108
Doral, Florida 33172
Tel: +1 305 395 6062
E Fax: +1 305 395 6262
Email: USoffice@auxein.com

MEXICO

Auxein México, S.A. de C.V.
Tepic 139 int 801, Colonia Roma Sur,
Alcaldía Cuauhtémoc, CDMX,
México, C.P. 06760
Tel: +521 55 7261 0318
Email: info@auxein.mx

INDIA

Auxein Medical Pvt. Ltd.
Plot No. 168-169-170, Phase-4,
Kundli Industrial Area,
HSIIDC, Sector-57, Sonapat - 131028, Haryana
Tel: +91 99106 43638 | Fax: +91 86077 70197
Email: info@auxein.com