






U-Spine color

				
Pink 8mm	Type II 10mm	Green 12mm	Blue 14mm	Brown 16mm

U-Spine








Material:

Wrought titanium 6-aluminium 4-vanadium alloy according to ISO 5832-3

The **U-Spine** implant is delivered in sterile packaging.


During the course of the clinical trial, the wings were modified slightly for ease of stacking two devices at adjacent levels. The holes in the wings were also removed. The modification was not the result of any clinical problems, safety issues or adverse events, product complaints, or surgeon requests from within or outside the United States. As this modification was minor, it did not affect the mechanical behavior of the device or the anticipated clinical outcome.

U-Spine color

				
Pink 8mm	Type II 10mm	Green 12mm	Blue 14mm	Brown 16mm

4-073-08TI

U-Spine, 8mm


4-073-10TI

U-Spine, 10mm


4-073-12TI


U-Spine, 12mm


4-073-14TI

U-Spine, 14mm


4-073-16TI

U-Spine, 16mm



7-098-01 U-Spine Template Inserter, 8mm



7-098-02 U-Spine Template Inserter, 10mm



7-098-03 U-Spine Template Inserter, 12mm




7-098-04 U-Spine Template Inserter, 14mm




7-098-05 U-Spine Template Inserter, 16mm



7-098-06 U-Spine Bending Plier




7-098-07 U-Spine Instrument and Implant box with lid




7-098-08 U-Spine Crimping Plier




7-098-09 U-Spine Mallet



7-098 U-Spine Instrument and Implant Set



7-098 U-Spine Instrument and Implant Set


Code	Set Consisting of	Qty.
4-073-08TI	U-Spine, 8mm	3
4-073-10TI	U-Spine, 10mm	3
4-073-12TI	U-Spine, 12mm	3
4-073-14TI	U-Spine, 14mm	3
4-073-16TI	U-Spine, 16mm	3
7-098-01	U-Spine Template Inserters, 8mm	1
7-098-02	U-Spine Template Inserters, 10mm	1
7-098-03	U-Spine Template Inserters, 12mm	1
7-098-04	U-Spine Template Inserters, 14mm	1
7-098-05	U-Spine Template Inserters, 16mm	1
7-098-06	U-Spine Bending Plier	1
7-098-07	U-Spine Instrument and Implant box with lid	1
7-098-08	U-Spine Crimping Plier	1
7-098-09	U-Spine Mallet	1