


VERTAUX-MIS SCREW

# VERTAUX-MIS

Screw Spine System



edition  
2023

 **Auxein**<sup>®</sup>





# about us

Auxein Medical is an integrated, research based, orthopaedic Implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in india and serves customers in over 75 countries worldwide.






## Our Achievements




## **Orthopaedic Implants**

# VERTAUX MIS SCREW SYSTEM

Colour Parameters as per Ø

				
Silver Ø 5.0	Green Ø 5.5	Blue Ø 6.0	Pink Ø 6.5	Brown Ø 7.0



### Integrated Long Arm

- More stable than extended sleeve
- More convenient for rod implantation & cap tightening

### Break-off Tail Tabs

- Protect the long arms from deformation during screw insertion

### Full Round


- Better strength


### Dual-core Screw Design

- Enlarged core at proximal screw and thickened thread
- Enhanced purchase with pedicle cortex

### Double-thread Design


- Faster screw insertion



**VERTAUX MIS Polyaxial Screw with Cap**


Code	Dia X Length
<b>4-109-01TI</b>	Ø5mm x 30mm
<b>4-109-02TI</b>	Ø5mm x 35mm
<b>4-109-03TI</b>	Ø5mm x 40mm
<b>4-109-04TI</b>	Ø5mm x 45mm
<b>4-110-01TI</b>	Ø5.5mm x 30mm
<b>4-110-02TI</b>	Ø5.5mm x 35mm
<b>4-110-03TI</b>	Ø5.5mm x 40mm
<b>4-110-04TI</b>	Ø5.5mm x 45mm
<b>4-110-05TI</b>	Ø5.5mm x 50mm
<b>4-111-01TI</b>	Ø6mm x 30mm
<b>4-111-02TI</b>	Ø6mm x 35mm
<b>4-111-03TI</b>	Ø6mm x 40mm
<b>4-111-04TI</b>	Ø6mm x 45mm
<b>4-111-05TI</b>	Ø6mm x 50mm
<b>4-111-06TI</b>	Ø6mm x 55mm


Code	Dia X Length
<b>4-112-01TI</b>	Ø6.5mm x 30mm
<b>4-112-02TI</b>	Ø6.5mm x 35mm
<b>4-112-03TI</b>	Ø6.5mm x 40mm
<b>4-112-04TI</b>	Ø6.5mm x 45mm
<b>4-112-05TI</b>	Ø6.5mm x 50mm
<b>4-112-06TI</b>	Ø6.5mm x 55mm
<b>4-112-07TI</b>	Ø6.5mm x 60mm
<b>4-113-01TI</b>	Ø7mm x 35mm
<b>4-113-02TI</b>	Ø7mm x 40mm
<b>4-113-03TI</b>	Ø7mm x 45mm
<b>4-113-04TI</b>	Ø7mm x 50mm
<b>4-113-05TI</b>	Ø7mm x 55mm
<b>4-113-06TI</b>	Ø7mm x 60mm

**VERTAUX MIS Polyaxial Screw with Cap**

**STERILE**

Code	Dia X Length
4-109-01TI-S	Ø5mm x 30mm
4-109-02TI-S	Ø5mm x 35mm
4-109-03TI-S	Ø5mm x 40mm
4-109-04TI-S	Ø5mm x 45mm
4-110-01TI-S	Ø5.5mm x 30mm
4-110-02TI-S	Ø5.5mm x 35mm
4-110-03TI-S	Ø5.5mm x 40mm
4-110-04TI-S	Ø5.5mm x 45mm
4-110-05TI-S	Ø5.5mm x 50mm
4-111-01TI-S	Ø6mm x 30mm
4-111-02TI-S	Ø6mm x 35mm
4-111-03TI-S	Ø6mm x 40mm
4-111-04TI-S	Ø6mm x 45mm
4-111-05TI-S	Ø6mm x 50mm
4-111-06TI-S	Ø6mm x 55mm

Code	Dia X Length
4-112-01TI-S	Ø6.5mm x 30mm
4-112-02TI-S	Ø6.5mm x 35mm
4-112-03TI-S	Ø6.5mm x 40mm
4-112-04TI-S	Ø6.5mm x 45mm
4-112-05TI-S	Ø6.5mm x 50mm
4-112-06TI-S	Ø6.5mm x 55mm
4-112-07TI-S	Ø6.5mm x 60mm
4-113-01TI-S	Ø7mm x 35mm
4-113-02TI-S	Ø7mm x 40mm
4-113-03TI-S	Ø7mm x 45mm
4-113-04TI-S	Ø7mm x 50mm
4-113-05TI-S	Ø7mm x 55mm
4-113-06TI-S	Ø7mm x 60mm



**VERTAUX MIS Monoaxial Screw with Cap**


Code	Dia X Length
4-117-01TI	Ø5mm x 30mm
4-117-02TI	Ø5mm x 35mm
4-117-03TI	Ø5mm x 40mm
4-117-04TI	Ø5mm x 45mm
4-118-01TI	Ø5.5mm x 30mm
4-118-02TI	Ø5.5mm x 35mm
4-118-03TI	Ø5.5mm x 40mm
4-118-04TI	Ø5.5mm x 45mm
4-118-05TI	Ø5.5mm x 50mm
4-119-01TI	Ø6mm x 30mm
4-119-02TI	Ø6mm x 35mm
4-119-03TI	Ø6mm x 40mm
4-119-04TI	Ø6mm x 45mm
4-119-05TI	Ø6mm x 50mm
4-119-06TI	Ø6mm x 55mm


Code	Dia X Length
4-120-01TI	Ø6.5mm x 30mm
4-120-02TI	Ø6.5mm x 35mm
4-120-03TI	Ø6.5mm x 40mm
4-120-04TI	Ø6.5mm x 45mm
4-120-05TI	Ø6.5mm x 50mm
4-120-06TI	Ø6.5mm x 55mm
4-120-07TI	Ø6.5mm x 60mm
4-121-01TI	Ø7mm x 35mm
4-121-02TI	Ø7mm x 40mm
4-121-03TI	Ø7mm x 45mm
4-121-04TI	Ø7mm x 50mm
4-121-05TI	Ø7mm x 55mm
4-121-06TI	Ø7mm x 60mm

**VERTAUX MIS Inner Screw Cap**


Code

**4-114-01TI**

**Star Head**

**VERTAUX MIS Monoaxial Screw with Cap**

**STERILE**

Code	Dia X Length
4-117-01TI-S	Ø5mm x 30mm
4-117-02TI-S	Ø5mm x 35mm
4-117-03TI-S	Ø5mm x 40mm
4-117-04TI-S	Ø5mm x 45mm
4-118-01TI-S	Ø5.5mm x 30mm
4-118-02TI-S	Ø5.5mm x 35mm
4-118-03TI-S	Ø5.5mm x 40mm
4-118-04TI-S	Ø5.5mm x 45mm
4-118-05TI-S	Ø5.5mm x 50mm
4-119-01TI-S	Ø6mm x 30mm
4-119-02TI-S	Ø6mm x 35mm
4-119-03TI-S	Ø6mm x 40mm
4-119-04TI-S	Ø6mm x 45mm
4-119-05TI-S	Ø6mm x 50mm
4-119-06TI-S	Ø6mm x 55mm

Code	Dia X Length
4-120-01TI-S	Ø6.5mm x 30mm
4-120-02TI-S	Ø6.5mm x 35mm
4-120-03TI-S	Ø6.5mm x 40mm
4-120-04TI-S	Ø6.5mm x 45mm
4-120-05TI-S	Ø6.5mm x 50mm
4-120-06TI-S	Ø6.5mm x 55mm
4-120-07TI-S	Ø6.5mm x 60mm
4-121-01TI-S	Ø7mm x 35mm
4-121-02TI-S	Ø7mm x 40mm
4-121-03TI-S	Ø7mm x 45mm
4-121-04TI-S	Ø7mm x 50mm
4-121-05TI-S	Ø7mm x 55mm
4-121-06TI-S	Ø7mm x 60mm

**VERTAUX MIS Inner Screw Cap**
**STERILE**

Code

**4-114-01TI**

**Star Head**

**VERTAUX MIS Pre Bend Rod Ø5.5mm**


Code	Length	Code	Length
4-122-01	30mm	4-122-50	275mm
4-122-02	35mm	4-122-51	280mm
4-122-03	40mm	4-122-52	285mm
4-122-04	45mm	4-122-53	290mm
4-122-05	50mm	4-122-54	295mm
4-122-06	55mm	4-122-55	300mm
4-122-07	60mm	4-122-56	305mm
4-122-08	65mm	4-122-57	310mm
4-122-09	70mm	4-122-58	315mm
4-122-10	75mm	4-122-59	320mm
4-122-11	80mm	4-122-60	325mm
4-122-12	85mm	4-122-61	330mm
4-122-13	90mm	4-122-62	335mm
4-122-14	95mm	4-122-63	340mm
4-122-15	100mm	4-122-64	345mm
4-122-16	105mm	4-122-65	350mm
4-122-17	110mm	4-122-66	355mm
4-122-18	115mm	4-122-67	360mm
4-122-19	120mm	4-122-68	365mm
4-122-20	125mm	4-122-69	370mm
4-122-21	130mm	4-122-70	375mm
4-122-22	135mm	4-122-71	380mm
4-122-23	140mm	4-122-72	385mm
4-122-24	145mm	4-122-73	390mm
4-122-25	150mm	4-122-74	395mm
4-122-26	155mm	4-122-75	400mm
4-122-27	160mm	4-122-76	405mm
4-122-28	165mm	4-122-77	410mm
4-122-29	170mm	4-122-78	415mm
4-122-30	175mm	4-122-79	420mm
4-122-31	180mm	4-122-80	425mm
4-122-32	185mm	4-122-81	430mm
4-122-33	190mm	4-122-82	435mm
4-122-34	195mm	4-122-83	440mm
4-122-35	200mm	4-122-84	445mm
4-122-36	205mm	4-122-85	450mm
4-122-37	210mm	4-122-86	455mm
4-122-38	215mm	4-122-87	460mm
4-122-39	220mm	4-122-88	465mm
4-122-40	225mm	4-122-89	470mm
4-122-41	230mm	4-122-90	475mm
4-122-42	235mm	4-122-91	480mm
4-122-43	240mm	4-122-92	485mm
4-122-44	245mm	4-122-93	490mm
4-122-45	250mm	4-122-94	495mm
4-122-46	255mm	4-122-95	500mm
4-122-47	260mm		
4-122-48	265mm		
4-122-49	270mm		

**VERTAUX MIS Pre Bend Rod Ø5.5mm**
**STERILE**

Code	Length	Code	Length
4-122-01-S	30mm	4-122-50-S	275mm
4-122-02-S	35mm	4-122-51-S	280mm
4-122-03-S	40mm	4-122-52-S	285mm
4-122-04-S	45mm	4-122-53-S	290mm
4-122-05-S	50mm	4-122-54-S	295mm
4-122-06-S	55mm	4-122-55-S	300mm
4-122-07-S	60mm	4-122-56-S	305mm
4-122-08-S	65mm	4-122-57-S	310mm
4-122-09-S	70mm	4-122-58-S	315mm
4-122-10-S	75mm	4-122-59-S	320mm
4-122-11-S	80mm	4-122-60-S	325mm
4-122-12-S	85mm	4-122-61-S	330mm
4-122-13-S	90mm	4-122-62-S	335mm
4-122-14-S	95mm	4-122-63-S	340mm
4-122-15-S	100mm	4-122-64-S	345mm
4-122-16-S	105mm	4-122-65-S	350mm
4-122-17-S	110mm	4-122-66-S	355mm
4-122-18-S	115mm	4-122-67-S	360mm
4-122-19-S	120mm	4-122-68-S	365mm
4-122-20-S	125mm	4-122-69-S	370mm
4-122-21-S	130mm	4-122-70-S	375mm
4-122-22-S	135mm	4-122-71-S	380mm
4-122-23-S	140mm	4-122-72-S	385mm
4-122-24-S	145mm	4-122-73-S	390mm
4-122-25-S	150mm	4-122-74-S	395mm
4-122-26-S	155mm	4-122-75-S	400mm
4-122-27-S	160mm	4-122-76-S	405mm
4-122-28-S	165mm	4-122-77-S	410mm
4-122-29-S	170mm	4-122-78-S	415mm
4-122-30-S	175mm	4-122-79-S	420mm
4-122-31-S	180mm	4-122-80-S	425mm
4-122-32-S	185mm	4-122-81-S	430mm
4-122-33-S	190mm	4-122-82-S	435mm
4-122-34-S	195mm	4-122-83-S	440mm
4-122-35-S	200mm	4-122-84-S	445mm
4-122-36-S	205mm	4-122-85-S	450mm
4-122-37-S	210mm	4-122-86-S	455mm
4-122-38-S	215mm	4-122-87-S	460mm
4-122-39-S	220mm	4-122-88-S	465mm
4-122-40-S	225mm	4-122-89-S	470mm
4-122-41-S	230mm	4-122-90-S	475mm
4-122-42-S	235mm	4-122-91-S	480mm
4-122-43-S	240mm	4-122-92-S	485mm
4-122-44-S	245mm	4-122-93-S	490mm
4-122-45-S	250mm	4-122-94-S	495mm
4-122-46-S	255mm	4-122-95-S	500mm
4-122-47-S	260mm		
4-122-48-S	265mm		
4-122-49-S	270mm		



**VERTAUX MIS Straight Rod Ø5.5mm**


Code	Length	Code	Length
4-130-01TI	30mm	4-130-50TI	275mm
4-130-02TI	35mm	4-130-51TI	280mm
4-130-03TI	40mm	4-130-52TI	285mm
4-130-04TI	45mm	4-130-53TI	290mm
4-130-05TI	50mm	4-130-54TI	295mm
4-130-06TI	55mm	4-130-55TI	300mm
4-130-07TI	60mm	4-130-56TI	305mm
4-130-08TI	65mm	4-130-57TI	310mm
4-130-09TI	70mm	4-130-58TI	315mm
4-130-10TI	75mm	4-130-59TI	320mm
4-130-11TI	80mm	4-130-60TI	325mm
4-130-12TI	85mm	4-130-61TI	330mm
4-130-13TI	90mm	4-130-62TI	335mm
4-130-14TI	95mm	4-130-63TI	340mm
4-130-15TI	100mm	4-130-64TI	345mm
4-130-16TI	105mm	4-130-65TI	350mm
4-130-17TI	110mm	4-130-66TI	355mm
4-130-18TI	115mm	4-130-67TI	360mm
4-130-19TI	120mm	4-130-68TI	365mm
4-130-20TI	125mm	4-130-69TI	370mm
4-130-21TI	130mm	4-130-70TI	375mm
4-130-22TI	135mm	4-130-71TI	380mm
4-130-23TI	140mm	4-130-72TI	385mm
4-130-24TI	145mm	4-130-73TI	390mm
4-130-25TI	150mm	4-130-74TI	395mm
4-130-26TI	155mm	4-130-75TI	400mm
4-130-27TI	160mm	4-130-76TI	405mm
4-130-28TI	165mm	4-130-77TI	410mm
4-130-29TI	170mm	4-130-78TI	415mm
4-130-30TI	175mm	4-130-79TI	420mm
4-130-31TI	180mm	4-130-80TI	425mm
4-130-32TI	185mm	4-130-81TI	430mm
4-130-33TI	190mm	4-130-82TI	435mm
4-130-34TI	195mm	4-130-83TI	440mm
4-130-35TI	200mm	4-130-84TI	445mm
4-130-36TI	205mm	4-130-85TI	450mm
4-130-37TI	210mm	4-130-86TI	455mm
4-130-38TI	215mm	4-130-87TI	460mm
4-130-39TI	220mm	4-130-88TI	465mm
4-130-40TI	225mm	4-130-89TI	470mm
4-130-41TI	230mm	4-130-90TI	475mm
4-130-42TI	235mm	4-130-91TI	480mm
4-130-43TI	240mm	4-130-92TI	485mm
4-130-44TI	245mm	4-130-93TI	490mm
4-130-45TI	250mm	4-130-94TI	495mm
4-130-46TI	255mm	4-130-95TI	500mm
4-130-47TI	260mm		
4-130-48TI	265mm		
4-130-49TI	270mm		




**VERTAUX MIS Straight Rod Ø5.5mm**
**STERILE**

Code	Length	Code	Length
4-130-01TI-S	30mm	4-130-50TI-S	275mm
4-130-02TI-S	35mm	4-130-51TI-S	280mm
4-130-03TI-S	40mm	4-130-52TI-S	285mm
4-130-04TI-S	45mm	4-130-53TI-S	290mm
4-130-05TI-S	50mm	4-130-54TI-S	295mm
4-130-06TI-S	55mm	4-130-55TI-S	300mm
4-130-07TI-S	60mm	4-130-56TI-S	305mm
4-130-08TI-S	65mm	4-130-57TI-S	310mm
4-130-09TI-S	70mm	4-130-58TI-S	315mm
4-130-10TI-S	75mm	4-130-59TI-S	320mm
4-130-11TI-S	80mm	4-130-60TI-S	325mm
4-130-12TI-S	85mm	4-130-61TI-S	330mm
4-130-13TI-S	90mm	4-130-62TI-S	335mm
4-130-14TI-S	95mm	4-130-63TI-S	340mm
4-130-15TI-S	100mm	4-130-64TI-S	345mm
4-130-16TI-S	105mm	4-130-65TI-S	350mm
4-130-17TI-S	110mm	4-130-66TI-S	355mm
4-130-18TI-S	115mm	4-130-67TI-S	360mm
4-130-19TI-S	120mm	4-130-68TI-S	365mm
4-130-20TI-S	125mm	4-130-69TI-S	370mm
4-130-21TI-S	130mm	4-130-70TI-S	375mm
4-130-22TI-S	135mm	4-130-71TI-S	380mm
4-130-23TI-S	140mm	4-130-72TI-S	385mm
4-130-24TI-S	145mm	4-130-73TI-S	390mm
4-130-25TI-S	150mm	4-130-74TI-S	395mm
4-130-26TI-S	155mm	4-130-75TI-S	400mm
4-130-27TI-S	160mm	4-130-76TI-S	405mm
4-130-28TI-S	165mm	4-130-77TI-S	410mm
4-130-29TI-S	170mm	4-130-78TI-S	415mm
4-130-30TI-S	175mm	4-130-79TI-S	420mm
4-130-31TI-S	180mm	4-130-80TI-S	425mm
4-130-32TI-S	185mm	4-130-81TI-S	430mm
4-130-33TI-S	190mm	4-130-82TI-S	435mm
4-130-34TI-S	195mm	4-130-83TI-S	440mm
4-130-35TI-S	200mm	4-130-84TI-S	445mm
4-130-36TI-S	205mm	4-130-85TI-S	450mm
4-130-37TI-S	210mm	4-130-86TI-S	455mm
4-130-38TI-S	215mm	4-130-87TI-S	460mm
4-130-39TI-S	220mm	4-130-88TI-S	465mm
4-130-40TI-S	225mm	4-130-89TI-S	470mm
4-130-41TI-S	230mm	4-130-90TI-S	475mm
4-130-42TI-S	235mm	4-130-91TI-S	480mm
4-130-43TI-S	240mm	4-130-92TI-S	485mm
4-130-44TI-S	245mm	4-130-93TI-S	490mm
4-130-45TI-S	250mm	4-130-94TI-S	495mm
4-130-46TI-S	255mm	4-130-95TI-S	500mm
4-130-47TI-S	260mm		
4-130-48TI-S	265mm		
4-130-49TI-S	270mm		



**7-145-01** Guide Wire,  $\varnothing 1.5 \times 500$  for Vertaux MIS Screw



**7-145-02** Inserted needle for Vertaux MIS Screw



**7-145-03** Multiple Level Dilator 5.3mm for Vertaux MIS Screw



**7-145-04** Multiple Level Dilator 9.4mm for Vertaux MIS Screw



**7-145-05** Multiple Level Dilator 12.8mm for Vertaux MIS Screw



**7-145-06** Multiple Level Dilator 16.8mm for Vertaux MIS Screw



**7-145-07** Multiple Level Dilator 20.8mm for Vertaux MIS Screw



**7-145-08** Cannulated Tap Ø1.5/Ø5.0 for Vertaux MIS Screw



**7-145-09** Cannulated Tap Ø1.5/Ø5.5 for Vertaux MIS Screw



**7-145-10** Cannulated Tap Ø1.5/Ø6.0 for Vertaux MIS Screw



**7-145-11** Cannulated Tap Ø1.5/Ø6.5 for Vertaux MIS Screw



**7-145-12** Cannulated Tap Ø1.5/Ø7.0 for Vertaux MIS Screw





**7-145-13** Cannulated Tap  $\varnothing 1.5/\varnothing 7.5$  for Vertaux MIS Screw



**7-145-14** Cannulated Quick Change Handle for Vertaux MIS Screw



**7-145-15** Awl for Vertaux MIS Screw



**7-145-16** Bone Reamer  $\varnothing 8$  for Vertaux MIS Screw



**7-145-17** Reduction Multi-axial Screwdriver for Vertaux MIS Screw



**7-145-18** Counter Torque for Vertaux MIS Screw



**7-145-19** Ratchet Handle for Vertaux MIS Screw



**7-145-20** Screwdriver for Vertaux MIS Screw



**7-145-21** Break Off Iron for Vertaux MIS Screw



**7-145-22** Plug Screwdriver (Self Hold with Inner Locking Pin) for Vertaux MIS Screw




**7-145-23** Screwdriver Shaft for Vertaux MIS Screw




**7-145-24** Torque-limiting Handle, 9Nm for Vertaux MIS Screw



**7-145-25** Rod Bender for Vertaux MIS Screw



**7-145-26** Rod Holder, Straight for Vertaux MIS Screw



**7-145-27** Rod Holder, Long for Vertaux MIS Screw



**7-145-28** Break Off Iron for Vertaux MIS Screw



**7-145-29** Reduction Sleeve for Vertaux MIS Screw



**7-145-30** Reduction Extension A for Vertaux MIS Screw



**7-145-31** Reduction Extension B for Vertaux MIS Screw




**7-145-32** Torx Screwdriver for Vertaux MIS Screw



**7-145-33** Distractor/Compressor for Vertaux MIS Screw



**7-145-34** Puller for Vertaux MIS Screw




**7-145-35** Screwdriver SW11 for Vertaux MIS Screw



**7-145-36** Rod Gauge for Vertaux MIS Screw



**7-145-37**      Template for Vertaux MIS Screw



**7-145-38**      Rod Inserter for Vertaux MIS Screw



**7-145-39**      Rod Pusher for Vertaux MIS Screw



**7-145-40**      Tissue Dissector for Vertaux MIS Screw




**7-145-41**      Instrument Trays for Vertaux MIS Screw



**7-145-42**      Container for Vertaux MIS Screw



**7-145 Vertaux MIS Instrument Set**



**7-145 Vertaux MIS Instrument Set**

Code	Set Consisting of	Qty.
7-145-01	Guide Wire, Ø1.5X500 for Vertaux MIS Screw	6
7-145-02	Inserted needle for Vertaux MIS Screw	2
7-145-03	Multiple Level Dilator 5.3mm for Vertaux MIS Screw	1
7-145-04	Multiple Level Dilator 9.4mm for Vertaux MIS Screw	1
7-145-05	Multiple Level Dilator 12.8mm for Vertaux MIS Screw	1
7-145-06	Multiple Level Dilator 16.8mm for Vertaux MIS Screw	1
7-145-07	Multiple Level Dilator 20.8mm for Vertaux MIS Screw	1
7-145-08	Cannulated Tap Ø1.5/Ø5.0 for Vertaux MIS Screw	1
7-145-09	Cannulated Tap Ø1.5/Ø5.5 for Vertaux MIS Screw	1
7-145-10	Cannulated Tap Ø1.5/Ø6.0 for Vertaux MIS Screw	1
7-145-11	Cannulated Tap Ø1.5/Ø6.5 for Vertaux MIS Screw	1
7-145-12	Cannulated Tap Ø1.5/Ø7.0 for Vertaux MIS Screw	1
7-145-13	Cannulated Tap Ø1.5/Ø7.5 for Vertaux MIS Screw	1
7-145-14	Cannulated Quick Change Handle for Vertaux MIS Screw	2
7-145-15	Awl for Vertaux MIS Screw	1
7-145-16	Bone Reamer Ø8 for Vertaux MIS Screw	1
7-145-17	Reduction Multi-axial Screwdriver for Vertaux MIS Screw	2
7-145-18	Counter Torque for Vertaux MIS Screw	1
7-145-19	Ratchet Handle for Vertaux MIS Screw	1
7-145-20	Screwdriver for Vertaux MIS Screw	1
7-145-21	Break Off Iron for Vertaux MIS Screw	1
7-145-22	Plug Screwdriver (Self Hold with Inner Locking Pin) for Vertaux MIS Screw	2
7-145-23	Screwdriver Shaft for Vertaux MIS Screw	2
7-145-24	Torque-limiting Handle, 9Nm for Vertaux MIS Screw	1
7-145-25	Rod Bender for Vertaux MIS Screw	1
7-145-26	Rod Holder, Straight for Vertaux MIS Screw	1
7-145-27	Rod Holder, Long for Vertaux MIS Screw	1
7-145-28	Break Off Iron for Vertaux MIS Screw	1
7-145-29	Reduction Sleeve for Vertaux MIS Screw	1
7-145-30	Reduction Extension A for Vertaux MIS Screw	2
7-145-31	Reduction Extension B for Vertaux MIS Screw	2
7-145-32	Torx Screwdriver for Vertaux MIS Screw	1
7-145-33	Distractor/Compressor for Vertaux MIS Screw	1
7-145-34	Puller for Vertaux MIS Screw	1
7-145-35	Screwdriver SW11 for Vertaux MIS Screw	1
7-145-36	Rod Gauge for Vertaux MIS Screw	1
7-145-37	Template for Vertaux MIS Screw	1
7-145-38	Rod Inserter for Vertaux MIS Screw	1
7-145-39	Rod Pusher for Vertaux MIS Screw	1
7-145-40	Tissue Dissector for Vertaux MIS Screw	1
7-145-41	Instrument Trays for Vertaux MIS Screw	4
7-145-42	Container for Vertaux MIS Screw	1



**USA**

Auxein Inc.  
1500 Nw 89th Court, Suite 107-108  
Doral, Florida 33172  
Tel: +1 305 395 6062  
E Fax: +1 305 395 6262  
Email: USoffice@auxein.com

**MEXICO**

Auxein México, S.A. de C.V.  
Tepic 139 int 801, Colonia Roma Sur,  
Alcaldía Cuauhtémoc, CDMX,  
México, C.P. 06760  
Tel: +521 55 7261 0318  
Email: info@auxein.mx

**INDIA**

Auxein Medical Pvt. Ltd.  
Plot No. 168-169-170, Phase-4,  
Kundli Industrial Area,  
HSIIDC, Sector-57, Sonapat - 131028, Haryana  
Tel: +91 99106 43638 | Fax: +91 86077 70197  
Email: info@auxein.com