


VERTAUX Basico Polyaxial Screw

VERTAUX Basico Polyaxial Screw System available in standard configuration is a top loading low profile posterior spinal fixation system designed to provide stabilization as an adjunct to fusion. VERTAUX Basico Polyaxial Screw System is a comprehensive system which consists of pedicle screws, rods, cross-links and locking caps. All of the components are available on a variety of sizes to more closely match the patient's anatomy.

Vertaux Basico Polyaxial Screw with Cap

Ø Dia	Length (mm)	Titanium
4.5mm	25mm	4-066-025TI
4.5mm	30mm	4-066-030TI
4.5mm	35mm	4-066-035TI
4.5mm	40mm	4-066-040TI
4.5mm	45mm	4-066-045TI
4.5mm	50mm	4-066-050TI
4.5mm	55mm	4-066-055TI
5.0mm	25mm	4-067-025TI
5.0mm	30mm	4-067-030TI
5.0mm	35mm	4-067-035TI
5.0mm	40mm	4-067-040TI
5.0mm	45mm	4-067-045TI
5.0mm	50mm	4-067-050TI
5.0mm	55mm	4-067-055TI
5.5mm	25mm	4-068-025TI
5.5mm	30mm	4-068-030TI
5.5mm	35mm	4-068-035TI
5.5mm	40mm	4-068-040TI
5.5mm	45mm	4-068-045TI
5.5mm	50mm	4-068-050TI
5.5mm	55mm	4-068-055TI
6.0mm	25mm	4-069-025TI
6.0mm	30mm	4-069-030TI
6.0mm	35mm	4-069-035TI
6.0mm	40mm	4-069-040TI
6.0mm	45mm	4-069-045TI
6.0mm	50mm	4-069-050TI
6.0mm	55mm	4-069-055TI
6.5mm	25mm	4-070-025TI
6.5mm	30mm	4-070-030TI
6.5mm	35mm	4-070-035TI
6.5mm	40mm	4-070-040TI
6.5mm	45mm	4-070-045TI
6.5mm	50mm	4-070-050TI
6.5mm	55mm	4-070-055TI
7.0mm	25mm	4-071-025TI
7.0mm	30mm	4-071-030TI
7.0mm	35mm	4-071-035TI
7.0mm	40mm	4-071-040TI
7.0mm	45mm	4-071-045TI
7.0mm	50mm	4-071-050TI
7.0mm	55mm	4-071-055TI



Vertaux Basico Polyaxial Screw Rod

Length (mm)	Titanium
400mm	4-072-400TI



VERTAUX - Pedicle Screws Cap, M9

Product Specification	Titanium
M9	SP-303



7-108-02

Vertaux Basico Reducer

**7-108-03**

Vertaux Basico Screwdriver

**7-108-04**

Vertaux Basico Screw Shaft with Quick Coupling

**7-108-05**



Implant Tray for Vertaux Basico Screw

**7-108-06**

Container with Lid for Vertaux Basico Screw



7-108-000 Vertaux Basico Screw Implant & Instrument Set



7-108-000 Vertaux Basico Screw Implant & Instrument Set

Code	Set Consisting of	Unit
4-066-030TI	Vertaux Basico Polyaxial Screw with Cap, Ø4.5 x 30mm	10
4-066-035TI	Vertaux Basico Polyaxial Screw with Cap, Ø4.5 x 35mm	10
4-066-040TI	Vertaux Basico Polyaxial Screw with Cap, Ø4.5 x 40mm	10
4-066-045TI	Vertaux Basico Polyaxial Screw with Cap, Ø4.5 x 45mm	10
4-066-050TI	Vertaux Basico Polyaxial Screw with Cap, Ø4.5 x 50mm	10
4-066-055TI	Vertaux Basico Polyaxial Screw with Cap, Ø4.5 x 55mm	10
4-068-030TI	Vertaux Basico Polyaxial Screw with Cap, Ø5.5 x 30mm	10
4-068-035TI	Vertaux Basico Polyaxial Screw with Cap, Ø5.5 x 35mm	10
4-068-040TI	Vertaux Basico Polyaxial Screw with Cap, Ø5.5 x 40mm	10
4-068-045TI	Vertaux Basico Polyaxial Screw with Cap, Ø5.5 x 45mm	10
4-068-050TI	Vertaux Basico Polyaxial Screw with Cap, Ø5.5 x 50mm	10
4-068-055TI	Vertaux Basico Polyaxial Screw with Cap, Ø5.5 x 55mm	10
4-070-030TI	Vertaux Basico Polyaxial Screw with Cap, Ø6.5 x 30mm	10
4-070-035TI	Vertaux Basico Polyaxial Screw with Cap, Ø6.5 x 35mm	10
4-070-040TI	Vertaux Basico Polyaxial Screw with Cap, Ø6.5 x 40mm	10
4-070-045TI	Vertaux Basico Polyaxial Screw with Cap, Ø6.5 x 45mm	10
4-070-050TI	Vertaux Basico Polyaxial Screw with Cap, Ø6.5 x 50mm	10
4-070-055TI	Vertaux Basico Polyaxial Screw with Cap, Ø6.5 x 55mm	10
4-072-400TI	Vertaux Basico Polyaxial Screw Rod, Length 400mm	5
SP-303	VERTAUX - Pedicle Screw Cap, M9	10
7-108-05	Implant Tray for Vertaux Basico Screw	1
7-108-02	Vertaux Basico Reducer	1
7-108-03	Vertaux Basico Screwdriver	2
7-108-04	Vertaux Basico Screw Shaft with Quick Coupling	1
7-108-06	Container with Lid for Vertaux Basico Screw	1

SP-061 Set Screw Holder - Long



SP-052 Crosslink Hook Holder



SP-053 Rod Bender



SP-063 Star Screwdriver with Quick Coupling - Long



SP-064 Pedicle Sounding Probe, Straight



SP-065 Pedicle Sounding Probe, Curved



SP-066 Breaker



SP-067 Marker(L)



SP-068 Marker(R)



SP-069 Marker Inserter



SP-070 In-Situ Bender (L)



SP-071 In-Situ Bender (R)




SP-054 Hex Rod Rotation Wrench



SP-055 Power Gripper



7-116-01 Rod Reduction Forcep - Frog Type



SP-056 Distractor for 5.5mm Pedicle Screw



SP-057 Compressor for 5.5mm Pedicle Screw



SP-073 Awl for 5.5mm Pedicle Screw



SP-074 Straight Probe**SP-075** Curved Probe**SP-077** Tap, Ø4.5mm**SP-1950** Tap, Ø5.0mm**SP-078** Tap, Ø5.5mm**SP-079** Tap, Ø6.0mm

SP-080 Tap, Ø6.5mm



SP-1951 Tap, Ø7.0mm



SP-081 Screwdriver, Hex 3.5mm, for 5.5mm Pedicle Screw



SP-082 Mono-Axial Pedicle Screwdriver



SP-083 Mono-Axial Reduction Pedicle Screwdriver



SP-084 Poly-Axial Pedicle Screwdriver



SP-085 Poly-Axial Reduction Pedicle Screwdriver



SP-086 T-Handle with Quick Coupling for 5.5mm Pedicle Screw



SP-087 Template 120mm



SP-1952 Template 450mm




SP-088 Rod Pusher



SP-089 Counter Torque



SP-090 Ratchet Handle for 5.5mm Pedicle Screw**SP-091** Pedicle Screw Measuring Device (White)**SP-092** Measuring Plate for Crosslink (Black)**SP-1953** Rod Holder, Curved**SP-093** Hex Rod Rotation Forcep**SP-094** Rod Push Forcep

SP-1493 Torque Limiter Handle, 10.0Nm



7-116-02 Ratchet Handle T-Type for 5.5mm Pedicle Screw




4-106-01 Screwdriver, Hex 3.0mm for Crosslink Connector Type-III







SP-468 Instrument Trays for Vertaux 5.5 Pedicle Screw Instrument Set



SP-059 Container with PPSU Lid for Vertaux 5.5 Pedicle Screw Instrument Set



SP-058 VERTAUX 5.5mm Pedicle Screw Instrument Set



SP-058 VERTAUX 5.5mm Pedicle Screw Instrument Set

Code	Set Consisting of	Qty.
SP-061	Set Screw Holder - Long	2
SP-052	Crosslink Hook Holder	1
SP-053	Rod Bender	1
SP-063	Star Screwdriver with Quick Coupling - Long	2
SP-064	Pedicle Sounding Probe, Straight	1
SP-065	Pedicle Sounding Probe, Curved	1
SP-066	Breaker	1
SP-067	Marker(L)	5
SP-068	Marker(R)	5
SP-069	Marker Inserter	1
SP-070	In-Situ Bender (L)	1
SP-071	In-Situ Bender (R)	1
SP-054	Hex Rod Rotation Wrench	1
SP-055	Power Gripper	1
7-116-01	Rod Reduction Forcep - Frog Type	1
SP-056	Distractor for 5.5mm Pedicle Screw	1
SP-057	Compressor for 5.5mm Pedicle Screw	1
SP-073	Awl for 5.5mm Pedicle Screw	1
SP-074	Straight Probe	1
SP-075	Curved Probe	1
SP-077	Tap, Ø4.5mm	1
SP-1950	Tap, Ø5.0mm	1
SP-078	Tap, Ø5.5mm	1
SP-079	Tap, Ø6.0mm	1
SP-080	Tap, Ø6.5mm	1
SP-1951	Tap, Ø7.0mm	1
SP-081	Screwdriver, Hex 3.5mm, for 5.5mm Pedicle Screw	1
SP-082	Mono-Axial Pedicle Screwdriver	1
SP-083	Mono-Axial Reduction Pedicle Screwdriver	1
SP-084	Poly-Axial Pedicle Screwdriver	1
SP-085	Poly-Axial Reduction Pedicle Screwdriver	1
SP-086	T-Handle with Quick Coupling for 5.5mm Pedicle Screw	1
SP-087	Template 120mm	1
SP-1952	Template 450mm	1
SP-088	Rod Pusher	1
SP-089	Counter Torque	1
SP-090	Ratchet Handle for 5.5mm Pedicle Screw	1
SP-091	Pedicle Screw Measuring Device (White)	1
SP-092	Measuring Plate for Crosslink (Black)	1
SP-1953	Rod Holder, Curved	1
SP-093	Hex Rod Rotation Forcep	1
SP-094	Rod Push Forcep	1
SP-1493	Torque Limiter Handle, 10.0Nm	1
7-116-02	Ratchet Handle T-Type for 5.5mm Pedicle Screw	1
4-106-01	Screwdriver, Hex 3.0mm for Crosslink Connector Type-III	1
SP-468	Instrument Trays for Vertaux 5.5 Pedicle Screw Instrument Set	4
SP-059	Container with PPSU Lid for Vertaux 5.5 Pedicle Screw Instrument Set	1