

# ASTERIUS


## Anterior Cervical Plate System

The Asterius(Anterior Cervical Plate) is used in anterior plating of the cervical spine(C2–T2) for the internal fixation in the treatment of instabilities associated with:


- Fractures/dislocations
- Degenerative diseases
- Tumours
- Partial or total Spondylectomy



## Asterius-Anterior Cervical Plate



Unlocked



Locked



**STERILE**

Code	Length	
SP-145	22.5mm	Level-1
SP-146	25mm	
SP-147	27.5mm	
SP-148	30mm	Level-2
SP-149	32.5mm	
SP-150	35mm	
SP-151	37.5mm	
SP-152	40mm	
SP-153	42.5mm	
SP-154	45mm	Level-3
SP-155	47.5mm	
SP-156	50mm	
SP-157	52.5mm	
SP-158	55mm	
SP-159	57.5mm	
SP-160	60mm	
SP-161	62.5mm	
SP-162	65mm	
SP-163	67.5mm	
SP-164	70mm	

Code	Length	
SP-145-S	22.5mm	Level-1
SP-146-S	25mm	
SP-147-S	27.5mm	
SP-148-S	30mm	Level-2
SP-149-S	32.5mm	
SP-150-S	35mm	
SP-151-S	37.5mm	
SP-152-S	40mm	
SP-153-S	42.5mm	
SP-154-S	45mm	Level-3
SP-155-S	47.5mm	
SP-156-S	50mm	
SP-157-S	52.5mm	
SP-158-S	55mm	
SP-159-S	57.5mm	
SP-160-S	60mm	
SP-161-S	62.5mm	
SP-162-S	65mm	
SP-163-S	67.5mm	
SP-164-S	70mm	

**Asterius-3.5mm Cancellous Screw, Self-Tapping, (Hex Head)**



Code	Length
SP-1789	12
SP-1790	14
SP-1791	16
SP-1792	18
SP-1793	20
SP-1794	22

**Asterius-4.0mm Cancellous Screw, Self-Tapping, (Hex Head)**



Code	Length
SP-1795	12
SP-1796	14
SP-1797	16
SP-1798	18
SP-1799	20
SP-1800	22

**Asterius-3.5mm Cancellous Screw, Self-Drilling, (Hex Head)**



Code	Length
SP-1801	12
SP-1802	14
SP-1803	16
SP-1804	18
SP-1805	20
SP-1806	22

**Asterius-4.0mm Cancellous Screw, Self-Drilling, (Hex Head)**



Code	Length
SP-1807	12
SP-1808	14
SP-1809	16
SP-1810	18
SP-1811	20
SP-1812	22

**Asterius-3.5mm Cancellous Screw,  
Self-Tapping, (Hex Head)**

**STERILE**

Code	Length
SP-1789-S	12
SP-1790-S	14
SP-1791-S	16
SP-1792-S	18
SP-1793-S	20
SP-1794-S	22

**Asterius-4.0mm Cancellous Screw,  
Self-Tapping, (Hex Head)**

**STERILE**

Code	Length
SP-1795-S	12
SP-1796-S	14
SP-1797-S	16
SP-1798-S	18
SP-1799-S	20
SP-1800-S	22

**Asterius-3.5mm Cancellous Screw,  
Self-Drilling, (Hex Head)**


**STERILE**

Code	Length
SP-1801-S	12
SP-1802-S	14
SP-1803-S	16
SP-1804-S	18
SP-1805-S	20
SP-1806-S	22

**Asterius-4.0mm Cancellous Screw,  
Self-Drilling, (Hex Head)**

**STERILE**

Code	Length
SP-1807-S	12
SP-1808-S	14
SP-1809-S	16
SP-1810-S	18
SP-1811-S	20
SP-1812-S	22

**SP-1609-000 Asterius- Anterior Cervical Plate Implant Set**


Code	Product Specification	Qty.
<b>SP-145</b>	Asterius -Anterior Cervical Plate, Length, 22.5mm Titanium	1
<b>SP-146</b>	Asterius -Anterior Cervical Plate, Length, 25.0mm Titanium	1
<b>SP-147</b>	Asterius -Anterior Cervical Plate, Length, 27.5mm Titanium	1
<b>SP-148</b>	Asterius -Anterior Cervical Plate, Length, 30.0mm Titanium	1
<b>SP-149</b>	Asterius -Anterior Cervical Plate, Length, 32.5mm Titanium	1
<b>SP-150</b>	Asterius -Anterior Cervical Plate, Length, 35.0mm Titanium	1
<b>SP-151</b>	Asterius -Anterior Cervical Plate, Length, 37.5mm Titanium	1
<b>SP-152</b>	Asterius -Anterior Cervical Plate, Length, 40.0mm Titanium	1
<b>SP-153</b>	Asterius -Anterior Cervical Plate, Length, 42.5mm Titanium	1
<b>SP-154</b>	Asterius -Anterior Cervical Plate, Length, 45.0mm Titanium	1
<b>SP-155</b>	Asterius -Anterior Cervical Plate, Length, 47.5mm Titanium	1
<b>SP-156</b>	Asterius -Anterior Cervical Plate, Length, 50.0mm Titanium	1
<b>SP-157</b>	Asterius -Anterior Cervical Plate, Length, 52.5mm Titanium	1
<b>SP-158</b>	Asterius -Anterior Cervical Plate, Length, 55.0mm Titanium	1
<b>SP-159</b>	Asterius -Anterior Cervical Plate, Length, 57.5mm Titanium	1
<b>SP-160</b>	Asterius -Anterior Cervical Plate, Length, 60.0mm Titanium	1
<b>SP-161</b>	Asterius -Anterior Cervical Plate, Length, 62.5mm Titanium	1
<b>SP-162</b>	Asterius -Anterior Cervical Plate, Length, 65.0mm Titanium	1
<b>SP-163</b>	Asterius -Anterior Cervical Plate, Length, 67.5mm Titanium	1
<b>SP-164</b>	Asterius -Anterior Cervical Plate, Length, 70mm Titanium	1

Code	Product Specification	Qty.
<b>SP-1790</b>	Asterius - 3.5mm Cancellous Screws, Self-Tapping, (Hex Head), Length 14mm	8
<b>SP-1791</b>	Asterius - 3.5mm Cancellous Screws, Self-Tapping, (Hex Head), Length 16mm	8
<b>SP-1792</b>	Asterius - 3.5mm Cancellous Screws, Self-Tapping, (Hex Head), Length 18mm	8
<b>SP-1796</b>	Asterius - 4.0mm Cancellous Screws, Self-Tapping, (Hex Head), Length 14mm	8
<b>SP-1797</b>	Asterius - 4.0mm Cancellous Screws, Self-Tapping, (Hex Head), Length 16mm	8
<b>SP-1798</b>	Asterius - 4.0mm Cancellous Screws, Self-Tapping, (Hex Head), Length 18mm	8
<b>OR</b>		
<b>SP-1802</b>	Asterius - 3.5mm Cancellous Screws, Self-Drilling, (Hex Head), Length 14mm	8
<b>SP-1803</b>	Asterius - 3.5mm Cancellous Screws, Self-Drilling, (Hex Head), Length 16mm	8
<b>SP-1804</b>	Asterius - 3.5mm Cancellous Screws, Self-Tapping, (Hex Head), Length 18mm	8
<b>SP-1808</b>	Asterius - 4.0mm Cancellous Screws, Self-Drilling, (Hex Head), Length 14mm	8
<b>SP-1809</b>	Asterius - 4.0mm Cancellous Screws, Self-Drilling, (Hex Head), Length 16mm	8
<b>SP-1810</b>	Asterius - 4.0mm Cancellous Screws, Self-Drilling, (Hex Head), Length 18mm	8
<b>4-028-01</b>	Tray for Asterius - Anterior Cervical Plate Implant Set	1



Hex Head

**SP-003**      Fixation Pin for Anterior Cervical Plate




**SP-004**      Plate Bender



**SP-005**      Screwdriver Sleeve



**SP-006**      Positioning Rod for Anterior Cervical Plate



**SP-007**      Distractor for Anterior Cervical Plate



**SP-008**      Awl, 2.0mm, for Anterior Cervical Plate



**SP-009** Drill Bit,  $\Phi$ 2.2mm, for Anterior Cervical Plate



**SP-010** Screwdriver, Hex 2.0mm



**SP-011** Locking Cap Pusher



**SP-012** Inserter for Positioning Rod



**SP-013** Multifunctional Holder



**SP-014** Inserter for Fixation Pin





**SP-015** Plate Holder



**SP-017** Mono Drilling Guide




**SP-016** Dual Drilling Guide



**SP-018** Tap Ø3.0mm



**SP-466** Mesh Tray for Asterius- Anterior Cervical Plate Instrument Set




**4-028-01** Implant Box for Asterius Anterior Cervical Plate System



**SP-002** Sterilization Container for Asterius- Anterior Cervical Plate Instrument Set



**SP-001 Asterius- Anterior Cervical Plate Instrument Set**


Code	Product Specification	Qty.
<b>SP-003</b>	Fixation Pin for Anterior Cervical Plate	2
<b>SP-004</b>	Plate Bender	1
<b>SP-005</b>	Screwdriver Sleeve	1
<b>SP-006</b>	Positioning Rod for Anterior Cervical Plate	2
<b>SP-007</b>	Distractor for Anterior Cervical Plate	1
<b>SP-008</b>	Awl, 2.0mm, for Anterior Cervical Plate	1
<b>SP-009</b>	Drill Bit, Ø2.2mm, for Anterior Cervical Plate	1
<b>SP-010</b>	Screwdriver, Hex 2.0mm	1
<b>SP-011</b>	Locking Cap Pusher	1
<b>SP-012</b>	Inserter for Positioning Rod	1
<b>SP-013</b>	Multifunctional Holder	1
<b>SP-014</b>	Inserter for Fixation Pin	1
<b>SP-015</b>	Plate Holder	1
<b>SP-017</b>	Mono Drilling Guide	1
<b>SP-016</b>	Dual Drilling Guide	1
<b>SP-018</b>	Tap Ø3.0mm	1
<b>SP-466</b>	Mesh Tray for Asterius- Anterior Cervical Plate Instrument Set	2
<b>4-028-01</b>	Implant Box for Asterius Anterior Cervical Plate System	1
<b>SP-002</b>	Sterilization Container for Asterius- Anterior Cervical Plate Instrument Set	1