


# AUXILOCK® AUXSUTURE ANCHOR

Shoulder Instability Repair



AUXILOCK® AuxSuture Anchor Provides numerous advantages over traditional anchors for shoulder instability repair procedures with its strong fixations, less bone removal and a push in techniques. AuxSuture anchor is available in diameter of 1.5mm and 1.8mm. It is available in single and double loaded options with various combination of BioBraid and suture tapes. The drill guide and 1.5mm drill bit are specially designed for AuxSuture Anchor while assuming an accurate hole placement for the anchor. The small 1.5mm footprint AuxSuture anchor provides placement and fixation advantages when bone real estate is limited, especially during revision surgeries. The anchor is also available with needles which are ideal for mini-open surgeries.

## Ordering Information

Code	Product Description	Suture/Tape Colour	Needle
6-034-01	AUXILOCK® 1.5mm AuxSuture Anchor With One #2 BioBraid: White/Blue		
6-034-02	AUXILOCK® 1.5mm AuxSuture Anchor With One 1.4mm Suture Tape: White/Blue		
6-034-03	AUXILOCK® 1.5mm AuxSuture Anchor With One #2 BioBraid: White/Blue, With Needles: MO-6		
6-035-01	AUXILOCK® 1.8mm AuxSuture Anchor With Two #2 BioBraid: White/Blue & White/Green		
6-035-02	AUXILOCK® 1.8mm AuxSuture Anchor With One #2 BioBraid and One 1.4mm Suture Tape: White/Black & White/Blue		
6-035-03	AUXILOCK® 1.8mm AuxSuture Anchor With Two #2 BioBraid: White/Blue & White/Green, With Needles: MO-6		
6-035-04	AUXILOCK® 1.8mm AuxSuture Anchor With One #2 BioBraid and One 1.4mm Suture Tape: White/Black & White/Blue, With Needles: MO-6		

## Recommended Instrumentation

Code	Product Description	Compatible Implants
7-013-03	V-Mouth Drill Guide	6-034-01, 6-034-02, 6-035-01, 6-035-02
7-013-04	8 Point Drill Guide	6-034-01, 6-034-02, 6-035-01, 6-035-02
7-013-05	V-Mouth Drill Guide - Curve	6-034-01, 6-034-02, 6-035-01, 6-035-02
7-013-06	Trocar Tip Obturator	6-034-01, 6-034-02, 6-035-01, 6-035-02
7-013-07	Trocar Tip Obturator - Curve	6-034-01, 6-034-02, 6-035-01, 6-035-02
7-013-20	Ø1.5mm X 300mm Drill Bit, Stainless Steel	6-034-01, 6-034-02, 6-034-03, 6-035-01, 6-035-02, 6-035-03, 6-035-04
7-013-21	Ø1.5mm X 300mm Drill Bit, Nitinol	6-034-01, 6-034-02, 6-035-01, 6-035-02
7-013-01	8 Point Slotted Drill Guide, Length 190mm	6-034-03, 6-035-03, 6-035-04
7-013-02	8 Point Drill Guide, Length 100mm	6-034-03, 6-035-03, 6-035-04

STERILE EO