

AUXILOCK® GFS BUTTON EXTENDER



GFS Button Extender is an extension device designed to provide extended surface area during cruciate graft reconstructions. It can be used in conjunction with the GFS Mini and GFS Ultimate Mini button providing an effective solution for cortical blowouts. It also helps overcome intraoperative complications seen in revision surgeries.

Ordering Information

Code	Product Description	Length	Width	Height	Slot Depth
6-032-01	AUXILOCK® GFS Button Extender, L: 18mm, W: 5.0mm	18mm	5mm	2.3mm	1.7mm

STERILE EO