

# AUXILOCK® GFS CONCAVE BUTTON



The GFS Concave Button is an attachable Button System that has revolutionized tibial fixation of ACL and PCL grafts. GFS Concave Button can be used on all graft types and attached to a variety of button configurations for fixation over sockets or full tunnels.

**The advantages of the GFS Concave Button implant include:**

- Maximum graft-to-bone contact improves incorporation and healing
- The ability to retension grafts after fixation and knee cycling
- Several different button options for sockets and full tunnels

## Ordering Information

Code	Product Description
6-050-01	AUXILOCK® GFS Concave Button 11mm with 4mm Collar
6-050-02	AUXILOCK® GFS Concave Button 14mm with 7mm Collar
6-050-03	AUXILOCK® GFS Concave Button 20mm with 9mm Collar
6-050-04	AUXILOCK® GFS Round Button without Collar, 14mm



Concave Button 20mm  
with 9mm Collar

Use with 9 mm-13 mm  
tunnels



Concave Button 14mm  
with 7mm Collar

Use with 7 mm-9 mm  
tunnels



Concave Button  
11mm with 4mm Collar

Use with 4 mm-7 mm  
tunnels



Round Button without  
Collar, 14mm

STERILE EO