


### AUXILOCK® GFS ULTIMATE MINI



GFS Ultimate Mini is an adjustable loop with an oblong shaped titanium button used for Cruciate Graft Reconstructions. GFS Ultimate provides a double locking mechanism which eliminates the need for knot tying.

**GFS Ultimate Mini provides three sutures (UHMWPE):**

- Adjustable Suture (#7, White/Blue): It allows the surgeon to maximize the amount of graft inside the femoral tunnel, thereby optimizing the healing process. It also enables calibration of the loop to its optimum size.
- Pulling Suture (#5, White): It is available to pull the graft inside the tunnel.
- Flipping Suture (#5, White/Black): It ensures the flipping of the button.

### Ordering Information

Code	Product Description	Length of Loop (mm)
6-006-01	AUXILOCK® GFS Ultimate Mini Button, Button L: 12mm, W: 3.9mm (Adjustable Loop)	Max: 70, Min: 12

### Recommended Instrumentation

Code	Product Description
7-013-57	Ø4.5mm Cannulated Headed Reamer

**STERILE EO**