


AUXILOCK® MENISCO INSIDE-OUT

Menisco Inside-out system is specifically designed for inside-out meniscal repair procedures features cannulas, needles and rasps.

Five pre-bent Single Lumen cannulas provide optimum access to all zones of the meniscus with funnel ends that make the loading of needles safe and easy.

Double-arm meniscal repair needles are thin, strong and flexible for easy passage. The needles are available with Pre-attached #2-0, #0 BioBraid Suture & 1.0mm Suture Tape. Meniscal repair needles are made of Stainless Steel and have a large eyelet at one end to permit easy threading of suture material.



Ordering Information

| Code | Product Description |
|----------|---|
| 6-047-01 | AUXILOCK® MENISCO Inside-Out Meniscal Repair Needle with #2-0 BioBraid: White/Blue, Length 36in, Needle Length 25cm |
| 6-047-02 | AUXILOCK® MENISCO Inside-Out Meniscal Repair Needle with #0 BioBraid: White/Blue, Length 36in, Needle Length 25cm |
| 6-047-03 | AUXILOCK® MENISCO Inside-Out Meniscal Repair Needle with 1.0mm Suture Tape: White/Blue, Length 36in, Needle Length 25cm |

Related Instrumentation

| Code | Product Description |
|----------|---|
| 7-156-01 | Single Lumen Cannula for Meniscal Repair, Straight |
| 7-156-02 | Single Lumen Cannula for Meniscal Repair, 10° Curve |
| 7-156-03 | Single Lumen Cannula for Meniscal Repair, 20° Curve |
| 7-156-04 | Single Lumen Cannula for Meniscal Repair, 30° Curve |
| 7-156-05 | Single Lumen Cannula for Meniscal Repair, 40° Curve |
| 7-156-06 | Bender for Cannula |
| 7-156-07 | Needle Pusher |
| 7-156-08 | Slotted Cannula |
| 7-156-09 | Knot Pusher for Meniscal Repair |
| 7-156-10 | Curved Probe for Meniscal Repair |
| 7-156-11 | Diamond Rasp, 45° for Meniscal Repair |
| 7-156-12 | Diamond Rasp, 90° for Meniscal Repair |
| 7-156-13 | Ring Curette for Meniscal Repair |
| 7-156-14 | Knot Pusher with Suture Cutter for Meniscal Repair |

STERILE EO