

AUXILOCK® NAKED BUTTON



GFS Mini Naked Button



GFS II Large Naked Button



Ultimate Mini Naked Button



Ultimate Large Naked Button

GFS Naked button is made up of titanium alloy. The system provides an excellent combination of strength and stiffness required for successful cortical fixation. The large naked buttons of GFS may eliminate the need for a stepped tunnel technique and are considered to be an excellent choice for revision surgeries.

Ordering Information

Code	Product Description	Length	Width	Height
6-007-11	AUXILOCK® GFS Mini Naked Button L: 12mm, W: 3.9mm	12mm	3.9mm	1.5mm
6-008-09	AUXILOCK® GFS II Large Naked Button L: 16.5mm, W: 4.4mm	16.5mm	4.4mm	2.8mm
6-006-03	AUXILOCK® Ultimate Mini Naked Button L: 12mm, W: 3.9mm	12mm	3.9mm	1.5mm
6-006-04	AUXILOCK® Ultimate Large Naked Button L: 16.5mm, W: 3.9mm	16.5mm	3.9mm	1.5mm

STERILE EO