


AUXILOCK® TITANIUM SCREW-IN SUTURE ANCHOR



Ordering Information

Code	Product Description	Suture/Tape Colour	Needle
6-005-16	AUXILOCK® 2.0mm Titanium Screw-In Suture Anchor With One #0 BioBraid: White/Blue		
6-005-11	AUXILOCK® 2.0mm Titanium Screw-In Suture Anchor With One #0 BioBraid: White/Blue, With Needles		
6-005-15	AUXILOCK® 2.5mm Titanium Screw-In Suture Anchor With One #1 BioBraid: White/Blue		
6-005-12	AUXILOCK® 2.5mm Titanium Screw-In Suture Anchor With One #1 BioBraid: White/Blue, With Needles: MO-6		

Recommended Instrumentation

Code	Product Description	Compatible Implants
7-139-01	Ø1.2mm X 100mm Drill Bit	6-005-11, 6-005-16
7-139-02	8 Point Slotted Drill Guide, Length 50mm	6-005-11, 6-005-12, 6-005-16, 6-005-15
7-139-66	Ø1.8mm X 100mm Drill Bit	6-005-15, 6-005-12

STERILE|EO