

HERBERT SCREW

# HERBERT SCREW

2.5mm | 3.0mm



edition  
2023





# about us

Auxein Medical is an integrated, research based, orthopaedic Implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in india and serves customers in over 75 countries worldwide.

## Our Achievements



# **Orthopaedic Implants**

## INDICATIONS


For fixation of bone fractures or for bone reconstruction.

### Examples Include:

- Fixation of Osteotomies for Hallux Valgus treatment (such as Scarf, Chevron, etc.).
- Mono or Bi-cortical Osteotomies in the Foot or Hand.
- Distal or Proximal Metatarsal or Metacarpal Osteotomies.
- Arthrodesis in Hand or Foot surgery.
- Fixation of small bone fragments, in long bones or small bones fractures.

1.0mm Guide wire

Screw	Compression
Dia 2.5mm	Bicortical compression
Dia 3.0mm	Bicortical compression



Self drill &  
Self tapping flutes



Star Head



- Totally Intra-Osseous Implants.
- Dual threads.
- Material: Titanium Alloy

## 2.5mm Cannulated Herbert Screw, Short Thread, (Star Head)



**STERILE**

Code	Length
HS25-110	10mm
HS25-112	12mm
HS25-114	14mm
HS25-116	16mm
HS25-118	18mm
HS25-120	20mm
HS25-122	22mm
HS25-124	24mm
HS25-126	26mm
HS25-128	28mm
HS25-130	30mm
HS25-132	32mm
HS25-134	34mm
HS25-136	36mm
HS25-138	38mm
HS25-140	40mm

Code	Length
HS25-110-S	10mm
HS25-112-S	12mm
HS25-114-S	14mm
HS25-116-S	16mm
HS25-118-S	18mm
HS25-120-S	20mm
HS25-122-S	22mm
HS25-124-S	24mm
HS25-126-S	26mm
HS25-128-S	28mm
HS25-130-S	30mm
HS25-132-S	32mm
HS25-134-S	34mm
HS25-136-S	36mm
HS25-138-S	38mm
HS25-140-S	40mm



## 2.5mm Cannulated Herbert Screw, Long Thread, (Star Head)



**STERILE**

Code	Length
HS25-216	16mm
HS25-218	18mm
HS25-220	20mm
HS25-222	22mm
HS25-224	24mm
HS25-226	26mm
HS25-228	28mm
HS25-230	30mm
HS25-232	32mm
HS25-234	34mm
HS25-236	36mm
HS25-238	38mm
HS25-240	40mm

Code	Length
HS25-216-S	16mm
HS25-218-S	18mm
HS25-220-S	20mm
HS25-222-S	22mm
HS25-224-S	24mm
HS25-226-S	26mm
HS25-228-S	28mm
HS25-230-S	30mm
HS25-232-S	32mm
HS25-234-S	34mm
HS25-236-S	36mm
HS25-238-S	38mm
HS25-240-S	40mm



## 3.0mm Cannulated Herbert Screw, Short Thread, (Star Head)



**STERILE**

Code	Length
HS30-110	10mm
HS30-112	12mm
HS30-114	14mm
HS30-116	16mm
HS30-118	18mm
HS30-120	20mm
HS30-122	22mm
HS30-124	24mm
HS30-126	26mm
HS30-128	28mm
HS30-130	30mm
HS30-132	32mm
HS30-134	34mm
HS30-136	36mm
HS30-138	38mm
HS30-140	40mm

Code	Length
HS30-110-S	10mm
HS30-112-S	12mm
HS30-114-S	14mm
HS30-116-S	16mm
HS30-118-S	18mm
HS30-120-S	20mm
HS30-122-S	22mm
HS30-124-S	24mm
HS30-126-S	26mm
HS30-128-S	28mm
HS30-130-S	30mm
HS30-132-S	32mm
HS30-134-S	34mm
HS30-136-S	36mm
HS30-138-S	38mm
HS30-140-S	40mm



## 3.0mm Cannulated Herbert Screw, Long Thread, (Star Head)



**STERILE**

Code	Length
HS30-216	16mm
HS30-218	18mm
HS30-220	20mm
HS30-222	22mm
HS30-224	24mm
HS30-226	26mm
HS30-228	28mm
HS30-230	30mm
HS30-232	32mm
HS30-234	34mm
HS30-236	36mm
HS30-238	38mm
HS30-240	40mm

Code	Length
HS30-216-S	16mm
HS30-218-S	18mm
HS30-220-S	20mm
HS30-222-S	22mm
HS30-224-S	24mm
HS30-226-S	26mm
HS30-228-S	28mm
HS30-230-S	30mm
HS30-232-S	32mm
HS30-234-S	34mm
HS30-236-S	36mm
HS30-238-S	38mm
HS30-240-S	40mm



**HS101-19** Screw Caddy for 2.5/3.0mm Cannulated Herbert Screw System



**2.5mm Cannulated Herbert Screw,  
Short Thread, (Star Head)**

Code	Length	Unit
HS25-110	10mm	2
HS25-112	12mm	2
HS25-114	14mm	2
HS25-116	16mm	2
HS25-118	18mm	2
HS25-120	20mm	2
HS25-122	22mm	2
HS25-124	24mm	2
HS25-126	26mm	2
HS25-128	28mm	2
HS25-130	30mm	2
HS25-132	32mm	2
HS25-134	34mm	2
HS25-136	36mm	2
HS25-138	38mm	2
HS25-140	40mm	2

**2.5mm Cannulated Herbert Screw,  
Long Thread, (Star Head)**

Code	Length	Unit
HS25-216	16mm	2
HS25-218	18mm	2
HS25-220	20mm	2
HS25-222	22mm	2
HS25-224	24mm	2
HS25-226	26mm	2
HS25-228	28mm	2
HS25-230	30mm	2
HS25-232	32mm	2
HS25-234	34mm	2
HS25-236	36mm	2
HS25-238	38mm	2
HS25-240	40mm	2

**3.0mm Cannulated Herbert Screw,  
Short Thread, (Star Head)**

Code	Length	Unit
HS30-110	10mm	2
HS30-112	12mm	2
HS30-114	14mm	2
HS30-116	16mm	2
HS30-118	18mm	2
HS30-120	20mm	2
HS30-122	22mm	2
HS30-124	24mm	2
HS30-126	26mm	2
HS30-128	28mm	2
HS30-130	30mm	2
HS30-132	32mm	2
HS30-134	34mm	2
HS30-136	36mm	2
HS30-138	38mm	2
HS30-140	40mm	2

**3.0mm Cannulated Herbert Screw,  
Long Thread, (Star Head)**

Code	Length	Unit
HS30-216	16mm	2
HS30-218	18mm	2
HS30-220	20mm	2
HS30-222	22mm	2
HS30-224	24mm	2
HS30-226	26mm	2
HS30-228	28mm	2
HS30-230	30mm	2
HS30-232	32mm	2
HS30-234	34mm	2
HS30-236	36mm	2
HS30-238	38mm	2
HS30-240	40mm	2

## Herbert Screw Instruments

**HS101-01** Screw Holding Forcep for Herbert Screw



**HS101-02** Holding Sleeve for 3.0mm Herbert Screw



**HS101-03** Holding Sleeve for 2.5mm Herbert Screw



**HS101-04** Holding Sleeve Handle for Herbert Screw



**HS101-05** Guide Sleeve, Ø4.9mm, for Herbert Screw



**HS101-06** Double Drill Guide, Ø2.0/1.1mm, for Herbert Screw



**HS101-07** Guide Wire, Ø0.8mm x Length 150mm



**HS101-08** Guide Wire, Ø1.0mm x Length 150mm



**HS101-09** Sharp Hook for Herbert Screw



**HS101-10** Depth Gauge for Herbert Screw



**HS101-11** Double Periosteal Elevator for Herbert Screw



**HS101-12** Cannulated Screwdriver Shaft, T8, for Herbert Screw



**HS101-13** Screwdriver Shaft, T8, for Herbert Screw



**HS101-14** Cannulated Drill Bit with Quick Coupling End, Ø2.0mm x Length 150mm, for Herbert Screw



**HS101-15** Cannulated Reamer Bit with Quick Coupling End, Ø4.8mm x Length 150mm, for Herbert Screw



**HS101-16** Handle With Quick Coupling for Herbert Screw




**HS101-18** Cleaning Stylet Ø1.0mm for Herbert Screw




**HS101-19** Screw Caddy for 2.5/3.0mm Cannulated Herbert Screw System




**HS101-21** Instrument Tray for Herbert Screw Instrument Set



**HS101-20** Container for Herbert Screw Instrument Set



HS101 Herbert Screw Instrument Set



**HS101 Herbert Screw Instrument Set**

Code	Set Consisting of	Units
<b>HS101-01</b>	Screw Holding Forcep for Herbert Screw	1
<b>HS101-02</b>	Holding Sleeve for 3.0mm Herbert Screw	1
<b>HS101-03</b>	Holding Sleeve for 2.5mm Herbert Screw	1
<b>HS101-04</b>	Holding Sleeve Handle for Herbert Screw	1
<b>HS101-05</b>	Guide Sleeve, Ø4.9mm, for Herbert Screw	1
<b>HS101-06</b>	Double Drill Guide, Ø2.0/1.1mm, for Herbert Screw	1
<b>HS101-07</b>	Guide Wire, Ø0.8mm x Length 150mm	4
<b>HS101-08</b>	Guide Wire, Ø1.0mm x Length 150mm	4
<b>HS101-09</b>	Sharp Hook for Herbert Screw	1
<b>HS101-10</b>	Depth Gauge for Herbert Screw	1
<b>HS101-11</b>	Double Periosteal Elevator for Herbert Screw	1
<b>HS101-12</b>	Cannulated Screwdriver Shaft, T8, for Herbert Screw	1
<b>HS101-13</b>	Screwdriver Shaft, T8, for Herbert Screw	1
<b>HS101-14</b>	Cannulated Drill Bit with Quick Coupling End, Ø2.0mm x Length 150mm, for Herbert Screw	1
<b>HS101-15</b>	Cannulated Reamer Bit with Quick Coupling End, Ø4.8mm x Length 150mm, for Herbert Screw	1
<b>HS101-16</b>	Handle With Quick Coupling for Herbert Screw	1
<b>HS101-18</b>	Cleaning Stylet Ø1.0mm for Herbert Screw	1
<b>HS101-19</b>	Screw Caddy for 2.5/3.0mm Cannulated Herbert Screw System	1
<b>HS101-21</b>	Instrument Tray for Herbert Screw Instrument Set	1
<b>HS101-20</b>	Container for Herbert Screw Instrument Set	1



**USA**

Auxein Inc.  
1500 Nw 89th Court, Suite 107-108  
Doral, Florida 33172  
Tel: +1 305 395 6062  
E Fax: +1 305 395 6262  
Email: USoffice@auxein.com

**MEXICO**

Auxein México, S.A. de C.V.  
Tepic 139 int 801, Colonia Roma Sur,  
Alcaldía Cuauhtémoc, CDMX,  
México, C.P. 06760  
Tel: +521 55 7261 0318  
Email: info@auxein.mx

**INDIA**

Auxein Medical Pvt. Ltd.  
Plot No. 168-169-170, Phase-4,  
Kundli Industrial Area,  
HSIIDC, Sector-57, Sonapat - 131028, Haryana  
Tel: +91 99106 43638 | Fax: +91 86077 70197  
Email: info@auxein.com