




AUXILOCK[®] AUXPIK SUTURE-BASED ANCHOR



AUXILOCK[®] AuxPik Suture-Based Anchor is made from UHMWPE suture anchor body with a PEEK OPTIMA eyelet tip. The AuxPik Suture-Based Anchor is designed for efficiency and ease of use. AuxPik Suture-Based Anchor is available in 1.8mm diameter with various combination of BioBraid and suture tapes. The anchor can be used for Bankart or SLAP repair surgeries. The anchor provides small foot print and subcortical fixation for anchor insertions. The drill guide and anchor driver combination are well designed to improve performance and reliability.

Ordering Information

| Code | Product Description | Suture/Tape Colour |
|----------|---|---|
| 6-001-01 | AUXILOCK [®] 1.8mm AuxPik Suture-Based Anchor With One #2 BioBraid: White/Blue |  |
| 6-001-08 | AUXILOCK [®] 1.8mm AuxPik Suture-Based Anchor With One 1.4mm Suture Tape: White/Blue |  |

Recommended Instrumentation

| Code | Product Description | Compatible Implants |
|----------|--------------------------|---------------------|
| 7-013-03 | V-Mouth Drill Guide | 6-001-01, 6-001-08 |
| 7-013-04 | 8 Point Drill Guide | 6-001-01, 6-001-08 |
| 7-013-06 | Trocar Tip Obturator | 6-001-01, 6-001-08 |
| 7-013-53 | Ø1.8mm X 325mm Drill Bit | 6-001-01, 6-001-08 |

STERILE EO