




## Surgical Technique

Naked Button for ACL Reconstruction




## AUXILOCK<sup>®</sup> NAKED BUTTON

GFS Naked button is made up of titanium alloy. The system provides an excellent combination of strength and stiffness required for successful cortical fixation. GFS Naked button range includes GFS Mini naked button, GFS II Large Naked button, GFS Ultimate Mini Naked button, GFS Ultimate Large naked button. The large naked buttons of GFS may eliminate the need for a stepped tunnel technique and are considered to be an excellent choice for revision surgeries.




Ideal for use over tibial sockets created with a Flykuter reamer, the GFS loops pass easily through small-diameter tunnels and allow attachment of GFS Naked buttons against the tibial cortex. GFS Naked buttons are available in multiple sizes and shapes and provide strong, reliable cortical fixation.


Ideal for full tunnels, the centering feature of these buttons maintains position over the tunnel and provides a better seal at the cortex than standard flat buttons. The concave surface countersinks sutures and knots. The 14 mm and 20mm buttons have slots for the GFS Ultimate Mini loop along with 2 holes for additional sutures.




**GFS Mini Naked Button**  
L: 12mm X 3.9mm



**GFS II Large Naked Button**  
L: 16.5mm X 4.4mm



**Ultimate Mini Naked Button**  
L: 12mm X 3.9mm



**Ultimate II Large Naked Button**  
L: 16.5mm X 3.9mm

## **USA**

### **Auxein Inc.**

1500 Nw 89th Court, Suite 107-108  
Doral, Florida 33172  
Tel: +1 305 395 6062  
E Fax: +1 305 395 6262  
Email: USoffice@auxein.com

## **MEXICO**

### **Auxein México, S.A. de C.V.**

Tepic 139 int 801, Colonia Roma Sur,  
Alcaldía Cuauhtémoc, CDMX,  
México, C.P. 06760  
Tel: +521 55 7261 0318  
Email: info@auxein.mx

## **INDIA**

### **Auxein Medical Pvt. Ltd.**

Plot No. 168-169-170, Phase-4,  
Kundli Industrial Area,  
HSIIDC, Sector-57, Sonapat - 131028, Haryana  
Tel: +91 99106 43638 | Fax: +91 86077 70197  
Email: info@auxein.com