








# AUXILOCK® BIOBRAID INFINITY LOOP NEEDLE



AUXILOCK® BioBraid Infinity Loop is a continuous loop of BioBraid with a straight needle. BioBraid Infinity Loop utilises the whipstitch technique for graft preparation and reduces the time spent. It also compresses the graft uniformly and improves strength.

## Ordering Information

Code	Product Description	Suture Colour	Needle
6-012-01	AUXILOCK® #2 BioBraid Infinity Loop: White, With Needle: Straight, 20in Loop Length, 24in Total Length		
6-012-02	AUXILOCK® #2 BioBraid Infinity Loop: White/Blue, With Needle: Straight, 20in Loop Length, 24in Total Length		
6-012-03	AUXILOCK® #5 BioBraid Infinity Loop: White/Blue, With Needle: Straight, 20in Loop Length, 24in Total Length		

**STERILE EO**