


AUXILOCK® BIOBRAID SUTURE TAPE



AUXILOCK® BioBraid Suture Tape made from UHMWPE. The suture tape design accommodates the tape in the centre with a suture on the either sides of the tape. It can be used in Knee, Shoulder, and Foot & Ankle surgeries including AC joint reconstruction, knotless double-row rotator cuff repair, patella fracture repair, MPFL Reconstruction etc. The broad footprint of the BioBraid Suture Tape is appropriate for repairs in degenerative cuff tissue where tissue pull-through may be a concern.

Ordering Information

Code	Product Description	Needle
6-013-01	AUXILOCK® 1.8mm BioBraid Suture Tape: White/Blue, 39in Total length	
6-013-02	AUXILOCK® 1.4mm BioBraid Suture Tape: White/Blue, 39in Total length	
6-013-03	AUXILOCK® 1.4mm BioBraid Suture Tape: White/Black, 39in Total length	
6-013-04	AUXILOCK® 1.6mm BioBraid Suture Tape: White/Blue, 39in Total length	
6-013-05	AUXILOCK® 2.0mm BioBraid Suture Tape: White/Blue, 39in Total length	
6-013-06	AUXILOCK® 1.0mm BioBraid Suture Tape: White/Blue, 39in Total length	