


AUXILOCK[®] DUWLOC KNOTLESS PEEK OPTIMA PUSH-IN ANCHOR



PEEK OPTIMA
Eyelet

PEEK OPTIMA
Anchor Body

DUWLOC delivers secure, knotless labral repair with precise tension control—minimizing surgical steps and maximizing glenoid preservation.

AUXILOCK[®] DUWLOC PEEK OPTIMA Knotless Push-in Anchor is engineered for efficient, secure labral repair without the complications of arthroscopic knots. It is available in the diameter of 2.8mm. Its innovative “suture-first” technique allows surgeons to pass sutures through the capsule or labrum and fine-tune tissue tension before anchor insertion. It comes with a PEEK OPTIMA Eyelet. The DUWLOC anchor reduces the risk of knot-related impingement, offering multiple stitch configurations and precise tension control for optimal outcomes.

Ordering Information

Code	Product Description
6-064-01	AUXILOCK [®] Duwloc 2.8mm x 15.5mm Knotless PEEK OPTIMA Push-In Anchor

Related Instrumentation

Code	Product Description
7-013-83	Ø2.5mm X 200mm Drill Bit

STERILE EO