


AUXILOCK[®] DUWLOC KNOTLESS PEEK OPTIMA PUSH-IN ANCHOR



DUWLOC delivers secure, knotless rotator cuff repair with precise tension control—minimizing surgical steps and maximizing glenoid preservation.

AUXILOCK[®] **DUWLOC** PEEK OPTIMA Knotless Push-in Anchor is engineered for efficient, secure rotator cuff repair without the complications of arthroscopic knots. It is available in the diameter of 3.5mm. Its innovative “suture-first” technique allows surgeons to pass sutures through the rotator cuff repair and fine-tune tissue tension before anchor insertion. It comes with a **PEEK OPTIMA Eyelet and Self-punching metal eyelet**. **DUWLOC** Self Punching anchor is designed to eliminate the need for pre-drilling, simplifying the procedure and saving operative time by eliminating knots, the **DUWLOC** anchor reduces the risk of knot-related impingement, offering multiple stitch configurations and precise tension control for optimal outcomes.

Ordering Information

Code	Product Description
6-064-03	AUXILOCK [®] Duwloc 3.5mm x 15.5mm Knotless PEEK OPTIMA Push-In Anchor
6-064-04	AUXILOCK [®] Duwloc SP 3.5mm x 18.5mm Knotless PEEK OPTIMA Push-In Anchor

Related Instrumentation

Code	Product Description
7-013-84	Ø3.2mm X 200mm Drill Bit

This is not for SP Anchor

STERILE EO