


## AUXILOCK<sup>®</sup> GFS ULTIMATE MINI BUTTON WITH BRACE



GFS ULTIMATE MINI BUTTON WITH BRACE implant delivers enhanced handling characteristics and superior resistance to graft abrasion and tissue pull-through compared to traditional round sutures. Designed with precision in mind, its adjustable-loop mechanism facilitates incremental retensioning of the graft construct even after cortical fixation, ensuring optimal tensioning and stability.

### GFS Ultimate Mini provides three sutures (UHMWPE):

- Adjustable Suture (#7, White/Blue): It allows the surgeon to maximize the amount of graft inside the femoral tunnel, thereby optimizing the healing process. It also enables the calibration of the loop to its optimum size.
- Pulling Suture (#5, White): It is available to pull the graft inside the tunnel.
- Flipping Suture (#5, White/Black): It ensures the flipping of the button.
- 1.4mm BioBraid Suture Tape, White/black): Offers better handling.

### Ordering Information

Code	Product Description	Length of Loop (mm)
6-006-05	AUXILOCK <sup>®</sup> GFS Ultimate Mini Button With Brace, Button L: 12mm, W: 3.9mm (Adjustable Loop)	Max: 100, Min: 12

STERILE EO