


LAPIDUS

Intramedullary Nailing System



about us

Auxein Medical is an integrated, research based, orthopaedic Implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in india and serves customers in over 75 countries worldwide.

Our Achievements





LAPIDUS

Intramedullary Nailing System

INTRODUCTION

The Lapidus Intramedullary Nail System is designed to improve on existing technology for the Lapidus Arthrodesis procedure by providing a structurally sound implant that minimizes hardware prominence, improves compression capability and helps to preserve the periosteum. This Nail can be used in primary arthrodesis or for revision Lapidus procedures. An extensive offering of sizes of Lapidus Nails are available to fit variations in patient anatomy and allow for use of a bone graft to restore the length of the first ray.





Lapidus Intramedullary Nail, (3 Holes)


Ø Dia	Length (mm)	Titanium (Left)	Titanium (Right)
5.5mm	38mm	1-065-38LTI	1-065-38RTI
5.5mm	40mm	1-065-40LTI	1-065-40RTI
5.5mm	42mm	1-065-42LTI	1-065-42RTI
5.5mm	44mm	1-065-44LTI	1-065-44RTI
5.5mm	46mm	1-065-46LTI	1-065-46RTI
5.5mm	48mm	1-065-48LTI	1-065-48RTI
5.5mm	50mm	1-065-50LTI	1-065-50RTI
5.5mm	52mm	1-065-52LTI	1-065-52RTI
5.5mm	54mm	1-065-54LTI	1-065-54RTI
5.5mm	56mm	1-065-56LTI	1-065-56RTI
5.5mm	58mm	1-065-58LTI	1-065-58RTI
5.5mm	60mm	1-065-60LTI	1-065-60RTI

STERILE

Ø Dia	Length (mm)	Titanium (Left)	Titanium (Right)
5.5mm	38mm	1-065-38LTI-S	1-065-38RTI-S
5.5mm	40mm	1-065-40LTI-S	1-065-40RTI-S
5.5mm	42mm	1-065-42LTI-S	1-065-42RTI-S
5.5mm	44mm	1-065-44LTI-S	1-065-44RTI-S
5.5mm	46mm	1-065-46LTI-S	1-065-46RTI-S
5.5mm	48mm	1-065-48LTI-S	1-065-48RTI-S
5.5mm	50mm	1-065-50LTI-S	1-065-50RTI-S
5.5mm	52mm	1-065-52LTI-S	1-065-52RTI-S
5.5mm	54mm	1-065-54LTI-S	1-065-54RTI-S
5.5mm	56mm	1-065-56LTI-S	1-065-56RTI-S
5.5mm	58mm	1-065-58LTI-S	1-065-58RTI-S
5.5mm	60mm	1-065-60LTI-S	1-065-60RTI-S

Lapidus Intramedullary Nail, (4 Holes)


Ø Dia	Length (mm)	Titanium (Left)	Titanium (Right)
5.5mm	38mm	1-066-38LTI	1-066-38RTI
5.5mm	40mm	1-066-40LTI	1-066-40RTI
5.5mm	42mm	1-066-42LTI	1-066-42RTI
5.5mm	44mm	1-066-44LTI	1-066-44RTI
5.5mm	46mm	1-066-46LTI	1-066-46RTI
5.5mm	48mm	1-066-48LTI	1-066-48RTI
5.5mm	50mm	1-066-50LTI	1-066-50RTI
5.5mm	52mm	1-066-52LTI	1-066-52RTI
5.5mm	54mm	1-066-54LTI	1-066-54RTI
5.5mm	56mm	1-066-56LTI	1-066-56RTI
5.5mm	58mm	1-066-58LTI	1-066-58RTI
5.5mm	60mm	1-066-60LTI	1-066-60RTI

STERILE

Ø Dia	Length (mm)	Titanium (Left)	Titanium (Right)
5.5mm	38mm	1-066-38LTI-S	1-066-38RTI-S
5.5mm	40mm	1-066-40LTI-S	1-066-40RTI-S
5.5mm	42mm	1-066-42LTI-S	1-066-42RTI-S
5.5mm	44mm	1-066-44LTI-S	1-066-44RTI-S
5.5mm	46mm	1-066-46LTI-S	1-066-46RTI-S
5.5mm	48mm	1-066-48LTI-S	1-066-48RTI-S
5.5mm	50mm	1-066-50LTI-S	1-066-50RTI-S
5.5mm	52mm	1-066-52LTI-S	1-066-52RTI-S
5.5mm	54mm	1-066-54LTI-S	1-066-54RTI-S
5.5mm	56mm	1-066-56LTI-S	1-066-56RTI-S
5.5mm	58mm	1-066-58LTI-S	1-066-58RTI-S
5.5mm	60mm	1-066-60LTI-S	1-066-60RTI-S

3.5mm Peg Screw for Lapidus Intramedullary Nail

Length	Titanium
10mm	1-070-010TI
12mm	1-070-012TI
14mm	1-070-014TI
16mm	1-070-016TI
18mm	1-070-018TI
20mm	1-070-020TI
22mm	1-070-022TI
24mm	1-070-024TI
26mm	1-070-026TI
28mm	1-070-028TI
30mm	1-070-030TI
32mm	1-070-032TI
34mm	1-070-034TI
36mm	1-070-036TI
38mm	1-070-038TI
40mm	1-070-040TI
42mm	1-070-042TI
44mm	1-070-044TI
46mm	1-070-046TI



Length	Titanium
10mm	1-070-010TI-S
12mm	1-070-012TI-S
14mm	1-070-014TI-S
16mm	1-070-016TI-S
18mm	1-070-018TI-S
20mm	1-070-020TI-S
22mm	1-070-022TI-S
24mm	1-070-024TI-S
26mm	1-070-026TI-S
28mm	1-070-028TI-S
30mm	1-070-030TI-S
32mm	1-070-032TI-S
34mm	1-070-034TI-S
36mm	1-070-036TI-S
38mm	1-070-038TI-S
40mm	1-070-040TI-S
42mm	1-070-042TI-S
44mm	1-070-044TI-S
46mm	1-070-046TI-S

STERILE


3.5mm Wise-Lock Screw, Self Tapping, (Star Head)

Dia	Length	Titanium
Ø3.5	10mm	TI-1473-010
Ø3.5	12mm	TI-1473-012
Ø3.5	14mm	TI-1473-014
Ø3.5	16mm	TI-1473-016
Ø3.5	18mm	TI-1473-018
Ø3.5	20mm	TI-1473-020
Ø3.5	22mm	TI-1473-022
Ø3.5	24mm	TI-1473-024
Ø3.5	26mm	TI-1473-026
Ø3.5	28mm	TI-1473-028
Ø3.5	30mm	TI-1473-030
Ø3.5	32mm	TI-1473-032
Ø3.5	34mm	TI-1473-034
Ø3.5	36mm	TI-1473-036
Ø3.5	38mm	TI-1473-038
Ø3.5	40mm	TI-1473-040
Ø3.5	42mm	TI-1473-042
Ø3.5	44mm	TI-1473-044
Ø3.5	46mm	TI-1473-046




Dia	Length	Titanium
Ø3.5	10mm	TI-1473-010-S
Ø3.5	12mm	TI-1473-012-S
Ø3.5	14mm	TI-1473-014-S
Ø3.5	16mm	TI-1473-016-S
Ø3.5	18mm	TI-1473-018-S
Ø3.5	20mm	TI-1473-020-S
Ø3.5	22mm	TI-1473-022-S
Ø3.5	24mm	TI-1473-024-S
Ø3.5	26mm	TI-1473-026-S
Ø3.5	28mm	TI-1473-028-S
Ø3.5	30mm	TI-1473-030-S
Ø3.5	32mm	TI-1473-032-S
Ø3.5	34mm	TI-1473-034-S
Ø3.5	36mm	TI-1473-036-S
Ø3.5	38mm	TI-1473-038-S
Ø3.5	40mm	TI-1473-040-S
Ø3.5	42mm	TI-1473-042-S
Ø3.5	44mm	TI-1473-044-S
Ø3.5	46mm	TI-1473-046-S


STERILE

7-164-01 Outrigger, Right for Lapidus Nail Instrument





7-164-22 Outrigger, Left for Lapidus Nail Instrument



7-164-02 Outrigger Slider, Right for Lapidus Nail Instrument



7-164-23 Outrigger Slider, Left for Lapidus Nail Instrument



7-164-03 Thumb Screw for Lapidus Nail Instrument



7-164-04 Sphere Wire for Lapidus Nail Instrument



7-164-05

Nail Positioning Guide for Lapidus Nail Instrument


7-164-06

Polyaxial Targeting Guide - Sleeve Insert for Lapidus Nail Instrument


7-164-07

Polyaxial Targeting Guide for Lapidus Nail Instrument


7-164-08

Ø2.75mm Drill Depth Gauge for Lapidus Nail Instrument


7-164-09

Ø2.75mm Drill-Pin for Lapidus Nail Instrument


7-164-10

Ø2.3mm K-wire for Lapidus Nail Instrument



7-164-11 Ø5.5mm Cannulated Nail Drill for Lapidus Nail Instrument



7-164-12 Screw Guide for Lapidus Nail Instrument



7-164-13 Drill-Pin Guide for Lapidus Nail Instrument



7-164-14 Obturator for Lapidus Nail Instrument



7-164-15 Depth Probe, 60.0mm for Lapidus Nail Instrument



7-164-16 Tissue Protector for Lapidus Nail Instrument



7-164-17 Solid Driver for Lapidus Nail Instrument



7-164-18 Locking Guide for Lapidus Nail Instrument



7-164-19 Slaphammer for Lapidus Nail Instrument



7-164-20 Subchondral Perforating Drill for Lapidus Nail Instrument



7-164-21 Joint Preparation Chisel for Lapidus Nail Instrument




7-164-24 Lapidus Nipper



7-164-25

Pin Distractor



7-164-26

Torque Indicating Driver



7-164-27

Lid for Intramedullary Lapidus Nail Instrument Set

7-164-28

Container for Intramedullary Lapidus Nail Instrument Set

9/386"41pvtco gf wnet{ "Ncr kf wu" Pckllpvtwo gpv/Ugv

Code	Description	Qty.
7-164-01	Outrigger, Right for Lapidus Nail Instrument	1
7-164-22	Outrigger, Left for Lapidus Nail Instrument	1
7-164-02	Outrigger Slider, Right for Lapidus Nail Instrument	6
7-164-23	Outrigger Slider, Left for Lapidus Nail Instrument	6
7-164-03	Thumb Screw for Lapidus Nail Instrument	1
7-164-04	Sphere Wire for Lapidus Nail Instrument	1
7-164-05	Nail Positioning Guide for Lapidus Nail Instrument	1
7-164-06	Polyaxial Targeting Guide - Sleeve Insert for Lapidus Nail Instrument	1
7-164-07	Polyaxial Targeting Guide for Lapidus Nail Instrument	1
7-164-08	Ø2.75mm Drill Depth Gauge for Lapidus Nail Instrument	1
7-164-09	Ø2.75mm Drill-Pin for Lapidus Nail Instrument	1
7-164-10	Ø2.3mm K-wire for Lapidus Nail Instrument	2
7-164-11	Ø5.5mm Cannulated Nail Drill for Lapidus Nail Instrument	1
7-164-12	Screw Guide for Lapidus Nail Instrument	2
7-164-13	Drill-Pin Guide for Lapidus Nail Instrument	2
7-164-14	Obturator for Lapidus Nail Instrument	2
7-164-15	Depth Probe, 60.0mm for Lapidus Nail Instrument	1
7-164-16	Tissue Protector for Lapidus Nail Instrument	1
7-164-17	Solid Driver for Lapidus Nail Instrument	1
7-164-18	Locking Guide for Lapidus Nail Instrument	1
7-164-19	Slaphammer for Lapidus Nail Instrument	1
7-164-20	Subchondral Perforating Drill for Lapidus Nail Instrument	1
7-164-21	Joint Preparation Chisel for Lapidus Nail Instrument	1
7-164-24	Lapidus Nipper	1
7-164-25	Pin Distractor	1
7-164-26	Torque Indicating Driver	1
7-164-27	Lid for Intramedullary Lapidus Nail Instrument Set	1
7-164-28	Container for Intramedullary Lapidus Nail Instrument Set	1



USA

Auxein Inc.
1500 Nw 89th Court, Suite 107-108
Doral, Florida 33172
Tel: +1 305 395 6062
E Fax: +1 305 395 6262
Email: USoffice@auxein.com

MEXICO

Auxein México, S.A. de C.V.
Tepic 139 int 801, Colonia Roma Sur,
Alcaldía Cuauhtémoc, CDMX,
México, C.P. 06760
Tel: +521 55 7261 0318
Email: info@auxein.mx

INDIA

Auxein Medical Pvt. Ltd.
Plot No. 168-169-170, Phase-4,
Kundli Industrial Area,
HSIIDC, Sector-57, Sonapat - 131028, Haryana
Tel: +91 99106 43638 | Fax: +91 86077 70197
Email: info@auxein.com