



The background is a solid light gray. It is decorated with several sets of thin, white, wavy lines that resemble topographical map contours. These lines are located in the top-left, top-right, and bottom-left corners, creating a sense of depth and movement.

Auxein

5.0mm WISE-LOCK Large Fragment System

NEXT-LEVEL CARE







Auxein®



WISE-LOCK SYSTEM

The content of this catalogue is intended for healthcare professionals only. All product information, illustrations, descriptions, and specifications provided here are for general informational purposes and are subject to change without notice. Auxein Medical Pvt. Ltd. reserves the right to improve or discontinue products at any time without prior notification.

This document does not constitute medical advice. Surgeons and medical practitioners must exercise their own professional judgment when using any medical devices or instruments, including those provided by Auxein.

INTENDED USE

The implants and instruments in this catalogue are intended for use by qualified orthopedic surgeons trained in trauma surgery and fracture-fixation procedures. Improper selection, misuse, or lack of experience may result in complications.

REGULATORY COMPLIANCE

US FDA 510(k)

CE Mark (Class IIb)

MDSAP (Medical Device Single Audit Program)

ISO 13485:2016 Quality Management System

CDSCO – Indian Medical Device Rules, 2017

EU MDR- European Union Medical Device Regulation, 2017

INTELLECTUAL PROPERTY

All trademarks, trade names, brand names, product names, and logos mentioned in this catalogue are the property of Auxein Medical Pvt. Ltd. Unauthorized use or reproduction is strictly prohibited.

Auxein



ABOUT US

Auxein Medical is an integrated, research based, orthopaedic implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavour to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in india and serves customers in over 75 countries worldwide.



OUR ACHIEVEMENTS



OVER THE YEARS



**5.0MM WISE-LOCK
LARGE FRAGMENT
FEMUR**




Indications

- Distal shaft fractures
- Supracondylar fractures
- Intra-articular fractures
- Periprosthetic fractures



4.5/5.0mm Wise-Lock Distal Femur Plate



Holes	Direction	TITANIUM	Stainless Steel
5	Left	TI-748.105L	748.105L
5	Right	TI-748.105R	748.105R
7	Left	TI-748.107L	748.107L
7	Right	TI-748.109R	748.107R
9	Left	TI-748.109L	748.109L
9	Right	TI-748.109R	748.109R
11	Left	TI-748.111L	748.111L
11	Right	TI-748.111R	748.111R
13	Left	TI-748.113L	748.113L
13	Right	TI-748.113R	748.113R
15	Left	TI-748.115L	748.115L
15	Right	TI-748.115R	748.115R

STERILE

Holes	Direction	TITANIUM	Stainless Steel
5	Left	TI-748.105L-S	748.105L-S
5	Right	TI-748.105R-S	748.105R-S
7	Left	TI-748.107L-S	748.107L-S
7	Right	TI-748.107R-S	748.107R-S
9	Left	TI-748.109L-S	748.109L-S
9	Right	TI-748.109R-S	748.109R-S
11	Left	TI-748.111L-S	748.111L-S
11	Right	TI-748.111R-S	748.111R-S
13	Left	TI-748.113L-S	748.113L-S
13	Right	TI-748.113R-S	748.113R-S
15	Left	TI-748.115L-S	748.115L-S
15	Right	TI-748.115R-S	748.115R-S



Indications

- Buttreassing of multifragmentary distal femur fractures
- Supracondylar fractures
- Intra-articular and extra-articular condylar fractures
- Malunions and nonunions of the distal femur
- Periprosthetic fractures
- Osteopenic bone



4.5/5.0mm Wise-Lock Condylar Femur Plate



Holes	Direction	TITANIUM	Stainless Steel
7	Left	TI-731.207L	731.207L
7	Right	TI-731.207R	731.207R
9	Left	TI-731.209L	731.209L
9	Right	TI-731.209R	731.209R
11	Left	TI-731.211L	731.211L
11	Right	TI-731.211R	731.211R
12	Left	TI-731.212L	731.212L
12	Right	TI-731.212R	731.212R
13	Left	TI-731.213L	731.213L
13	Right	TI-731.213R	731.213R
15	Left	TI-731.215L	731.215L
15	Right	TI-731.215R	731.215R

STERILE

Holes	Direction	TITANIUM	Stainless Steel
7	Left	TI-731.207L-S	731.207L-S
7	Right	TI-731.207R-S	731.207R-S
9	Left	TI-731.209L-S	731.209L-S
9	Right	TI-731.209R-S	731.209R-S
11	Left	TI-731.211L-S	731.211L-S
11	Right	TI-731.211R-S	731.211R-S
12	Left	TI-731.212L-S	731.212L-S
12	Right	TI-731.212R-S	731.212R-S
13	Left	TI-731.213L-S	731.213L-S
13	Right	TI-731.213R-S	731.213R-S
15	Left	TI-731.215L-S	731.215L-S
15	Right	TI-731.215R-S	731.215R-S



Indications

Wise-Lock Proximal Femur Plate 4.5/5.0 is intended for fractures of the femur including:

- Fractures of the trochanteric region, trochanteric simple, cervicotrochanteric, trochanterodiaphyseal, multifragmentary pertrochanteric, intertrochanteric, trochanteric reversed or transverse or with additional fracture of the medial cortex
- Fractures of the proximal end of the femur combined with ipsilateral shaft fractures
- Metastatic fracture of the proximal femur
- Osteotomies of the proximal femur
- Also for use in fixation of osteopenic



4.5/5.0mm Wise-Lock Proximal Femur Plate

Left

Right




Holes	Direction	TITANIUM	Stainless Steel
2	Left	TI-798.002L	798.002L
2	Right	TI-798.002R	798.002R
3	Left	TI-798.003L	798.003L
3	Right	TI-798.003R	798.003R
4	Left	TI-798.004L	798.004L
4	Right	TI-798.004R	798.004R
6	Left	TI-798.006L	798.006L
6	Right	TI-798.006R	798.006R
8	Left	TI-798.008L	798.008L
8	Right	TI-798.008R	798.008R
10	Left	TI-798.010L	798.010L
10	Right	TI-798.010R	798.010R
12	Left	TI-798.012L	798.012L
12	Right	TI-798.012R	798.012R
14	Left	TI-798.014L	798.014L
14	Right	TI-798.014R	798.014L
16	Left	TI-798.016L	798.016L
16	Right	TI-798.016R	798.016R

STERILE

Holes	Direction	TITANIUM	Stainless Steel
2	Left	TI-798.002L-S	798.002L-S
2	Right	TI-798.002R-S	798.002R-S
3	Left	TI-798.003L-S	798.003L-S
3	Right	TI-798.003R-S	798.003R-S
4	Left	TI-798.004L-S	798.004L-S
4	Right	TI-798.004R-S	798.004R-S
6	Left	TI-798.006L-S	798.006L-S
6	Right	TI-798.006R-S	798.006R-S
8	Left	TI-798.008L-S	798.008L-S
8	Right	TI-798.008R-S	798.008R-S
10	Left	TI-798.010L-S	798.010L-S
10	Right	TI-798.010R-S	798.010R-S
12	Left	TI-798.012L-S	798.012L-S
12	Right	TI-798.012R-S	798.012R-S
14	Left	TI-798.014L-S	798.014L-S
14	Right	TI-798.014R-S	798.014L-S
16	Left	TI-798.016L-S	798.016L-S
16	Right	TI-798.016R-S	798.016R-S



4.5/5.0mm Wise-Lock Proximal Femur Lateral Plate (Type-II)



Holes	Direction	TITANIUM	Stainless Steel
4	Left	TI-1128.04L	1128.04L
4	Right	TI-1128.04R	1128.04R
5	Left	TI-1128.05L	1128.05L
5	Right	TI-1128.05R	1128.05R
6	Left	TI-1128.06L	1128.06L
6	Right	TI-1128.06R	1128.06R
7	Left	TI-1128.07L	1128.07L
7	Right	TI-1128.07R	1128.07R
8	Left	TI-1128.08L	1128.08L
8	Right	TI-1128.08R	1128.08R
9	Left	TI-1128.09L	1128.09L
9	Right	TI-1128.09R	1128.09R
10	Left	TI-1128.10L	1128.10L
10	Right	TI-1128.10R	1128.10R
12	Left	TI-1128.12L	1128.12L
12	Right	TI-1128.12R	1128.12R
14	Left	TI-1128.14L	1128.14L
14	Right	TI-1128.14R	1128.14R
16	Left	TI-1128.16L	1128.16L
16	Right	TI-1128.16R	1128.16R

STERILE

Holes	Direction	TITANIUM	Stainless Steel
4	Left	TI-1128.04L-S	1128.04L-S
4	Right	TI-1128.04R-S	1128.04R-S
5	Left	TI-1128.05L-S	1128.05L-S
5	Right	TI-1128.05R-S	1128.05R-S
6	Left	TI-1128.06L-S	1128.06L-S
6	Right	TI-1128.06R-S	1128.06R-S
7	Left	TI-1128.07L-S	1128.07L-S
7	Right	TI-1128.07R-S	1128.07R-S
8	Left	TI-1128.08L-S	1128.08L-S
8	Right	TI-1128.08R-S	1128.08R-S
9	Left	TI-1128.09L-S	1128.09L-S
9	Right	TI-1128.09R-S	1128.09R-S
10	Left	TI-1128.10L-S	1128.10L-S
10	Right	TI-1128.10R-S	1128.10R-S
12	Left	TI-1128.12L-S	1128.12L-S
12	Right	TI-1128.12R-S	1128.12R-S
14	Left	TI-1128.14L-S	1128.14L-S
14	Right	TI-1128.14R-S	1128.14R-S
16	Left	TI-1128.16L-S	1128.16L-S
16	Right	TI-1128.16R-S	1128.16R-S




Indications

Wise-Lock Proximal Femur Plate 4.5/5.0 is intended for fractures of the femur including:


- Fractures of the trochanteric region, trochanteric simple, cervicotrochanteric, trochanterodiaphyseal, multifragmentary pertrochanteric, intertrochanteric, introchanteric reversed or transverse or with additional fracture of the medial cortex
- Fractures of the proximal end of the femur combined with ipsilateral shaft fractures
- Metastatic fracture of the proximal femur
- Osteotomies of the proximal femur
- Also for use in fixation of osteopenic bone and fixation of nonunions or malunions



4.5/5.0mm Wise-Lock Proximal Femur Hook Plate



Compression and Locking Hole



Holes	TITANIUM	Stainless Steel
2	TI-579.002	579.002
4	TI-579.004	579.004
6	TI-579.006	579.006
8	TI-579.008	579.008
10	TI-579.010	579.010
12	TI-579.012	579.012
14	TI-579.014	579.014
16	TI-579.016	579.016
18	TI-579.018	579.018



Holes	TITANIUM	Stainless Steel
2	TI-579.002-S	579.002-S
4	TI-579.004-S	579.004-S
6	TI-579.006-S	579.006-S
8	TI-579.008-S	579.008-S
10	TI-579.010-S	579.010-S
12	TI-579.012-S	579.012-S
14	TI-579.014-S	579.014-S
16	TI-579.016-S	579.016-S
18	TI-579.018-S	579.018-S


The background of the entire page is a solid light gray. Overlaid on this are numerous thin, white, wavy lines that flow from the top left towards the bottom right, creating a sense of movement and depth. These lines are more densely packed in some areas and more spread out in others, resembling a stylized topographical map or a fluid motion graphic.

INSTRUMENTS



7-041-01

DIRECT MEASURING DEVICE FOR 7.3MM SCREWS - PROXIMAL FEMUR PLATE



7-041-02

3.5MM CANNULATED HEXAGONAL SCREWDRIVER FOR 7.3MM CANNULATED SCREWS - PROXIMAL FEMUR PLATE



7-041-03

WIRE GUIDE 5.0MM, FOR GUIDE WIRE Ø2.0MM - PROXIMAL FEMUR PLATE




7-041-04

WIRE GUIDE 7.3MM, FOR GUIDE WIRE Ø2.0MM - PROXIMAL FEMUR PLATE



7-041-05

THREADED DRILL GUIDE 7.3MM FOR Ø5.0MM DRILL BIT - PROXIMAL FEMUR PLATE



7-041-06

THREADED DRILL GUIDE 5.0MM FOR Ø4.3MM DRILL BIT - PROXIMAL FEMUR PLATE





PGW-2.0

GUIDE WIRE Ø2.0MM X LENGTH 250MM




PDB-5.0

CANNULATED DRILL BIT Ø5.0MM X LENGTH 230MM




7-041-07

SCREW CADDY FOR PROXIMAL FEMUR PLATE SYSTEM




7-041-10

IMPLANT TRAY FOR PROXIMAL FEMUR PLATE SYSTEM



7-041-09

INSTRUMENT TRAY FOR PROXIMAL FEMUR PLATE INSTRUMENT SET



7-041-08

STERILIZATION CONTAINER FOR PROXIMAL FEMUR PLATE INSTRUMENT SET





7-041	Proximal Femur Plate Instrument Set
-------	-------------------------------------

CODE	SET CONSISTING OF	UNIT
7-041-01	Direct Measuring Device for 7.3mm Screws - Proximal Femur Plate	1
7-041-02	3.5mm Cannulated hexagonal Screwdriver for 7.3mm Cannulated Screws - Proximal Femur Plate	1
7-041-03	Wire Guide 5.0mm, for Guide Wire Ø2.0mm - Proximal Femur Plate	1
7-041-04	Wire Guide 7.3mm, for Guide Wire Ø2.0mm - Proximal Femur Plate	1
7-041-05	Threaded Drill Guide 7.3mm for Ø5.0mm Drill Bit - Proximal Femur Plate	1
7-041-06	Threaded Drill Guide 5.0mm for Ø4.3mm Drill Bit - Proximal Femur Plate	1
PGW-2.0	Guide Wire Ø2.0mm x Length 250mm	5
PDB-5.0	Cannulated Drill Bit Ø5.0mm x Length 230mm	2
7-041-07	Screw Caddy for Proximal Femur Plate System	1
7-041-10	Implant Tray for Proximal Femur Plate System	2
7-041-09	Instrument Tray for Proximal Femur Plate Instrument Set	1
7-041-08	Sterilization Container for Proximal Femur Plate Instrument Set	1

The background features a series of thin, light blue wavy lines that flow from the top left towards the bottom right, creating a sense of movement and depth. The lines are more densely packed in some areas, creating a mesh-like effect.

**5.0MM WISE-LOCK
LARGE FRAGMENT
TIBIA**




Indications


- Proximal shaft fractures
- Metaphyseal fractures
- Intra-articular fractures
- Periprosthetic fractures




4.5/5.0mm Wise-Lock Proximal Lateral Tibia Plate



Round Locking Hole



Compression and Locking Hole



STERILE

Holes	Direction	TITANIUM	Stainless Steel
5	Left	TI-733.205L	733.205L
5	Right	TI-733.205R	733.205R
7	Left	TI-733.207L	733.207L
7	Right	TI-733.207R	733.207R
9	Left	TI-733.209L	733.209L
9	Right	TI-733.209R	733.209R
11	Left	TI-733.211L	733.211L
11	Right	TI-733.211R	733.211R
13	Left	TI-733.213L	733.213L
13	Right	TI-733.213R	733.213R

Holes	Direction	TITANIUM	Stainless Steel
5	Left	TI-733.205L-S	733.205L-S
5	Right	TI-733.205R-S	733.205R-S
7	Left	TI-733.207L-S	733.207L-S
7	Right	TI-733.207R-S	733.207R-S
9	Left	TI-733.209L-S	733.209L-S
9	Right	TI-733.209R-S	733.209R-S
11	Left	TI-733.211L-S	733.211L-S
11	Right	TI-733.211R-S	733.211R-S
13	Left	TI-733.213L-S	733.213L-S
13	Right	TI-733.213R-S	733.213R-S




Indications

Wise-Lock Medial Proximal Tibia Plates 4.5/5.0 are intended to buttress metaphyseal fractures of the medial tibia plateau, split-type fractures of the medial tibia plateau, medial split fractures with associated depressions and split or depression fractures of the medial tibia plateau. The plates may also be used for fixation of the proximal quarter (lateral and medial) of the tibia, as well as segmental fractures of the proximal tibia. The Wise-Lock Medial Proximal Tibia Plates 4.5/5.0 may also be used for fixation of nonunions and malunions of the medial proximal tibia and tibia shaft, as well as opening and closing wedge tibia osteotomies.



4.5/5.0mm Wise-Lock Medial Proximal Tibia Plate



STERILE

Holes	Direction	TITANIUM	Stainless Steel
4	Left	TI-747.204L	747.204L
4	Right	TI-747.204R	747.204R
6	Left	TI-747.206L	747.206L
6	Right	TI-747.206R	747.206R
8	Left	TI-747.208L	747.208L
8	Right	TI-747.208R	747.208R
10	Left	TI-747.210L	747.210L
10	Right	TI-747.210R	747.210R
12	Left	TI-747.212L	747.212L
12	Right	TI-747.212R	747.212R
14	Left	TI-747.214L	747.214L
14	Right	TI-747.214R	747.214R
16	Left	TI-747.216L	747.216L
16	Right	TI-747.216R	747.216R

Holes	Direction	TITANIUM	Stainless Steel
4	Left	TI-747.204L-S	747.204L-S
4	Right	TI-747.204R-S	747.204R-S
6	Left	TI-747.206L-S	747.206L-S
6	Right	TI-747.206R-S	747.206R-S
8	Left	TI-747.208L-S	747.208L-S
8	Right	TI-747.208R-S	747.208R-S
10	Left	TI-747.210L-S	747.210L-S
10	Right	TI-747.210R-S	747.210R-S
12	Left	TI-747.212L-S	747.212L-S
12	Right	TI-747.212R-S	747.212R-S
14	Left	TI-747.214L-S	747.214L-S
14	Right	TI-747.214R-S	747.214R-S
16	Left	TI-747.216L-S	747.216L-S
16	Right	TI-747.216R-S	747.216R-S




Indications

Wise-Lock 4.5mm Proximal Tibia Plate is intended for treatment of osteopenic bone, tibia osteotomies, nonunions, malunions and fractures of the proximal tibia including:


- Simple, comminuted fractures
- Lateral wedge fractures
- Depression medial wedge fractures
- Bicondylar combination of lateral wedge and depression fractures
- Periprosthetic fractures
- Proximal fractures with associated shaft fractures



4.5/5.0mm Wise-Lock Proximal Tibia Plate



Round Locking Hole



Compression and Locking Hole



STERILE

Holes	Direction	TITANIUM	Stainless Steel
4	Left	TI-749.104L	749.104L
4	Right	TI-749.104R	749.104R
6	Left	TI-749.106L	749.106L
6	Right	TI-749.106R	749.106R
8	Left	TI-749.108L	749.108L
8	Right	TI-749.108R	749.108R
10	Left	TI-749.110L	749.110L
10	Right	TI-749.110R	749.110R
12	Left	TI-749.112L	749.112L
12	Right	TI-749.112R	749.112R
14	Left	TI-749.114L	749.114L
14	Right	TI-749.114R	749.114L

Holes	Direction	TITANIUM	Stainless Steel
4	Left	TI-749.104L-S	749.104L-S
4	Right	TI-749.104R-S	749.104R-S
6	Left	TI-749.106L-S	749.106L-S
6	Right	TI-749.106R-S	749.106R-S
8	Left	TI-749.108L-S	749.108L-S
8	Right	TI-749.108R-S	749.108R-S
10	Left	TI-749.110L-S	749.110L-S
10	Right	TI-749.110R-S	749.110R-S
12	Left	TI-749.112L-S	749.112L-S
12	Right	TI-749.112R-S	749.112R-S
14	Left	TI-749.114L-S	749.114L-S
14	Right	TI-749.114R-S	749.114L-S

The background features a series of thin, light blue wavy lines that flow from the top left towards the bottom right, creating a sense of movement and depth. The lines are more densely packed in some areas, creating a mesh-like effect.

**5.0MM WISE-LOCK
LARGE FRAGMENT
GENERAL**



4.5/5.0mm Wise-Lock Reconstruction Plate



Compression and
Locking Hole




Holes	TITANIUM	Stainless Steel
3	TI-728.203	728.203
4	TI-728.204	728.204
5	TI-728.205	728.205
6	TI-728.206	728.206
7	TI-728.207	728.207
8	TI-728.208	728.208
9	TI-728.209	728.209
10	TI-728.210	728.210
11	TI-728.211	728.211
12	TI-728.212	728.212
13	TI-728.213	728.213
14	TI-728.214	728.214
15	TI-728.215	728.215
16	TI-728.216	728.216



Holes	TITANIUM	Stainless Steel
3	TI-728.203-S	728.203-S
4	TI-728.204-S	728.204-S
5	TI-728.205-S	728.205-S
6	TI-728.206-S	728.206-S
7	TI-728.207-S	728.207-S
8	TI-728.208-S	728.208-S
9	TI-728.209-S	728.209-S
10	TI-728.210-S	728.210-S
11	TI-728.211-S	728.211-S
12	TI-728.212-S	728.212-S
13	TI-728.213-S	728.213-S
14	TI-728.214-S	728.214-S
15	TI-728.215-S	728.215-S
16	TI-728.216-S	728.216-S



4.5/5.0mm Wise-Lock Narrow Dynamic Compression Plate with LC under cuts



STERILE

Holes	TITANIUM	Stainless Steel
2	TI-723.202	723.202
3	TI-723.203	723.203
4	TI-723.204	723.204
5	TI-723.205	723.205
6	TI-723.206	723.206
7	TI-723.207	723.207
8	TI-723.208	723.208
9	TI-723.209	723.209
10	TI-723.210	723.210
11	TI-723.211	723.211
12	TI-723.212	723.212
14	TI-723.214	723.214
16	TI-723.216	723.216
18	TI-723.218	723.218
20	TI-723.220	723.220
22	TI-723.222	723.222
24	TI-723.224	723.224

Holes	TITANIUM	Stainless Steel
2	TI-723.202	723.202
3	TI-723.203	723.203
4	TI-723.204	723.204
5	TI-723.205	723.205
6	TI-723.206	723.206
7	TI-723.207	723.207
8	TI-723.208	723.208
9	TI-723.209	723.209
10	TI-723.210	723.210
11	TI-723.211	723.211
12	TI-723.212	723.212
14	TI-723.214	723.214
16	TI-723.216	723.216
18	TI-723.218	723.218
20	TI-723.220	723.220
22	TI-723.222	723.222
24	TI-723.224	723.224



4.5/5.0mm Wise-Lock Broad Dynamic Compression Plate with LC under cuts



Compression and Locking Hole




Holes	TITANIUM	Stainless Steel
5	TI-725.205	725.205
6	TI-725.206	725.206
7	TI-725.207	725.207
8	TI-725.208	725.208
9	TI-725.209	725.209
10	TI-725.210	725.210
11	TI-725.211	725.211
12	TI-725.212	725.212
13	TI-725.213	725.213
14	TI-725.214	725.214
15	TI-725.215	725.215
16	TI-725.216	725.216
18	TI-725.218	725.218
20	TI-725.220	725.220
22	TI-725.222	725.222
24	TI-725.224	725.224



Holes	TITANIUM	Stainless Steel
5	TI-725.205-S	725.205-S
6	TI-725.206-S	725.206-S
7	TI-725.207-S	725.207-S
8	TI-725.208-S	725.208-S
9	TI-725.209-S	725.209-S
10	TI-725.210-S	725.210-S
11	TI-725.211-S	725.211-S
12	TI-725.212-S	725.212-S
13	TI-725.213-S	725.213-S
14	TI-725.214-S	725.214-S
15	TI-725.215-S	725.215-S
16	TI-725.216-S	725.216-S
18	TI-725.218-S	725.218-S
20	TI-725.220-S	725.220-S
22	TI-725.222-S	725.222-S
24	TI-725.224-S	725.224-S



4.5/5.0mm Wise-Lock Curved Broad Plate




Holes	TITANIUM	Stainless Steel
12	TI-750.112	750.112
13	TI-750.113	750.113
14	TI-750.114	750.114
15	TI-750.115	750.115
16	TI-750.116	750.116
17	TI-750.117	750.117
18	TI-750.118	750.118



Holes	TITANIUM	Stainless Steel
12	TI-750.112-S	750.112-S
13	TI-750.113-S	750.113-S
14	TI-750.114-S	750.114-S
15	TI-750.115-S	750.115-S
16	TI-750.116-S	750.116-S
17	TI-750.117-S	750.117-S
18	TI-750.118-S	750.118-S



4.5/5.0mm Wise-Lock "L" Buttress Plate




Holes	Direction	TITANIUM	Stainless Steel
3	Left	TI-706.203L	706.203L
3	Right	TI-706.203R	706.203R
4	Left	TI-706.204L	706.204L
4	Right	TI-706.204R	706.204R
5	Left	TI-706.205L	706.205L
5	Right	TI-706.205R	706.205R
6	Left	TI-706.206L	706.206L
6	Right	TI-706.206R	706.206R
7	Left	TI-706.207L	706.207L
7	Right	TI-706.207R	706.207R
8	Left	TI-706.208L	706.208L
8	Right	TI-706.208R	706.208R
9	Left	TI-706.209L	706.209L
9	Right	TI-706.209R	706.209R
10	Left	TI-706.210L	706.210L
10	Right	TI-706.210R	706.210R
12	Left	TI-706.212L	706.212L
12	Right	TI-706.212R	706.212R

STERILE

Holes	Direction	TITANIUM	Stainless Steel
3	Left	TI-706.203L-S	706.203L-S
3	Right	TI-706.203R-S	706.203R-S
4	Left	TI-706.204L-S	706.204L-S
4	Right	TI-706.204R-S	706.204R-S
5	Left	TI-706.205L-S	706.205L-S
5	Right	TI-706.205R-S	706.205R-S
6	Left	TI-706.206L-S	706.206L-S
6	Right	TI-706.206R-S	706.206R-S
7	Left	TI-706.207L-S	706.207L-S
7	Right	TI-706.207R-S	706.207R-S
8	Left	TI-706.208L-S	706.208L-S
8	Right	TI-706.208R-S	706.208R-S
9	Left	TI-706.209L-S	706.209L-S
9	Right	TI-706.209R-S	706.209R-S
10	Left	TI-706.210L-S	706.210L-S
10	Right	TI-706.210R-S	706.210R-S
12	Left	TI-706.212L-S	706.212L-S
12	Right	TI-706.212R-S	706.212R



4.5/5.0mm Wise-Lock "T" Buttress Plate



Round Locking Hole



Compression and Locking Hole



Holes	TITANIUM	Stainless Steel
4	TI-714.204	714.204
5	TI-714.205	714.205
6	TI-714.206	714.206
7	TI-714.207	714.207
8	TI-714.208	714.208
9	TI-714.209	714.209
10	TI-714.210	714.210
12	TI-714.212	714.212

STERILE

Holes	TITANIUM	Stainless Steel
4	TI-714.204-S	714.204-S
5	TI-714.205-S	714.205-S
6	TI-714.206-S	714.206-S
7	TI-714.207-S	714.207-S
8	TI-714.208-S	714.208-S
9	TI-714.209-S	714.209-S
10	TI-714.210-S	714.210-S
12	TI-714.212 -S	714.212-S



4.5/5.0mm Wise-Lock "T" Plate



Compression and Locking Hole



Holes	TITANIUM	Stainless Steel
3	TI-709.203	709.203
4	TI-709.204	709.204
5	TI-709.205	709.205
6	TI-709.206	709.206
7	TI-709.207	709.207
8	TI-709.208	709.208
9	TI-709.209	709.209
10	TI-709.210	709.210
12	TI-709.212	709.212



Holes	TITANIUM	Stainless Steel
3	TI-709.203-S	709.203-S
4	TI-709.204-S	709.204-S
5	TI-709.205-S	709.205-S
6	TI-709.206-S	709.206-S
7	TI-709.207-S	709.207-S
8	TI-709.208-S	709.208-S
9	TI-709.209-S	709.209-S
10	TI-709.210-S	709.210-S
12	TI-709.212-S	709.212-S



SCREWS



4.5mm Cortical Screw, Self-Tapping, (Hex Head)



STERILE

Length	Stainless Steel	Titanium
* 12mm	105.212	TI-105.212
* 14mm	105.214	TI-105.214
* 16mm	105.216	TI-105.216
* 18mm	105.218	TI-105.218
20mm	105.220	TI-105.220
22mm	105.222	TI-105.222
24mm	105.224	TI-105.224
26mm	105.226	TI-105.226
28mm	105.228	TI-105.228
30mm	105.230	TI-105.230
32mm	105.232	TI-105.232
34mm	105.234	TI-105.234
36mm	105.236	TI-105.236
38mm	105.238	TI-105.238
40mm	105.240	TI-105.240
42mm	105.242	TI-105.242
44mm	105.244	TI-105.244
46mm	105.246	TI-105.246
48mm	105.248	TI-105.248
50mm	105.250	TI-105.250
52mm	105.252	TI-105.252
54mm	105.254	TI-105.254
* 55mm	105.255	TI-105.255
56mm	105.256	TI-105.256
58mm	105.258	TI-105.258
60mm	105.260	TI-105.260
65mm	105.265	TI-105.265
70mm	105.270	TI-105.270
80mm	105.280	TI-105.280
* 85mm	105.285	TI-105.285
* 90mm	105.290	TI-105.290

* Sizes not available in Screw Caddy

Length	Stainless Steel	Titanium
* 12mm	105.212-S	TI-105.212-S
* 14mm	105.214-S	TI-105.214-S
* 16mm	105.216-S	TI-105.216-S
* 18mm	105.218-S	TI-105.218-S
20mm	105.220-S	TI-105.220-S
22mm	105.222-S	TI-105.222-S
24mm	105.224-S	TI-105.224-S
26mm	105.226-S	TI-105.226-S
28mm	105.228-S	TI-105.228-S
30mm	105.230-S	TI-105.230-S
32mm	105.232-S	TI-105.232-S
34mm	105.234-S	TI-105.234-S
36mm	105.236-S	TI-105.236-S
38mm	105.238-S	TI-105.238-S
40mm	105.240-S	TI-105.240-S
42mm	105.242-S	TI-105.242-S
44mm	105.244-S	TI-105.244-S
46mm	105.246-S	TI-105.246-S
48mm	105.248-S	TI-105.248-S
50mm	105.250-S	TI-105.250-S
52mm	105.252-S	TI-105.252-S
54mm	105.254-S	TI-105.254-S
* 55mm	105.255-S	TI-105.255-S
56mm	105.256-S	TI-105.256-S
58mm	105.258-S	TI-105.258-S
60mm	105.260-S	TI-105.260-S
65mm	105.265-S	TI-105.265-S
70mm	105.270-S	TI-105.270-S
80mm	105.280-S	TI-105.280-S
* 85mm	105.285-S	TI-105.285-S
* 90mm	105.290-S	TI-105.290-S

* Sizes not available in Screw Caddy






5.0mm Wise-Lock Screw, Self-Tapping, (Hex Head)



Length	Stainless Steel	Titanium
12mm	119.012	TI-119.012
14mm	119.014	TI-119.014
16mm	119.016	TI-119.016
18mm	119.018	TI-119.018
20mm	119.020	TI-119.020
22mm	119.022	TI-119.022
24mm	119.024	TI-119.024
26mm	119.026	TI-119.026
28mm	119.028	TI-119.028
30mm	119.030	TI-119.030
32mm	119.032	TI-119.032
34mm	119.034	TI-119.034
36mm	119.036	TI-119.036
38mm	119.038	TI-119.038
40mm	119.040	TI-119.040
42mm	119.042	TI-119.042
44mm	119.044	TI-119.044
46mm	119.046	TI-119.046
48mm	119.048	TI-119.048
50mm	119.050	TI-119.050
52mm	119.052	TI-119.052
54mm	119.054	TI-119.054
56mm	119.056	TI-119.056
58mm	119.058	TI-119.058
60mm	119.060	TI-119.060
65mm	119.065	TI-119.065
70mm	119.070	TI-119.070
75mm	119.075	TI-119.075
80mm	119.080	TI-119.080
85mm	119.085	TI-119.085
90mm	119.090	TI-119.090

STERILE

Length	Stainless Steel	Titanium
12mm	119.012-S	TI-119.012-S
14mm	119.014-S	TI-119.014-S
16mm	119.016-S	TI-119.016-S
18mm	119.018-S	TI-119.018-S
20mm	119.020-S	TI-119.020-S
22mm	119.022-S	TI-119.022-S
24mm	119.024-S	TI-119.024-S
26mm	119.026-S	TI-119.026-S
28mm	119.028-S	TI-119.028-S
30mm	119.030-S	TI-119.030-S
32mm	119.032-S	TI-119.032-S
34mm	119.034-S	TI-119.034-S
36mm	119.036-S	TI-119.036-S
38mm	119.038-S	TI-119.038-S
40mm	119.040-S	TI-119.040-S
42mm	119.042-S	TI-119.042-S
44mm	119.044-S	TI-119.044-S
46mm	119.046-S	TI-119.046-S
48mm	119.048-S	TI-119.048-S
50mm	119.050-S	TI-119.050-S
52mm	119.052-S	TI-119.052-S
54mm	119.054-S	TI-119.054-S
56mm	119.056-S	TI-119.056-S
58mm	119.058-S	TI-119.058-S
60mm	119.060-S	TI-119.060-S
65mm	119.065-S	TI-119.065-S
70mm	119.070-S	TI-119.070-S
75mm	119.075-S	TI-119.075-S
80mm	119.080-S	TI-119.080-S
85mm	119.085-S	TI-119.085-S
90mm	119.090-S	TI-119.090-S






5.0mm Wise-Lock Screw, Self-Tapping & Self-Drilling, (Hex Head)



Length	Stainless Steel	Titanium
12mm	119.212	TI-119.212
14mm	119.214	TI-119.214
16mm	119.216	TI-119.216
18mm	119.218	TI-119.218
20mm	119.220	TI-119.220
22mm	119.222	TI-119.222
24mm	119.224	TI-119.224
26mm	119.226	TI-119.226
28mm	119.228	TI-119.228
30mm	119.230	TI-119.230
32mm	119.232	TI-119.232
34mm	119.234	TI-119.234
36mm	119.236	TI-119.236
38mm	119.238	TI-119.238
40mm	119.240	TI-119.240
42mm	119.242	TI-119.242
44mm	119.244	TI-119.244
46mm	119.246	TI-119.246
48mm	119.248	TI-119.248
50mm	119.250	TI-119.250
55mm	119.255	TI-119.255
60mm	119.260	TI-119.260
65mm	119.265	TI-119.265
70mm	119.270	TI-119.270
75mm	119.275	TI-119.275
80mm	119.280	TI-119.280
85mm	119.285	TI-119.285
90mm	119.290	TI-119.290

STERILE

Length	Stainless Steel	Titanium
12mm	119.212-S	TI-119.212-S
14mm	119.214-S	TI-119.214-S
16mm	119.216-S	TI-119.216-S
18mm	119.218-S	TI-119.218-S
20mm	119.220-S	TI-119.220-S
22mm	119.222-S	TI-119.222-S
24mm	119.224-S	TI-119.224-S
26mm	119.226-S	TI-119.226-S
28mm	119.228-S	TI-119.228-S
30mm	119.230-S	TI-119.230-S
32mm	119.232-S	TI-119.232-S
34mm	119.234-S	TI-119.234-S
36mm	119.236-S	TI-119.236-S
38mm	119.238-S	TI-119.238-S
40mm	119.240-S	TI-119.240-S
42mm	119.242-S	TI-119.242-S
44mm	119.244-S	TI-119.244-S
46mm	119.246-S	TI-119.246-S
48mm	119.248-S	TI-119.248-S
50mm	119.250-S	TI-119.250-S
55mm	119.255-S	TI-119.255-S
60mm	119.260-S	TI-119.260-S
65mm	119.265-S	TI-119.265-S
70mm	119.270-S	TI-119.270-S
75mm	119.275-S	TI-119.275-S
80mm	119.280-S	TI-119.280-S
85mm	119.285-S	TI-119.285-S
90mm	119.290-S	TI-119.290-S






5.0mm Wise-Lock Cancellous Screw, Full Thread, Self-Tapping, (Hex Head)



Length	Stainless Steel	Titanium
20mm	1360-20	TI-1360-20
22mm	1360-22	TI-1360-22
24mm	1360-24	TI-1360-24
26mm	1360-26	TI-1360-26
28mm	1360-28	TI-1360-28
30mm	1360-30	TI-1360-30
32mm	1360-32	TI-1360-32
34mm	1360-34	TI-1360-34
36mm	1360-36	TI-1360-36
38mm	1360-38	TI-1360-38
40mm	1360-40	TI-1360-40
42mm	1360-42	TI-1360-42
44mm	1360-44	TI-1360-44
46mm	1360-46	TI-1360-46
48mm	1360-48	TI-1360-48
50mm	1360-50	TI-1360-50
55mm	1360-55	TI-1360-55
60mm	1360-60	TI-1360-60
65mm	1360-65	TI-1360-65
70mm	1360-70	TI-1360-70
75mm	1360-75	TI-1360-75
80mm	1360-80	TI-1360-80
85mm	1360-85	TI-1360-85
90mm	1360-90	TI-1360-90
95mm	1360-95	TI-1360-95
100mm	1360-100	TI-1360-100
105mm	1360-105	TI-1360-105
110mm	1360-110	TI-1360-110

STERILE

Length	Stainless Steel	Titanium
20mm	1360-20-S	TI-1360-20-S
22mm	1360-22-S	TI-1360-22-S
24mm	1360-24-S	TI-1360-24-S
26mm	1360-26-S	TI-1360-26-S
28mm	1360-28-S	TI-1360-28-S
30mm	1360-30-S	TI-1360-30-S
32mm	1360-32-S	TI-1360-32-S
34mm	1360-34-S	TI-1360-34-S
36mm	1360-36-S	TI-1360-36-S
38mm	1360-38-S	TI-1360-38-S
40mm	1360-40-S	TI-1360-40-S
42mm	1360-42-S	TI-1360-42-S
44mm	1360-44-S	TI-1360-44-S
46mm	1360-46-S	TI-1360-46-S
48mm	1360-48-S	TI-1360-48-S
50mm	1360-50-S	TI-1360-50-S
55mm	1360-55-S	TI-1360-55-S
60mm	1360-60-S	TI-1360-60-S
65mm	1360-65-S	TI-1360-65-S
70mm	1360-70-S	TI-1360-70-S
75mm	1360-75-S	TI-1360-75-S
80mm	1360-80-S	TI-1360-80-S
85mm	1360-85-S	TI-1360-85-S
90mm	1360-90-S	TI-1360-90-S
95mm	1360-95-S	TI-1360-95-S
100mm	1360-100-S	TI-1360-100-S
105mm	1360-105-S	TI-1360-105-S
110mm	1360-110-S	TI-1360-110-S





5.0mm Wise-Lock Cancellous Screw, Short Thread, Self-Tapping, (Hex Head)



Length	Stainless Steel	Titanium
20mm	1256-20	TI-1256-20
22mm	1256-22	TI-1256-22
24mm	1256-24	TI-1256-24
26mm	1256-26	TI-1256-26
28mm	1256-28	TI-1256-28
30mm	1256-30	TI-1256-30
32mm	1256-32	TI-1256-32
34mm	1256-34	TI-1256-34
36mm	1256-36	TI-1256-36
38mm	1256-38	TI-1256-38
40mm	1256-40	TI-1256-40
42mm	1256-42	TI-1256-42
44mm	1256-44	TI-1256-44
46mm	1256-46	TI-1256-46
48mm	1256-48	TI-1256-48
50mm	1256-50	TI-1256-50
55mm	1256-55	TI-1256-55
60mm	1256-60	TI-1256-60
65mm	1256-65	TI-1256-65
70mm	1256-70	TI-1256-70
75mm	1256-75	TI-1256-75
80mm	1256-80	TI-1256-80
85mm	1256-85	TI-1256-85
90mm	1256-90	TI-1256-90
95mm	1256-95	TI-1256-95
100mm	1256-100	TI-1256-100
105mm	1256-105	TI-1256-105
110mm	1256-110	TI-1256-110

STERILE

Length	Stainless Steel	Titanium
20mm	1256-20-S	TI-1256-20-S
22mm	1256-22-S	TI-1256-22-S
24mm	1256-24-S	TI-1256-24-S
26mm	1256-26-S	TI-1256-26-S
28mm	1256-28-S	TI-1256-28-S
30mm	1256-30-S	TI-1256-30-S
32mm	1256-32-S	TI-1256-32-S
34mm	1256-34-S	TI-1256-34-S
36mm	1256-36-S	TI-1256-36-S
38mm	1256-38-S	TI-1256-38-S
40mm	1256-40-S	TI-1256-40-S
42mm	1256-42-S	TI-1256-42-S
44mm	1256-44-S	TI-1256-44-S
46mm	1256-46-S	TI-1256-46-S
48mm	1256-48-S	TI-1256-48-S
50mm	1256-50-S	TI-1256-50-S
55mm	1256-55-S	TI-1256-55-S
60mm	1256-60-S	TI-1256-60-S
65mm	1256-65-S	TI-1256-65-S
70mm	1256-70-S	TI-1256-70-S
75mm	1256-75-S	TI-1256-75-S
80mm	1256-80-S	TI-1256-80-S
85mm	1256-85-S	TI-1256-85-S
90mm	1256-90-S	TI-1256-90-S
95mm	1256-95-S	TI-1256-95-S
100mm	1256-100-S	TI-1256-100-S
105mm	1256-105-S	TI-1256-105-S
110mm	1256-110-S	TI-1256-110-S






5.0mm Wise-Lock Cannulated Screw, Full Thread, Self-Tapping



Length	Stainless Steel	Titanium
20mm	156.020	TI-156.020
25mm	156.025	TI-156.025
30mm	156.030	TI-156.030
35mm	156.035	TI-156.035
40mm	156.040	TI-156.040
45mm	156.045	TI-156.045
50mm	156.050	TI-156.050
55mm	156.055	TI-156.055
60mm	156.060	TI-156.060
65mm	156.065	TI-156.065
70mm	156.070	TI-156.070
75mm	156.075	TI-156.075
80mm	156.080	TI-156.080
85mm	156.085	TI-156.085
90mm	156.090	TI-156.090
95mm	156.095	TI-156.095
100mm	156.100	TI-156.100
105mm	156.105	TI-156.105
110mm	156.110	TI-156.110
115mm	156.115	TI-156.115
120mm	156.120	TI-156.120

STERILE

Length	Stainless Steel	Titanium
20mm	156.020-S	TI-156.020-S
25mm	156.025-S	TI-156.025-S
30mm	156.030-S	TI-156.030-S
35mm	156.035-S	TI-156.035-S
40mm	156.040-S	TI-156.040-S
45mm	156.045-S	TI-156.045-S
50mm	156.050-S	TI-156.050-S
55mm	156.055-S	TI-156.055-S
60mm	156.060-S	TI-156.060-S
65mm	156.065-S	TI-156.065-S
70mm	156.070-S	TI-156.070-S
75mm	156.075-S	TI-156.075-S
80mm	156.080-S	TI-156.080-S
85mm	156.085-S	TI-156.085-S
90mm	156.090-S	TI-156.090-S
95mm	156.095-S	TI-156.095-S
100mm	156.100-S	TI-156.100-S
105mm	156.105-S	TI-156.105-S
110mm	156.110-S	TI-156.110-S
115mm	156.115-S	TI-156.115-S
120mm	156.120-S	TI-156.120-S





5.0mm Wise-Lock Cannulated Screw, Partial Thread, Self-Tapping



Length	Stainless Steel	Titanium
30mm	1269.30	TI-1269.30
35mm	1269.35	TI-1269.35
40mm	1269.40	TI-1269.40
45mm	1269.45	TI-1269.45
50mm	1269.50	TI-1269.50
55mm	1269.55	TI-1269.55
60mm	1269.60	TI-1269.60
65mm	1269.65	TI-1269.65
70mm	1269.70	TI-1269.70
75mm	1269.75	TI-1269.75
80mm	1269.80	TI-1269.80
85mm	1269.85	TI-1269.85
90mm	1269.90	TI-1269.90
95mm	1269.95	TI-1269.95
100mm	1269.100	TI-1269.100

STERILE

Length	Stainless Steel	Titanium
30mm	1269.30-S	TI-1269.30-S
35mm	1269.35-S	TI-1269.35-S
40mm	1269.40-S	TI-1269.40-S
45mm	1269.45-S	TI-1269.45-S
50mm	1269.50-S	TI-1269.50-S
55mm	1269.55-S	TI-1269.55-S
60mm	1269.60-S	TI-1269.60-S
65mm	1269.65-S	TI-1269.65-S
70mm	1269.70-S	TI-1269.70-S
75mm	1269.75-S	TI-1269.75-S
80mm	1269.80-S	TI-1269.80-S
85mm	1269.85-S	TI-1269.85-S
90mm	1269.90-S	TI-1269.90-S
95mm	1269.95-S	TI-1269.95-S
100mm	1269.100-S	TI-1269.100-S





6.5mm Cancellous Screw, 16mm Thread



Length	Stainless Steel	Titanium
25mm	111.025	TI-111.025
30mm	111.030	TI-111.030
35mm	111.035	TI-111.035
40mm	111.040	TI-111.040
45mm	111.045	TI-111.045
50mm	111.050	TI-111.050
55mm	111.055	TI-111.055
60mm	111.060	TI-111.060
65mm	111.065	TI-111.065
70mm	111.070	TI-111.070
75mm	111.075	TI-111.075
80mm	111.080	TI-111.080
85mm	111.085	TI-111.085
90mm	111.090	TI-111.090
95mm	111.095	TI-111.095
100mm	111.100	TI-111.100
* 105mm	111.105	TI-111.105
* 110mm	111.110	TI-111.110
* 115mm	111.115	TI-111.115
* 120mm	111.120	TI-111.120

* Sizes not available in Screw Caddy



Length	Stainless Steel	Titanium
25mm	111.025-S	TI-111.025-S
30mm	111.030-S	TI-111.030-S
35mm	111.035-S	TI-111.035-S
40mm	111.040-S	TI-111.040-S
45mm	111.045-S	TI-111.045-S
50mm	111.050-S	TI-111.050-S
55mm	111.055-S	TI-111.055-S
60mm	111.060-S	TI-111.060-S
65mm	111.065-S	TI-111.065-S
70mm	111.070-S	TI-111.070-S
75mm	111.075-S	TI-111.075-S
80mm	111.080-S	TI-111.080-S
85mm	111.085-S	TI-111.085-S
90mm	111.090-S	TI-111.090-S
95mm	111.095-S	TI-111.095-S
100mm	111.100-S	TI-111.100-S
* 105mm	111.105-S	TI-111.105-S
* 110mm	111.110-S	TI-111.110-S
* 115mm	111.115-S	TI-111.115-S
* 120mm	111.120-S	TI-111.120-S

* Sizes not available in Screw Caddy





6.5mm Cancellous Screw, 32mm Thread



Length	Stainless Steel	Titanium
35mm	112.035	TI-112.035
40mm	112.040	TI-112.040
45mm	112.045	TI-112.045
50mm	112.050	TI-112.050
55mm	112.055	TI-112.055
60mm	112.060	TI-112.060
65mm	112.065	TI-112.065
70mm	112.070	TI-112.070
75mm	112.075	TI-112.075
80mm	112.080	TI-112.080
85mm	112.085	TI-112.085
90mm	112.090	TI-112.090
95mm	112.095	TI-112.095
100mm	112.100	TI-112.100
* 105mm	112.105	TI-112.105
* 110mm	112.110	TI-112.110
* 115mm	112.115	TI-112.115
* 120mm	112.120	TI-112.120

* Sizes not available in Screw Caddy



Length	Stainless Steel	Titanium
35mm	112.035-S	TI-112.035-S
40mm	112.040-S	TI-112.040-S
45mm	112.045-S	TI-112.045-S
50mm	112.050-S	TI-112.050-S
55mm	112.055-S	TI-112.055-S
60mm	112.060-S	TI-112.060-S
65mm	112.065-S	TI-112.065-S
70mm	112.070-S	TI-112.070-S
75mm	112.075-S	TI-112.075-S
80mm	112.080-S	TI-112.080-S
85mm	112.085-S	TI-112.085-S
90mm	112.090-S	TI-112.090-S
95mm	112.095-S	TI-112.095-S
100mm	112.100-S	TI-112.100-S
* 105mm	112.105-S	TI-112.105-S
* 110mm	112.110-S	TI-112.110-S
* 115mm	112.115-S	TI-112.115-S
* 120mm	112.120-S	TI-112.120-S

* Sizes not available in Screw Caddy





6.5mm Cancellous Screw, Full Thread




Length	Stainless Steel	Titanium
35mm	113.035	TI-113.035
40mm	113.040	TI-113.040
45mm	113.045	TI-113.045
50mm	113.050	TI-113.050
55mm	113.055	TI-113.055
60mm	113.060	TI-113.060
65mm	113.065	TI-113.065
70mm	113.070	TI-113.070
75mm	113.075	TI-113.075
80mm	113.080	TI-113.080
85mm	113.085	TI-113.085
90mm	113.090	TI-113.090
95mm	113.095	TI-113.095
100mm	113.100	TI-113.100
* 105mm	113.105	TI-113.105
* 110mm	113.110	TI-113.110
* 115mm	113.115	TI-113.115
* 120mm	113.120	TI-113.120

* Sizes not available in Screw Caddy



Length	Stainless Steel	Titanium
35mm	113.035-S	TI-113.035-S
40mm	113.040-S	TI-113.040-S
45mm	113.045-S	TI-113.045-S
50mm	113.050-S	TI-113.050-S
55mm	113.055-S	TI-113.055-S
60mm	113.060-S	TI-113.060-S
65mm	113.065-S	TI-113.065-S
70mm	113.070-S	TI-113.070-S
75mm	113.075-S	TI-113.075-S
80mm	113.080-S	TI-113.080-S
85mm	113.085-S	TI-113.085-S
90mm	113.090-S	TI-113.090-S
95mm	113.095-S	TI-113.095-S
100mm	113.100-S	TI-113.100-S
* 105mm	113.105-S	TI-113.105-S
* 110mm	113.110-S	TI-113.110-S
* 115mm	113.115-S	TI-113.115-S
* 120mm	113.120-S	TI-113.120-S

* Sizes not available in Screw Caddy






7.3mm Wise-Lock Cannulated Screw, Self-Tapping Partial Thread



Length	Stainless Steel	Titanium
20mm	157.020	TI-157.020
25mm	157.025	TI-157.025
30mm	157.030	TI-157.030
35mm	157.035	TI-157.035
40mm	157.040	TI-157.040
45mm	157.045	TI-157.045
50mm	157.050	TI-157.050
55mm	157.055	TI-157.055
60mm	157.060	TI-157.060
65mm	157.065	TI-157.065
70mm	157.070	TI-157.070
75mm	157.075	TI-157.075
80mm	157.080	TI-157.080
85mm	157.085	TI-157.085
90mm	157.090	TI-157.090
95mm	157.095	TI-157.095
100mm	157.100	TI-157.100
105mm	157.105	TI-157.105
110mm	157.110	TI-157.110
115mm	157.115	TI-157.115
120mm	157.120	TI-157.120

STERILE

Length	Stainless Steel	Titanium
20mm	157.020-S	TI-157.020-S
25mm	157.025-S	TI-157.025-S
30mm	157.030-S	TI-157.030-S
35mm	157.035-S	TI-157.035-S
40mm	157.040-S	TI-157.040-S
45mm	157.045-S	TI-157.045-S
50mm	157.050-S	TI-157.050-S
55mm	157.055-S	TI-157.055-S
60mm	157.060-S	TI-157.060-S
65mm	157.065-S	TI-157.065-S
70mm	157.070-S	TI-157.070-S
75mm	157.075-S	TI-157.075-S
80mm	157.080-S	TI-157.080-S
85mm	157.085-S	TI-157.085-S
90mm	157.090-S	TI-157.090-S
95mm	157.095-S	TI-157.095-S
100mm	157.100-S	TI-157.100-S
105mm	157.105-S	TI-157.105-S
110mm	157.110-S	TI-157.110-S
115mm	157.115-S	TI-157.115-S
120mm	157.120-S	TI-157.120-S






7.3mm Wise-Lock Cannulated Screw, Self Tapping, Full Thread



Length	Stainless Steel	Titanium
20mm	158.020	TI-158.020
25mm	158.025	TI-158.025
30mm	158.030	TI-158.030
35mm	158.035	TI-157.035
40mm	158.040	TI-158.040
45mm	158.045	TI-158.045
50mm	158.050	TI-158.050
55mm	158.055	TI-158.055
60mm	158.060	TI-158.060
65mm	158.065	TI-158.065
70mm	158.070	TI-158.070
75mm	158.075	TI-158.075
80mm	158.080	TI-158.080
85mm	158.085	TI-158.085
90mm	158.090	TI-158.090
95mm	158.095	TI-158.095
100mm	158.100	TI-158.100
105mm	158.105	TI-158.105
110mm	158.110	TI-158.110
115mm	158.115	TI-158.115
120mm	158.120	TI-158.120

STERILE

Length	Stainless Steel	Titanium
20mm	158.020-S	TI-158.020-S
25mm	158.025-S	TI-158.025-S
30mm	158.030-S	TI-158.030-S
35mm	158.035-S	TI-158.035-S
40mm	158.040-S	TI-158.040-S
45mm	158.045-S	TI-158.045-S
50mm	158.050-S	TI-158.050-S
55mm	158.055-S	TI-158.055-S
60mm	158.060-S	TI-158.060-S
65mm	158.065-S	TI-158.065-S
70mm	158.070-S	TI-158.070-S
75mm	158.075-S	TI-158.075-S
80mm	158.080-S	TI-158.080-S
85mm	158.085-S	TI-158.085-S
90mm	158.090-S	TI-158.090-S
95mm	158.095-S	TI-158.095-S
100mm	158.100-S	TI-158.100-S
105mm	158.105-S	TI-158.105-S
110mm	158.110-S	TI-158.110-S
115mm	158.115-S	TI-158.115-S
120mm	158.120-S	TI-158.120-S





TI-007NW	Screw Caddy for 5.0mm Wise-Lock Large Fragment System
----------	---





4.5mm Cortical Screw, Self-Tapping, (Hex Head)

Length	Stainless Steel	Titanium	Unit
20mm	105.220	TI-105.220	3
22mm	105.222	TI-105.222	3
24mm	105.224	TI-105.224	3
26mm	105.226	TI-105.226	3
28mm	105.228	TI-105.228	3
30mm	105.230	TI-105.230	3
32mm	105.232	TI-105.232	3
34mm	105.234	TI-105.234	3
36mm	105.236	TI-105.236	3
38mm	105.238	TI-105.238	3
40mm	105.240	TI-105.240	3
42mm	105.242	TI-105.242	3
44mm	105.244	TI-105.244	3
46mm	105.246	TI-105.246	3
48mm	105.248	TI-105.248	3
50mm	105.250	TI-105.250	3
52mm	105.252	TI-105.252	3
54mm	105.254	TI-105.254	3
56mm	105.256	TI-105.256	3
58mm	105.258	TI-105.258	3
60mm	105.260	TI-105.260	3
65mm	105.265	TI-105.265	3
70mm	105.270	TI-105.270	3
80mm	105.280	TI-105.280	3

5.0mm Wise-Lock Screw, Self-Tapping, (Hex Head)

Length	Stainless Steel	Titanium	Unit
12mm	119.012	TI-119.012	5
14mm	119.014	TI-119.014	5
16mm	119.016	TI-119.016	5
18mm	119.018	TI-119.018	5
20mm	119.020	TI-119.020	5
22mm	119.022	TI-119.022	5
24mm	119.024	TI-119.024	5
26mm	119.026	TI-119.026	5
28mm	119.028	TI-119.028	5
30mm	119.030	TI-119.030	5
32mm	119.032	TI-119.032	5
34mm	119.034	TI-119.034	5
36mm	119.036	TI-119.036	5
38mm	119.038	TI-119.038	5
40mm	119.040	TI-119.040	5
42mm	119.042	TI-119.042	5
44mm	119.044	TI-119.044	5
46mm	119.046	TI-119.046	5
48mm	119.048	TI-119.048	5
50mm	119.050	TI-119.050	5
52mm	119.052	TI-119.052	5
54mm	119.054	TI-119.054	5
56mm	119.056	TI-119.056	5
58mm	119.058	TI-119.058	5
60mm	119.060	TI-119.060	5
65mm	119.065	TI-119.065	5
70mm	119.070	TI-119.070	5
75mm	119.075	TI-119.075	5
80mm	119.080	TI-119.080	5
85mm	119.085	TI-119.085	5
90mm	119.090	TI-119.090	5



5.0mm Wise-Lock Cancellous Screw, Full Thread, Self-Tapping, (Hex Head)

Length	Stainless Steel	Titanium	Unit
20mm	1360-20	TI-1360-20	2
22mm	1360-22	TI-1360-22	2
24mm	1360-24	TI-1360-24	2
26mm	1360-26	TI-1360-26	2
28mm	1360-28	TI-1360-28	2
30mm	1360-30	TI-1360-30	2
32mm	1360-32	TI-1360-32	2
34mm	1360-34	TI-1360-34	2
36mm	1360-36	TI-1360-36	2
38mm	1360-38	TI-1360-38	2
40mm	1360-40	TI-1360-40	2
42mm	1360-42	TI-1360-42	2
44mm	1360-44	TI-1360-44	2
46mm	1360-46	TI-1360-46	2
48mm	1360-48	TI-1360-48	2
50mm	1360-50	TI-1360-50	2
55mm	1360-55	TI-1360-55	2
60mm	1360-60	TI-1360-60	2
65mm	1360-65	TI-1360-65	2
70mm	1360-70	TI-1360-70	2
75mm	1360-75	TI-1360-75	2
80mm	1360-80	TI-1360-80	2
85mm	1360-85	TI-1360-85	2
90mm	1360-90	TI-1360-90	2
95mm	1360-95	TI-1360-95	2
100mm	1360-100	TI-1360-100	2
105mm	1360-105	TI-1360-105	2
110mm	1360-110	TI-1360-110	2

5.0mm Wise-Lock Cancellous Screw, Short Thread, Self-Tapping, (Hex Head)

Length	Stainless Steel	Titanium	Unit
20mm	1256-20	TI-1256-20	2
22mm	1256-22	TI-1256-22	2
24mm	1256-24	TI-1256-24	2
26mm	1256-26	TI-1256-26	2
28mm	1256-28	TI-1256-28	2
30mm	1256-30	TI-1256-30	2
32mm	1256-32	TI-1256-32	2
34mm	1256-34	TI-1256-34	2
36mm	1256-36	TI-1256-36	2
38mm	1256-38	TI-1256-38	2
40mm	1256-40	TI-1256-40	2
42mm	1256-42	TI-1256-42	2
44mm	1256-44	TI-1256-44	2
46mm	1256-46	TI-1256-46	2
48mm	1256-48	TI-1256-48	2
50mm	1256-50	TI-1256-50	2
55mm	1256-55	TI-1256-55	2
60mm	1256-60	TI-1256-60	2
65mm	1256-65	TI-1256-65	2
70mm	1256-70	TI-1256-70	2
75mm	1256-75	TI-1256-75	2
80mm	1256-80	TI-1256-80	2
85mm	1256-85	TI-1256-85	2
90mm	1256-90	TI-1256-90	2
95mm	1256-95	TI-1256-95	2
100mm	1256-100	TI-1256-100	2
105mm	1256-105	TI-1256-105	2
110mm	1256-110	TI-1256-110	2



6.5mm Cancellous Screw,
16mm Thread

Length	Stainless Steel	Titanium	Unit
25mm	111.025	TI-111.025	1
30mm	111.030	TI-111.030	1
35mm	111.035	TI-111.035	1
40mm	111.040	TI-111.040	1
45mm	111.045	TI-111.045	1
50mm	111.050	TI-111.050	1
55mm	111.055	TI-111.055	1
60mm	111.060	TI-111.060	1
65mm	111.065	TI-111.065	1
70mm	111.070	TI-111.070	1
75mm	111.075	TI-111.075	1
80mm	111.080	TI-111.080	1
85mm	111.085	TI-111.085	1
90mm	111.090	TI-111.090	1
95mm	111.095	TI-111.095	1
100mm	111.100	TI-111.100	1

6.5mm Cancellous Screw,
32mm Thread

Length	Stainless Steel	Titanium	Unit
35mm	112.035	TI-112.035	1
40mm	112.040	TI-112.040	1
45mm	112.045	TI-112.045	1
50mm	112.050	TI-112.050	1
55mm	112.055	TI-112.055	1
60mm	112.060	TI-112.060	1
65mm	112.065	TI-112.065	1
70mm	112.070	TI-112.070	1
75mm	112.075	TI-112.075	1
80mm	112.080	TI-112.080	1
85mm	112.085	TI-112.085	1
90mm	112.090	TI-112.090	1
95mm	112.095	TI-112.095	1
100mm	112.100	TI-112.100	1



6.5mm Cancellous Screw, Full Thread

Length	Stainless Steel	Titanium	Unit
35	113.035	TI-113.035	1
40	113.040	TI-113.040	1
45	113.045	TI-113.045	1
50	113.050	TI-113.050	1
55	113.055	TI-113.055	1
60	113.060	TI-113.060	1
65	113.065	TI-113.065	1
70	113.070	TI-113.070	1
75	113.075	TI-113.075	1
80	113.080	TI-113.080	1
85	113.085	TI-113.085	1
90	113.090	TI-113.090	1
95	113.095	TI-113.095	1
100	113.100	TI-113.100	1

5.0mm Wise-Lock Screw, Self-Tapping & Self- rilling, (Hex Head)

Length	Stainless Steel	Titanium	Unit
12	119.212	TI-119.212	5
14	119.214	TI-119.214	5
16	119.216	TI-119.216	5
18	119.218	TI-119.218	5
20	119.220	TI-119.220	5
22	119.222	TI-119.222	5
24	119.224	TI-119.224	5
26	119.226	TI-119.226	5
28	119.228	TI-119.228	5
30	119.230	TI-119.230	5
32	119.232	TI-119.232	5
34	119.234	TI-119.234	5
36	119.236	TI-119.236	5
38	119.238	TI-119.238	5
40	119.240	TI-119.240	5
42	119.242	TI-119.242	5
44	119.244	TI-119.244	5
46	119.246	TI-119.246	5
48	119.248	TI-119.248	5
50	119.250	TI-119.250	5
60	119.260	TI-119.260	5
65	119.265	TI-119.265	5
70	119.270	TI-119.270	5
75	119.275	TI-119.275	5
80	119.280	TI-119.280	5
85	119.285	TI-119.285	5
90	119.290	TI-119.290	5



INSTRUMENTS

**2100-3.2-145**

DRILL BIT WITH QUICK COUPLING END, Ø3.2MM X LENGTH 145MM, FOR LARGE FRAGMENT

**2100-4.5-145**

DRILL BIT WITH QUICK COUPLING END, Ø4.5MM X LENGTH 145MM, FOR LARGE FRAGMENT

**2103-07**

BONE TAP QUICK COUPLING FOR Ø4.5MM CORTICAL SCREWS, FOR LARGE FRAGMENT

**2103-09**

BONE TAP QUICK COUPLING FOR Ø6.5MM CANCELLOUS SCREWS, FOR LARGE FRAGMENT

**2103-4.3-225**

DRILL BIT QUICK COUPLING WITH STOPPER Ø4.3MM X LENGTH 225MM FOR LARGE FRAGMENT

**2106-2.0**

GUIDE SLEEVE FOR Ø2.0MM K. WIRES - LARGE FRAGMENT





3400-02

T-HANDLE COUNTERSINK FOR Ø4.5/6.5MM SCREWS, FOR LARGE FRAGMENT



3420-02

INSERT DRILL SLEEVE, Ø4.5/3.2MM, FOR LARGE FRAGMENT



3443-38-5.0

EXTRACTION SCREW, Ø5.0MM, FOR LARGE FRAGMENT



3408-04

HEXAGONAL SCREWDRIVER SHAFT - 3.5MM TIP, QUICK COUPLING, SHORT



2186-3.5

HSS DRILL BIT, Ø3.5MM, FOR LARGE FRAGMENT



3767-4.5

DOUBLE DRILL GUIDE, Ø4.5/3.2MM





1472-058

SELF-CENTERING DOUBLE DRILL GUIDE, Ø4.5/3.2MM, FOR LARGE FRAGMENT



3767-6.5

DOUBLE DRILL GUIDE, Ø6.5/3.2MM



BT-LF-S

BENDING TEMPLATE, SMALL, FOR LARGE FRAGMENT-LOCKING



BT-LF-M

BENDING TEMPLATE, MEDIUM, FOR WISE-LOCK LARGE FRAGMENT



BT-LF-L

BENDING TEMPLATE, LARGE, FOR WISE-LOCK LARGE FRAGMENT



7-051-05

THREADED DRILL SLEEVE FOR DRILL BIT Ø4.3 FOR WISE-LOCK LARGE FRAGMENT FOR T-HANDLE





GW-2.0-230

GUIDE WIRE WITH THREADED TIP, Ø2.0MM X THREAD LENGTH 10MM X LENGTH 230MM



3766-003

ALLEN KEY, HEX 3.0MM FOR LARGE FRAGMENT



1472-052

PULLING DRILL BIT, Ø4.0MM, FOR LARGE FRAGMENT



3443-400

DEPTH GAUGE MEASURING UPTO 110MM



3406-04

SCREWDRIVER, HEX 3.5MM



3406-04S

HOLDING SLEEVE FOR SCREWDRIVER, HEX 3.5MM, FOR LOCKING SCREW





3445-3.5

T-HANDLE SCREWDRIVER, HEX 3.5MM, SELF-RETAINING FOR LARGE FRAGMENT




1472-056

SCREWDRIVER HOLDING SLEEVE FOR CORTICAL AND CANCELLOUS SCREWS



3408-04L

HEXAGONAL SCREWDRIVER SHAFT WITH QUICK COUPLING - 3.5MM TIP, LONG



7-051-04

TREPHINE FOR LARGE FRAGMENT



3402-000

T-HANDLE WITH QUICK COUPLING



1472-054

QUICK COUPLING SHAFT






TQ-3.5

TORQUE LIMITING ATTACHMENT, 4.0NM (FOR LARGE FRAGMENT)



1472-064

TORQUE SCREWDRIVER HANDLE (FOR LARGE FRAGMENT)



3409-02L

BENDING IRON, LEFT FOR LARGE FRAGMENT




3409-02R

BENDING IRON, RIGHT FOR LARGE FRAGMENT



1472-032

ARTICULATED TENSION DEVICE



1472-034

ARTICULATED TENSION DEVICE KEY





2150-1008

PERIOSTEAL ELEVATOR, ROUND, 8MM



2150-1016

PERIOSTEAL ELEVATOR, FLAT, 16MM



2107-1240

REDUCTION FORCEP, POINTED, RATCHET LOCK, 240MM, FOR LARGE FRAGMENT




2106-1240

REDUCTION FORCEP, SERRATED, SPEED LOCK, 240MM, FOR LARGE FRAGMENT



2105-000

SELF-CENTERING BONE HOLDING FORCEP, SPEED LOCK



1472-070

HOHMANN RETRACTOR, 44MM, FOR LARGE FRAGMENT






2146-016

HOHMANN RETRACTOR, 16MM, FOR LARGE FRAGMENT



7-050-05

T- HANDLE FOR THREADED DRILL SLEEVE




7-051-03

SCREW HOLDING FORCEP FOR WISE-LOCK LARGE FRAGMENT




TI-007NW

SCREW CADDY FOR 5.0MM WISE-LOCK LARGE FRAGMENT SYSTEM



7-051-02

INSTRUMENT TRAYS FOR 5.0MM WISE-LOCK LARGE FRAGMENT INSTRUMENT SET






SC-004

CONTAINER FOR 5.0MM WISE-LOCK LARGE FRAGMENT INSTRUMENT SET





2303-000	Large Fragment Wise-Lock Instrument Set
----------	---





2303-000

Large Fragment Wise-Lock Instrument Set

CODE	SET CONSISTING OF	UNIT
2100-3.2-145	Drill Bit with Quick Coupling End, Ø3.2mm x Length 145mm, for Large Fragment	1
2100-4.5-145	Drill Bit with Quick Coupling End, Ø4.5mm x Length 145mm, for Large Fragment	1
2103-07	Bone Tap Quick Coupling for Ø4.5mm Cortical Screws, for Large Fragment	1
2103-09	Bone Tap Quick Coupling for Ø6.5mm Cancellous Screws, for Large Fragment	1
2103-4.3-225	Drill Bit Quick Coupling with Stopper Ø4.3mm x Length 225mm for Large Fragment	2
2106-2.0	Guide Sleeve for Ø2.0mm K. Wires - Large Fragment	1
3400-02	T-Handle Countersink for Ø4.5/6.5mm Screws, for Large Fragment	1
3420-02	Insert Drill Sleeve, Ø4.5/3.2mm, for Large Fragment	1
3443-38-5.0	Extraction Screw, Ø5.0mm for Large Fragment	1
3408-04	Hexagonal Screwdriver Shaft - 3.5mm Tip, Quick Coupling, Short	1
2186-3.5	HSS Drill Bit, Ø3.5mm, for Large Fragment	1
3767-4.5	Double Drill Guide, Ø4.5/3.2mm	1
1472-058	Self-Centering Double Drill Guide, Ø4.5/3.2mm, for Large Fragment	1
3767-6.5	Double Drill Guide, Ø6.5/3.2mm	1
BT-LF-S	Bending Template, Small, for Large Fragment-Locking	1
BT-LF-M	Bending Template, Medium, for Wise-Lock Large Fragment	1
BT-LF-L	Bending Template, Large, for Wise-Lock Large Fragment	1
7-051-05	Threaded Drill Sleeve for Drill Bit Ø4.3 for Wise-lock Large Fragment for T-Handle	3
GW-2.0-230	Guide Wire with Threaded Tip, Ø2.0mm x Thread Length 10mm x Length 230mm	3
3766-003	Allen Key, Hex 3.0mm for Large Fragment	1
1472-052	Pulling Drill Bit, Ø4.0mm, for Large Fragment	1
3443-400	Depth Gauge measuring upto 110mm	1
3406-04	Screwdriver, Hex 3.5mm	1
3406-04S	Holding Sleeve for Screwdriver, Hex 3.5mm, for Locking Screw	1
3445-3.5	T-Handle Screwdriver, Hex 3.5mm, Self-Retaining for Large Fragment	1
1472-056	Screwdriver Holding Sleeve for Cortical and Cancellous Screws	1
3408-04L	Hexagonal Screwdriver Shaft with Quick Coupling - 3.5mm Tip, Long	1
7-051-04	Trephine for Large Fragment	1
3402-000	T-Handle with Quick Coupling	1
1472-054	Quick Coupling Shaft	1
TQ-3.5	Torque Limiting Attachment, 4.0Nm, (for Large Fragment)	1
1472-064	Torque Screwdriver Handle , (for Large Fragment)	1
3409-02L	Bending Iron, Left for Large Fragment	1
3409-02R	Bending Iron, Right for Large Fragment	1
1472-032	Articulated Tension Device	1
1472-034	Articulated Tension Device Key	1
2150-1008	Periosteal Elevator, Round, 8mm	1
2150-1016	Periosteal Elevator, Flat, 16mm	1
2107-1240	Reduction Forcep, Pointed, Ratchet Lock, 240mm, for Large Fragment	1
2106-1240	Reduction Forcep, Serrated, Speed Lock, 240mm, for Large Fragment	1
2105-000	Self-Centering Bone Holding Forcep, Speed Lock	2
1472-070	Hohmann Retractor, 44mm, for Large Fragment	1
2146-016	Hohmann Retractor, 16mm, for Large Fragment	1
7-050-05	T- Handle for Threaded Drill Sleeve	1
7-051-03	Screw Holding Forcep for Wise-Lock Large Fragment	1
TI-007NW	Screw Caddy for 5.0mm Wise-Lock Large Fragment System	1
7-051-02	Instrument Trays for 5.0mm Wise-Lock Large Fragment Instrument Set	3
SC-004	Container for 5.0mm Wise-Lock Large Fragment Instrument Set	1

