


AUXILOCK® AUXPIK SUTURE- BASED ANCHOR

SHOULDER INSTABILITY REPAIR



AUXILOCK® AuxPik Suture-Based Anchor is made from UHMWPE suture anchor body with a PEEK OPTIMA eyelet tip. The AuxPik Suture-Based Anchor is designed for efficiency and ease of use. AuxPik Suture-Based Anchor is available in 1.8mm diameter with a BioBraid and suture tape. The anchor can be used for Bankart or SLAP repair surgeries. The anchor provides small foot print and subcortical fixation for anchor insertions. The drill guide and anchor driver combination are well designed to improve performance and reliability.

ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION	SUTURE/TAPE COLOUR
6-001-01	AUXILOCK® 1.8mm AuxPik Suture-Based Anchor With One #2 BioBraid: White/Blue	
6-001-08	AUXILOCK® 1.8mm AuxPik Suture-Based Anchor With One 1.4mm Suture Tape: White/Blue	

RECOMMENDED INSTRUMENTATION

CODE	PRODUCT DESCRIPTION	COMPATIBLE IMPLANTS
7-013-03	V-Mouth Drill Guide	6-001-01, 6-001-08
7-013-04	8 Point Drill Guide	6-001-01, 6-001-08
7-013-06	Trocar Tip Obturator	6-001-01, 6-001-08
7-013-53	Ø1.8mm X 325mm Drill Bit	6-001-01, 6-001-08