




AUXILOCK® BIOBRAID INFINITY LOOP NEEDLE



AUXILOCK® BioBraid Infinity Loop is a continuous loop of BioBraid with a straight needle. BioBraid Infinity Loop utilises the whipstitch technique on graft preparation and reduces the time spent on graft preparation. It also compresses the graft uniformly and improves strength.

ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION	SUTURE COLOUR	NEEDLE
6-012-01	AUXILOCK® #2 BioBraid Infinity Loop: White, With Needle: Straight, 20in Loop Length, 24in Total Length		
6-012-02	AUXILOCK® #2 BioBraid Infinity Loop: White/Blue, With Needle: Straight, 20in Loop Length, 24in Total Length		
6-012-03	AUXILOCK® #5 BioBraid Infinity Loop: White/Blue, With Needle: Straight, 20in Loop Length, 24in Total Length		

Note: Each box consist of 10 Pieces.