




AUXILOCK® BIOBRAID SUTURE TAPE



AUXILOCK® BioBraid Suture Tape made from UHMWPE The suture tape design accommodates the tape in the centre with a suture on the either sides of the tape.

It can be used in Knee, Shoulder, and Foot & Ankle surgeries including AC joint reconstruction, knotless double-row rotator cuff repair, patella fracture repair, MPFL Reconstruction etc. The broad footprint of the BioBraid Suture Tape is appropriate for repairs in degenerative cuff tissue where tissue pull-through may be a concern.

ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION	SUTURE COLOUR
6-013-01	AUXILOCK® 1.8mm BioBraid Suture Tape: White/Blue, 39in Total length	
6-013-02	AUXILOCK® 1.4mm BioBraid Suture Tape: White/Blue, 39in Total length	
6-013-03	AUXILOCK® 1.4mm BioBraid Suture Tape: White/Black, 39in Total length	
6-013-05	AUXILOCK® 2.0mm BioBraid Suture Tape: White/Blue, 39in Total length	
6-013-06	AUXILOCK® 1.0mm BioBraid Suture Tape: White/Blue, 39in Total length	
6-013-07	AUXILOCK® 1.4mm BioBraid Suture Tape: Solid Blue, 39in Total length	

Note: Each box consist of 10 Pieces.