


SHOULDER INSTABILITY REPAIR

AUXILOCK®
DUWLOC KNOTLESS
PEEK OPTIMA
PUSH-IN ANCHOR



DUWLOC delivers secure, knotless rotator cuff repair with precise tension control—minimizing surgical steps and maximizing glenoid preservation.

AUXILOCK® DUWLOC PEEK OPTIMA Knotless Push-in Anchor is engineered for efficient, secure bankart and SLAP repair without the complications of arthroscopic knots. It is available in the diameter of 2.8mm & 3.5mm Its innovative “suture-first” technique allows surgeons to pass sutures through the labrum repair and fine-tune tissue tension before anchor insertion. It comes with a PEEK OPTIMA Eyelet a two piece design allowing for tension adjustment before locking. The DUWLOC anchor reduces the risk of knot-related impingement, offering multiple stitch configurations and precise tension control for optimal outcomes.

ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION
6-064-01	AUXILOCK® Duwloc 2.8mm x 15.5mm Knotless PEEK OPTIMA Push-In Anchor
6-064-03	AUXILOCK® Duwloc 3.5mm x 20.0mm Knotless PEEK OPTIMA Push-In Anchor

RECOMMENDED INSTRUMENTATION

CODE	PRODUCT DESCRIPTION	COMPATIBLE IMPLANTS
7-013-84	Ø3.2mm X 200mm Drill Bit	6-064-03
7-013-01	8 Point Slotted Drill Guide, Length 190mm	6-064-01, 6-064-03
7-013-02	8 Point Drill Guide, Length 100mm	6-064-01, 6-064-03
7-013-83	Ø2.5mm X 200mm Drill Bit	6-064-01