





AUXILOCK® GFS BUTTON EXTENDER




EXTENDER BONE MODEL



HOLES TO ACCOMMODATE
PULLING & FLIPPING SUTURES

GFS MINI BUTTON &
GFS ULTIMATE MINI BUTTON



SLOTTED DESIGN

BUTTON
EXTENDER

GFS Button Extender is an extension device designed to provide extended surface area during cruciate graft reconstructions. It can be used in conjunction with the GFS Mini and GFS Ultimate Mini button providing an effective solution for cortical blowouts. It also helps overcome intraoperative complications seen in revision surgeries.

ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION	LENGTH	WIDTH	HEIGHT	SLOT DEPTH
6-032-01	AUXILOCK® GFS Button Extender, L: 18mm, W: 5.0mm	18mm	5mm	2.3mm	1.7mm