



AUXILOCK® GFS CONCAVE BUTTON



CONCAVE BUTTON
11MM WITH 4MM COLLAR

USE WITH 4 MM-7 MM
TUNNELS



CONCAVE BUTTON 14MM
WITH 7MM COLLAR

USE WITH 7 MM-9 MM
TUNNELS



CONCAVE BUTTON 20MM
WITH 9MM COLLAR

USE WITH 9 MM-13 MM
TUNNELS



ROUND BUTTON WITHOUT
COLLAR, 14MM

The GFS Concave Button is an attachable button system with the GFS No-Button adjustable loop that has revolutionized tibial fixation for ACL and PCL grafts. GFS Concave Button can be used on all graft types and attached to a variety of button configurations for fixation over sockets or tunnels.

THE ADVANTAGES OF THE GFS CONCAVE BUTTON IMPLANT INCLUDE:

- Maximum graft-to-bone contact improves incorporation and healing.
- The ability to retention grafts after fixation and knee cycling.
- Several different button options for sockets and tunnels.

ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION
6-050-01	AUXILOCK® GFS Concave Button 11mm with 4mm Collar
6-050-02	AUXILOCK® GFS Concave Button 14mm with 7mm Collar
6-050-03	AUXILOCK® GFS Concave Button 20mm with 9mm Collar
6-050-04	AUXILOCK® GFS Round Button without Collar, 14mm