




AUXILOCK® GFS NO-BUTTON ADJUSTABLE LOOP



GFS No-Button is an adjustable loop without a button, used for the tibial fixation in all-inside ACL & RCL reconstructions. The adjustable Loop is made of UHMWPE (Ultra High Molecular Weight Polyethylene).

Adjustable Suture (#7, White/Blue): It allows the surgeon to maximize the amount of graft inside the tibial tunnel, thereby optimizing the healing process. It also enables the calibration of the loop to its optimum size.

ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION	LENGTH OF LOOP
6-049-01	AUXILOCK® GFS No-Button Adjustable Loop	Package loop length: 70 mm Min: 12 mm