




AUXILOCK® MENISCO INSIDE-OUT



Menisco Inside-out system is specifically designed for inside-out meniscal repair procedures features cannulas, needles and rasps.

Five pre-bent Single Lumen cannulas provide optimum access to all zones of the meniscus with funnel ends that make the loading of needles safe and easy.

Double-arm meniscal repair needles are thin, strong and flexible for easy passage. The needles are available with Pre-attached #2-0, #0 BioBraid Suture & 1.0mm Suture Tape. Meniscal repair needles are made of Stainless Steel and have a large eyelet at one end to permit easy threading of suture material.

ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION
6-047-01	AUXILOCK® MENISCO Inside-Out Meniscal Repair Needle with #2-0 BioBraid: White/Blue, Length 36in, Needle Length 25cm
6-047-02	AUXILOCK® MENISCO Inside-Out Meniscal Repair Needle with #0 BioBraid: White/Blue, Length 36in, Needle Length 25cm
6-047-03	AUXILOCK® MENISCO Inside-Out Meniscal Repair Needle with 1.0mm Suture Tape: White/Blue, Length 36in, Needle Length 25cm

RELATED INSTRUMENTATION

CODE	PRODUCT DESCRIPTION
7-156-01	Single Lumen Cannula for Meniscal Repair, Straight
7-156-02	Single Lumen Cannula for Meniscal Repair, 10° Curve
7-156-03	Single Lumen Cannula for Meniscal Repair, 20° Curve
7-156-04	Single Lumen Cannula for Meniscal Repair, 30° Curve
7-156-05	Single Lumen Cannula for Meniscal Repair, 40° Curve
7-156-06	Bender for Cannula
7-156-07	Needle Pusher
7-156-08	Slotted Cannula
7-156-09	Knot Pusher for Meniscal Repair
7-156-10	Curved Probe for Meniscal Repair
7-156-11	Diamond Rasp, 45° for Meniscal Repair
7-156-12	Diamond Rasp, 90° for Meniscal Repair
7-156-13	Ring Curette for Meniscal Repair
7-156-14	Knot Pusher with Suture Cutter for Meniscal Repair
7-156-15	Container for Meniscal Repair Instrument Set
7-156-16	Lid for Container, Meniscal Repair Instrument Set