



AUXILOCK® NAKED BUTTON



GFS MINI NAKED BUTTON



ULTIMATE MINI NAKED BUTTON



GFS II LARGE NAKED BUTTON



ULTIMATE LARGE NAKED BUTTON

GFS Naked button is made up of titanium alloy. The system provides an excellent combination of strength and stiffness required for successful cortical fixation.

ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION
6-007-11	AUXILOCK GFS Mini Naked Button L: 12mm, W: 3.9mm
6-008-09	AUXILOCK GFS II Large Naked Button L: 16.5mm, W: 4.4mm
6-006-03	AUXILOCK Ultimate Mini Naked Button L: 12mm, W: 3.9mm
6-006-04	AUXILOCK Ultimate Large Naked Button L: 16.5mm, W: 3.9mm