



AUXILOCK® PEEK CF INTERFERENCE SCREW




The PEEK CF Interference Screw is made of a carbon fiber reinforced PEEK OPTIMA from Invibio, USA. The polymer in PEEK CF enhances the physical strength by nearly twice as much as natural PEEK.

Its mechanical properties are much closer to the cortical bone than the natural PEEK, PLLA or Titanium. The PEEK CF Interference Screw has a fully threaded design that provides strong mechanical fixation for both bone tendon bone (BTB) and soft tissue grafts. It is renowned for its radiolucent properties, less imaging artifact and for being MRI safe. PEEK CF interference screw offers the revision ability of an absorbable screw.

ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION
6-009-01	AUXILOCK® 7mm x 20mm PEEK CF Interference Screw
6-009-02	AUXILOCK® 7mm x 25mm PEEK CF Interference Screw
6-009-03	AUXILOCK® 7mm x 30mm PEEK CF Interference Screw
6-009-04	AUXILOCK® 8mm x 20mm PEEK CF Interference Screw
6-009-05	AUXILOCK® 8mm x 25mm PEEK CF Interference Screw
6-009-06	AUXILOCK® 8mm x 30mm PEEK CF Interference Screw
6-009-07	AUXILOCK® 8mm x 35mm PEEK CF Interference Screw
6-009-08	AUXILOCK® 9mm x 20mm PEEK CF Interference Screw
6-009-09	AUXILOCK® 9mm x 25mm PEEK CF Interference Screw
6-009-10	AUXILOCK® 9mm x 30mm PEEK CF Interference Screw
6-009-11	AUXILOCK® 9mm x 35mm PEEK CF Interference Screw
6-009-12	AUXILOCK® 10mm x 25mm PEEK CF Interference Screw
6-009-13	AUXILOCK® 10mm x 30mm PEEK CF Interference Screw
6-009-14	AUXILOCK® 10mm x 35mm PEEK CF Interference Screw
6-009-15	AUXILOCK® 11mm x 25mm PEEK CF Interference Screw
6-009-16	AUXILOCK® 11mm x 30mm PEEK CF Interference Screw
6-009-17	AUXILOCK® 11mm x 35mm PEEK CF Interference Screw
6-009-18	AUXILOCK® 12mm x 25mm PEEK CF Interference Screw



ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION
6-009-19	AUXILOCK® 12mm x 30mm PEEK CF Interference Screw
6-009-20	AUXILOCK® 12mm x 35mm PEEK CF Interference Screw

RECOMMENDED INSTRUMENTATION

CODE	PRODUCT DESCRIPTION	COMPATIBLE IMPLANTS
7-014-43	Nitinol Guide Wire Ø1.2mm x Length 300mm (For Interference Screws)	6-009-01, 6-009-02, 6-009-03, 6-009-04, 6-009-05, 6-009-06, 6-009-07, 6-009-08,
7-014-04	Screwdriver For Polymer Interference Screw	6-009-09, 6-009-10, 6-009-11, 6-009-12,
7-014-96	Screwdriver Shaft For Polymer Interference Screw, Hudson Coupling	6-009-13, 6-009-14, 6-009-15, 6-009-16,
7-014-42	Ratchet Handle For Hudson Coupling Screwdriver Shaft	6-009-17, 6-009-18, 6-009-19, 6-009-20