




AUXILOCK® PEEK OPTIMA INTERFERENCE SCREW




The PEEK OPTIMA (From InVivo, USA) interference screw is made up of PEEK (Poly-Ether- Ether-Ketone). PEEK is a thermoplastic material and has radiolucent properties for easy artifact free monitoring and assessment of the healing site with X-ray, CT or MRI. It also offers the revision ability of an absorbable screw.

PEEK OPTIMA Interference Screw has a fully threaded design. It provides strong mechanical fixation for both Bone-Tendon-Bone (BTB) and soft tissue grafts. PEEK OPTIMA is considered to be extremely strong, durable, and highly resistant to creep & fatigue. The bonelike modulus help minimize stress shielding and also stimulates bone healing.

ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION
6-014-01	AUXILOCK® 7mm x 20mm PEEK OPTIMA Interference Screw
6-014-02	AUXILOCK® 7mm x 25mm PEEK OPTIMA Interference Screw
6-014-03	AUXILOCK® 7mm x 30mm PEEK OPTIMA Interference Screw
6-014-04	AUXILOCK® 8mm x 20mm PEEK OPTIMA Interference Screw
6-014-05	AUXILOCK® 8mm x 25mm PEEK OPTIMA Interference Screw
6-014-06	AUXILOCK® 8mm x 30mm PEEK OPTIMA Interference Screw
6-014-07	AUXILOCK® 8mm x 35mm PEEK OPTIMA Interference Screw
6-014-08	AUXILOCK® 9mm x 20mm PEEK OPTIMA Interference Screw
6-014-09	AUXILOCK® 9mm x 25mm PEEK OPTIMA Interference Screw
6-014-10	AUXILOCK® 9mm x 30mm PEEK OPTIMA Interference Screw
6-014-11	AUXILOCK® 9mm x 35mm PEEK OPTIMA Interference Screw
6-014-12	AUXILOCK® 10mm x 25mm PEEK OPTIMA Interference Screw
6-014-13	AUXILOCK® 10mm x 30mm PEEK OPTIMA Interference Screw
6-014-14	AUXILOCK® 10mm x 35mm PEEK OPTIMA Interference Screw
6-014-15	AUXILOCK® 11mm x 25mm PEEK OPTIMA Interference Screw
6-014-16	AUXILOCK® 11mm x 30mm PEEK OPTIMA Interference Screw
6-014-17	AUXILOCK® 11mm x 35mm PEEK OPTIMA Interference Screw
6-014-18	AUXILOCK® 12mm x 25mm PEEK OPTIMA Interference Screw



ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION
6-014-19	AUXILOCK® 12mm x 30mm PEEK OPTIMA Interference Screw
6-014-20	AUXILOCK® 12mm x 35mm PEEK OPTIMA Interference Screw

RECOMMENDED INSTRUMENTATION

CODE	PRODUCT DESCRIPTION	COMPATIBLE IMPLANTS
7-014-43	Nitinol Guide Wire Ø1.2mm x Length 300mm (For Interference Screws)	6-014-01, 6-014-02, 6-014-03, 6-014-04, 6-014-05, 6-014-06, 6-014-07, 6-014-08,
7-014-04	Screwdriver For Polymer Interference Screw	6-014-09, 6-014-10, 6-014-11, 6-014-12,
7-014-96	Screwdriver Shaft For Polymer Interference Screw, Hudson Coupling	6-014-13, 6-014-14, 6-014-15, 6-014-16,
7-014-42	Ratchet Handle For Hudson Coupling Screwdriver Shaft	6-014-17, 6-014-18, 6-014-19, 6-014-20