




AUXILOCK® PEEK OPTIMA SCREW-IN SUTURE ANCHOR

ROTATOR CUFF REPAIR




AUXILOCK® PEEK OPTIMA Screw-In Suture anchor is a fully threaded suture anchor featuring dual threads to maximize cortical and cancellous fixation. It is available in diameter of 4.5, 5.5 and 6.5mm with double and triple loaded options.

It has a flat tip to facilitate the insertion. The anchor is particularly suitable for repairing rotator cuff and associated pathologies. The anchors are also available with needles which are ideal for mini-open rotator cuff repair procedures.

ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION	SUTURE/TAPE COLOUR
6-016-05	AUXILOCK® 4.5mm PEEK OPTIMA Screw-In Suture Anchor With Two #2 BioBraid: White/Blue & White/Black	
6-016-07	AUXILOCK® 5.5mm PEEK OPTIMA Screw-In Suture Anchor With Two #2 BioBraid: White/Blue & White/Black	
6-016-08	AUXILOCK® 6.5mm PEEK OPTIMA Screw-In Suture Anchor With Two #2 BioBraid: White/Blue & White/Black	

6-016-13	AUXILOCK® 4.5mm PEEK OPTIMA Screw-In Suture Anchor With Three #2 BioBraid: White, White/Blue & White/Black	
6-016-14	AUXILOCK® 5.5mm PEEK OPTIMA Screw-In Suture Anchor With Three #2 BioBraid: White, White/Blue & White/Black	
6-016-15	AUXILOCK® 4.5mm PEEK OPTIMA Screw-In Suture Anchor With One #2 BioBraid and One 1.4mm Suture Tape: White/Blue & White/Black	
6-016-16	AUXILOCK® 5.5mm PEEK OPTIMA Screw-In Suture Anchor With One #2 BioBraid and One 1.4mm Suture Tape: White/Blue & White/Black	
6-016-17	AUXILOCK® 6.5mm PEEK OPTIMA Screw-In Suture Anchor With One #2 BioBraid and One 1.4mm Suture Tape: White/Blue & White/Black	
6-016-18	AUXILOCK® 4.5mm PEEK OPTIMA Screw-In Suture Anchor With Two 1.4mm Suture Tape: White/Blue & White/Black	
6-016-19	AUXILOCK® 5.5mm PEEK OPTIMA Screw-In Suture Anchor With Two 1.4mm Suture Tape: White/Blue & White/Black	
6-016-20	AUXILOCK® 6.5mm PEEK OPTIMA Screw-In Suture Anchor With Two 1.4mm Suture Tape: White/Blue & White/Black	



ROTATOR CUFF REPAIR

ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION	SUTURE/TAPE COLOUR	NEEDLE
6-016-01	AUXILOCK® 4.5mm PEEK OP-TIMA Screw-In Suture Anchor With Two #2 BioBraid: White & White/Blue, With Needles: MO-6	 	
6-016-03	AUXILOCK® 5.5mm PEEK OP-TIMA Screw-In Suture Anchor With Two #2 BioBraid: White & White/Blue, With Needles: MO-6	 	
6-016-04	AUXILOCK® 6.5mm PEEK OP-TIMA Screw-In Suture Anchor With Two #2 BioBraid: White & White/Blue, With Needles: MO-6	 	

RECOMMENDED INSTRUMENTATION

CODE	PRODUCT DESCRIPTION	BONE DENSITY	COMPATIBLE IMPLANTS
7-013-25	Awl For 4.5/5.5mm PEEK Screw-In Anchor	Soft Bone	6-016-05, 6-016-07, 6-016-13, 6-016-14, 6-016-15, 6-016-16, 6-016-18, 6-016-19, 6-016-01, 6-016-03
7-013-26	Awl For 6.5mm PEEK Screw-In Anchor		6-016-08, 6-016-17, 6-016-20, 6-016-04
7-013-27	Awl/ Tap For 4.5/5.5mm PEEK Screw-In Anchor	Normal Bone	6-016-05, 6-016-07, 6-016-13, 6-016-14, 6-016-15, 6-016-16, 6-016-18, 6-016-19, 6-016-01, 6-016-03
7-013-28	Awl/ Tap For 6.5mm PEEK Screw-In Anchor		6-016-08, 6-016-17, 6-016-20, 6-016-04
7-013-29	Tap For 4.5/5.5mm PEEK Screw-In Anchor	Hard Bone	6-016-05, 6-016-07, 6-016-13, 6-016-14, 6-016-15, 6-016-16, 6-016-18, 6-016-19, 6-016-01, 6-016-03
7-013-30	Tap For 6.5mm PEEK Screw-In Anchor		6-016-08, 6-016-17, 6-016-20, 6-016-04